

# GC Analysis of Derivatized Chlorinated Acetic Acids

Anila I Khan, Thermo Fisher Scientific, Runcorn, Cheshire, UK

## Key Words

Derivatization reagent, extractive alkylation, pentafluorobenzyl bromide (PFBBR), chloroacetic acid (CA), dichloroacetic acid (DCA), trichloroacetic acid (TCA)

## Abstract

Chloroacetic acid (CA), Dichloroacetic acid (DCA) and Trichloroacetic acid (TCA) are difficult to analyze by GC due to highly polar acidic groups and therefore must be derivatised to increase their detectability. This application note demonstrates derivatization is achieved using the alkylation reagent pentafluorobenzyl bromide (PFBBR) to form their corresponding fluorinated derivatives. Separation of the derivatized chlorinated acetic acids with highly symmetrical peak shapes was achieved using a 5% phenylmethylpolysiloxane phase column.

## Introduction

Chloroacetic acid (CA), Dichloroacetic acid (DCA) and Trichloroacetic acid (TCA) play important roles in chemical and pharmaceutical industries. CA is actively used in industrial processes such as in the production of dyes, insecticides and active pharmaceutical ingredients<sup>1</sup>. DCA has been used in cancer treatment<sup>2</sup> and TCA has been used in herbicides<sup>3</sup>.

Trace levels of low molecular weight chlorinated acetic acids TCA, DCA and CA can be difficult to analyze by GC due to their highly polar acidic nature. Derivatization of chlorinated acetic acids improves the volatility and peak shapes on a chromatographic column.

Pentafluorobenzyl bromide (PFBBR) is commonly used for converting carboxylic acids, as well as mercaptans, phenols and sulphonamides to their respective fluorinated derivatives. In this case chlorinated acetic acid derivatives were obtained via an extractive alkylation procedure using the alkylation reagent PFBBR (Figures 1 and 2). PFBBR is manufactured to meet the exact needs for sensitive derivatization reactions. This procedure simultaneously combines extraction and derivatisation without the need for purification. In addition, the reaction time is fast and can take less than 30 min to complete.



## Experimental Details

### Sample Preparation

10 mg each of TCA, DA and CA were weighed into a 50 mL volumetric flask and made up to volume with Dichloromethane. A 1 mL aliquot was taken and placed into a Thermo Scientific Reacti-Vial containing a Reacti-Vial magnetic stirrer. To the Reacti-Vial, 1 mL of aqueous 0.1 M tetrabutylammonium hydrogen sulphate and 1 mL of 0.2 M NaOH were added. This was followed by the addition of 20  $\mu$ L of PFBBR. The vials were capped and placed in the Reacti-Therm Sample Incubation System and stirred for 30 minutes at 60 °C. The final sample was then transferred to a 2 mL autosampler vial and 1  $\mu$ L was injected into the GC/FID.

<b>Reagents</b>	<b>Part Number</b>
Thermo Scientific alkylation reagent PFBBR 5g	TS-58220
Fisher Scientific tetrabutylammonium hydrogen sulphate	16838-0250
Fisher Scientific sodium hydroxide	BPE359-500

<b>Sample Handling Equipment</b>	<b>Part Number</b>
Thermo Scientific Reacti-Therm III Heating/Stirring Module	TS-18823
Thermo Scientific Reacti-Vap III Evaporator	TS-18826
Thermo Scientific Reacti-Block Q-1 (Holds 8 × 10 mL Reacti-Vials)	TS-18814
Thermo Scientific Reacti-Vial clear glass reaction vials 10 mL	TS-13225
Thermo Scientific 2 mL amber vial and screw tops	60180-565

<b>Separation Conditions</b>	<b>Part Number</b>	
Instrumentation:	Thermo Scientific TRACE GC Ultra	
Column:	TRACE TR-5 30 m × 0.25 mm × 0.25 µm	260E142P
Thermo Scientific BTO 17 mm septa		31303211
5 mm ID Focus Split Liner, 105 mm long		453T1905
Graphite liner seal		29033406
10 µL, 50 mm needle length gauge 25 Syringe		36500525
Graphite ferrules to fit 0.1-0.25 mm ID columns		29053488
Carrier gas:	Helium	
Split flow:	60 mL/min	
Column flow:	1.2 mL/min, constant flow	
Split ratio:	50:1	
Oven temperature:	40 °C (1 min), 10 °C/min, 300 °C (5 min)	
Injector type:	Split/Splitless	
Injector mode:	Split	
Injector temperature:	240 °C	
Detector type:	FID	
Detector temperature:	280 °C	
Detector Hydrogen flow:	35 mL/min	
Detector Air flow:	350 mL/min	
Detector Nitrogen flow:	30 mL/min	
Thermo Scientific TriPlus Autosampler		
Injection Volume:	1 µL	

<b>Data Processing</b>	
Software:	Thermo Scientific XCalibur™




Figure 1: Reaction products of pentafluorobenzyl derivative of TCA<sup>4</sup>




Figure 2 Pentafluoro derivative of (a) DCA and (b) CA

## Results

Two peaks were observed for derivatised TCA-PFB; the product structures are shown in Figure 1. The reaction itself does not comprise of two products but one product 1,1,1-trichloro-2-pentafluorophenylethane decomposing to 1,1-dichloro-2-pentafluorophenylethane<sup>4</sup>. Each of the derivatized chlorinated peaks was identified by individually derivatizing the TCA, DCA and CA. Figure 1 shows the derivatization of TCA and Figure 2 shows the structures of the derivatized form of DCA and CA. Figure 3 shows the chromatograms of the corresponding esters of the derivatized TCA, DCA and CA.

## Conclusion

PFBBR is an ideal derivatization reagent for increasing the detectability of low molecular weight chlorinated acetic acids. This enabled enhanced separation and detection with a TRACE TR-5 GC column.




Figure 3: Chromatogram of PFB-derivatized acetic acids

Peak Number	Compounds	t <sub>R</sub> (min)
1	Derivatized TCA	10.2
2	Derivatized TCA	10.9
3	PFB-by Product	11.3
4	Derivatized CA	13.3
5	Derivatized DCA	16.7

## References

- Günter Koenig, Elmar Lohmar, Norbert Rupprich "Chloroacetic Acids" in Ullmann's Encyclopedia of Industrial Chemistry 2002, Wiley-VCH, Weinheim.
- Matthew G. Vander Heiden, Targeting Cell Metabolism in Cancer Patients, *Sci Transl Med* 12 May 2010 2:31ed1
- Monaco T.J, Weller S.C, Ashton F. M, Weed science: principles and practices. 2002, Technology and Engineering.
- Sinkkonen S, Kolehmainen E, Paasivirta J, Hamalainen S, Lahtipera M, *Journal of chromatography A*, 1995, 718, 391-396.
- ACD labs software to draw chemical structures.
- Thermo Scientific reagents, solvents and accessories brochure (Ref: BR20535\_E 06/12S). Available upon request.

## [thermoscientific.com/chromatography](http://thermoscientific.com/chromatography)

©2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**USA and Canada** +1 800 332 3331  
**Australia** 1300 735 292 (free call domestic)  
**China** 800 810 5118 (free call domestic)  
 400 650 5118  
**France** +33 (0)1 60 92 48 34  
**Germany** +49 (0) 2423 9431 20 or 21  
**India** +91 22 6742 9494  
 +91 27 1766 2352

**Japan** 0120 753 670 (free call domestic)  
 0120 753 671 (fax)  
**Korea** +82 2 3420 8600  
**United Kingdom** +44 (0) 1928 534 110  
**New Zealand** 0800 933 966 (free call domestic)  
**Singapore** +65 6289 1190  
**All Other Enquiries** +44 (0) 1928 534 050

### Technical Support

For advice and support, please visit our website:  
[www.thermoscientific.com/chromexpert](http://www.thermoscientific.com/chromexpert)

**Thermo**  
 SCIENTIFIC

A Thermo Fisher Scientific Brand