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

Positively charged metal ions can attract negatively charged biomolecules through **chelation**

Iron can also cause **oxidation** of some proteins

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## Corrosion Resistance

MP35N alloy possesses excellent resistance to sulfidation, high temperature oxidation, hydrogen embrittlement, saline solutions and most mineral acids.

Nitric Acid	Good	Sulfuric Acid	Good
Phosphoric Acid	Good	Acetic Acid	Excellent
Sodium Hydroxide	Good	Salt Spray (NaCl)	Excellent
Sea Water	Excellent	Sour Oil/Gas	Excellent
Humidity	Excellent		

## Type Analysis

Carbon (Maximum)	0.03 %	Manganese (Maximum)	0.15 %
Phosphorus (Maximum)	0.015 %	Sulfur (Maximum)	0.010 %
Silicon (Maximum)	0.15 %	Chromium	19.00 to 21.00 %
Nickel	33.00 to 37.00 %	Molybdenum	9.00 to 10.50 %
Cobalt	Balance	Titanium (Maximum)	1.00 %
Boron	0.010 %	Iron (Maximum)	1.00 %

## Galvanic Series of Metals and Alloys Ladder

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Technical Overview

## Highly Confident Peptide Mapping of Protein Digests Using Agilent LC/Q-TOFs

### Authors

Stephen Madden,  
Crystal Cody, and  
Jungkap Park  
Agilent Technologies, Inc.  
Santa Clara, California, USA

### Introduction

The biopharmaceutical industry is adding an increasing number of protein therapeutics to its pipelines. With a five-year compound annual growth rate (CAGR) of 7.3 %, it is predicted to attain sales of \$248.7 billion worldwide by 2020<sup>1</sup>. The fastest growing segment with a five-year CAGR of 16.6 % is monoclonal antibodies (mAbs), which have been approved for a range of conditions such as oncology and autoimmune diseases. During the development of protein therapeutics, monitoring is required for a broad set of Critical Quality Attributes (CQAs), among which is confirming the protein sequence and post-translational modifications (PTMs).

Typically, a bottoms-up approach is taken in which the protein is digested with an enzyme to create peptides in the range 5 to 60 amino acids. The peptide mixture is separated by liquid chromatography then detected in a mass spectrometer (MS), which generates mass spectra containing information about them. Given the richness of MS data, there are challenges in finding spectra that have peptide information, then extracting the sequence and PTMs accurately and efficiently. A manual approach by inspection requires significant time and expertise.

The Agilent MassHunter BioConfirm B.09.00 has made a step forward in this peptide mapping workflow for finding peptides and visualizing the results. A software algorithm uses MS/MS data to narrow the list of possible peptides, saving time versus manual inspection. The results are displayed in the Sequence Coverage Map, which allows the user to review multiple samples (for example, from different digest enzymes) at the same time. A new mode in Agilent MassHunter Acquisition software on the Agilent 6545XT AdvanceBio LC/Q-TOF called Iterative MS/MS focuses on low-concentration PTMs, and enables them to be viewed in the Sequence Coverage Map.

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
File View Find and Identify Method Sequence Configuration Help

Sample Table: NIST mAb Digest\_250 ng-uL\_01.d

Results Workflow

Confirmation Status	File Name	Saved Results Method	Last Run Method	Workflow	Condition	Sequence / Mass	Modification	Enzyme
Partially confirmed	NIST mAb Digest_250 ng-uL_01.d	BioConfirmProteinDigest-NISTmAb 2.m	Protein Digest	reduced	NISTmAb	mAb     Protein Digest (Reduced + Alkylated)	Trypsin,Lys	

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Biomolecule Fragment Spectrum

**Anotované MS/MS spektrum**

Sequence Coverage Map: NISTmAb (Protein Digest) (64.50%)

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