

GLYCOWORKS *Rapi*Fluor-MS QUICK START PROTOCOL

STREAMLINED PROTOCOL – 24 SAMPLE (3 x 8 FORMAT)

- Set heat blocks to at least 90 °C and 50 °C.
- Protocol is based on 1.5 mg/mL glycoprotein starting concentration.



STEP 1: Rapid Deglycosylation

1. Reconstitute 1 vial of the Intact mAb Mass Check Standard (1 mg/vial) in 670 μ L 18.2 M Ω water to create a 1.5 mg/mL solution.

Note: For glycoproteins with a formulation buffer containing nucleophiles or anionic reagents (e.g., His, Gly, Tris, PO₄²⁻), if a Glycan C₁₈ AX Column is applied for LC separation, it is highly recommended to desalt the sample with water prior to Step 1.

2. Prepare 3% (w/v) *Rapi*Gest™ SF by dissolving 3 mg of *Rapi*Gest SF Surfactant in 60 μ L of Rapid Buffer and 40 μ L of water, vortex.
3. Dilute PNGase F enzyme (35 μ L) with 220 μ L water for a total of 255 μ L.
4. Add 10 μ L of 1.5 mg/mL glycoprotein into the provided tube.
5. Add 10 μ L of buffered 3% (w/v) *Rapi*Gest SF solution to above tube, aspirate to mix.
6. Heat at least to 90 °C for 3 minutes.
7. Cool at room temperature for 3 minutes.
8. Add 10 μ L Rapid PNGase F and aspirate to mix.
9. Incubate at 50 °C for 5 minutes.
10. Cool at room temperature for 3 minutes.

STEP 2: Rapid Labeling of Glycosylamines

1. Add 110 μ L of anhydrous DMF directly to one vial of 9 mg of *Rapi*Fluor-MS™ Reagent. Mix to solubilize.
2. Add 10 μ L of the *Rapi*Fluor-MS solution to the deglycosylation mixture and aspirate to mix.
3. Allow the labeling to proceed at room temperature for 5 minutes.
4. Dilute the reaction with 360 μ L of acetonitrile (ACN) and aspirate to mix.

STEP 3: HILIC Cleanup of Labeled Glycosylamines

1. Set up a GlycoWorks™ HILIC μ Elution Plate and add in shims or spacer and waste tray.
2. Condition wells by adding 200 μ L of water per well.
3. Equilibrate wells by adding 200 μ L 85% ACN.
4. Load ACN-diluted samples (~400 μ L).
5. Wash wells with two (2) 600 μ L volumes of 1% formic acid, 90% ACN.
6. Replace waste tray with sample collection tray loaded with 600 μ L tubes.
7. Elute glycans with three (3) 30 μ L volumes of SPE Elution Buffer into 600 μ L tapered bottom inserts.
8. Dilute SPE eluate with 310 μ L of the GlycoWorks SPE Diluent (DMF/ACN). Aspirate to mix.
Note: For a Glycan C₁₈ AX separation sample, either skip dilution Step 8 or dilute with 310 μ L of water.
9. Cap the tubes with pre-slit cap mats.

► For the complete Care and Use Manual, visit [waters.com](https://www.waters.com) and search [715004903EN](#).

► For more details on this method, download Application Notes [720005506EN](#) and [720007038EN](#).