

TheReporter

Reprinted from Volume 15, No. 6, 1996

T296016

© 1999 Sigma-Aldrich Co.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA • Sigma-Aldrich de Argentina, S.A. • Av. Pueyrredon 2446/50 • Piso 5-B • Buenos Aires 1119
AUSTRALIA • Sigma-Aldrich Pty. Ltd. • Unit #2, 14 Anella Avenue • Castle Hill NSW 2154
AUSTRIA • Sigma-Aldrich Handels GmbH • Hebbelplatz 7 • A-1110 Wien
BELGIUM • Sigma-Aldrich N.V./S.A. • K. Cardijnplein 8 • B-2880 Bornem
BRAZIL • Sigma-Aldrich Quimica Brasil Ltda. • Rua Sabara, 566-Conj. 53 • 01239-010 São Paulo, SP
CANADA • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC • Sigma-Aldrich s.r.o. • Pobrezni 46 • 186 21 Praha 8
DENMARK • Sigma-Aldrich Denmark A/S • Vejlegaardsvej 65B • DK-2665 Vallengbaek Strand
FINLAND • Sigma-Aldrich Finland/YA-Kemia Oy • Teerisuonkuja 4 • FIN-00700 Helsinki
FRANCE • Sigma-Aldrich Chimie • Chromatographie Supelco • L'Isle d'Abeau Chesnes - B.P. 701 • 38297 Saint-Quentin Fallavier Cedex
GERMANY • Sigma-Aldrich Chemie GmbH • Geschäftsbereich Supelco • Grünwalder Weg 30 • D-82041 Deisenhofen
GREECE • Sigma-Aldrich (o.m.) Ltd. • 72 Argonafton Str. • 16346 Ilioupoli, Athens
HUNGARY • Sigma-Aldrich Kft. • Nagy Diófa u. 7., IV fl. • H-1067 Budapest
INDIA • Sigma-Aldrich Co. • Survey No. 31/1, Sitharamapalaya • Mahadevapura P.O. • Bangalore 560 048
IRELAND • Sigma-Aldrich Ireland Ltd. • Airton Road • Tallaght • Dublin 24
ISRAEL • Sigma Israel Chemicals Ltd. • Park Rabin • Rohovot 76100
ITALY • Sigma-Aldrich s.r.l. • Via Gallarate, 154 • 20151 Milano
JAPAN • Sigma-Aldrich Japan K.K. • Division Supelco • JL Nihonbashi Building • 1-10-15 Nihonbashi Horidome-cho, Chuo-ku • Tokyo 103
KOREA • Sigma-Aldrich Korea • Samhan Camus Annex, 10th Floor • 17-26 Yoido-dong, Yungdeungpo-ku • Seoul
MALAYSIA • Sigma-Aldrich (M) Sdn. Bhd. • 9-2, Jalan 2/128, Taman Gembira • Off Jalan Kuchai Lama • 58200 Kuala Lumpur • Selangor
MEXICO • Sigma-Aldrich Quimica S.A. de C.V. • Calle 6 North No. 107 • Parque Industrial Toluca 2000 • 50200 Toluca
NETHERLANDS • Sigma-Aldrich Chemie BV • Postbus 27 • 3330 AA Zwijndrecht
NORWAY • Sigma-Aldrich Norway • Sandakerveien 102 • N-0483 Oslo
POLAND • Sigma-Aldrich Sp. z o.o. • Szelagowska 30 • 61-626 Poznań
PORTUGAL • Sigma-Aldrich Quimica, S.A. • P.O. Box 131 • Sintra 2710
RUSSIA • Sigma-Aldrich Russia • TOO Techmedbiochem • Makarenko Str. 2/21 • Building 1, Flat 22 • Moscow 103062
SINGAPORE • Sigma-Aldrich Pte. Ltd. • 102E Pasir Panjang Road • #08-01 Citilink Warehouse • Singapore 118529
SOUTH AFRICA • Sigma-Aldrich (pty) Ltd. • CNR Kelly & Ackerman Streets • Southern Life Industrial Park Unit • Unit 16/17 • Jet Park 1459
SPAIN • Sigma-Aldrich Quimica, S.A. • Apt. Correos 161 • 28100 Alcobendas, Madrid
SWEDEN • Sigma-Aldrich Sweden AB • Solkraftsvägen 14C • 135 70 Stockholm
SWITZERLAND • Supelco Switzerland • Industriestrasse 25 • P.O. Box 260 • CH-9471 Buchs
UNITED KINGDOM • Sigma-Aldrich Company Ltd. • Supelco UK • Fancy Road, Poole • Dorset BH12 4QH
UNITED STATES • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email:supelco@sial.com

H

This article is archived from a past issue of The Supelco Reporter. Information in the article was appropriate at the time of publication, but product specifications, catalog numbers, and availability may have changed over time.

If you have questions about applying methodology described in this article to a current application, please contact our technical service chemists.

 **SUPELCO**

ISO 9001
REGISTERED

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

New Pharmaceutical and Food & Beverage Applications for SPME/HPLC

G. Gora-Maslak
Supelco, Bellefonte, PA USA

Solid phase microextraction (SPME) has been established as a practical alternative for sample preparation for gas chromatography. A new SPME/HPLC interface now allows SPME to be used in analyses of many weakly volatile or thermally labile compounds, such as pharmaceutical compounds. New SPME/HPLC applications include monitoring tricyclic antidepressant drugs in serum and antioxidants and preservatives in foods and beverages.

Because analytes are rapidly extracted from a sample matrix, with virtually no solvent consumption, solid phase microextraction (SPME)* has become established as a practical approach to sample preparation for gas chromatography. SPME saves preparation time and solvent purchase and disposal costs, and can improve the limits of detection. Now, our SPME/HPLC interface enables analysts to use SPME in analyses of numerous weakly volatile or thermally labile analytes (1). New applications for SPME/HPLC are summarized here.

Tricyclic Antidepressant Drugs in Serum

In analyses of drugs in blood, urine, or other biological fluids sample preparation usually consists of removing the analytes through liquid-liquid extraction, solid phase extraction, or other techniques. These methods can involve excessive preparation time, and require expensive organic solvents. SPME eliminates these drawbacks.

Investigators from the departments of legal medicine at Showa University School of Medicine (Tokyo) and Hamamatsu University School of Medicine (Hamamatsu) developed a headspace SPME/capillary GC method for extracting the most volatile tricyclic antidepressants (amitriptyline, chlorimipramine, imipramine, trimipramine) from urine (2). In developing their method, these analysts extracted four additional tricyclic antidepressants (carpipramine, clocapramine, desipramine, lofepramine) from urine samples, but these compounds decomposed in the GC and eluted as multiple peaks. Such heat-sensitive analytes should be analyzed by combining SPME extraction with HPLC analysis. Our SPME/HPLC interface makes this combination possible.

We have developed a procedure for SPME/HPLC analysis of three nonvolatile tricyclic antidepressants, desipramine, nordoxepin, and nortriptyline, in serum. We add methanol to the serum sample to denature the serum proteins and release protein-bound drugs, then extract the drugs by immersion SPME, using a specially prepared 40 μ m polydimethylsiloxane/divinylbenzene SPME fiber. Extraction and analytical conditions are summarized in Figure A. The data summarized in Table 1 show that the extraction is

Figure A. Tricyclic Antidepressant Drugs in Serum

Sample: 0.8mL thawed serum containing 100ng/mL each analyte, 0.16mL methanol (20% of serum volume) added to denature proteins/release protein-bound drugs
SPME Fiber: 40 μ m PDMS/DVB-HPLC
Cat. No.: 57317
Extraction: immersion, 30 min (rapid stirring)
Desorption: dynamic, methanol:acetonitrile, 60:40, 10 min
Column: SUPELCOSIL™ LC-PCN, 25cm x 4.6mm ID, 5 μ m particles (with Supelguard™ LC-PCN guard column)
Cat. No.: 58378
Mobile Phase: acetonitrile:methanol:0.01M K₂HPO₄ (pH 7.0 with phosphoric acid), 60:15:25
Flow Rate: 2mL/min
Temp.: 40°C
Det.: UV, 215nm




Table 1. Consistent Recovery of Tricyclic Antidepressants by SPME/HPLC

Analyte	Precision [^]		Calibration Curve ^{^^}	
	Mean \pm Std. Dev.	%RSD	Mean \pm Std. Dev.	%RSD
Nordoxepin	0.68 \pm 0.06	10.0	0.74 \pm 0.07	9.9
Nortriptyline	1.51 \pm 0.14	9.1	1.52 \pm 0.06	4.1
Desipramine	1.86 \pm 0.09	5.2	1.90 \pm 0.09	4.9

[^]100ng/mL each analyte in serum, n = 8 extractions (4 fibers, 2 extractions each).

^{^^}50–400ng/mL each analyte in serum.

Internal standard: protriptyline. Conditions listed in Figure A.

consistent and sensitive, and the analytes are efficiently delivered to the HPLC column.

In addition to ensuring speed and economy in extractions of pharmaceuticals from biological fluids, SPME can be more effective than liquid-liquid extractions in separating analytes of interest from potentially interfering matrix components or coadministered

drugs (2). Antioxidants and Preservatives in Foods and Beverages


Antioxidants and preservatives maintain the quality of food products during storage. These additives must be monitored to ensure that they are within specified limits. Currently, BHA (2- and 3- tert-butyl-4-hydroxyanisole) and BHT (3,5-di-tert-butyl-4-hydroxytoluene) in oils and fats are analyzed through a process that includes tedious liquid extraction/concentration (3). SPME/HPLC offers a considerable savings of time and solvents, while providing very satisfactory results, as shown by analyses of BHA and BHT extracted from olive oil and soft drink powder (Figure B).

The procedure for monitoring benzoic acid in foods includes analyte derivatization, as well as liquid extraction and concentration (4). SPME/HPLC eliminates all of these steps. Note in Figure C that caffeine, as well as benzoic acid, is extracted by the 40 μ m PDMS/DVB-HPLC fiber.

These examples show SPME now offers a faster, better approach to sample preparation for HPLC. Alternatively, SPME may be the key to an analysis that, until now, has been difficult or impossible to perform. In either case, the combination of SPME with HPLC should prove valuable in many environmental, pharmacological, and food and beverage applications.

Figure B. Antioxidants in Foods

Sample: 3mL olive oil containing 1mg/g each analyte
0.6g soft drink powder in 3mL water
SPME Fiber: **40 μ m PDMS/DVB-HPLC**
Cat. No.: **57317**
Extraction: immersion, 15 min (olive oil) or 20 min (powder)(rapid stirring)
Desorption: static, acetonitrile:methanol, 50:50, 5 min (olive oil) or 3 min (powder)
Column: **SUPELCO SILLC-18, 15cm x 4.6mm ID, 5 μ m particles**
Cat. No.: **58230**
Mobile Phase: A = 5% acetic acid in water
B = acetonitrile:methanol, 1:1
30% B for 2 min, then to 100% B over 10 min, hold 10 min
Flow Rate: 2mL/min
Det.: UV, 280nm




796-0494, 0495

Trademarks

Rheodyne – Rheodyne, Inc.
SUPELCO SILL, Supelguard – Sigma-Aldrich Co.
Valco – Valco Instruments Co., Inc.

Figure C. Benzoic Acid in Diet Cola

Sample: cola in water, 1:10 (pH to 2.2 with acetic acid), 3mL
SPME Fiber: **40 μ m PDMS/DVB-HPLC**
Cat. No.: **57317**
Extraction: immersion, 25 min (rapid stirring)
Desorption: static, mobile phase, 1 min
Column: **SUPELCO SILLC-18, 15cm x 4.6mm ID, 3 μ m particles**
Cat. No.: **58985**
Mobile Phase: acetonitrile:0.02M sodium acetate (pH 4.3 with acetic acid), 20:80
Flow Rate: 1.5mL/min
Det.: UV, 254nm



796-0496

Ordering Information:

Description	Cat. No.
SPME Holder**	
For manual sampling	57330-U
For HPLC	57331
SPME/HPLC Interface	
with Valco® valve	57350-U
with Rheodyne® valve	57353
SPME/HPLC Fiber Assembly	
40 μ m PDMS/DVB-HPLC	57317
HPLC Columns	
SUPELCO SILLC-PCN	
25cm x 4.6mm ID, 5 μ m particles	58378
SUPELCO SILLC-18	
15cm x 4.6mm ID, 5 μ m particles	58230-U
SUPELCO SILLC-18	
15cm x 4.6mm ID, 3 μ m particles	58985

References

- For additional information about the HPLC/SPME interface request Product Specification 496049.
 - Supelco Bulletin 901 (available on request).
 - AOAC Official Method 983.15, *AOAC Official Methods of Analysis*, Supplement March 1995.
 - AOAC Official Method 983.16, *AOAC Official Methods of Analysis*, Supplement March 1995.
- References 3 and 4 not available from Supelco (obtain from Association of Official Analytical Chemists, Arlington, VA USA).

For descriptions of our SPME fibers and accessories, request publication 413019.

For chemical standards for tricyclic antidepressants, refer to the Supelco catalog. For BHA, BHT, and other antioxidant/preservative standards, request our Food & Beverage Analysis Guide.

* Technology licensed exclusively to Supelco. US patent #5,691,206; European patent #0523092.

** Initially you must order both holder and fiber assembly. Holder is reusable indefinitely. Cat. No. 5-7331 also is used with Varian 8100/8200 AutoSampler (requires Varian SPME upgrade kit, available from Varian).