



Application Note 2003/1



Stir Bar Sorptive Extraction (Twister™)
RTL-CGC-MS.
A Versatile Method to Monitor more
than 400 Pesticides in Different
Matrices (Water, Beverages, Fruits,
Vegetables, Baby Food)

Pat Sandra, Bart Tienpont, Frank David
*Research Institute for Chromatography, Kennedypark
20, 8500 Kortrijk, Belgium*

Keywords

Pesticides, Multi-Residue Analysis, SBSE,
TDU, CGC-MS

Introduction

The analysis of pesticides in solid samples (fruits, vegetables, baby food) is discussed in detail because this is much more challenging than aqueous samples (water, beverages). For beverages e.g. wine see Pat Sandra, Bart Tienpont, Joeri Vercammen, Andreas Tredoux, Tom Sandra, Frank David, *J. Chromatogr. A*, 928 (2001) 117. Pesticides from vegetables, fruits and baby food are extracted with methanol and an aliquot is diluted with water for SBSE. Series of 98 stir bars are automatically analyzed by application of a new thermal desorption unit (TDU) developed by Gerstel, in combination with retention time locked

capillary gas chromatography-mass spectrometry (RTL-capillary GC-MS) analysis. Elucidation and identification is finally performed automatically using the Agilent RTL library for pesticides. Accurate quantitation is performed by standard addition or isotope dilution.

Experimental

Sample Preparation

Sample: 15 g vegetable, fruit, baby food
← 30 mL methanol
Ultraturrax: 5 min + Ultrasonic bath: 15 min
1 mL extract
+ 10 mL water
SBSE (Twister 10 mm L x 0.5 mm d _r): 60 min
TD-RTL-CGC-MS

Instrumental setup



Figure 1. Sample rack of MPS2-TDU with a capacity of 98 positions; Twisters are put in a clean empty glass liner and capped with a special tube head.



Figure 2. MPS-2-TDU-CGC-MS setup; The MPS-2 robot inserts the liners into the TDU for thermal desorption.

Instrumental

TDU	Gerstel TDU with Gerstel MPS-2 XYZ robot 40°C (1 min) – 60°C/min – 280°C (5 min)
PTV	Gerstel CIS-4 -150°C – 600°C/min – 280°C (2 min)
GC	Agilent 6890 (Agilent Technologies)
Oven	70°C (2 min) – 25°C/min – 150°C – 3°C/min – 200°C – 8°C/min – 300°C
Column	30 m L x 250 µm I.D., 0.25 µm d _r HP5-MS (Agilent Technologies)
Carrier	Helium, constant pressure eluting p,p'DDT @ 26.98 min
MSD	Agilent 5973N (Agilent Technologies) Full scan @ m/z 40-500
RTL library	Agilent RTLpest.L

Applications

SBSE of pesticides

The recovery of SBSE for pesticides can be calculated from their log K_{ow} using the [SBSE recovery calculator](#) package (Application Note 2003/2 of this website). Tables 1 to 4 show alphabetical lists of GC amenable pesticides from the Agilent RTL pesticide library that can be analyzed by SBSE. RTLocked retention times (RT), four MSD qualifier ions (Tion, Q1, Q2, Q3), the CAS number, log K_{ow} and SBSE recoveries for Twisters of 10 mm L x 0.5 mm df (Twister 1) and 20 mm L x 1.0 mm df (Twister 2) are given. The tables comprise four separate lists of pesticides that are grouped in function of SBSE recoveries for Twister 1 between 10-24% (Table 1), 25-49% (Table 2), 50-74% (Table 3) and 75-100% (Table 4). Recoveries are calculated for 10 mL water samples. From all 440 pesticides, 367 compounds are in Table 4 and show very high enrichments in the PDMS coated stir bar. More polar compounds in Table 1 to 3 can also be recovered at minimum 10% and allows their detection at sub- $\mu\text{g/L}$ level.

These recovery values are calculated for water samples and at equilibrium conditions. In the case of solid samples, an extract in methanol is made and diluted with water and SBSE is performed during 60 min. Matrix interferences (10% methanol, sample matrix) and sampling at non-equilibrium conditions decrease the recovery percentages slightly. However, Tables 1 to 4 give a good approximation of the feasibility of SBSE analysis for the listed pesticides.

Thermally sensitive pesticides are annotated with an asterisk in the tables and possible partial degradation can occur during TD-CGC-MS analysis. Thermal desorption using the TDU is however equal or better than classical split/splitless injection. However, it was demonstrated that thermal desorption under the prescribed conditions occurs under mild conditions and thermo-sensitive pesticides like captan can be recovered for more than 70%.

RTL result screener

The pesticides are automatically identified using the Agilent Result Screener software in combination with the Agilent RTL library for pesticides (RTLPest.I, Agilent Technologies). An example of the detection of permethrin II in grapes is shown in Figure 3 and shows the extracted ion chromatograms at the selective target and qualifier ions of the pesticide (m/z 163, 165, 183, 184) and the mass spectrum at peak apex.




Figure 3. Result screener window of the automatic detection of permethrin II in grapes by SBSE-TD-RTL-CGC-MS; **A:** extracted ion chromatograms at m/z 183, 163, 165, 184; **B:** recorded mass spectrum; **C:** identification results.

Besides permethrin II, permethrin I, vinclozolin and procymidone were detected (Figure 4).




Figure 4. Extracted ion chromatograms at m/z 183, 197 and 283 for vinclozolin, procymidone, permethrin I and permethrin II.

An other example is shown for the detection of benalaxyl, vinclozolin, tolylfluand and procymidone in apple (Figure 5 and 6).




Figure 5. Result screener window of the automatic detection of benalaxyl in apple by SBSE-TD-RTL-CGC-MS; **A:** extracted ion chromatograms at m/z 148, 91, 204, 206; **B:** recorded mass spectrum; **C:** identification results.




Figure 6. Extracted ion chromatograms at m/z 137, 148, 169, 283 for vinclozolin, tolylfluand, procymidone and benalaxyl.

Quantitative analysis

Quantitation of pesticides can be performed into two manners: first, external calibration is performed by spiking a pesticide mixture in an aliquot of methanol (1 mL) before dilution with water and SBSE-TD-CGC-MS analysis. Pat Sandra *et al.* demonstrated negligible matrix effects in the case of vegetable and fruit samples and this permits semi-quantitative measurements of pesticides and differentiation of positives from negatives for in depth quantitation.

Standard addition calibration or isotope dilution are the only solutions for accurate MRM analysis. The paper in reference [1] describes the quantification procedures in detail and illustrates that the required limits of quantitation (LOQ) even for baby food can easily be reached.

References

- [1] P. Sandra, B. Tienpont, F. David, *accepted for publication in J. Chromatogr. A.*

Table 1. Pesticides that result in SBSE recoveries (for 10 mL water) between 10 and 24% for Twister 1 (containing 24 μ L PDMS) and 38 and 62% for Twister 2 (containing 116 μ L PDMS).

Compound Name	Cpd	RT	Tlon	Q1	Q2	Q3	CAS	log Kow	Twister 1	Twister 2
3-Chloroaniline	16	5.45	127	129	65	92	108-42-9	1,72	11	40
4-Chloroaniline	17	5.48	127	129	65	92	106-47-8	1,72	11	40
Allidochlor	22	6.18	41	56	138	132	93-71-0	1,79	13	44
Aminocarb	148	13.58	151	150	136	134	2032-59-9	1,9	16	50
Ancymidol	413	25.05	228	107	121	215	12771-68-5	1,99	19	55
Bentazone	292	19.71	198	119	161	92	25057-89-0	1,67	10	37
Bromacil	260	18.35	205	207	206	42	314-40-9	1,68	10	38
Crotoxyphos	355	22.21	127	105	193	104	7700-17-6	1,89	16	49
Demeton-S	122	12.63	88	60	170	89	8065-48-3	2,09	23	61
Diethyl dithiobis(thionoformate) (EXD)	123	12.65	76	89	94	60	502-55-6	2,06	22	59
Ethiofencarb	196	15.61	107	168	77	57	29973-13-5	2,04	21	58
Ethiolate	15	5.41	100	72	161	44	2941-55-1	2,04	21	58
Famphur	447	26.64	218	125	217	93	52-85-7	1,99	19	55
Fenthion sulfoxide	422	25.54	278	125	109	279	-	1,92	17	51
Metalaxyl	238	17.34	206	45	160	249	57837-19-1	1,7	11	39
Metolcarb	46	7.99	108	107	77	79	1129-41-5	1,72	11	40
N-Methyl-N-1-naphthyl acetamide	146	13.56	157	156	56	199	-	2,1	23	61
Paraoxon	239	17.34	109	149	139	275	311-45-5	1,97	18	54
Propoxur	76	10.35	110	152	57	111	114-26-1	1,9	16	50
Pyracarbolid	302	20.02	125	97	217	55	24691-76-7	2,08	22	60
Pyroquilon	154	13.79	173	130	172	144	57369-32-1	1,8	13	44
Terbacil	178	14.77	161	160	117	163	5902-51-2	1,75	12	41
Thiabendazole	320	20.94	201	174	202	175	148-79-8	2	19	56
Thionazin	74	10.26	96	97	107	143	297-97-2	1,86	15	48
Tycor (SMY 1500)	245	17.58	212	200	213	61	64529-56-2	1,98	19	55

Table 2. Pesticides that result in SBSE recoveries (for 10 mL water) between 25 and 49% for Twister 1 (containing 24 µL PDMS) and 63 and 83% for Twister 2 (containing 116 µL PDMS); (*) thermally sensitive pesticides

Compound Name	Cpd	RT	Tlon	Q1	Q2	Q3	CAS	log Kow	Twister 1	Twister 2
2,3,5-Trimethacarb	96	11.20	136	121	91	135	2655-15-4	2,56	47	82
2,4-Dimethylaniline	11	5.19	121	120	106	77	1300-73-8	2,17	26	65
2,6-Dimethylaniline	12	5.20	121	106	120	77	87-62-7	2,17	26	65
3,4-Dichloroaniline	37	7.66	161	163	99	90	95-76-1	2,37	36	75
3,5-Dichloroaniline	34	7.32	161	163	90	99	626-43-7	2,37	36	75
4,6-Dinitro-o-cresol (DNOC)	78	10.36	198	105	121	51	534-52-1	2,27	31	70
Azinphos-methyl	514	29.65	160	132	77	105	86-50-0	2,53	45	81
Bendiocarb	102	11.54	151	126	166	57	22781-23-3	2,55	46	82
Carbaryl	225	16.81	144	115	116	145	63-25-2	2,35	35	74
Carbofuran	131	13.03	164	149	131	123	1563-66-2	2,3	32	72
Chlorotoluron*	261	18.37	132	72	167	212	15545-48-9	2,58	48	83
Chlozolinate	336	21.38	188	259	186	187	72391-46-9	2,22	28	68
Cyanazine	282	19.35	212	213	214	68	21725-46-2	2,51	44	80
Cyanophos	152	13.77	109	243	125	79	2636-26-2	2,48	42	79
Dicamba	91	10.95	173	175	220	222	1918-00-9	2,14	25	63
Dichlormid	28	6.79	41	172	124	56	37764-25-3	2,28	31	71
Dicyclopentadiene	4	4.00	66	132	67	65	542-92-7	2,25	30	69
Dimethachlor	208	16.18	134	197	132	210	50563-36-5	2,33	34	73
Drazoxolon	297	19.91	127	125	44	129	5707-69-7	2,44	40	78
Fensulfothion	421	25.56	292	293	156	140	115-90-2	2,35	35	74
Fluometuron*	89	10.81	72	187	232	168	2164-17-2	2,35	35	74
Flutriafol	369	23.18	123	219	164	83	76674-21-0	2,52	44	81
Fuberidazole	216	16.55	155	156	129	92	3878-19-1	2,37	36	75
Hexazinone	464	27.39	171	83	128	71	51235-04-2	2,15	25	64
Isoprocarb	58	9.09	121	136	91	77	2631-40-5	2,37	36	75
Malathion	270	18.80	173	127	125	93	121-75-5	2,29	32	71
Malathion-o-analog	229	16.87	127	99	109	125	121-75-5	2,29	32	71
Mecarbam	344	21.74	131	97	159	125	2595-54-2	2,29	32	71
Metazachlor	317	20.83	133	81	132	134	67129-08-2	2,38	37	75
Metobromuron*	192	15.51	61	197	199	90	3060-89-7	2,51	44	80
Monolinuron*	134	13.10	61	153	127	125	1746-81-2	2,26	30	70
N,N-Diethyl-m-toluamide	69	9.75	119	190	91	191	134-62-3	2,26	30	70
Norflurazon	456	26.95	303	145	102	305	27314-13-2	2,19	27	66
Phosmet	488	28.50	160	161	77	93	732-11-6	2,48	42	79
p-Nitrotoluene	19	5.57	91	137	65	107	99-99-0	2,36	35	74
Procymidone	350	21.96	96	283	285	67	32809-16-8	2,59	48	83
Propachlor	77	10.36	120	176	77	93	1918-16-7	2,42	39	77
Propetamphos	158	13.91	138	194	236	222	31218-83-4	2,51	44	80
Prothoate	213	16.46	286	288	125	109	2275-18-5	2,17	26	65
Simazine	129	12.91	201	186	44	173	122-34-9	2,4	38	76
Thiofanox	100	11.47	57	68	115	56	39196-18-4	2,16	26	65
Tricyclazole	379	23.61	189	162	161	135	41814-78-2	2,48	42	79

Table 3. Pesticides that result in SBSE recoveries (for 10 mL water) between 50 and 74% for Twister 1 (containing 24 μ L PDMS) and 84 and 92% for Twister 2 (containing 116 μ L PDMS); (*) thermally sensitive pesticides

Compound Name	Cpd	RT	Tlon	Q1	Q2	Q3	CAS	log Kow	Twister 1	Twister 2
2,3,5-Trimethylphenyl methyl carbamate	97	11.23	136	121	91	135	12407-86-2	2,81	61	89
2,6-Dichlorobenzonitrile	24	6.75	171	173	136	100	1194-65-6	2,83	62	89
2-Hydroxyestradiol	555	33.81	288	289	229	176	50-27-1	2,81	61	89
3,4,5-Trimethacarb	127	12.78	121	136	135	57	2686-99-9	2,81	61	89
Atrazine	136	13.16	200	215	202	58	1912-24-9	2,82	61	89
Azaconazole	402	24.55	217	219	173	175	60207-31-0	2,73	56	87
Benazolin-ethyl	321	20.96	170	198	271	172	25059-80-7	2,86	63	90
Benfuresate	203	15.97	163	256	121	164	68505-69-1	2,8	60	89
Captan*	330	21.23	79	80	151	77	133-06-2	2,74	57	87
Chlorbufam	133	13.06	153	127	53	155	1967-16-4	3,02	72	93
Chlordimeform	94	11.20	196	181	117	152	6164-98-3	2,89	65	91
Chlormefos	39	7.73	121	97	154	234	24934-91-6	3,04	72	93
Chlorthiamid	205	16.04	171	173	170	172	1918-13-4	2,96	69	92
Clomazone	138	13.25	125	204	127	205	81777-89-1	2,86	63	90
Cycluron*	149	13.63	72	89	127	45	2163-69-1	2,84	62	90
Desmetryn	204	15.99	213	198	171	58	1014-69-3	2,82	61	89
Dichlofluanid	262	18.41	123	167	224	226	1085-98-9	2,72	56	87
Dichlone	166	14.29	191	226	228	163	117-80-6	2,65	52	85
Dichlorprop	126	12.69	162	164	44	234	120-36-5	3,03	72	93
Dicloran	120	12.56	206	176	178	208	99-30-9	2,76	58	88
Diphenamid	307	20.20	167	72	165	239	957-51-7	2,86	63	90
Endosulfan lactone	309	20.56	277	321	239	237	-	2,91	66	91
EPTC	27	6.80	128	43	86	132	759-94-4	3,02	72	93
Ethofumesate	258	18.31	207	161	286	137	26225-79-6	2,89	65	91
Etrimfos	184	15.16	292	181	153	125	38260-54-7	2,94	68	92
Fenfuram	174	14.62	109	201	43	110	24691-80-3	2,62	50	84
Fenobucarb	75	10.28	121	150	91	122	3766-81-2	2,86	63	90
Folpet	338	21.59	104	76	260	147	133-07-3	2,84	62	90
Furalaxyl	349	21.92	95	242	152	146	57646-30-7	2,7	55	86
Iprodione*	485	28.39	187	314	189	244	36734-19-7	2,85	63	90
Isoprothiolane	388	23.88	118	162	189	290	50512-35-1	2,79	60	89
Isoproturon*	231	17.09	146	72	206	128	34123-59-6	2,84	62	90
Linuron*	254	18.19	61	187	189	124	330-55-2	2,91	66	91
Mefenacet	518	29.97	192	77	120	106	73250-68-7	2,8	60	89
Mefluidide	319	20.82	135	268	310	136	53780-34-0	2,72	56	87
Methiocarb	250	18.05	168	153	109	57	2032-65-7	2,87	64	90
Methoprotryne	406	24.71	256	213	226	240	841-06-5	3,04	72	93
Methyl parathion	217	16.59	263	109	125	79	298-00-0	2,75	57	88
Molinat	57	9.11	126	55	187	83	2212-67-1	2,91	66	91
Nitralin	480	28.18	316	274	300	317	4726-14-1	2,92	67	91
Oryzalin	528	31.25	317	275	58	258	19044-88-3	2,73	56	87
Oxabetrinil	188	15.31	73	103	45	77	74782-23-3	2,78	59	88
Pindone	165	14.28	173	174	146	230	83-26-1	2,87	64	90

Probenazole	259	18.33	103	130	76	159	27605-76-1	2,89	65	91
Propanil	206	16.12	161	163	217	219	709-98-8	2,88	65	91
Propham	42	7.91	93	179	119	137	122-42-9	2,66	52	85
Pyridinitril	340	21.64	273	275	237	202	1086-02-8	2,95	68	92
Quinalphos	341	21.65	146	157	156	118	13593-03-8	3,04	72	93
Simetryn	227	16.82	213	170	155	198	1014-70-6	2,9	66	91
Swep	132	13.05	187	189	124	219	1918-18-9	3,04	72	93
Thiometon	116	12.34	88	125	89	93	640-15-3	2,88	65	91
Triadimefon	283	19.39	57	208	85	210	43121-43-3	2,94	68	92
Triadimenol	342	21.67	112	168	128	70	55219-65-3	2,95	68	92
Triazophos	442	26.46	161	162	172	77	24017-47-8	2,92	67	91
Vinclozolin	219	16.63	212	285	198	187	50471-44-8	3,03	72	93

Table 4. Pesticides that result in SBSE recoveries (for 10 mL water) between 75 and 100 % for Twister 1 (containing 24 µL PDMS) and 93 and 100% for Twister 2 (containing 116 µL PDMS); (*) thermally sensitive pesticides

Compound Name	Cpd	RT	Tlon	Q1	Q2	Q3	CAS	log Kow	Twister 1	Twister 2
2-(Octylthio)ethanol	56	8.98	145	61	190	159	3547-33-9	3,38	85	97
2,4-D sec-butyl ester	176	14.73	57	175	177	162	-	4,3	98	100
4,4'-Dichlorobenzophenone	276	19.20	139	111	141	250	90-98-2	4,44	99	100
5,7-Dihydroxy-4'-methoxyisoflavone	551	33.41	284	132	283	285	552-59-0	3,41	86	97
Acetochlor	215	16.54	146	162	59	223	34256-82-1	3,37	85	97
Alachlor	230	17.03	160	188	45	146	15972-60-8	3,37	85	97
Aldrin*	265	18.53	263	66	265	261	309-00-2	6,75	100	100
Ametryn	232	17.11	227	212	170	185	834-12-8	3,32	83	96
Amitraz	519	30.14	121	162	132	147	33089-61-1	5,55	100	100
Anilazine	322	20.95	239	241	178	143	101-05-3	3,64	91	98
Atraton	125	12.70	196	211	169	58	1610-17-9	3,15	77	95
Azinphos-ethyl	522	30.65	132	160	77	105	2642-71-9	3,51	89	98
Aziprotryne	173	14.58	199	184	139	68	4658-28-0	3,27	82	96
Azobenzene	83	10.59	77	182	105	51	103-33-3	4,11	97	99
Barban	392	24.07	153	127	155	125	101-27-9	3,41	86	97
Benalaxyl	448	26.75	148	91	206	204	71626-11-4	3,69	92	98
Benfluralin	104	11.73	292	264	276	293	1861-40-1	5,31	100	100
Benodanil	434	25.95	231	323	203	76	15310-01-7	3,87	95	99
Benthiocarb	267	18.58	100	72	125	257	28249-77-6	3,9	95	99
Benzophenone	84	10.67	105	182	77	51	119-61-9	3,15	77	95
Benzoylprop ethyl	484	28.33	105	77	292	106	22212-55-1	4,27	98	100
BHC alpha isomer	113	12.08	181	219	183	217	58-89-9	4,26	98	100
BHC beta isomer	137	13.20	219	181	183	217	93-72-1	3,68	92	98
BHC delta isomer	169	14.54	181	219	183	217	93-72-1	3,68	92	98
Bifenox	506	29.18	341	343	311	189	42576-02-3	4,15	97	99
Bifenthrin	500	28.84	181	165	166	182	82657-04-3	8,15	100	100
Binapacryl	416	25.18	83	55	82	84	485-31-4	4,49	99	100
Bioallethrin	347	21.80	123	136	79	107	584-79-2	5,52	100	100
Bioallethrin S-cyclopentenyl isomer	345	21.75	123	136	79	107	28434-00-6	5,52	100	100
Bioresmethrin	477	27.98	123	171	128	143	28434-01-7	7,11	100	100
Bitertanol I	527	31.21	170	57	168	171	55179-31-2	4,07	97	99
Bitertanol II	529	31.34	170	57	168	171	55179-31-2	4,07	97	99
Bromobutide	211	16.25	119	120	118	232	74712-19-9	4,33	98	100
Bromocyclen	185	15.19	357	359	237	361	1715-40-8	6,13	100	100
Bromophos	304	20.08	331	329	125	333	2104-96-3	5,11	100	100
Bromophos-ethyl	360	22.53	359	303	357	301	4824-78-6	6,09	100	100
Bromopropylate	494	28.62	341	339	343	183	18181-80-1	4,9	99	100
Bromoxynil	101	11.55	277	275	279	88	1689-84-5	3,39	85	97
Bromoxynil octanoic acid ester	452	26.83	127	57	128	88	1689-99-2	5,86	100	100
Buprofezin	404	24.58	105	106	104	172	69327-76-0	4,3	98	100
Butachlor	371	23.22	176	160	188	146	23184-66-9	4,84	99	100
Butamifos	378	23.59	286	200	96	202	36335-67-8	4,77	99	100
Butralin	308	20.23	266	224	267	220	33629-47-9	5,15	100	100

Butylate	36	7.61	57	146	156	174	2008-41-5	3,85	94	99
Captafol*	468	27.62	79	80	77	151	2425-06-1	3,42	86	97
Carbophenothion	446	26.64	157	121	153	97	786-19-6	5,19	100	100
Carbosulfan	497	28.66	164	160	149	118	55285-14-8	5,57	100	100
Chinomethionat	348	21.90	206	234	116	148	2439-01-2	3,37	85	97
Chloranocryl	222	16.74	69	125	229	231	2164-09-2	3,29	82	96
Chlorbenside	346	21.82	125	127	268	270	103-17-3	5,59	100	100
Chlorbromuron*	323	20.99	61	233	231	124	13360-45-7	3,15	77	95
Chlordecone	436	26.02	272	274	270	237	143-50-0	4,91	99	100
Chlorfenethol	326	21.06	251	139	253	178	80-06-8	4,45	99	100
Chlorfenprop-methyl	70	9.96	125	165	196	197	14437-17-3	3,4	86	97
Chlorfenson	372	23.32	175	111	177	302	80-33-1	4,21	97	100
Chlorfenvinphos	337	21.56	267	323	269	325	470-90-6	4,15	97	99
Chlorflurecol-methyl ester	353	22.08	215	152	217	216	-	3,27	82	96
Chlornitrofen	444	26.53	317	319	287	236	1836-77-7	4,96	100	100
Chlorobenzilate	419	25.40	251	139	253	111	510-15-6	3,99	96	99
Chloroneb	50	8.68	191	193	206	208	2675-77-6	3,44	87	97
Chloropropylate	420	25.42	251	253	139	111	5836-10-2	4,41	98	100
Chlorothalonil	179	14.78	266	264	268	270	1897-45-6	3,66	92	98
Chlorpropham	92	11.05	127	213	153	154	101-21-3	3,3	83	96
Chlorpyrifos	278	19.23	197	199	314	97	2921-88-2	4,66	99	100
Chlorpyrifos Methyl	218	16.59	286	288	125	290	2921-88-2	4,66	99	100
Chlorthal-dimethyl	286	19.43	301	299	303	332	1861-32-1	4,24	98	100
Chlorthion	293	19.75	109	125	297	267	500-28-7	3,39	85	97
Chlorthiophos	438	26.13	269	97	325	271	60238-56-4	5,8	100	100
Chlorthiophos sulfone	508	29.37	301	357	97	303	-	3,77	93	99
Chlorthiophos sulfoxide	491	28.56	97	341	269	285	-	3,64	91	98
cis-Chlordane	365	22.83	373	375	377	371	57-74-9	6,26	100	100
Coumaphos	534	31.67	362	226	109	210	56-72-4	4,47	99	100
Crufomate	298	19.90	256	169	182	276	299-86-5	3,3	83	96
Cyanofenphos	451	26.82	157	169	141	185	13067-93-1	4,2	97	100
Cycloate	85	10.76	83	154	72	84	1134-23-2	3,81	94	99
Cyfluthrin I	539	32.22	163	206	165	227	68359-37-5	5,74	100	100
Cyfluthrin II	540	32.36	163	206	165	227	68359-37-5	5,74	100	100
Cyfluthrin III	541	32.48	163	206	165	227	68359-37-5	5,74	100	100
Cyfluthrin IV	542	32.54	163	206	227	199	68359-37-5	5,74	100	100
Cyhalothrin I (lambda)	520	30.37	181	197	208	209	68085-85-8	6,85	100	100
Cymoxanil	80	10.45	44	70	111	167	57966-95-7	4,24	98	100
Cypermethrin I	543	32.69	181	163	165	77	52315-07-8	6,38	100	100
Cypermethrin II	544	32.84	181	163	165	209	52315-07-8	6,38	100	100
Cypermethrin III	546	32.97	163	181	165	209	52315-07-8	6,38	100	100
Cypermethrin IV	547	33.02	163	181	165	209	52315-07-8	6,38	100	100
Cyprazine	209	16.18	212	227	170	214	22936-86-3	3,12	76	94
d-(cis-trans)-Phenothrin-I	510	29.45	123	183	81	184	26002-80-2	7,54	100	100
d-(cis-trans)-Phenothrin-II	511	29.61	123	183	81	184	26002-80-2	7,54	100	100
Deltamethrin*	566	36.00	181	253	251	255	52918-63-5	6,18	100	100
Desbromo-bromobutide	108	11.85	119	120	118	233	-	3,74	93	99
Desmedipham	118	12.37	181	135	109	122	13684-56-5	3,22	80	95
Dialifos	526	30.83	208	173	210	76	10311-84-9	4,14	97	99
Di-allate I	110	11.96	86	234	128	236	2303-16-4	4,08	97	99

Di-allate II	115	12.29	86	234	236	128	2303-16-4	4,08	97	99
Diazinon	168	14.47	179	137	152	199	333-41-5	3,86	95	99
Dicapthon	284	19.41	262	125	263	79	2463-84-5	3,39	85	97
Dichlofenthion	207	16.17	279	223	97	162	97-17-6	5,2	100	100
Dichlorophen	486	28.37	128	141	268	270	97-23-4	4,34	98	100
Diclobutrazol	398	24.42	270	272	159	82	75736-33-3	4,01	96	99
Diclofop methyl	469	27.68	253	340	255	342	51338-27-3	4,54	99	100
Dieldrin	387	23.87	79	263	277	279	60-57-1	5,45	100	100
Diethyl ethyl	381	23.60	188	162	238	160	38727-55-8	3,62	91	98
Diethofencarb	274	19.12	151	207	150	267	87130-20-9	3,29	82	96
Difenoconazol I	563	35.09	323	265	325	267	119446-68-3	5,2	100	100
Difenoconazol II	564	35.23	323	265	325	267	119446-68-3	5,2	100	100
Diflufenican	472	27.79	266	394	267	246	83164-33-4	3,53	89	98
Dimethametryn	327	21.11	212	213	255	71	22936-75-0	4,22	98	100
Dimethylvinphos(z)	275	19.15	295	297	109	299	2274-67-1	3,16	78	95
Diniconazole	423	25.55	268	270	70	269	83657-24-3	3,92	95	99
Dinitramine	180	14.81	305	307	261	216	29091-05-2	3,96	96	99
Dinobuton	351	21.98	45	44	211	207	973-21-7	3,94	95	99
Dinocap I	474	27.85	69	44	41	159	39300-45-3	5,98	100	100
Dinocap II	482	28.22	69	44	68	70	39300-45-3	5,98	100	100
Dinocap III	495	28.60	69	70	197	77	39300-45-3	5,98	100	100
Dinocap IV	504	28.93	69	44	207	41	39300-45-3	5,98	100	100
Dinoseb	171	14.54	211	163	147	117	88-85-7	3,67	92	98
Dinoseb acetate	241	17.53	60	240	211	205	2813-95-8	3,17	78	95
Dinoterb	163	14.07	225	177	131	77	1420-07-1	3,64	91	98
Dinoterb acetate	248	17.91	225	77	177	240	3204-27-1	3,13	76	94
Dioxathion	537	31.84	97	125	271	153	78-34-2	3,45	87	97
Diphacinone	557	34.26	173	168	167	165	82-66-6	4,85	99	100
Diphenylamine	82	10.52	169	168	170	83	122-39-4	3,29	82	96
Dipropetryn	268	18.72	255	240	184	222	4147-51-7	4,22	98	100
Disulfoton	170	14.55	88	89	97	142	298-04-4	3,86	95	99
Ditalimfos	370	23.17	130	148	299	209	5131-24-8	3,74	93	99
Dithiopyr	252	18.10	354	306	286	237	97886-45-8	4,46	99	100
Diuron*	334	21.37	72	187	124	189	470-90-6	4,15	97	99
Dodemorph I	299	19.92	154	141	155	281	1593-77-7	5,7	100	100
Dodemorph II	316	20.80	154	141	155	238	1593-77-7	5,7	100	100
Edifenphos	449	26.79	109	173	110	310	17109-49-8	3,61	91	98
Endosulfan (alpha isomer)	362	22.64	241	195	239	237	115-29-7	3,5	88	98
Endosulfan (beta isomer)	414	25.16	195	237	207	241	115-29-7	3,5	88	98
Endosulfan ether	189	15.35	69	241	239	277	-	4,15	97	99
Endosulfan sulfate	450	26.76	272	274	229	237	1031-07-8	3,64	91	98
Endrin	408	24.76	317	263	315	319	72-20-8	5,45	100	100
Endrin aldehyde	431	25.91	67	345	250	347	7421-93-4	4,8	99	100
Endrin ketone	481	28.23	317	67	315	319	53494-70-5	4,99	100	100
EPN	496	28.66	157	169	141	185	2104-64-5	4,47	99	100
Epoxiconazole	478	28.02	192	194	138	165	133855-98-8	3,47	88	97
Erbon	411	24.89	169	171	97	135	136-25-4	5,63	100	100
Esfenvalerate	560	34.69	125	167	181	209	51630-58-1	6,76	100	100
Esprocarb	256	18.25	91	222	71	162	85785-20-2	4,58	99	100
Etaconazole	427	25.76	173	245	55	175	60207-93-4	4,13	97	99

Ethalfuralin	98	11.28	276	316	292	333	55283-68-6	5,23	100	100
Ethion	435	26.00	231	153	97	125	563-12-2	4,75	99	100
Ethoprophos	86	10.74	158	97	126	139	13194-48-4	3,14	77	95
Ethoxyquin	128	12.83	202	174	203	145	91-53-2	3,87	95	99
Etridiazole	44	7.96	211	183	213	185	2593-15-9	3,6	91	98
Fenarimol	521	30.40	139	219	107	251	60168-88-9	3,62	91	98
Fenazaflor	439	26.26	77	374	234	376	14255-88-0	4,54	99	100
Fenbuconazole	538	32.23	129	198	125	127	114369-43-6	4,23	98	100
Fenchlorphos	236	17.33	285	287	125	289	299-84-3	4,86	99	100
Fenitrothion	251	18.07	277	125	109	260	122-14-5	3,3	83	96
Fenoprop	228	16.80	196	198	97	200	93-72-1	3,68	92	98
Fenoxycarb	498	28.70	255	186	77	185	72490-01-8	4,24	98	100
Fenpropathrin	505	28.99	97	181	125	265	39515-41-8	5,62	100	100
Fenpropimorph	279	19.27	128	129	303	117	67564-91-4	5,5	100	100
Fenson	290	19.66	77	141	268	51	80-38-6	3,57	90	98
Fenthion	273	19.12	278	125	109	169	55-38-9	4,08	97	99
Fenvalerate I	558	34.27	167	125	181	152	51630-58-1	6,76	100	100
Fenvalerate II	559	34.68	125	167	181	169	51630-58-1	6,76	100	100
Flamprop-isopropyl	429	25.82	105	77	276	106	52756-22-6	4,24	98	100
Flamprop-methyl	403	24.54	105	77	276	106	52756-25-9	3,33	84	96
Fluazifop-p-butyl	417	25.32	282	254	383	255	69806-50-4	5,34	100	100
Flubenzimine	397	24.36	135	186	77	416	37893-02-0	6,66	100	100
Fluchloralin	175	14.64	306	326	264	63	33245-39-5	5,07	100	100
Flucythrinate I	548	33.12	199	157	44	207	70124-77-5	6,56	100	100
Flucythrinate II	552	33.46	199	157	44	181	70124-77-5	6,56	100	100
Flumetralin	373	23.34	143	145	157	404	62924-70-3	6,09	100	100
Fluorodifen	382	23.73	190	126	160	298	15457-05-3	3,81	94	99
Fluotrimazole	479	28.04	165	311	312	379	31251-03-3	5,76	100	100
Fluridone	554	33.77	328	329	330	154	2051-24-3	10,2	100	100
Flurochloridone I	300	19.96	174	311	187	145	61213-25-0	3,45	87	97
Flurochloridone II	305	20.10	174	311	187	145	31213-25-0	3,45	87	97
Fluroxypyr-1-methylheptyl ester	470	27.73	57	71	209	181	81406-37-3	4,82	99	100
Flusilazole	405	24.59	233	206	234	315	85509-19-9	4,89	99	100
Flutolanil	385	23.84	173	145	281	323	66332-96-5	4,65	99	100
Fluvalinate-tau-I	561	34.72	250	252	209	181	69409-94-5	6,81	100	100
Fluvalinate-tau-II	562	34.85	250	252	209	181	69409-94-5	6,81	100	100
Fonofos	157	13.89	109	246	137	110	944-22-9	4,02	96	99
Furathiocarb	512	29.61	163	164	149	76	65907-30-4	4,43	98	100
Furmecyclox	197	15.64	123	124	139	221	60568-05-0	4,38	98	100
Heptachlor	224	16.80	272	100	274	270	76-44-8	5,86	100	100
Heptachlor epoxide	313	20.72	353	355	351	357	1024-57-3	4,56	99	100
Heptachlor exo-epoxide isomer B	312	20.72	353	355	351	357	1024-57-3	4,56	99	100
Hexabromobenzene	487	28.45	551	549	553	232	87-82-1	7,33	100	100
Hexachlorobenzene	117	12.38	284	286	282	288	118-74-1	5,86	100	100
Hexachlorophene	550	33.20	196	198	211	210	70-30-4	6,92	100	100
Hexaconazole	376	23.52	83	214	216	82	79983-71-4	3,66	92	98
Hexestrol	457	26.98	135	107	134	136	84-16-2	5,6	100	100
Imazalil	384	23.79	41	215	173	217	35554-44-0	4,1	97	99
Ioxynil	247	17.86	371	117	88	89	1689-83-4	3,94	95	99
Iprobenfos	190	15.35	91	204	123	122	26087-47-8	3,57	90	98

Isazophos	182	15.05	161	119	162	97	42509-80-8	3,26	81	96
Isobenzan	285	19.44	311	103	313	309	297-78-9	4,51	99	100
Isobornyl thiocynoacetate	186	15.24	121	93	95	136	115-31-1	3,75	93	99
Isodrin	303	20.03	193	195	263	66	465-73-6	6,75	100	100
Isofenphos	339	21.62	213	58	121	255	25311-71-1	4,65	99	100
Isomethiozin	281	19.35	225	198	184	199	57052-04-7	3,99	96	99
Isopropalin	314	20.70	280	238	281	264	33820-53-0	5,8	100	100
Isoxaben	525	30.80	165	68	96	163	82558-50-7	3,98	96	99
Isoxathion	412	24.98	105	77	177	313	18854-01-8	3,9	95	99
Jodfenphos	377	23.58	377	379	125	109	18181-70-9	5,39	100	100
Kinoprene	272	19.09	149	79	105	71	42588-37-4	6,45	100	100
Lenacil	453	26.88	153	154	110	136	2164-08-1	3,09	75	94
Leptophos	516	29.73	171	377	375	155	21609-90-5	6,34	100	100
Leptophos oxon	490	28.52	155	361	359	77	21609-90-5	6,34	100	100
Lindane	143	13.46	181	183	219	111	58-89-9	4,26	98	100
Mepronil	440	26.29	119	91	269	120	55814-41-0	4,24	98	100
Methfuroxam	287	19.47	137	229	138	67	28730-17-8	3,71	92	98
Methoprene I	318	20.90	73	111	81	109	40596-69-8	6,34	100	100
Methoprene II	357	22.33	73	111	109	107	40596-69-8	6,34	100	100
Methoxychlor	502	28.86	227	228	152	113	72-43-5	5,67	100	100
Methyldymron	333	21.37	107	146	106	77	42609-73-4	3,58	90	98
Metolachlor	271	18.91	162	238	240	146	51218-45-2	3,24	81	96
Mirex	517	29.84	272	274	270	237	2385-85-5	7,01	100	100
Monalide	191	15.47	85	127	197	239	7287-36-7	4,09	97	99
Myclobutanil	399	24.44	179	150	181	82	88671-89-0	3,5	88	98
Napropamide	374	23.44	72	128	271	100	15299-99-7	3,33	84	96
Nitrapyrin	43	7.93	194	196	198	160	1929-82-4	3,35	84	97
Nitrofen	410	24.86	283	285	202	253	1836-75-5	4,32	98	100
Nitrothal-isopropyl	296	19.87	236	194	212	254	10552-74-6	3,3	83	96
Nuarimol	466	27.45	235	107	203	139	63284-71-9	3,17	78	95
o,p'-DDD	396	24.35	235	237	165	199	53-19-0	5,87	100	100
o,p'-DDE	359	22.50	246	248	318	316	3424-82-6	6	100	100
o,p'-DDT	428	25.77	235	237	165	236	50-29-3	6,79	100	100
Octachlorostyrene	310	20.61	308	380	310	343	29082-74-4	7,46	100	100
o-Phenylphenol	51	8.78	170	169	141	115	90-43-7	3,28	82	96
Oxadiazon	400	24.42	175	177	258	260	19666-30-9	4,81	99	100
Oxychlorane	315	20.79	115	185	187	149	27304-13-8	5,48	100	100
Oxyfluorfen	407	24.73	252	302	331	361	42874-03-3	5,21	100	100
p,p'-DDD	425	25.69	235	237	165	236	72-54-8	5,87	100	100
p,p'-DDE	391	24.02	246	318	316	248	72-55-9	6	100	100
p,p'-DDT	458	26.98	235	237	165	236	50-29-3	6,79	100	100
Paclobutrazol	361	22.55	236	125	238	167	76738-62-0	3,36	85	97
Parathion	280	19.27	291	109	97	139	56-38-2	3,73	93	99
p-Dichlorobenzene	3	3.88	146	148	111	75	106-46-7	3,28	82	96
Pebulate	45	7.97	128	57	72	161	1114-71-2	3,51	89	98
Penconazole	325	21.03	248	159	161	250	66246-88-6	4,67	99	100
Pendimethalin	324	20.99	252	253	281	162	40487-42-1	4,82	99	100
Pentachloroaniline	193	15.54	265	267	263	269	527-20-8	4,3	98	100
Pentachloroanisole	121	12.58	280	265	237	282	1825-21-4	5,3	100	100
Pentachlorobenzene	55	8.95	250	252	248	254	608-93-5	5,22	100	100

Pentachloronitrobenzene	151	13.68	237	249	295	214	82-68-8	5,03	100	100
Pentachlorophenol	139	13.26	266	268	264	165	87-86-5	4,74	99	100
Pentanochlor	255	18.23	141	71	143	239	2307-68-8	4,18	97	99
Permethrin I	530	31.37	183	163	165	184	52645-53-1	7,43	100	100
Permethrin II	532	31.55	183	163	165	184	52645-53-1	7,43	100	100
Perthane	415	25.15	223	224	165	178	72-56-0	6,66	100	100
Phenamiphos	380	23.59	303	154	288	217	22224-92-6	3,29	82	96
Phenkapton	503	28.89	121	153	97	191	2275-14-1	5,84	100	100
Phenthoate	343	21.71	274	125	121	93	2597-03-7	3,47	88	97
Phorate	111	11.96	75	121	260	97	298-02-2	3,37	85	97
Phosalone	515	29.68	182	121	184	367	2310-17-0	4,29	98	100
Piperalin	463	27.24	112	314	316	173	3478-94-2	5,08	100	100
Piperonyl butoxide	475	27.90	176	177	149	178	51-03-6	4,29	98	100
Piperophos	501	28.88	122	140	320	97	24151-93-7	4,23	98	100
Pirimiphos-ethyl	311	20.65	333	318	304	168	23505-41-1	4,42	98	100
Pirimiphos-methyl	257	18.31	290	276	305	233	29232-93-7	3,44	87	97
Plifenat	220	16.64	43	217	175	219	21757-82-4	5,14	100	100
Pretilachlor	394	24.14	162	238	176	202	51218-49-6	4,29	98	100
Prochloraz	535	31.73	180	70	308	310	67747-09-5	4,13	97	99
Profenofos	389	23.90	208	139	206	339	41198-08-7	4,82	99	100
Profluralin	162	14.04	55	318	330	264	26399-36-0	5,62	100	100
Promecarb	109	11.92	135	150	91	136	2631-37-0	3,18	78	95
Prometon	130	12.99	210	225	168	183	1610-18-0	3,57	90	98
Prometryn	240	17.34	241	184	226	105	7287-19-6	3,73	93	99
Propargite	471	27.73	135	150	231	64	2312-35-8	5,57	100	100
Propazine	141	13.36	214	229	172	58	139-40-2	3,24	81	96
Propiconazole-I	454	26.94	173	69	259	175	60207-90-1	4,13	97	99
Propiconazole-II	462	27.15	173	69	259	175	60207-90-1	4,13	97	99
Propyzamide	159	13.95	173	175	145	255	23950-58-5	3,57	90	98
Prothiofos	383	23.75	309	267	162	113	34643-46-4	5,69	100	100
Pyrazophos	523	30.72	221	232	237	373	13457-18-6	3,53	89	98
Pyrazoxyfen	565	35.52	105	77	173	51	71561-11-0	4,97	100	100
Pyributicarb	483	28.35	165	108	181	93	88678-67-5	5,34	100	100
Pyridaben	531	31.52	147	117	148	132	96489-71-3	5,47	100	100
Pyridaphenthion	489	28.52	340	77	199	97	119-12-0	3,66	92	98
Pyridate	556	33.97	207	205	71	206	55512-33-9	5,73	100	100
Pyrifenox I	332	21.21	171	173	262	100	88283-41-4	4,2	97	100
Pyrifenox II	363	22.62	171	173	262	92	88283-41-4	4,2	97	100
Pyrimethanil	164	14.13	198	199	200	77	53112-28-0	3,19	79	95
Quizalofop-ethyl	545	32.91	299	372	163	301	76578-14-8	4,35	98	100
Resmethrin	476	27.98	123	171	128	143	10453-86-8	7,11	100	100
S,S,S-Tributylphosphorotrithioate (tribufos)	393	24.10	169	202	170	147	78-48-8	5,75	100	100
Sebuthylazine	187	15.26	200	202	214	229	7286-69-3	3,31	83	96
Secbumeton	177	14.74	196	169	210	225	26259-45-0	3,64	91	98
Sulfotep	106	11.83	322	202	97	238	3689-24-5	3,98	96	99
Sulprofos	441	26.34	322	156	140	139	35400-43-2	5,65	100	100
Tamoxifen	536	31.75	58	72	371	253	10540-29-1	6,3	100	100
TCMTB	368	23.10	180	181	136	238	21564-17-0	3,12	76	94
Tebuconazole	465	27.43	125	250	70	83	107534-96-3	3,89	95	99

Tebutam	107	11.83	91	57	190	233	35256-85-0	3,49	88	97
Tecnazene	73	10.24	203	215	261	201	117-18-0	4,39	98	100
Temephos	567	40.74	466	125	93	109	3383-96-8	6,17	100	100
Terbucarb	221	16.69	205	220	206	57	1918-11-2	5,28	100	100
Terbufos	153	13.80	231	57	103	153	13071-79-9	4,24	98	100
Terbumeton	144	13.44	210	169	225	154	33693-04-8	3,6	91	98
Terbutylazine	155	13.81	214	173	216	229	5915-41-3	3,27	82	96
Terbutryne	249	17.97	226	185	241	170	886-50-0	3,77	93	99
Tetrachlorvinphos	366	22.95	329	331	109	333	22248-79-9	3,81	94	99
Tetradifon	509	29.38	159	111	229	227	116-29-0	5,18	100	100
Tetramethrin I	493	28.63	164	123	79	107	7696-12-0	5,54	100	100
Tetramethrin II	499	28.86	164	123	165	79	7696-12-0	5,54	100	100
Tetrapropyl thiodiphosphate	277	19.21	211	210	253	115	3244-90-4	5,94	100	100
Tetrasul	437	26.07	252	254	324	322	2227-13-6	6,87	100	100
Thenylchlor	467	27.49	127	288	141	287	96491-05-3	4,21	97	100
Tolclofos-methyl	226	16.81	265	267	125	266	57018-04-9	4,77	99	100
Tolyfluanid	331	21.25	137	238	106	63	731-27-1	3,27	82	96
Tri-allate	181	14.97	86	268	270	128	2303-17-5	4,57	99	100
Tributyl phosphorotrithioite	328	21.18	209	153	298	97	150-50-5	7,67	100	100
Trichloronate	295	19.84	109	297	299	271	327-98-0	5,86	100	100
Tridiphane	235	17.25	173	187	189	175	58138-08-2	5,18	100	100
Trietazine	156	13.86	200	229	214	186	1912-26-1	3,44	87	97
Trifluralin	103	11.64	306	264	307	290	1582-09-8	5,31	100	100
Tryclopyrbutoxyethyl	409	24.82	57	56	85	210	-	3,94	95	99
Vernolate	40	7.81	128	43	86	146	1929-77-7	3,51	89	98
Uniconizole-P	390	24.02	234	236	70	235	83657-22-1	3,28	82	96

Information, descriptions and specifications in this publication are subject to change without notice.

