



# Halogenated hydrocarbons, aromatics

Separation of volatile halogenated hydrocarbons and some aromatics on a fused silica capillary column

## Application Note

Environmental

### Authors

Agilent Technologies, Inc.

### Introduction

Gas chromatography using an Agilent CP-Wax 57 CB column separates ten volatile halogenated hydrocarbons and some aromatics in 20 minutes.




**Agilent Technologies**

## Conditions

Technique : GC-capillary  
Column : Agilent CP-Wax 57 CB, 0.32 mm fused silica WCOT  
CP-Wax 57 CB (1.2  $\mu\text{m}$ ) (Part no. CP97773)  
Temperature : 60 °C  $\rightarrow$  120 °C, 1 °C/min  
Carrier Gas :  $\text{N}_2$ , 200 kPa (2.0 bar, 29 psi), 16 cm/s  
Injector : on-column  
T = 200 °C  
Detector : FID,  $32 \times 10^{-12}$  Afs  
T = 200 °C  
Sample Size : 30  $\mu\text{L}$  headspace

## Peak identification

1. tetrachloromethane
2. dichloromethane
3. benzene
4. trichloroethene
5. trichloromethane
6. toluene
7. ethylbenzene
8. p-xylene
9. m-xylene
10. o-xylene



[www.agilent.com/chem](http://www.agilent.com/chem)

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A00068



**Agilent Technologies**