



Sulfur compounds

Analysis of sulfur compounds in natural gas

Application Note

Energy & Fuels

Authors

Agilent Technologies, Inc.

Introduction

Gas chromatography using an Agilent CP-Sil 5 CB UltiMetal column separates sulfur compounds in natural gas in eight minutes.




Agilent Technologies

Conditions

Technique : GC-wide-bore
Column : Agilent CP-Sil 5 CB Ultimetal, 0.53 mm x 50 m,
(df = 5 μ m) (Part no. CP6670)
Temperature : 35 $^{\circ}$ C
Carrier Gas : He, 150 kPa (1.5 bar, 21 psi)
Injector : Valve
Detector : Sievers SCD 355
Sample Size : 10 μ L
Concentration Range : \pm 6 ppm/component

Peak identification

1. hydrogen sulfide
2. carbonyl sulfide
3. methanethiol (methyl mercaptan)
4. ethanethiol (ethyl mercaptan)
5. carbon disulfide



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01296



Agilent Technologies