



Drugs

Analysis of drugs of abuse (underivatized) standard mixture

Application Note

Forensic Toxicology

Authors

Agilent Technologies, Inc.

Introduction

Gas chromatography using an Agilent CP-Sil 8 CB column separates 18 underivatized drugs of abuse in a standard mixture in 15 minutes.

As a retention gap, a high temperature stable, thin film coated piece of a nonpolar fused silica column was used. This resulted in a better peak shape for the basic compounds and a longer lifetime of the precolumn under these injection conditions.




Agilent Technologies

Conditions

Technique : GC-capillary
Column : Agilent CP-Sil 8 CB, 0.32 mm x 25 m, fused silica
WCOT CP-Sil 8 CB (df = 0.25 µm) (Part no. 7452)
Precolumn : Agilent CP-SimDist, 0.53 mm x 2 m, fused silica
WCOT (df = 0.1 µm) (Part no. CP7541)
(for 10 m column)
Temperature : 75 °C (1 min) → 200 °C, 20 °C/min;
200 °C → 280 °C, 15 °C/min; 280 °C (3 min)
Carrier Gas : He, 80 kPa (0.8 bar, 11 psi)
Injector : on-column, T = 75 °C
Detector : NPD, T = 300 °C
Sample Size : 1.0 µL
Concentration Range : 0-10 mg/L
Solvent Sample : hexane

Peak identification

1. α -phenylethylamine (a-MBMA)
2. amphetamine
3. α,α -dimethylphenethylamine
(phentermine)
4. methamphetamine
5. ephedrine
6. MDA (a-methyl-3,4-methylene-
dioxyphenethylamine)
7. MDMA (3,4-methylenedioxyamphetamine)
8. MDEA (3,4-methylenedioxyethylamphetamine)
9. MBDB (2-methylamino-1-(3,4-methylene-
dioxyphenyl)butane)
10. caffeine
11. lidocaine
12. chirald (internal standard)
13. methadone
14. cocaine
15. codeine
16. morphine
17. artefacts
18. heroin



www.agilent.com/chem

For Forensic Use only.

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01390



Agilent Technologies