



## **Sulfur simdist**

# Boiling point analysis of sulfur compounds in hydrogenated naphtha

## Application Note

Energy & Fuels

### **Authors**

Agilent Technologies, Inc.

### **Introduction**

Boiling point analysis of sulfur compounds in hydrogenated naphtha by gas chromatography is achieved using an Agilent CP-Sil 5 CB column in ten minutes.




**Agilent Technologies**

## Conditions

Technique : GC-capillary  
Column : Agilent CP-Sil 5 CB UltiMetal, 0.53 mm x 10 m  
(df = 5  $\mu$ m) (Part no. CP6666)  
Temperature : 35 °C (1 min)  $\rightarrow$  15 °C/min  $\rightarrow$  200 °C (2 min)  
Carrier Gas : Helium, 15 mL/min, constant flow  
Injector : PTV, 70 °C Hold 0 Rate 15 °C/min Final 200 °C  
Detector : SCLD  
Sample Size : 0.5  $\mu$ L  
Calculations : Envantage Simdist- 2000 with Sulfur module  
Total sulfur = 240 ppm

Ppm-Sulfur w/v

10



10 min

[www.agilent.com/chem](http://www.agilent.com/chem)

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01812



**Agilent Technologies**