



Steroids as HFBA esters

Determination of steroids in urine

Application Note

BioPharma

Authors

Agilent Technologies, Inc.

Introduction

The selectivity of the Agilent VF-17ms allows a separation between the interference and the methylandrostandiol. Deconvolution via MS is not possible because of similar mass.




Agilent Technologies

Conditions

Technique : GC/MS
Column : Agilent VF-17ms, 0.25 mm x 30 m fused silica (film thickness = 0.25 μ m) (Part no. CP8982)
Temperature : 100 °C, (1 min) \rightarrow 175 °C, 15 °C/min \rightarrow 250 °C, 2.5 °C/min, \rightarrow 290 °C, 25 °C/min
Carrier Gas : He, 1.0 mL/min
Injection Technique : Splitless, Initial time : 1 min
Injection Temperature : 250 °C
Injection Volume : 1 μ L
Detector : MS Transferline: 300 °C Source: 250 °C
SIM used on specific masses
Concentration : Standard, 1 μ g/L
Sample Preparation : 2 mL Urine, decoupled; SPE (C18 and Amino); HPLC pre-purification via diol;
Derivatization : with heptafluor butyric acid anhydride
Courtesy : Drs. R. Schilt and P. Boshuis, TNO Quality of Life, Dep. Analytical Services, The Netherlands

Peak identification

1. fluoxymesterone
2. methylandrostandiol
3. interference = mass separated
4. α -nortestosterone
5. β -boldenone
6. β -nortestosterone
7. α -boldenone
8. β -estradiol
9. ethynylestradiol



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A02390



Agilent Technologies