

Now sold under the  
Thermo Scientific brand

**Thermo**  
SCIENTIFIC

## Anions and Organic Acids in Wood Extracts

### **SUMMARY**

The aim of this application was to analyze inorganic anions (sulfite, sulfate, and thiosulfate) and anions of organic acids (formate and acetate) in waste liquor from wood extracts. The sample matrix and the contaminant (ligninsulfonates) were removed from the components of interest using an IonPac® NG1 column as an additional guard column. The periodic regeneration of the IonPac NG1 column after each sample was easily integrated into the analytical setup using an additional 10-port valve. With the configuration shown here, it was possible to have an IonPac NG1 column in line with the analytical column, while a second column was rinsed in counterflow direction with acetonitrile/water, followed by a conditioning step using sodium hydroxide (see appendix).

The anions and organic acids were successfully separated using an IonPac AS11-HC column with a hydroxide gradient.

### **INSTRUMENTATION**

#### **ICS-3000 Ion Chromatography System**

DP Dual Pump

EG Eluent Generator Module

DC Detector/Chromatography Module

ED Electrochemical Detector

AS Autosampler

Chromeleon® 6.8 Chromatography Data System

### **ANALYTICAL CONDITIONS**

Column:	IonPac NG1 (4 × 35 mm) IonPac AG11-HC (2 × 50 mm) + AS11-HC (2 × 250 mm)												
Eluent Source:	EGC II KOH cartridge												
Gradient:	<table><thead><tr><th>Time (min)</th><th>Concentration KOH (mMol/L)</th></tr></thead><tbody><tr><td>-7.0</td><td>1</td></tr><tr><td>9.0</td><td>1</td></tr><tr><td>17</td><td>15</td></tr><tr><td>25</td><td>20</td></tr><tr><td>37</td><td>70</td></tr></tbody></table>	Time (min)	Concentration KOH (mMol/L)	-7.0	1	9.0	1	17	15	25	20	37	70
Time (min)	Concentration KOH (mMol/L)												
-7.0	1												
9.0	1												
17	15												
25	20												
37	70												
Flow:	0.38 mL/min												
System Pressure:	2576 psi (17.76 MPa)												
Detection:	Suppressed Conductivity												
Background Signal:	< 0.8 µS/S												
Suppressor:	ASRS® 2 mm (Recycle Mode)												
Suppressor Current:	66 mA												
Sample Preparation:	Dilution in ultrapure water and then centrifugation (10 min at 14,500 rpm).												
Sample 1:	Dilution 1 to 1250 for the acetate, sulfate, sulfite, and thiosulfate; dilution 1 to 100 for the formate.												
Sample 2:	Dilution 1 to 1000 for the acetate, formate, sulfite, and sulfate; dilution 1 to 100 for the thiosulfate.												
Temperature:	30 °C												
Injection Volume:	10 µL												
Chemicals Used:	Ultrapure Water (18.2 MΩ-cm)												

## RESULTS




Figure 1. Chromatogram of sample 1 at a dilution of 1 to 1250.




Figure 3. Chromatogram of sample 2 at a dilution of 1 to 1000.




Figure 2. Magnification of the chromatogram of sample 1 at a dilution of 1 to 1250.




Figure 4. Magnification of the chromatogram of sample 2 at a dilution of 1 to 1250.

## APPENDIX

### IonPac NG1 Regeneration Conditions

Eluent: (A) 5 mM NaOH  
(B) Water/acetonitrile 20/80

Flow: 0.3 mL/min

Chemicals Used: Ultrapure Water (18.2 M $\Omega$ -cm),  
Sodium hydroxide (Fluka 72064,  
50-52% solution),  
Acetonitrile (VWR 100030)

Gradient:	Time (min)	% A	% B
	-7.0	100	0
	-2.0	0	100
	15	0	100
	17	100	0
	37	100	0




Figure 5. IonPac NG1 regeneration.

ASRS, Chromeleon, and IonPac are registered trademarks of Dionex Corporation.

Passion. Power. Productivity.



#### Dionex Corporation

1228 Titan Way  
P.O. Box 3603  
Sunnyvale, CA  
94088-3603  
(408) 737-0700

#### North America

U.S./Canada (847) 295-7500

#### South America

Brazil (55) 11 3731 5140

#### Europe

Austria (43) 1 616 51 25 Benelux (31) 20 683 9768 (32) 3 353 4294  
Denmark (45) 36 36 90 90 France (33) 1 39 30 01 10 Germany (49) 6126 991 0  
Ireland (353) 1 644 0064 Italy (39) 02 51 62 1267 Sweden (46) 8 473 3380  
Switzerland (41) 62 205 9966 United Kingdom (44) 1276 691722

#### Asia Pacific

Australia (61) 2 9420 5233 China (852) 2428 3282 India (91) 22 2764 2735  
Japan (81) 6 6885 1213 Korea (82) 2 2653 2580 Singapore (65) 6289 1190  
Taiwan (886) 2 8751 6655

[www.dionex.com](http://www.dionex.com)

LPN 2275 PDF 08/16  
©2016 Dionex Corporation