

Comprehensive tap water analysis inclusively oxygenated halides with TitrIC 2

Branch: water, waste water, environmental protection

Keywords

TitriC 2 / 861 / 855 / 712 / 800 / Aquatrode plus / 6.0257.000 / Dosino / tiamo / 6.1006.630 / Metrosep SUPP 7 – 250 / 6.1010.230 / Metrosep C 2 – 250 / 2.145.0320 / 853 / MCS / branch 2

Summary

The determination of anions / cations in tap water was done on the **Metrosep SUPP 7 – 250 / Metrosep C 2 – 250** using carbonate / tartaric acid – dipicolinic acid eluent with chemical suppression. **TitriC 2** enables a fully automated and parallel analysis of anions, cations, titration (p- and m-value), and direct measurements (temperature, pH, conductivity). In this application the determination of oxygenated halides is shown as well. For a faster analysis of the standard anions please refer to **AB 302 e**.

Sample

Sample 1: Tap water Herisau

Sample 2: Tap water Rondelle

Reagents

Eluent: **3.6 mmol/L sodium carbonate** in ultra pure water (resistivity >18 MΩ)

7 mM tartaric acid, 1.5 mM dipicolinic acid in ultra pure water

Suppressor solutions:
50 mmol/L sulfuric acid
ultra pure water

Titration: 0.1 mol/L HCl for p-/m-value determination
1 mol/L HNO₃ for acidification

Standards [ppm] in ultra pure water

Anions:

| Analyte / Level | 1 | 2 | 3 | 4 | 5 |
|-----------------|-------|------|------|------|-----|
| fluoride | 0.03 | 0.1 | 0.2 | 0.4 | 0.6 |
| chlorite | 0.005 | 0.01 | 0.02 | 0.05 | 0.1 |
| bromate | 0.005 | 0.01 | 0.02 | 0.05 | 0.1 |
| chloride | 3 | 6 | 10 | 20 | 30 |
| nitrite | 0.03 | 0.1 | 0.2 | 0.4 | 0.6 |
| bromide | 0.03 | 0.1 | 0.2 | 0.4 | 0.6 |
| chlorate | 0.005 | 0.01 | 0.02 | 0.05 | 0.1 |
| nitrate | 3 | 6 | 10 | 20 | 30 |
| phosphate | 0.06 | 0.2 | 0.4 | 0.8 | 1.2 |
| sulfate | 3 | 6 | 10 | 20 | 30 |

Cations:

| Analyte / Level | 1 | 2 | 3 | 4 | 5 |
|-----------------|------|-----|----|------|-----|
| sodium | 1 | 2 | 4 | 6 | 8 |
| potassium | 0.5 | 1 | 2 | 3 | 4 |
| calcium | 12.5 | 25 | 50 | 75 | 100 |
| magnesium | 3.75 | 7.5 | 15 | 22.5 | 30 |

Apparatus and Accessories

- delivered with TitriC 2 package:

| | |
|-------------------------------------|------------|
| 861 Adv. Compact IC with seq. supp. | 2.861.0020 |
| 861 Adv. Compact IC | 2.861.0010 |
| 855 Robotic Titrosampler | 2.855.0020 |
| 712 Conductometer | 2.712.0010 |
| Cond. measuring cell Pt 1000 | 6.0912.110 |
| Conductivity standard | 6.2301.060 |
| Aquatrode plus | 6.0257.000 |
| 3 Dosions | 2.800.0010 |
| tiamo 1.1 full | 6.6056.112 |
| USB Converter Edgeport/4 | 2.145.0320 |
| 802 Rod Stirrer | 2.802.0020 |

- used optional accessories for TitriC 2:

| | |
|------------------------------------|------------|
| Metrosep A SUPP 7 – 250 | 6.1006.630 |
| Metrosep RP Guard | 6.1011.020 |
| Metrosep C 2 – 250 | 6.1010.230 |
| Metrosep C 2 Guard | 6.1010.200 |
| Sample rack 59 x 120 mL | 6.2041.840 |
| Sample beaker 120 mL (x100) | 6.1459.300 |
| Column oven for 861 | 2.861.0500 |
| 853 MCS CO ₂ Suppressor | 2.853.0010 |



Sample Preparation

The water samples were injected directly.

Analysis

The standards and the sample were injected two 800 dosinos (10 and 50 mL Dosing unit) with automatic acidification for the cation run prior to injection. Loop volume anions: **20 µl**, cations: **10 µl**. The two (with the TitriC 2 package delivered) *tiamo* methods *TitriC 2 – complete run* and *TitriC 2 – IC calibration anion & cation* were used for all types of analysis mentioned.

Calculation (IC)

Automatic integration with IC Net 2.3 software using peak area.

Parameters

IC - Anions

STARTUP HARDWARE:

```

RECORDER METHOD ASupp7_250-TitrIC.mtw
RECORDER DATA Data acquisition [Cond]
861 Adv. Compact Unit version 4
861 Adv. Compact Polarity +
861 Adv. Compact Supp. autostep yes
861 Adv. Compact Autostep with Fill
861 Adv. Compact Peristaltic pump yes
861 Adv. Compact Flow 0.70 mL/min
861 Adv. Compact Pmax 15.0 MPa
861 Adv. Compact Pmin 0.0 MPa
861 Adv. Compact FullScale 50 uS/cm
861 Adv. Compact Remote 00001000
Column Thermostat Temperature 45.0 °C
    
```

START WITH DETERMINATION

```

0.0 Column Thermostat START
0.20 861 Adv. Compact Valve Inject
28.0 861 Adv. Compact Valve Fill
    
```

START WITH INJECT:

```

0.0 RECORDER START
    
```

SET pH 4,3(HCl) – Titration (m-value)

Control parameters

```

EP1 at pH 4.3
Titration rate user
Dynamics pH 1.00
Max. rate 10.00 mL/min
Min rate 15.00 µL/min
Stop criterion drift
Stop drift 30 µL/min
End point 2 off
    
```

SET pH 8,2 (HCl) – Titration (p-/m-value)

Control parameters:

```

EP1 at pH 8.2
Titration rate user
Dynamics pH 1.00
Max. rate 10.00 mL/min
Min rate 15.00 µL/min
Stop criterion drift
Stop drift 30 µL/min
End point 2 on
EP2 at pH 4.3
Dynamics pH 1.00
Max. rate 10.00 mL/min
Min rate 15.00 µL/min
Stop criterion drift
Stop drift 30 µL/min
    
```

IC - Cations

STARTUP HARDWARE:

```

RECORDER METHOD C2_250-TitrIC.mtw
RECORDER DATA Data acquisition [Cond]
861 Adv. Compact Unit version 1
861 Adv. Compact Polarity -
861 Adv. Compact Supp. autostep no
861 Adv. Compact Autostep with Fill
861 Adv. Compact Flow 1.0 mL/min
861 Adv. Compact Pmax 25.0 MPa
861 Adv. Compact Pmin 0.0 MPa
861 Adv. Compact FullScale 1000uS/cm
861 Adv. Compact Remote 00000000
    
```

START WITH DETERMINATION

```

0.20 861 Adv. Compact Valve Inject
25.2 861 Adv. Compact Valve Fill
    
```

START WITH INJECT:

```

0.0 RECORDER START
    
```

MEAS pH

Measuring parameters:

```

Measurement with drift control on
Signal drift 10.0 mV/min
Min. waiting time 15 s
Max. waiting time 52 s
Measuring interval 2.0 s
Stop measured value pH off
Measurement without drift control off
Temperature 25.0 °C
    
```

MEAS conductivity

Measuring parameters:

```

Measurement frequency auto
Measuring time 35 s
Measuring interval 2.0 s
Stop measured value off mS/cm
Temperature 20.0 °C
    
```

Results

TitriC-Report



| | |
|-----------------------|----------------------------|
| Print date | 4/11/2006 |
| General data | |
| ID | Tap water Herisau |
| Record date | 4/05/2006 |
| Record time | 14:35:14 |
| User comment | |
| System comment | recalculated, within limit |
| Titration data | |
| pH | 7.88 |
| Cond. [μ S/cm] | 588 |
| Temp. [$^{\circ}$ C] | 22.79 |
| m value [mmol/L] | 5.59 |
| p value [mmol/L] | 0 |

| | |
|-----------------------------|--------|
| Anions | |
| F [mg/L] | 0.056 |
| ClO ₂ [mg/L] | 0 |
| BrO ₃ [mg/L] | 0 |
| Cl [mg/L] | 13.11 |
| NO ₂ [mg/L] | 0 |
| Br [mg/L] | 0.004 |
| ClO ₃ [mg/L] | 0.004 |
| NO ₃ [mg/L] | 10.136 |
| PO ₄ [mg/L] | 0 |
| SO ₄ [mg/L] | 5.271 |
| Anion balance [mEq/L] | 6.07 |
| Cations | |
| Na [mg/L] | 6.21 |
| K [mg/L] | 0 |
| Mg [mg/L] | 17.49 |
| Ca [mg/L] | 86.502 |
| Cation balance [mEq/L] | 6.03 |
| ionic difference [mEq/L] | -0.04 |
| ionic difference [%] | -0.35 |

TitriC-Report: Sample 1

TitriC-Report



| | |
|-----------------------|----------------------------|
| Print date | 4/11/2006 |
| General data | |
| ID | Quelle Rondelle |
| Record date | 4/05/2006 |
| Record time | 14:36:17 |
| User comment | |
| System comment | recalculated, within limit |
| Titration data | |
| pH | 7.82 |
| Cond. [μ S/cm] | 469 |
| Temp. [$^{\circ}$ C] | 24.01 |
| m value [mmol/L] | 4.12 |
| p value [mmol/L] | 0 |

| | |
|-----------------------------|--------|
| Anions | |
| F [mg/L] | 0.054 |
| ClO ₂ [mg/L] | 0 |
| BrO ₃ [mg/L] | 0 |
| Cl [mg/L] | 14.235 |
| NO ₂ [mg/L] | 0 |
| Br [mg/L] | 0.005 |
| ClO ₃ [mg/L] | 0 |
| NO ₃ [mg/L] | 9.895 |
| PO ₄ [mg/L] | 0.009 |
| SO ₄ [mg/L] | 5.694 |
| Anion balance [mEq/L] | 4.66 |
| Cations | |
| Na [mg/L] | 6.559 |
| K [mg/L] | 1.292 |
| Mg [mg/L] | 13.6 |
| Ca [mg/L] | 65.188 |
| Cation balance [mEq/L] | 4.69 |
| ionic difference [mEq/L] | 0.03 |
| ionic difference [%] | 0.31 |

TitriC-Report: Sample 2

Sample 1 - Tap water Herisau:

Reproducibility: Measurement of 10 tap water samples:

IC:

| anions | fluoride | chlorite | bromate | chloride | nitrite | bromide | chlorate | nitrate | phosphate | sulfate |
|--------------------|--------------|--------------|--------------|---------------|--------------|--------------|--------------|---------------|--------------|--------------|
| Average (N=10) | 0.060 | 0.000 | 0.000 | 13.066 | 0.000 | 0.005 | 0.009 | 10.149 | 0.005 | 5.285 |
| Standard Deviation | 0.003 | 0.000 | 0.000 | 0.054 | 0.000 | 0.000 | 0.001 | 0.050 | 0.001 | 0.041 |
| % RSD (N=10) | 4.43 | - | - | 0.41 | - | 1.11 | 3.87 | 0.49 | 15.31 | 0.78 |

| cations | sodium | potassium | calcium | magnesium |
|--------------------|--------|-----------|---------|-----------|
| Average (N=10) | 6.228 | 1.729 | 84.911 | 17.424 |
| Standard Deviation | 0.104 | 0.098 | 0.896 | 0.260 |
| % RSD (N=10) | 1.67 | 5.72 | 1.06 | 1.51 |

Titration / ionic balance:

| | pH | Cond. [µS/cm] | Temp. [°C] | m value [mmol/L] | p value [mmol/L] | sum anions [meq/L] | sum cations [meq/L] | ionic difference [%] |
|--------------------|--------------|----------------|---------------|------------------|------------------|--------------------|---------------------|----------------------|
| Average (N=10) | 7.772 | 585.900 | 22.358 | 5.572 | 0.000 | 6.002 | 5.992 | -0.082 |
| Standard Deviation | 0.075 | 2.601 | 0.235 | 0.020 | 0.000 | 0.057 | 0.048 | - |
| % RSD (N=10) | 0.97 | 0.44 | 1.05 | 0.37 | - | 0.95 | 0.80 | - |

Recovery: Spiked sample Tap water Herisau:

| Anions [ppm] | fluoride | chlorite | bromate | chloride | nitrite | bromide | chlorate | nitrate | phosphate | sulfate |
|-------------------|-------------|--------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|
| sample (av. N=10) | 0.060 | 0.000 | 0.000 | 13.066 | 0.000 | 0.005 | 0.009 | 10.149 | 0.005 | 5.285 |
| spiked (av. N=4) | 0.084 | 0.006 | 0.005 | 15.917 | 0.023 | 0.029 | 0.014 | 13.044 | 0.050 | 8.088 |
| % RSD spiked | 1.39 | 5.64 | 5.08 | 0.29 | 1.02 | 1.73 | 4.85 | 0.11 | 2.66 | 0.12 |
| difference | 0.023 | 0.006 | 0.005 | 2.851 | 0.023 | 0.024 | 0.005 | 2.895 | 0.045 | 2.803 |
| theoretical spike | 0.025 | 0.005 | 0.005 | 3 | 0.025 | 0.025 | 0.005 | 3 | 0.05 | 3 |
| recovery % | 92.7 | 113.5 | 99.0 | 95.0 | 92.3 | 96.1 | 92.5 | 96.5 | 89.8 | 93.4 |

| Cations [ppm] | sodium | potassium | calcium | magnesium |
|-------------------|--------------|--------------|-------------|-------------|
| sample (av. N=10) | 6.228 | 1.729 | 84.911 | 17.424 |
| spiked (av. N=4) | 7.229 | 2.242 | 96.101 | 20.740 |
| difference | 1.001 | 0.513 | 11.190 | 3.316 |
| theoretical spike | 1 | 0.5 | 12.5 | 3.75 |
| recovery % | 100.1 | 102.6 | 89.5 | 88.4 |

Sample 2 - Tap water Rondelle:

IC:

| Anions [ppm] | fluoride | chlorite | bromate | chloride | nitrite | bromide | chlorate | nitrate | phosphate | sulfate |
|----------------|--------------|--------------|--------------|---------------|--------------|--------------|--------------|--------------|--------------|--------------|
| Average (N=10) | 0.054 | 0.000 | 0.000 | 14.235 | 0.000 | 0.005 | 0.000 | 9.895 | 0.009 | 5.694 |

| Cations [ppm] | sodium | potassium | calcium | magnesium |
|----------------|--------|-----------|---------|-----------|
| Average (N=10) | 6.559 | 1.292 | 65.188 | 13.600 |

Titration / ionic balance:

| | pH | Cond. [µS/cm] | Temp. [°C] | m value [mmol/L] | p value [mmol/L] | sum anions [meq/L] | sum cations [meq/L] | ionic difference [%] |
|----------------|-------------|---------------|--------------|------------------|------------------|--------------------|---------------------|----------------------|
| Average (N=10) | 7.82 | 469 | 24.01 | 4.12 | 0.00 | 4.66 | 4.69 | 0.31 |

Comments

This Application Bulletin is based on AW CH6-0864-042006.

Overview of the determined parameters:

- direct measurement:
conductivity, pH, temperature
- titrations:
p- and m-value (= acidic capacity)
- IC:
fluoride, chlorite, bromate, chloride, nitrite,
chlorate, bromide, nitrate, sulfate
sodium, potassium, calcium, magnesium

TitriC is absolutely flexible and can be adapted to determine the desired variables.

The calibration results are good which means that the filling of a loop with dosinos is working properly.

No sample preparation is necessary – even the acidification for the cation-IC-run is done automatically by a Dosino.

The recovery as well as the reproducibility for phosphate in tap water is poor while the calibration was good. There seems to be species in the sample which lower the phosphate concentrations. This problem (probably related to the presence of humic acids in the sample) is known and also appears with direct injection.

With this setup the oxohalogenides can be determined down to 5ppb. This could be improved with implementing an Advanced Modular system.

The ionic balance shown in the TitrIC-Report is calculated by comparison of total equivalent of anions (the carbonate concentration is calculated from the p- and m-values) and cations from the IC results:

$$\text{Ionic d. [\%]} = 100 * \left(\frac{\text{sum cations [meq/L]} - \text{sum anions [meq/L]}}{\text{sum cations [meq/L]} + \text{sum anions [meq/L]}} \right)$$

Example: If 5 ppm sulfate is found, this corresponds to 5 mg/L → 0.052 mmol/L ($M(\text{SO}_4) = 96 \text{ mg/mmol}$) → -0.104 meq/L (charge SO_4^{2-}).

For a more detailed information about the system setup of TitrIC 2 consult the document **AB 286 - Installation Instructions for TitrIC 2**.

For Application works of TitrIC 1 please have a look at AB 288 e.

The Aquatrode should be stored in the storage solution of Metrohm (6.2323.000).

Appendix


Calibration curves, reproducibility, chromatograms with peak tables and titration curves

Date

Thursday, 2006-04-05

Name

A. Rumi; IC Marketing,
Metrohm Ltd.; CH-9101 Herisau




TitrIC 2: System setup

Calibration

Anions:


fluoride: RSD: 0.81%

$$Q = -0.00020 \cdot A^2 + 0.034 \cdot A + 0.0038$$




nitrite: RSD: 3.350 %, corr. coeff.: 0.99938

$$Q = 0.07824 \cdot A + 0.008926$$




chlorite: RSD: 1.306 %, corr. coeff.: 0.999943

$$Q = 0.1362A + 0.000226$$




bromide: RSD: 1.804 %, corr. coeff.: 0.99984

$$Q = 0.1419 \cdot A + 0.004460$$




bromate: RSD: 0.496 % corr. coeff.: 0.999992

$$Q = 0.2546 \cdot A - 0.0006363$$




chlorate: RSD: 3.211 %, corr. coeff.: 0.99966

$$Q = 0.143708 \cdot A + 0.000734$$




chloride: RSD: 0.05%

$$Q = -7.96e-07 \cdot A^2 + 0.0428 \cdot A + 0.209$$




nitrate: RSD: 0.09%


$$Q = -6.998e-06 \cdot A^2 + 0.078 \cdot A + 0.390$$



phosphate: RSD: 2.550 %, corr. coeff.: 0.99968
 $Q = 0.18538 \cdot A + 0.0105642$




sulfate: RSD: 0.09%
 $Q = -2.759e-06 \cdot A^2 + 0.060 \cdot A + 0.392$




Cations:


sodium: RSD: 2.38%
 $Q = 0.1434 \cdot A + 0.0674$




calcium: RSD: 1.30%
 $Q = 0.1903 \cdot A + 0.1929$




potassium: RSD: 2.77%
 $Q = 0.5325 \cdot A - 0.3066$



magnesium: RSD: 0.87%
 $Q = 0.0775 \cdot A + 0.355$




Chromatograms with peak tables - Anions



Standard 1

File: q3211205
Volume: 20.0 µl
Dilution: 1.00
Amount: 1.00


| No | Retention min | Height uS/cm | Area uS/cm*sec | Conc. mg/L | Name |
|----|---------------|--------------|----------------|------------|-----------|
| 1 | 6.63 | 0.08 | 0.801 | 0.0309 | fluoride |
| 5 | 8.99 | 0.00 | 0.033 | 0.0050 | chlorite |
| 6 | 9.48 | 0.00 | 0.022 | 0.0050 | bromate |
| 7 | 10.24 | 5.92 | 65.089 | 2.9945 | chloride |
| 9 | 12.48 | 0.03 | 0.341 | 0.0291 | nitrite |
| 10 | 15.69 | 0.01 | 0.199 | 0.0295 | bromide |
| 11 | 16.21 | 0.00 | 0.018 | 0.0057 | chlorate |
| 12 | 18.17 | 1.77 | 33.358 | 3.0096 | nitrate |
| 13 | 24.43 | 0.01 | 0.336 | 0.0669 | phosphate |
| 14 | 27.11 | 1.59 | 43.355 | 3.0100 | sulfate |



Standard 5


File: q3211105
Volume: 20.0 µl
Dilution: 1.00
Amount: 1.00

| No | Retention min | Height uS/cm | Area uS/cm*sec | Conc. mg/L | Name |
|----|---------------|--------------|----------------|------------|-----------|
| 1 | 6.63 | 2.08 | 19.906 | 0.6001 | fluoride |
| 4 | 8.98 | 0.07 | 0.730 | 0.1000 | chlorite |
| 5 | 9.49 | 0.03 | 0.395 | 0.1000 | bromate |
| 6 | 10.29 | 58.94 | 704.574 | 30.0012 | chloride |
| 7 | 12.48 | 0.58 | 7.637 | 0.6002 | nitrite |
| 8 | 15.67 | 0.27 | 4.221 | 0.6005 | bromide |
| 9 | 16.22 | 0.03 | 0.378 | 0.0999 | chlorate |
| 10 | 17.98 | 18.67 | 389.425 | 29.9971 | nitrate |
| 11 | 24.41 | 0.24 | 6.477 | 1.2010 | phosphate |
| 12 | 26.99 | 17.81 | 500.850 | 29.9971 | sulfate |



Sample 1
 File: q3220119
 Volume: 20.0 µl
 Dilution: 1.00
 Amount: 1.00


| No | Retention min | Height uS/cm | Area uS/cm*sec | Conc. mg/L | Name |
|----|---------------|--------------|----------------|------------|----------|
| 1 | 6.63 | 0.15 | 1.743 | 0.0625 | fluoride |
| 4 | 10.26 | 23.94 | 303.280 | 13.1299 | chloride |
| 7 | 15.73 | 0.00 | 0.033 | 0.0048 | bromide |
| 8 | 16.27 | 0.00 | 0.031 | 0.0099 | chlorate |
| 9 | 18.14 | 6.58 | 125.507 | 10.1635 | nitrate |
| 10 | 27.11 | 2.99 | 82.005 | 5.3346 | sulfate |



Sample 2
 File: q3221750
 Volume: 20.0 µl
 Dilution: 1.00
 Amount: 1.00


| No | Retention min | Height uS/cm | Area uS/cm*sec | Conc. mg/L | Name |
|----|---------------|--------------|----------------|------------|-----------|
| 1 | 6.62 | 0.13 | 1.339 | 0.0490 | fluoride |
| 2 | 10.25 | 27.45 | 329.370 | 14.2346 | chloride |
| 3 | 11.76 | 0.01 | 0.089 | 0.0000 | |
| 4 | 15.72 | 0.00 | 0.042 | 0.0054 | bromide |
| 5 | 18.12 | 6.39 | 122.016 | 9.8952 | nitrate |
| 6 | 24.39 | 0.00 | 0.079 | 0.0092 | phosphate |
| 7 | 27.09 | 3.23 | 88.351 | 5.7155 | sulfate |

Cations:



Sample 1
 File: q3220120
 Volume: 10.0 µl
 Dilution: 1.00
 Amount: 1.00

| No | Retention min | Height uS/cm | Area uS/cm*sec | Conc. mg/L | Name |
|----|---------------|--------------|----------------|------------|-----------|
| 1 | 6.56 | 3.31 | 42.731 | 6.1884 | sodium |
| 2 | 9.79 | 0.14 | 3.155 | 1.5719 | potassium |
| 3 | 10.92 | 9.97 | 450.169 | 85.2894 | calcium |
| 4 | 15.86 | 4.16 | 225.558 | 17.7705 | magnesium |



Sample 2
 File: q3221751
 Volume: 10.0 µl
 Dilution: 1.00
 Amount: 1.00

| No | Retention min | Height uS/cm | Area uS/cm*sec | Conc. mg/L | Name |
|----|---------------|--------------|----------------|------------|-----------|
| 1 | 6.53 | 3.73 | 45.260 | 6.5590 | sodium |
| 2 | 9.79 | 0.11 | 3.002 | 1.2923 | potassium |
| 3 | 11.10 | 8.24 | 341.400 | 65.1877 | calcium |
| 4 | 15.98 | 3.54 | 170.864 | 13.5999 | magnesium |

Titration curves – Titration of *m*-value (to pH 4.3)

