


# GCxGC-TOFMS Utilized as a Broad-Spectrum Analysis for Endocrine Disruptor Compounds in Urban and Rural Watersheds

Allen Zhang, John Heim, and Joe Binkley | LECO Corporation, St. Joseph, Michigan USA



Endocrine-disrupting compounds (EDCs) encompass a variety of chemical classes, including drugs, pesticides, polymer additives, coatings materials, personal consumer products, industrial by-products and pollutants. There is worldwide concern over long-term environmental exposure to EDCs leading to serious health effects including a range of reproductive problems such as reduced fertility, male and female reproductive abnormalities, skewed male/female sex ratios, brain and behavior problems, impaired immune functions, and various cancers. This research presents a robust, broad range analysis for the detection of EDCs in impacted natural waters using Comprehensive Two-Dimensional Gas Chromatography-Timeof-Flight Mass Spectrometry (GCxGC-TOFMS). GCxGC facilitates enhanced detection, chromatographic resolution, and peak capacity while TOFMS allows the fast acquisition (up to 500 spectra per second) necessary to successfully acquire the data density needed to fully characterize low levels of targeted and untargeted compounds in complex samples. A reference standard of 108 known endocrinedisrupting compounds was prepared and analyzed. Methods for solid-phase extraction and GCxGC-TOFMS analysis were developed. GCxGC-TOFMS analysis was conducted on multiple water samples from a rural and urban Midwestern U.S watershed. Extraction of 1.0 liter water samples was conducted using Supel-Select HLB SPE cartridges (Supelco Analytical, Sigma-Aldrich) designed to recover a wide range of analytes. Subsequent analysis was conducted by GCxGC-TOFMS and the data was processed using a 152 component reference standard. The results are reported as targeted and untargeted analytes found. This research study presents a practical, robust, sensitive, and reliable method for the detection of EDCs in urban

Introduction

# **Experimental Methods**

### **Reference Standard Preparation**

- An endocrine disruptor reference stock standard at  $1 \text{ ng}/\mu\text{L}$  was prepared in acetone with EPA Method standards purchased from Restek Corp.
- The commercial standards purchased were Method 8270 Megamix, Method 527 pesticide mix # 1, Method 551.1 pesticide/herbicide mix, plus Bisphenol A.

### **Solid-Phase Extraction Procedure**

- Adjust 1 liter water sample to pH 2 with 37% HCl.
- Condition SPE Supel™ Select HLB, 500 mg cartridge with 5 mL HPLC Water/5% Methanol, then 5 mL Acetone, followed by 5 mL HPLC water.
- Load 1 Liter of water using the Supelco Visiprep Vacuum Manifold (Supelco Analytical, Sigma –Aldrich) slowly into the SPE cartridge.
- Dry SPE tube with vacuum for approximately 15 minutes.
- Elute slowly, 3 mL of acetone/5% methanol into a 20 mL clean glass test tube
- Elute slowly, 3 mL of dichloromethane into the same 20 mL glass test tube.
- Speedvac to dryness for approximately 2 hours.
- Reconstitute dried residue in 500  $\mu$ L of acetone, vortex, and pipet into
- Inject 1 μL sample for GCxGC-TOFMS analysis.

## **Method Development**

- The solid-phase extraction (SPE) procedure was developed using 1 liter of HPLC water spiked with the EDC reference standard at 50ppb.
- The GCxGC-TOFMS method was developed using the EDC reference standard injected at 5 ng on-column.

## Methods

### **GCxGC-TOFMS** Analysis Parameters

- Gas Chromatograph: Agilent 7890 equipped with a LECO dual stage, quad jet thermal modulator, and a GERSTEL MPS2 autosampler
- GC Primary Column: 30 m x 0.25 mm id. x 0.25  $\mu$ m film thickness Rxi-5SilMS
- GC Secondary Column: 1.0 m x 0.18 mm id. x 0.18  $\mu$ m film thickness Rxi-17Sil-MS (Restek Corp.)
- Carrier Gas: Helium set @ 1.5 mL/min
- **Injection Mode: Splitless**
- Injection Volume: 1 μL
- Inlet Temperature: 250°C
- Primary Column Temperature Program: Initial temperature 75°C for 1.0 min
- Secondary Column Temperature Program: Initial temperature 80°C for 1.0 min ramped @ 6.0°C/min to 305°C held for 10 min
- GCxGC Modulator Temperature Offset: 20°C
- Transfer Line Temperature: 280°C
- Total Run Time: 48.50 min

### **GCxGC** Variable Modulation

Variable modulation shown in Table 1 below allows the user to change and modify optimize the separation, thereby maximizing the peak capacity in the

Modu	lation Timing:	For 1D GC set second dimension time to 0					
#	Start	End	Modulation Period (s)	Hot Pulse Time	Cool Time Between Stages		
1×	Start of Run	614 s	2.00	0.50	0.50		
2	614 s	1838 ន	3.00	0.60	0.90		
3	1838 ន	End of Run	4.00	0.70	1.30		

## Mass Spectrometer: Pegasus® 4D GCxGC-TOFMS **Analysis Parameters**

- Mass Range: 35 800 υ
- Acquisition Rate: 200 spectra/s
- Ion Source Temperature: 230°C
- Detector Voltage: 1700 V
- Acquisition Delay: 150 sec.

# Using The Reference Standard Feature

A Reference Method is built in ChromaTOF® software from a user-created standard that is applied and compared to a sample. The purpose of a Reference Method is to determine the component differences between a sample and a reference standard within user-defined limits of retention time, peak area, and spectral match. In this study, untargeted components found in the water extracts were added to the original 108-component reference standard. The final Reference used in data analysis contained 152 compounds. The Reference was applied as part of the data processing method. The processed sample peak table displays each compound from the reference in the Type column as either a "Match", found but "Out of Tolerance" by percent, or "Not Found".

Table 2. An example of a peak table is shown below for the EDC Reference used for this analysis. The columns show the Compound Name, Type of Match, the Match similarity score, retention time, peak area, unique mass, signal to noise, and the library used to identify the mass spectral peak.

Name	Туре	Match	R.T. (s)	Area	UniqueMass	S/N	Library
Bisphenol A	Match	895	1724 , 2.600	740588	213	2979.2	mainlib
Ionol 2 * (ANTIOXIDANT- TO PREVENT GUMMING IN FUELS)	Out of Tolerance	918	1085 , 1.450	733755	219	5824.1	mainlib
Lindane/ EDC STD 8S rt dev 600 SIM	Not Found	600					mainlib

## GCxGC-TOFMS Results





Figure 1. The two-dimensional contour plot chromatogram above shows the 108-component endocrine disruptor Reference standard used to develop the SPE extraction procedure and GCxGC-TOFMS method. The on-column concentration for each component is 5 nanograms. The red hatched unprocessed classification region represents the unwanted solvent and erroneous background excluded from

# Calibration Linearity & Targeted GCXGC-TOFMS Quantitation




**Figure 2.** The six-point calibration curves for 1.4-dichlorobenzene and 1-methylnaphthalene are shown above in Figure 2. Calibration linearity for both components above show calculated Pearson's r values greater than 0.999. Targeted analysis was utilized in the data processing method with calibration curves for components contained within the endocrine disruptor standard. The peak table above shows the quantitative results calculated for two selected analytes detected in a 1 liter water sample.

aphthalene, 1-methyl- (CAS) 0.06 ng 791, 1.830 732 7533.1 605.24 141 70.566 Wiley9

0.08 ng 396,1.345 941 169194 12113 146 364.46 Wiley9

# GCxGC Increased Peak Capacity & NIST Searchable TOFMS



47 4-Chlorophenyl phenyl ether

59 | Ibuprofen-M (HO-) -H2O P329

3,5-DITERT-BUTYLBENZALDEHYDE

4 Gemfibrozil \* Medication to lower lipid levels

63 Myristic acid P412

Caffeine (CAS)

79 | Heptachlor epoxide

2 Carbazole

4-Bromophenyl phenyl ether \*USED AS A (PAST) FLAME RETARDANT

7-Acetyl-6-ethyl-1,1,4,4-tetramethyltetralin \* synthetic musk

Benzenesulfonanilide \* production of synthetic dyes, photochemicals, disinfectants,

Figure 3. The example above shows the two-dimensional Contour Plot for two peaks that are resolved in the second dimension by 105 milliseconds vet coeluted completely in the 1<sup>st</sup> dimension separation at 1764 s. The expanded peak capacity and advantage of a two-dimensional separation is illustrated by this example of two chemically similar pesticides that would be coeluted and difficult to characterize by one-dimensional chromatography. Notice the high mass spectral match scores by the NIST 08 library searches of 84% and 90.5% for both compounds.

## **EDCs & Pollutants Detected**

water samples prepared with solid-phase extraction and analyzed by GCxGC-TOFMS are shown below. A total of 102 chemical compounds were detected from six Midwestern watershed. The compounds detected matched the Reference standard with at least a 60% library match similarity. Endocrine-disrupting compounds (EDCs) encompass a variety of chemical classes, including

UNTARGETED AND TARGETED ANALYTES DETECTED	1ST DIMENSION RETENTION TIME (s)	COUNT # OF TIMES DETECTED	CHEMICAL TYPE		UNTARGETED AND TARGETED ANALYTES DETECTED	1ST DIMENSION RETENTION TIME (s)	COUNT # OF TIMES DETECTED	CHEMICAL TYPE
Phenol (CAS)	354	23	EDC					EDC, POP, INSECTICIDE, bann
Bis(2-chloroethyl) ether	364	1	EDC	80	Oxychlordane	1613	2	1988
	331	_		81	Bioallethrin	1625	4	EDC, PESTICIDE
1,4-Dichlorobenzene	392	21	EDC	82	FLUORANTHENE	1628	23	EDC, PAH
1,3-Dichlorobenzene	394	9	EDC	83	Lauryl acrylate * for manufacturing polymers used in hairstyling	1640	6	PCP
1,2-DICHLOROBENZENE	418	3	EDC	84	Naproxen * Pharmaceutical	1655	1	PCP, PHARMACEUTICAL
Benzylalcohol P194	420	19	EDC	85	Triclosan * antibacterial, antifungal agent used in toothpaste	1661	7	PCP, PHARMACEUTICAL
2-Methylphenol	440	22	EDC	86	Pyrene	1679	23	EDC, PAH
3-Methylphenol	458	24	EDC	87	Butyl citrate * used as a plasticizer, antifoam agent	1721	5	EDC, INDUSTRIAL CHEMICA
Nitrobenzene	480	7	EDC	88	Bisphenol A * Used in making plastics	1721	13	EDC, INDUSTRIAL CHEMICA
Isophorone * Solvent	532	26	EDC, INDUSTRIAL CHEMICAL					EDC, PESTICIDE, banned in EU
2-Nitrophenol	548	11	EDC	89	Nitrofen	1769	2	US
2,4-Dimethylphenol	566	11	EDC	90	Endrin P1112	1772	2	EDC, INSECTICIDE, banned 2
2,4 Dichlorophenol	600	7	EDC	91	TCPP _ Tris(1,3-dichloroisopropyl)phosphate * (PBDE) flame retardant REPLACEMENT	1838	16	EDC, INDUSTRIAL CHEMICA
Benzene, 1,2,4-trichloro- (CAS)	608	9	EDC					EDC, INDUSTRIAL CHEMICA
Naphthalene	623	25	EDC, PAH	92	2H-1-Benzopyran-2-one, 7-(diethylamino)-4-methyl- * (Optical bleach in textile industry)	1854	4	ED & LIEDDIAIDE
D-Glucitol, 1,4:3,6-dianhydro-2,5-di-O-methyl- * used in Personal Care Products;			PCP	93	Hexazinone P513	1882	1	EDC, HERBICIDE
cosmetics	731	7	rer	0.4	Dif- which	4050	-	EDC, INSECTICIDE, class C
Naphthalene, 2-methyl-	770	20	EDC, PAH	94		1958	5	carcinogen
Naphthalene, 1-methyl-	791	22	EDC, PAH	95	Benzo[a]anthracene	1962	4	EDC, PAH
1,3-Isobenzofurandione (CAS) * RUBBER RETARDER, CURING AGENT	794	14	EDC, INDUSTRIAL CHEMICAL	96	Chrysene (CAS)	1970	8	EDC, PAH EDC, INSECTICIDE U.S. BANN
Propofol	842	6	PCP, PHARMACEUTICAL	97	Methoxychlor	1974	,	2003
2,4,6-Trichlorophenol	842	6	EDC	37	IMETHOXYCHIOI	1974		EDC, INSECTICIDE U.S. BANN
2,4,5-Trichlorophenol	845	6	EDC	98	Mirex	2062	2	1976
à DAMASCONE	863	12	FOOD, FLAVOR, FRAGRANCE	99	Benzo[b or k]fluoranthene	2198	15	EDC, PAH
2-Chloronaphthalene	878	7	EDC, PAH	100	Benzo[a]pyrene (CAS)	2262	3	EDC, PAH
			EDC	101	Fenvalerate isomer-1 P1241	2302	1	EDC, INSECTICIDE
Phenol, 4-chloro-3,5-dimethyl- (CAS)	884	6	-50		Fenvalerate isomer-2 P1242	2322	2	EDC, INSECTICIDE
			FOOD, FLAVOR, FRAGRANCE	101			_	220, 111020110122
á-Patchoulene	887	7						
Naphthalene, 1,7-dimethyl-	911	7	EDC, PAH					
1H-3a,7-Methanoazulene, 2,3,4,7,8,8a-hexahydro-3,6,8,8-tetramethyl-, [3R-			EDC		Conclusion	nc		
(3à,3aá,7á,8aà)]-	929	5			Conclusion			
Acenaphthylene (CAS)	971	14	EDC, PAH					
Acenaphthene	1010	21	EDC, PAH		conclusion, comprehensive two-dimensional gas chromat			
BUTYL HYDROXY TOLUENE * BHT ANTIOXIDANT	1025	23	EDC, FOOD, FLAVOR, FRAGRANCE		pectrometry was utilized in this research to detect endocrine			
Tributyl phosphate * SOLVENT & PLASTICIZER	1031	19	EDC, INDUSTRIAL CHEMICAL		amples obtained from urban and rural point sources along			
Benzeneacetic acid, ethyl ester * flavor fragrance	1037	25	EDC, FOOD, FLAVOR, FRAGRANCE		xtraction method was developed using a hydrophilic modified			
à-N-METHYL IONONE	1043	6	EDC, FOOD, FLAVOR, FRAGRANCE		f compounds from aqueous samples. An optimized GCxGC m			
TRIPROPYLENE GLYCOL 5	1046	21	EDC, INDUSTRIAL CHEMICAL		on-polar and mid-polarity column set. The GCxGC method			
Dibenzofuran	1052	14	EDC		ptimization of the available peak capacity and chromatogr			
Lilial	1055	12	PCP, FOOD, FLAVOR, FRAGRANCE		imension separation. A TOFMS method was created which			
1H-Benzotriazole, 4-methyl- * CORROSION INHIBITOR	1070	4	EDC, INDUSTRIAL CHEMICAL	m	nass spectral information, True Signal Deconvolution $^{ exttt{ iny e}}$ , c	ınd fast (	acquisition	rates ideal for tl
1H-Benzotriazole, 5-methyl- *RETROCURE G USED IN PREVULCANIZATION IN RUBBER MANUFATURE	1076	4	EDC, INDUSTRIAL CHEMICAL		haracterization of EDCs and other contaminants in water. hass spectral data as illustrated in Figure 3.	TOFMS ge	nerates N	IST library searchab
Ionol 2 * (ANTIOXIDANT- TO PREVENT GUMMING IN FUELS)	1085	11	EDC, INDUSTRIAL CHEMICAL	• • •	,			
2-tert-Butylhydroquinone * TBHQ, FOOD PRESERVATIVE ANTIOXIDANT	1088	2	FOOD, FLAVOR, FRAGRANCE	Th	hirteen solid-phase extractions were conducted on 1 liter water	er sample	s obtained	from 6 different run
Dodecanamide, N,N-bis(2-hydroxyethyl)- •FOAM STABILIZER IN HOUSEHOLD DETERGENTS AND SHAMPOOS	1091	22	PCP	рі	nd urban point sources along a Midwestern watershed. GCx rocessing that utilized the "Reference" feature in ChromaTOF	software.	The final R	eference developed
á N METHYL IONONE	1109	6	PCP, FOOD, FLAVOR, FRAGRANCE		etect EDCs contained 152 compounds which was applied			
TAIN WEITHE TONOINE					argeted analysis was conducted using a six point calibration			

1139 5 EDC, FLAME RETARDAN

1160 5 PCP, PHARMACEUTICAL

1214 6 EDC, HERBICIDE

1298 14 EDC, HERBICIDE, banned by EU

1313 24 FOOD, FLAVOR, FRAGRANCE

1334 6 PCP, PHARMACEUTICAL

1406 10 PCP, FOOD, FLAVOR, PHARMA

1415 24 PCP, FOOD, FLAVOR, FRAGRANC

1415 12 EDC, INDUSTRIAL CHEMICAL

1490 11 PCP, PHARMACEUTICAL

1349 26 EDC, INDUSTRIAL CHEMICA

EDC, PHARMACEUTICAL,

DC, INDUSTRIAL CHEMICAL

EDC, FUNGICIDE, POP, banned

EDC, HERBICIDE, banned EU, US st

PCP, PHARMACEUTICAL

EDC, INDUSTRIAL CHEMICAL

PCP, PHARMACEUTICAL EDC, POP, INSECTICIDE, US limits

10 PCP, FOOD, FLAVOR, FRAGRANC

Targeted analysis was conducted using a six point calibration developed in the 50 ppb to 25 ppm range for selected EDC standard compounds illustrated in Figure 2. Calibration linearity of greater than 99.9% and quantitative examples are shown for selected analytes. Results of this research detected 102 chemicals in aqueous extractions from 6 different sampling sites. The detected compounds matched the reference standard with at least a 60% library match similarity. Furthermore, results show that 81% of the 102 detected chemicals were found at least 5 times.

The research presented in this study emphasizes the need for instrumentation that will detect, quantify and characterize sources of long-term environmental exposure to EDCs that can lead to ecological destruction and serious health effects. The application of GCxGC-TOFMS for this work presents an excellent instrumental option for the detection of targeted and untargeted pollutants in impacted natural waterways. The data presented illustrates the advantages and benefits of GCxGC-TOFMS to provide a robust analysis as well as a data mining strategy using the "Reference" feature and quantitative calibration available in ChromaTOF software to detect a broad range of chemical contaminants in impacted natural waterways.

For further information regarding this study, contact the authors at john heim@leco.com.

## References

] "Water Analysis: Emerging Contaminants and Current Issues", Susan D. Richardson, Anal. Chem. 2007, 79, 4295-4324

[2] "Determining Pharmaceuticals and personal Care Products in Water by Automated Solid Phase Extraction", Michael Ebitson, Horizon Technology, Inc., AN061-101101

[3] "Broad Range Analysis of Endocrine Disruptors in Pharmaceuticals Usina Gas Chromatography and Liquid Chromatography Tandem Mass Spectrometry", Rebecca A. Trenholm, Brett J. Vanderford, Janie C. Holady, David J. Rexing, Shane A. Snyder, Chemosphere 65 (2006) 1990-1998