

Extracting More Analytes from Complex Samples using Agilent's GC/TQ and GC/Q-TOF

GC/TQ (Triple Quad):

Chin-Kai (Kai) Meng, PhD

Senior Application Chemist

GC/Q-TOF (Quad-Time-of-Flight):

Terry L. Sheehan, PhD

GC/MS Marketing Manager



GO.



 Agilent Technologies




Extracting More Analytes from Complex Samples using Agilent's Pesticide GC/TQ Analyzer

-- **“Ready to Go” Analyzer based
on the most comprehensive
MRM database**

Ready to use Fast Track Analyzers

... no on-site break-in needed, ready for injections

At Factory




Plug and Play

At User's Lab



Quick installation and start-up





Why Do We Need Backflush?

- With the selectivity of MS/MS, **users cannot see “dirty matrix”**
- MS/MS users want the LOWEST detection limits – **inject more** with large volume injection
- Many late eluting peaks are not “chromatographically ideal” and **leave a residue** throughout the column
- Heavy matrix contaminates the source faster --- **performance is LOST!**


**Backflush ensures the highest data quality with GC/MS/MS!
The trade-off is 10-20% decrease in sensitivity.**



10% Fish Oil In Acetone: Retention Time Shifts Eliminated With Backflushing




Milk Extract - Using Bakeout To Remove High Boilers



Milk Extract

Once correct backflushing time is chosen, it will work for ALL high boilers thereafter, because heaviest elute first



Backflush: Many Advantages for GC/MS(/MS) Analysis of Complex Samples ('Dirty Matrices')

- Provides more consistent GC retention times
- Provides better, more consistent MS spectra through sequence
 - Reduces chemical noise that due to small carryover of matrix
 - Higher quality quantitation without increase in interfering ions
- Reduces contamination for the source
- Reduces analysis time
- Increases lifetime of analytical column

Benefits of Agilent GC/MS/MS Pesticide Analyzer

Retention Time Locking

- **No need to update the time segment RTs** after a column maintenance

Multimode inlet (MMI)

- injector adds flexibility by including standard, cold split/splitless, solvent vent (**LVI**) capabilities.

Capillary Flow Technology (CFT) and backflush

- Shorter analysis time, more consistent retention times and spectra, longer column life, and less frequent source cleaning -- **improve uptime.**

MS/MS MRM Database

- **Optimized and flexible MRM database** of hundreds of compounds

Pre-config. and factory setup analyzer

- **Factory setup and checked out on pesticide mixture** - ready to generate results on Day One



But...

How does the Analyzer work with the list of target pesticides in my lab?

Agilent's New Comprehensive MRM database with Extensive Flexibility

❖ Contains 7000 optimized MRMs for 1000 pesticides

-- over 3000 injections on \$70,000 worth of **chemical standards**

❖ Extensive flexibility allows method optimization

- average of 7 MRM transitions with relative intensity for each compound
 - provides alternatives to **avoid matrix interference**
- compound classification, CAS number etc. in excel format
 - allows **easy searching and sorting** for method customization
- three chromatographic methods (constant flow or pressure) with Retention Time and Retention Index
 - allows maximum freedom to follow **customer's workflow**
- absolute intensity for each MRM transition
 - allows **semi-quantitation** without standards

The Flexibility: 7 transitions; classifications; 3 RTs and RIs

MRM Database.xlsx - Microsoft Excel

Database has RTs (and RIs) to be used with three GC methods (CP, CF, and CF-screening).

1	Common Name	CAS # (format	Molecular Formula	Molecular Weight (average)	Molecular Weight (mono-isotopic)	CAS # (format 2	Classification	Classification	RT - 0502 screening (40.5 min), RTL = 18.111 min	RI - 0502 screenin	RT - 0501 (41.67 min), RTL = 16.593 min	RI - 0501	RT - 0502 (19.75 min), RTL = 9.143 min	RI - 0502
2	1 Acephate	30560-19-1	C4H10NO3PS	183.2	183.0119004	30560191	insecticide	Organophosphorus	6.234		7.610		5.651	
3	1 Acephate	30560-19-1	C4H10NO3PS	183.2	183.0119004	30560191	insecticide	Organophosphorus	6.234		7.610		5.651	
4	1 Acephate	30560-19-1	C4H10NO3PS	183.2	183.0119004	30560191	insecticide	Organophosphorus	6.234		7.610		5.651	
5	1 Acephate	30560-19-1	C4H10NO3PS	183.2	183.0119004	30560191	insecticide	Organophosphorus	6.234		7.610		5.651	
6	1 Acephate	30560-19-1	C4H10NO3PS	183.2	183.0119004	30560191	insecticide	Organophosphorus	6.234		7.610		5.651	
7	2 Etridiazole (Terrazole	2593-15-9	C5H5Cl3N2OS	259.2	259.0119004	2593159	fungicide	Thiadiazole	9.265		7.928		5.843	
8	2 Etridiazole (Terrazole	2593-15-9	C5H5Cl3N2OS	259.2	259.0119004	2593159	fungicide	Thiadiazole	9.265		7.928		5.843	
9	2 Etridiazole (Terrazole	2593-15-9	C5H5Cl3N2OS	259.2	259.0119004	2593159	fungicide	Thiadiazole	9.265		7.928		5.843	
10	2 Etridiazole (Terrazole	2593-15-9	C5H5Cl3N2OS	259.2	259.0119004	2593159	fungicide	Thiadiazole	9.265		7.928		5.843	
11	2 Etridiazole (Terrazole	2593-15-9	C5H5Cl3N2OS	259.2	259.0119004	2593159	fungicide	Thiadiazole	9.265		7.928		5.843	
12	2 Etridiazole (Terrazole	2593-15-9	C5H5Cl3N2OS	259.2	259.0119004	2593159	fungicide	Thiadiazole	9.265		7.928		5.843	
13	2 Etridiazole (Terrazole	2593-15-9	C5H5Cl3N2OS	259.2	259.0119004	2593159	fungicide	Thiadiazole	9.265		7.928		5.843	
14	3 Methabenzthiazuron	18691-97-9	C10H11N3OS	221.0	221.0622827	18691979	herbicide	Urea	13.461		11.176		7.269	
15	3 Methabenzthiazuron	18691-97-9	C10H11N3OS	221.0	221.0622827	18691979	herbicide	Urea	13.461		11.176		7.269	
16	3 Methabenzthiazuron	18691-97-9	C10H11N3OS	221.0	221.0622827	18691979	herbicide	Urea	13.461		11.176		7.269	
17	3 Methabenzthiazuron	18691-97-9	C10H11N3OS	221.0	221.0622827	18691979	herbicide	Urea	13.461		11.176		7.269	
18	3 Methabenzthiazuron	18691-97-9	C10H11N3OS	221.0	221.0622827	18691979	herbicide	Urea	13.461		11.176		7.269	
19	3 Methabenzthiazuron	18691-97-9	C10H11N3OS	221.0	221.0622827	18691979	herbicide	Urea	13.461		11.176		7.269	
20	3 Methabenzthiazuron	18691-97-9	C10H11N3OS	221.0	221.0622827	18691979	herbicide	Urea	13.461		11.176		7.269	
21	3 Methabenzthiazuron	18691-97-9	C10H11N3OS	221.0	221.0622827	18691979	herbicide	Urea	13.461		11.176		7.269	
22	4 Ethoxyquin	91-53-2	C14H19NO	217.3	217.1466642	91532	fungicide, plant	Quinoline	15.234		12.806		7.795	
23	4 Ethoxyquin	91-53-2	C14H19NO	217.3	217.1466642	91532	fungicide, plant	Quinoline	15.234		12.806		7.795	
24	4 Ethoxyquin	91-53-2	C14H19NO	217.3	217.1466642	91532	fungicide, plant	Quinoline	15.234		12.806		7.795	
25	4 Ethoxyquin	91-53-2	C14H19NO	217.3	217.1466642	91532	fungicide, plant	Quinoline	15.234		12.806		7.795	
26	4 Ethoxyquin	91-53-2	C14H19NO	217.3	217.1466642	91532	fungicide, plant	Quinoline	15.234		12.806		7.795	
27	4 Ethoxyquin	91-53-2	C14H19NO	217.3	217.1466642	91532	fungicide, plant	Quinoline	15.234		12.806		7.795	
28	5 Dicloran (Dichloran)	99-30-9	C6H4Cl2N2O2	207.0	205.9649828	99309	fungicide	Substituted benzen	14.775		12.536		7.798	
29	5 Dicloran (Dichloran)	99-30-9	C6H4Cl2N2O2	207.0	205.9649828	99309	fungicide	Substituted benzen	14.775		12.536		7.798	
30	5 Dicloran (Dichloran)	99-30-9	C6H4Cl2N2O2	207.0	205.9649828	99309	fungicide	Substituted benzen	14.775		12.536		7.798	

Average and exact Molecular Weight

Each pesticide is classified in two categories

The Flexibility: Excel format, relative and absolute Transition intensity

MassHunter Format

The absolute and relative intensities of transitions

(Color Scales): Red denotes strong intensity and blue denotes weak intensity among ALL transitions.

	P	Q	R	S	T	U	V	W	X	Y	Z	AA	AB	AC	AD	AE
				Precursor	MS1	Product	MS2	Dwell Time	CE (V)	Intensity	Relative Intensity	Quant	Chinese Name	China GB Method	Japanese Name	Note
1	RI - 0502	Common Name	ISTD?													
2		Acephate	FALSE	136.0	Wide	94.0	Wide	20	10	130	100%	Q0	乙酰甲胺磷	F	アセフエート	
3		Acephate	FALSE	142.1	Wide	96.0	Wide	20	10	30	21%	Q1	乙酰甲胺磷	F	アセフエート	
4		Acephate	FALSE	95.0	Wide	78.9	Wide	20	10	20	19%	Q2	乙酰甲胺磷	F	アセフエート	
5		Acephate	FALSE	95.0	Wide	79.9	Wide	20	10	20	17%	Q3	乙酰甲胺磷	F	アセフエート	
6		Acephate	FALSE	142.1	Wide	65.0	Wide	20	25	20	16%	Q4	乙酰甲胺磷	F	アセフエート	
7		Etridiazole (Terrazole, Echlomezol)	FALSE	183.0	Wide	140.0	Wide	20	15	2440	100%	Q0	土菌灵	A	エトリジアゾール	
8		Etridiazole (Terrazole, Echlomezol)	FALSE	211.1	Wide	183.0	Wide	20	15	2150	88%	Q1	土菌灵	A	エトリジアゾール	
9		Etridiazole (Terrazole, Echlomezol)	FALSE	185.0	Wide	142.1	Wide	20	15	1680	69%	Q2	土菌灵	A	エトリジアゾール	
10		Etridiazole (Terrazole, Echlomezol)	FALSE	211.1	Wide	140.0	Wide	20	15	1590	65%	Q3	土菌灵	A	エトリジアゾール	
11		Etridiazole (Terrazole, Echlomezol)	FALSE	213.1	Wide	185.0	Wide	20	15	1460	60%	Q4	土菌灵	A	エトリジアゾール	
12		Etridiazole (Terrazole, Echlomezol)	FALSE	213.1	Wide	142.1	Wide	20	15	1080	44%	Q5	土菌灵	A	エトリジアゾール	
13		Etridiazole (Terrazole, Echlomezol)	FALSE	183.0	Wide	108.0	Wide	20	45	500	20%	Q6	土菌灵	A	エトリジアゾール	
14		Methabenzthiazuron	FALSE	164.0	Wide	136.0	Wide	20	5	310	100%	Q0	甲基苯噻隆	D	メタベンズチアズロン	
15		Methabenzthiazuron	FALSE	163.1	Wide	136.0	Wide	20	15	190	60%	Q1	甲基苯噻隆	D	メタベンズチアズロン	
16		Methabenzthiazuron	FALSE	134.9	Wide	90.9	Wide	20	15	150	49%	Q2	甲基苯噻隆	D	メタベンズチアズロン	
17		Methabenzthiazuron	FALSE	134.9	Wide	108.0	Wide	20	15	100	32%	Q3	甲基苯噻隆	D	メタベンズチアズロン	
18		Methabenzthiazuron	FALSE	135.9	Wide	109.0	Wide	20	25	90	31%	Q4	甲基苯噻隆	D	メタベンズチアズロン	
19		Methabenzthiazuron	FALSE	135.9	Wide	64.9	Wide	20	35	80	25%	Q5	甲基苯噻隆	D	メタベンズチアズロン	
20		Methabenzthiazuron	FALSE	163.1	Wide	109.0	Wide	20	15	80	24%	Q6	甲基苯噻隆	D	メタベンズチアズロン	
21		Methabenzthiazuron	FALSE	164.0	Wide	108.0	Wide	20	30	50	16%	Q7	甲基苯噻隆	D	メタベンズチアズロン	
22		Ethoxyquin	FALSE	202.1	Wide	174.0	Wide	20	15	2890	100%	Q0	乙氧喹啉			
23		Ethoxyquin	FALSE	202.1	Wide	145.1	Wide	20	25	360	12%	Q1	乙氧喹啉			
24		Ethoxyquin	FALSE	203.0	Wide	175.0	Wide	20	15	360	12%	Q2	乙氧喹啉			
25		Ethoxyquin	FALSE	217.0	Wide	202.0	Wide	20	10	360	12%	Q3	乙氧喹啉			
26		Ethoxyquin	FALSE	174.0	Wide	146.1	Wide	20	10	310	11%	Q4	乙氧喹啉			
27		Ethoxyquin	FALSE	202.1	Wide	159.0	Wide	20	30	260	9%	Q5	乙氧喹啉			
28		Dicloran (Dichloran)	FALSE	206.0	Wide	176.0	Wide	20	15	2480	100%	Q0	氯硝胺	B	ジクロラン	
29		Dicloran (Dichloran)	FALSE	207.9	Wide	178.0	Wide	20	15	1560	63%	Q1	氯硝胺	B	ジクロラン	
30		Dicloran (Dichloran)	FALSE	124.0	Wide	73.1	Wide	20	15	1410	57%	Q2	氯硝胺	B	ジクロラン	

One Quant and several Qualifications for each compound

Click on the sorting button to show/hide Quant and Qualifier Ions

The screenshot shows an Excel spreadsheet titled 'MRM Database.xlsx'. The data table has columns for various parameters including 'Common Name', 'ISTD?', 'Precursor I_c', 'MS1 Resolution', 'Product I_c', 'MS2 Resolution', 'Dwell Time', 'CE (V)', 'Intensity Scale within the Database', 'Transition Relative Intensity', 'Quant (Q0)/Qual', 'China GB Method', and 'Japanese Name'. A sorting menu is open over the 'Quant' column, with a red arrow pointing to the 'Sort A to Z' button. The menu also includes options for 'Sort Z to A', 'Sort by Color', 'Clear Filter From "Quant (Q0)/Qual"', 'Filter by Color', and 'Text Filters'. The 'Text Filters' section is expanded, showing a list of checkboxes for Q0 through Q7, with Q0, Q1, and Q2 checked. The spreadsheet data includes entries for 'Acephate', 'Etridiazole (Terrazole, Echlomezol)', 'Methabenzthiazuron', and 'Ethoxyquin'.

Common Name	ISTD?	Precursor I _c	MS1 Resolution	Product I _c	MS2 Resolution	Dwell Time (ms)	CE (V)	Intensity Scale within the Database	Transition Relative Intensity	Quant (Q0)/Qual	China GB Method	Japanese Name
Acephate	FALSE	136.0	Wide	94.0	Wide	20	10				F	アセフエート
Etridiazole (Terrazole, Echlomezol)	FALSE	183.0	Wide	140.0	Wide	20	15				A	エトリジアゾール
Methabenzthiazuron	FALSE	164.0	Wide	136.0	Wide	20	5				D	メタベンズチアズロン
Ethoxyquin	FALSE	202.1	Wide	174.0	Wide	20	15	2890	100%	Q0		
Dicloran (Dichloran)	FALSE	206.0	Wide	176.0	Wide	20	15	2480	100%	Q0	B	ジクロラン

Use the sorting function to quickly select a Quant (Q0) and top three Qualifier ions (to build a method)!

MRM Database.xlsx - Microsoft Excel

	O	P	Q	R	S	T	U	V	W	X	Y	Z	AA	AB	AC	AD	AE
1	RI - 0502	Common Name	ISTD?	Precursor I _c	MS1 Resolution	Product I _c	MS2 Resolution	Dwell Time (ms)	CE (V)	Intensity Scale within the Database	Transition Relative Intensity (%)	Quant (Q0) / Qualifier	Chinese Name	China GB Method Group	Japanese Name	Note	
2		Acephate	FALSE	136.0	Wide	94.0	Wide	20	10	130	100%	Q0	乙酰甲胺磷	F	アセフエート		
3		Acephate	FALSE	142.1	Wide	96.0	Wide	20	10	30	21%	Q1	乙酰甲胺磷	F	アセフエート		
4		Acephate	FALSE	95.0	Wide	78.9	Wide	20	10	20	19%	Q2	乙酰甲胺磷	F	アセフエート		
5		Acephate	FALSE	95.0	Wide	79.9	Wide	20	10	20	17%	Q3	乙酰甲胺磷	F	アセフエート		
7		Etridiazole (Terrazole, Echlomezol)	FALSE	183.0	Wide	140.0	Wide	20	15	2440	100%	Q0	土菌灵	A	エトリジアゾール		
8		Etridiazole (Terrazole, Echlomezol)	FALSE	211.1	Wide	183.0	Wide	20	15	2150	88%	Q1	土菌灵	A	エトリジアゾール		
9		Etridiazole (Terrazole, Echlomezol)	FALSE	185.0	Wide	142.1	Wide	20	15	1680	69%	Q2	土菌灵	A	エトリジアゾール		
10		Etridiazole (Terrazole, Echlomezol)	FALSE	211.1	Wide	140.0	Wide	20	15	1590	65%	Q3	土菌灵	A	エトリジアゾール		
14		Methabenzthiazuron	FALSE	164.0	Wide	136.0	Wide	20	5	310	100%	Q0	甲基苯噻隆	D	メタベンズチアズロン		
15		Methabenzthiazuron	FALSE	163.1	Wide	136.0	Wide	20	15	190	60%	Q1	甲基苯噻隆	D	メタベンズチアズロン		
16		Methabenzthiazuron	FALSE	134.9	Wide	90.9	Wide	20	15	150	49%	Q2	甲基苯噻隆	D	メタベンズチアズロン		
17		Methabenzthiazuron	FALSE	134.9	Wide	108.0	Wide	20	15	100	32%	Q3	甲基苯噻隆	D	メタベンズチアズロン		
22		Ethoxyquin	FALSE	202.1	Wide	174.0	Wide	20	15	2890	100%	Q0	乙氧喹啉				
23		Ethoxyquin	FALSE	202.1	Wide	145.1	Wide	20	25	360	12%	Q1	乙氧喹啉				
24		Ethoxyquin	FALSE	203.0	Wide	175.0	Wide	20	15	360	12%	Q2	乙氧喹啉				
25		Ethoxyquin	FALSE	217.0	Wide	202.0	Wide	20	10	360	12%	Q3	乙氧喹啉				
28		Dicloran (Dichloran)	FALSE	206.0	Wide	176.0	Wide	20	15	2480	100%	Q0	氯硝胺	B	ジクロラン		
29		Dicloran (Dichloran)	FALSE	207.9	Wide	178.0	Wide	20	15	1560	63%	Q1	氯硝胺	B	ジクロラン		
30		Dicloran (Dichloran)	FALSE	124.0	Wide	73.1	Wide	20	15	1410	57%	Q2	氯硝胺	B	ジクロラン		
31		Dicloran (Dichloran)	FALSE	176.1	Wide	148.0	Wide	20	15	1100	44%	Q3	氯硝胺	B	ジクロラン		

Ready 20 of 278 records found


Click on the sorting button to select/show the compounds to be added into the acquisition method

The screenshot shows an Excel spreadsheet with a table of compounds. The table has columns for various parameters: Common Name, Precursor, MS1 Resolution, Product, MS2 Resolution, Dwell Time, CE (V), Intensity Scale within the Database, Transition Relative Intensity, Quant (Q0)/Quant, Chinese Name, China GB Method, Japanese Name, and Notes. A filter menu is open over the 'Common Name' column, showing options like 'Sort A to Z', 'Sort Z to A', and 'Filter by Color'. A 'Text Filters' dialog box is also open, showing a list of compounds with checkboxes for selection. Red arrows point to the filter button and the filter menu. A green box at the bottom contains the following text:


It is easy to add a column to associate each compound with a lab method or study to allow a quick sort to build an acq. method.

	O	P	Q	R	S	T	U	V	W	X	Y	Z	AA	AB	AC	AD	AE
1	RI - 0502	Common Name	ISTD?	Precursor	MS1 Resolution	Product	MS2 Resolution	Dwell Time (ms)	CE (V)	Intensity Scale within the Database	Transition Relative Intensity	Quant (Q0)/Quant	Chinese Name	China GB Method	Japanese Name	Notes	
2		Acephate	FALSE	136.0	Wide	94.0	Wide	20	10	130	100%	Q0	乙酰甲胺磷	F	アセフエート		
3		Acephate	FALSE	142.1	Wide	96.0	Wide	20	10	30	21%	Q1	乙酰甲胺磷	F	アセフエート		
4		Acephate	FALSE	95.0	Wide	78.9	Wide	20	10	20	19%	Q2	乙酰甲胺磷	F	アセフエート		
5		Acephate	FALSE	95.0	Wide	79.9	Wide	20	10	20	17%	Q3	乙酰甲胺磷	F	アセフエート		
7		Etridiazuron	FALSE	183.0	Wide	140.0	Wide	20	15	2440	100%	Q0	土菌灵	A	エトリジアゾール		
8		Etridiazuron	FALSE	211.1	Wide	183.0	Wide	20	15	2150	88%	Q1	土菌灵	A	エトリジアゾール		
9		Etridiazuron	FALSE	185.0	Wide	142.1	Wide	20	15	1680	69%	Q2	土菌灵	A	エトリジアゾール		
10		Etridiazuron	FALSE	211.1	Wide	140.0	Wide	20	15	1590	65%	Q3	土菌灵	A	エトリジアゾール		
14		Methidathion	FALSE	164.0	Wide	136.0	Wide	20	5	310	100%	Q0	甲基苯噻隆	D	メタベンズチアズロン		
15		Methidathion	FALSE	163.1	Wide	136.0	Wide	20	15	190	60%	Q1	甲基苯噻隆	D	メタベンズチアズロン		
16		Methidathion	FALSE	134.9	Wide	90.9	Wide	20	15	150	49%	Q2	甲基苯噻隆	D	メタベンズチアズロン		
17		Methidathion	FALSE	134.9	Wide	108.0	Wide	20	15	100	32%	Q3	甲基苯噻隆	D	メタベンズチアズロン		
22		Ethoxyquin	FALSE	202.1	Wide	174.0	Wide	20	15	2890	100%	Q0	乙氧喹啉				
23		Ethoxyquin	FALSE	202.1	Wide	145.1	Wide	20	25	360	12%	Q1	乙氧喹啉				
24		Ethoxyquin	FALSE	203.0	Wide	175.0	Wide	20	15	360	12%	Q2	乙氧喹啉				
25		Ethoxyquin	FALSE	217.0	Wide	202.0	Wide	20	10	360	12%	Q3	乙氧喹啉				
28		Dicloran (Dichloran)	FALSE	206.0	Wide	176.0	Wide	20	15	2480	100%	Q0	氯硝胺	B	ジクロラン		
29		Dicloran (Dichloran)	FALSE	207.9	Wide	178.0	Wide	20	15	1560	63%	Q1	氯硝胺	B	ジクロラン		
30		Dicloran (Dichloran)	FALSE	124.0	Wide	73.1	Wide	20	15	1410	57%	Q2	氯硝胺	B	ジクロラン		
31		Dicloran (Dichloran)	FALSE	176.1	Wide	148.0	Wide	20	15	1100	44%	Q3	氯硝胺	B	ジクロラン		

Why Do We Need More than 2 MRM Transitions?



MRM Transitions are not Universal, Should Choose them According to Matrices



Summary: Pesticide GC/MS/MS Analyzer

The Pesticide GC/MS/MS Analyzer is **tested as a unit** in the factory to ensure a quick and successful installation

New **flexible and comprehensive** Pesticide MRM Database expands the target compound list to 1000 to meet users' need

The MRM Database

- allows users to build acquisition methods without acquiring expensive or hard-to-get pesticide standards (saves time and money)
- applies to either constant flow or constant pressure method (adapts to user's preference or analytical method)
- has multiple MRM transitions (average 7) for each compound (helps to provide alternatives to work around matrix interference)
- shows relative intensity of each MRM transition (facilitates transition selection and acquisition method creation)

ASMS 2011 Introduction: 7200 High Resolution, Accurate Mass GC/Q-TOF



**Many ASMS Conferees Asked:
How Would a GC/Q-TOF Benefit My Lab?**

**Other ASMS Conferees Asked:
Why Q-TOF and Not Just TOF?**

**For This e-Seminar, You Might Ask:
How Can a Q-TOF Extract More from
Complex Samples?**

Fundamental Benefits (Agilent 7200 Q-TOF)

- High resolution ($> 10K$, typically $> 13K$ FWHM)
 - Increased detector selectivity (few interferences)

TOF Accurate Mass to Eliminate Matrix Interferants

Okra QuEChERS Extract

Matrix interferant ion (b-Tocopherol)


150.06839 Da

Analyte Indoxacarb ion (100pg)


150.01195 Da (fragment ion)

$\Delta m = 0.0564$ Da.

'MSD' resolution



TOF high resolution



TOF Accurate Mass to Eliminate Matrix Interferants

Okra QuEChERS Extract

Matrix interferant ion (b-Tocopherol)


150.06839 Da

Analyte Indoxacarb ion (100pg)


150.01195 Da (fragment ion)

$\Delta m = 0.0564$ Da.

'MSD' resolution



TOF high resolution




If even more selectivity is needed, option of MS/MS MS/MS with high resolution and accurate mass!



Agilent Technologies

Resolving Power




Resolving Power:

$$R = m/z / \text{FWHM}$$

Mass Accuracy: $\Delta m/z = dm/m/z * 10^6$,
parts per million (ppm)


PFTBA mass 614

C12F24N=613.964203



Agilent Technologies

Resolving Power




Resolving Power:

$$R = m/z / \text{FWHM}$$

Mass Accuracy: $\Delta m/z = dm/mz * 10^6$,
parts per million (ppm)


PFTBA mass 614

C12F24N=613.964203




Spectral Presentation (Tune File)

Most users think in “centroid”, but the MS operates in “profile”




Selectivity for Isobaric Ions



Accurate Mass Makes Mass Defect Important

Type	Element	Symbol	Integer Mass	Exact Mass	Abundance	X+1 Factor	X+2 Factor	Mass Defect
X	Hydrogen	H	1	1.0078	99.99			0.0078
		D or ² H	2	2.0141	0.01			0.0141
X+1	Carbon	¹² C	12	12	98.91			0
		¹³ C	13	13.0034	1.1	1.1n _C	0.0060n _C ²	0.0034
X+1	Nitrogen	¹⁴ N	14	14.0031	99.6			0.0031
		¹⁵ N	15	15.0001	0.4	0.37n _N		0.0001
X+2	Oxygen	¹⁶ O	16	15.9949	99.76			-0.0051
		¹⁷ O	17	16.9991	0.04	0.04n _O		-0.0009
		¹⁸ O	18	17.9992	0.2		0.20n _O	-0.0008
X	Fluorine	F	19	18.9984	100			-0.0016
X+2	Silicon	²⁸ Si	28	27.9769	92.2			-0.0231
		²⁹ Si	29	28.9765	4.7	5.1n _{Si}		-0.0235
		³⁰ Si	30	29.9738	3.1		3.4n _{Si}	-0.0262
X	Phosphorus	P	31	30.9738	100			-0.0262
X+2	Sulfur	³² S	32	31.9721	95.02			-0.0279
		³³ S	33	32.9715	0.76	0.8n _S		-0.0285
		³⁴ S	34	33.9679	4.22		4.4n _S	-0.0321
X+2	Chlorine	³⁵ Cl	35	34.9689	75.77			-0.0311
		³⁷ Cl	37	36.9659	24.23	32.5n _{Cl}		-0.0341
X+2	Bromine	⁷⁹ Br	79	78.9183	50.5			-0.0817
		⁸¹ Br	81	80.9163	49.5		98.0n _{Br}	-0.0837
X	Iodine	I	127	126.9045	100			-0.0955


Selectivity for Isobaric Ions




Fundamental Benefits (Agilent 7200 Q-TOF)

- High resolution ($> 10K$, typically $> 13K$ FWHM)
 - Increased detector selectivity (few interferences)
- Accurate mass measurements (low to sub-ppm)
 - < 5 ppm in TOF (typically < 2 ppm)
 - < 10 ppm in Q-TOF (typically < 5 ppm)
 - Valuable qualitative information about each ion

Many possible formulas with an MSD or IT But only a few with TOF



Octafluoronaphthalene (CAS 313-72-4)



mass uncertainty		# of Possible Formulas
ppm	amu	
1000	0.3	7657
368	0.1	4050
100	0.03	1223
37	0.01	466
10	0.003	120
4	0.001	43
1	0.0003	11
0.4	0.0001	5
0.1	0.00003	2

Formulas made of:
C, H, N, O, F, & Cl

Accurate mass reduces risk of investing effort on the wrong molecule

Fragment ion mass with high mass accuracy

Helps unambiguously identify corresponding formula and hence exact mass of a fragment

Most probable fragment ion

	Molecular Ion		Fragment Ions							
	m/z	Formula	Δ ppm	1 Formula	Δ ppm	2 Formula	Δ ppm	3 Formula	Δ ppm	4 Formula
Chlorpyrifos-methyl	320.8944	C7H7Cl3NO3PS	-0.7	C7 H7 Cl [37Cl] N O3 P S	0.0	C7 H7 Cl2 N O3 P S	0.0	C2 H6 O2 P S		
			215.4	C6 H2 Cl3 N O2 P S						
Dichlorvos	219.9454	C4H7Cl2O4P	1.6	C4H7ClO4P	-0.9	C2H6O3P	3.7	C4 H7 [37Cl] O4 P		
						336.5	C3 H2 Cl2 O3 P			
Endosulfan sulfate	419.8112	C9 H6 Cl6 O4 S	-2.1	C9 H6 Cl4 [37Cl] O4 S	-0.7	C5 Cl5 [37Cl]	0.0	C5 Cl4 [37Cl]2		
						-220.9	C8 H3 Cl5			
Propachlor	211.0758	C11 H14 Cl N O	-1.0	C10 H11 Cl N O	1.1	C11 H14 N O	1.8	C8 H8 Cl N O	-3.9	C6 H5
								771.6	C2 H2 Cl O	
Fluazifop-p-butyl	383.1339	C19 H20 F3 N O4	-2.2	C19 H20 F2 N O4	-1.1	C14 H11 F3 N O2	-1.3	C12 H7 F3 N O	2.1	C6 H3 F3 N
								-494.2	C7 H14 O3	
Triazophos	313.0645	C12 H16 N3 O3 P S	-1.4	C10 H12 N3 O3 P S	-2.7	C8 H8 N3 O3 P S	-0.6	C8 H7 N3 O	-1.9	C8 H8 N3 O
			-47.7	C11 H14 N2 O3 P S	-54.1	C9 H10 N2 O3 P S				

Examples from building accurate mass pesticide library

NIST 2011 MS Interpreter with Accurate Mass

COCAINE.MSP - MS Interpreter

File Edit View Options Help

m/z=182 (1/)

Formula Calculator


m = 182 C17H21NO4

Calculate Options Parent = 303
Loss = 121

18 Ions	O+E	RDB	Mass	C	H	N	O
C7H20NO4	+	0	182.13923	7	20	1	4
C8H8NO4	Even	5.5-6.5	182.04533	8	8	1	4
C9H10O4	Odd	5	182.05791	9	10	0	4
C9H12NO3	Even	4.5-5.5	182.08172	9	12	1	3
C10NO3	Even	11.5-12.5	181.98782	10	0	1	3
C10H14O3	Odd	4	182.09429	10	14	0	3
C10H16NO2	Even	3.5-4.5	182.11810	10	16	1	2
C11H2O3	Odd	11	182.00039	11	2	0	3
C11H4NO2	Even	10.5-11.5	182.02420	11	4	1	2
C11H18O2	Odd	3	182.13068	11	18	0	2
C11H20NO	Even	2.5-3.5	182.15449	11	20	1	1


Structure, Maximum Dissociation = 148

m/z	mass	formula	loss	type	rate	abund
182 (1/3)	182.118104	C10H16NO2	C7H5O2	dissociation	-65	716



Mass Spectrum for C17H21NO4; MW = 303; CAS = 50-36-2

Cocaine




L-Click/R-Click => Next/Prev fragment

Fundamental Benefits (Agilent 7200 Q-TOF)

- High resolution ($> 10K$, typically $> 13K$ FWHM)
 - Increased detector selectivity (few interferences)
- Accurate mass measurements (low to sub-ppm)
 - < 5 ppm in TOF (typically < 2 ppm)
 - < 10 ppm in Q-TOF (typically < 5 ppm)
 - Valuable qualitative information about each ion
- Structural elucidation with Accurate Mass MS/MS studies
 - High sensitivity tool to complement NMR

The Problem – Confirm Most Likely Structure

Kava Extract - Compound "B", $C_{16}H_{14}O_4$
(Rings + Double Bonds = 10)




Experimental measurements


m/z (experimental)	Formula	Error (ppm)	Score
- H			
269.0802	$C_{16}H_{13}O_4$	2.2	80.7
- C_6H_5			
193.0494	$C_{10}H_9O_4$	0.6	96.7
- $CH=CH-C_6H_5$			
167.0334	$C_8H_7O_4$	3.0	N/A
- $CH_2=CH-C_6H_5$			
166.0259	$C_8H_6O_4$	0.6	N/A
- CO			
138.0310	$C_7H_6O_3$	1.1	98.1
- CO			
110.0359	$C_6H_6O_2$	3.0	N/A
- CH_3			
95.0127	$C_5H_3O_2$	0.9	99.5

For the 5 candidate structures, only one fit the losses identified by CID experiments on multiple precursor ions

Problem – confirm most likely structure



Problem – confirm most likely structure




Mass at 138 is consistent with loss of **COCH**=CH-C₆H₅ (131.04969)

or

C₂H₄CH=CH-C₆H₅ (131.086075) from 269.08020.

Problem – confirm most likely structure



Mass at 138 is consistent with loss of **COCH=CH-C₆H₅** (131.04969)

or

C₂H₂CH=CH-C₆H₅ (131.086075) from 269.0802.

However, measured value of 269.0802 - 138.0310 = 131.04920 is consistent only with **COCH=CH-C₆H₅**.

Fundamental Benefits (Agilent 7200 Q-TOF)

- High resolution ($> 10K$, typically $> 13K$ FWHM)
 - Increased detector selectivity (few interferences)
- Accurate mass measurements (low to sub-ppm)
 - < 5 ppm in TOF (typically < 2 ppm)
 - < 10 ppm in Q-TOF (typically < 5 ppm)
 - Valuable qualitative information about each ion
- Structural elucidation with Accurate Mass MS/MS studies
 - High sensitivity tool to complement NMR
- “Fast”, full spectra acquisition with excellent sensitivity

What about TOF *SPEED*?


TOF always collects full mass range
Q-TOF always display full product ion spectrum

- **Acquisition Rate:** transients (pulses) /second
 - 10,000 transients/second
- Sum of transients = **Spectral Rate:**
 - Typical max rate: 25-200 spectra/sec (**Hz**) to disk
 - Usable rate is limited by signal level (ion count)
- New analysis opportunities for GC/MS:
 - High Throughput: ~20 Hz
 - Ultra high resolution GC: ~ 40Hz
 - GCxGC: ~50-200 Hz




“Speed” Enhances Deconvolution


Deconvolution Requires Time Offset



High Data Rate = Better Deconvolution



Time **Slow data rate will not pick each peak apex**



Time **Fast data rate will allow deconvolution of closely eluting peaks**



Successful GC/MS Applications Require
the Correct Combination of:

GC Resolution
+
MS Resolving Power



Agilent Technologies

Successful GC/MS Applications Require
the Correct Combination of:

GC Resolution

+

MS Resolving Power

+

Δ Mass Defect (**Mother Nature**)

+

Relative Ion Intensities (**Sample**)


There is more to learn about these new apps



Agilent Technologies


GC or MS: Relative intensity affects the result

Profile view of mass peaks



Centroid view of mass peaks obscures this fact

What Are the Application Drivers for Q-TOF?



**Thank You
For Joining Us Today**

**Time For
Q&A**