

# Utilization of Hydrogen Carrier Gas on a High Resolution GC-TOFMS System: An Application Compendium

Joe Binkley; David Alonso | LECO Corporation, St. Joseph, MI USA

## Background

Both the costs associated with a dwindling helium supply and the need for high sample throughput have fueled the desire to develop fast gas chromatography methods using hydrogen as a carrier gas. This poster demonstrates the ability to utilize hydrogen carrier gas on a high resolution GC-TOFMS system, the Pegasus® GC-HRT. Methods and data from multiple application markets, including specialty chemicals, forensics, and metabolomics, will be displayed. Various polymer additives were used to represent the specialty chemical market, drugs of abuse were used to demonstrate capabilities in the forensic market, and a derivatized metabolite mixture designed for metabolomics research was used for the metabolomics market.

GC-HRT methods utilizing hydrogen carrier gas were developed for analysis of representative specialty chemical, forensic, and metabolomic markets.

- Polymer additives were dissolved in chloroform prior to analysis.
- Drugs of abuse were dissolved in organic solvents and analyzed underivatized.
- Metabolites were derivatized prior to GC-HRT analysis using an optimized two-step procedure: 1) Treatment with methoxylamine hydrochloride and 2) MSTFA.



Helium




Figure 1. Van Deemter curves for He, N<sub>2</sub>, and H<sub>2</sub>. OPGV = Optimum practical gas velocity.

## Instrumentation



Figure 2. Pegasus GC-HRT High Resolution TOFMS.

### Pegasus GC-HRT Capabilities

- High Resolution (Up to R=50,000)
- Mass Accuracy (<1 ppm)
- Fast Acquisition Rates (Up to 200 spectra/second)
- High Resolution Deconvolution™ (HRD™)
- H<sub>2</sub> Carrier Gas Compatible

## Experimental

### GC-HRT Methods (Polymer Additives)

**GC Conditions**  
 Column: Rxi-1ms 20 m x 0.18 mm x 0.18 micron film (Restek)  
 Carrier: Hydrogen @ 0.8 mL/min  
 Injection: 1 µL, split 20:1 @ 300°C  
 Oven: 50°C (0.5 min) to 300°C @ 50°C/min  
 Tr. Line: 300°C

### HRT Conditions

Resolution: High (25k)  
 Mass Range: 40-700 m/z  
 Ionization: EI

### GC-HRT Methods (Drugs of Abuse)

**GC Conditions**  
 Column: Rxi-1ms 20 m x 0.18 mm x 0.18 micron film (Restek)  
 Carrier: Hydrogen @ 0.8 mL/min  
 Injection: 1 µL, split 50:1 @ 280°C  
 Oven: 40°C (2 min) to 300°C @ 30°C/min  
 Tr. Line: 300°C

### HRT Conditions

Resolution: High (25k)  
 Mass Range: 40-400 m/z  
 Ionization: EI

### GC-HRT Methods (Metabolomics)

**GC Conditions**  
 Column: Rxi-1ms 20 m x 0.18 mm x 0.18 micron film (Restek)  
 Carrier: Hydrogen @ 0.8 mL/min  
 Injection: 1 µL, split 10:1 @ 250°C  
 Oven: 80°C to 300°C @ 40°C/min  
 Tr. Line: 300°C

### HRT Conditions

Resolution: High (25k)  
 Mass Range: 50-600 m/z  
 Ionization: EI




Figure 3. AIC Showing Representative Polymer Additives.

Table I. Tabular Data for Polymer Additives

Compound	Retention (min)	Abundance	Observed (m/z)	Library (m/z)	Mass Accuracy (ppm)	Formula
Butylated Hydroxytoluene	813	230,182	230.182	1.2	-0.43	C <sub>19</sub> H <sub>22</sub> O
Dibutyl phthalate	875	278,153	NA	NA	NA	C <sub>24</sub> H <sub>38</sub> O <sub>4</sub>
Irganox 168	688	646,039	646.039	0.32	0.32	C <sub>24</sub> H <sub>38</sub> O <sub>2</sub>
Irganox 1076	566	529,803	529.803	0.49	0.49	C <sub>24</sub> H <sub>38</sub> O <sub>2</sub>




Figure 4. Deconvoluted (Peak True) and Library Hit Spectra for Butylated Hydroxytoluene (A) and Dibutyl Phthalate (B).




Figure 5. Example of Accurate Mass Data Used to Correctly Identify an Incorrect NIST Library Match.




Figure 6. Example of Accurate Mass Data Used to Confirm a NIST Library Match.

## Results




Figure 7. AIC Showing Representative Drugs of Abuse.

Table II. Tabular Data for Drugs of Abuse

Drug	Retention (min)	Abundance	Observed (m/z)	Library (m/z)	Mass Accuracy (ppm)	Formula
6-APB	682	375,092	375,092	-0.59	-0.59	C <sub>12</sub> H <sub>16</sub> N <sub>2</sub> O
Methylene	780	207,089	NA	NA	NA	C <sub>11</sub> H <sub>15</sub> N <sub>3</sub>
5-IT	719	174,155	174,155	2.26	2.26	C <sub>11</sub> H <sub>15</sub> N <sub>3</sub>
N,N-DMT	566	188,130	188,130	-0.27	-0.27	C <sub>11</sub> H <sub>15</sub> N <sub>3</sub>
2C-B	874	299,022	299,022	0.98	0.98	C <sub>12</sub> H <sub>16</sub> N <sub>2</sub> O
Acetyl norfentanyl	779	218,144	218,146	1.2	1.2	C <sub>17</sub> H <sub>21</sub> N <sub>3</sub> O
PV9	809	273,287	NA	NA	NA	C <sub>18</sub> H <sub>23</sub> N <sub>3</sub> O
5-methoxy DALT	763	270,172	270,172	-2.3	-2.3	C <sub>18</sub> H <sub>23</sub> N <sub>3</sub> O
AH 7921	761	328,104	NA	NA	NA	C <sub>18</sub> H <sub>23</sub> N <sub>3</sub> O
25I-NBOMe	799	427,069	NA	NA	NA	C <sub>19</sub> H <sub>25</sub> N <sub>3</sub> O
3OC-NBOMe	762	348,256	NA	NA	NA	C <sub>19</sub> H <sub>25</sub> N <sub>3</sub> O
AB-CHMINACA	384	NA	NA	NA	NA	C <sub>19</sub> H <sub>25</sub> N <sub>3</sub> O
THU 2201	737	360,162	360,162	-1.5	-1.5	C <sub>19</sub> H <sub>25</sub> N <sub>3</sub> O
AM201 benzimidazole analog	857	360,162	360,162	-3.34	-3.34	C <sub>19</sub> H <sub>25</sub> N <sub>3</sub> O




Figure 8. Deconvoluted (Peak True) and Library Spectra with Accurate Mass Utilized for Confirmation.




Figure 9. Deconvoluted (Peak True) and Library Spectra with Accurate Mass Utilized for Confirmation.




Figure 10. AIC Showing Representative Metabolites.

Table III. Tabular Data for TMS Derivatized Metabolites

Metabolite	Retention (min)	Abundance	Observed (m/z)	Library (m/z)	Mass Accuracy (ppm)	Formula
L-Proline, trimethylsilyl ester	771	187,1023	187,1023	0.034	0.034	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>3</sub> Si <sub>3</sub>
Serine, 2TMS	761	249,1211	NA	NA	NA	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>3</sub> Si <sub>3</sub>
Silanol, trimethyl-, phosphate (E:1)	898	334,0983	334,0983	1.172	1.172	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>3</sub> Si <sub>3</sub>
Glycine, 3TMS	826	291,1506	291,1498	-0.588	-0.588	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>3</sub> Si <sub>3</sub>
L-Serine, 2TMS	833	321,1606	NA	NA	NA	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>3</sub> Si <sub>3</sub>
meso-Erythritol, 4TMS	864	430,2147	NA	NA	NA	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>3</sub> Si <sub>3</sub>
L-Proline, 2TMS	747	273,1211	NA	NA	NA	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>3</sub> Si <sub>3</sub>
Ribitol, 5TMS	896	512,3036	NA	NA	NA	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>3</sub> Si <sub>3</sub>
Putrescine, 4TMS	899	378,2781	378,2613	-3.971	-3.971	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>3</sub> Si <sub>3</sub>
D-Fructofuranose, 5TMS	836	543,2047	NA	NA	NA	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>3</sub> Si <sub>3</sub>
D(-)-Fructofuranose, 5TMS	781	543,2047	NA	NA	NA	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>3</sub> Si <sub>3</sub>
OTIC ACID, 4TMS	848	480,1846	NA	NA	NA	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>3</sub> Si <sub>3</sub>
3-β-Norbornene, 8TMS	726	558,3188	NA	NA	NA	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>3</sub> Si <sub>3</sub>
L-Tryptophan, 3TMS	726	420,2075	NA	NA	NA	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>3</sub> Si <sub>3</sub>
Nonadecanoic acid, TMS	740	370,1206	NA	NA	NA	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>3</sub> Si <sub>3</sub>
Sucrose, 8TMS	806	918,4188	NA	NA	NA	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>3</sub> Si <sub>3</sub>




Figure 11. Deconvoluted (Peak True) and Library Spectra with Accurate Mass Utilized for Confirmation.




Figure 12. Deconvoluted (Peak True) and Library Spectra for Sucrose. Accurate Mass Measurement of m/z 361.16861 Confirms Chemical Formula of C<sub>14</sub>H<sub>27</sub>O<sub>11</sub>Si<sub>3</sub> Which is Part of the Sucrose Structure.