

Extraction of Polyaromatic Hydrocarbons (PAHs) from Natural Waters Using ISOLUTE® PAH Solid Phase Extraction Columns




Figure 1. Structure of Pyrene

Introduction

This application note describes the extraction of polyaromatic hydrocarbons (PAHs) from natural waters with a high concentration of humic acids, which can be problematic during analysis.

The ISOLUTE® PAH column contains C18 sorbent for retention of the PAHs, plus a layer of aminopropyl (NH₂) sorbent which acts as a scavenger and prevents humic interferences from contaminating the extract.

Extracts are suitable for analysis by GC or LC based techniques.

Analytes

Polyaromatic Hydrocarbons

Sample Preparation Procedure

Format:

ISOLUTE® PAH 1.5 g/6 mL Part number 927-0150-C



Sample Pre-treatment:

To a 1L water sample, add 10% isopropanol (v/v), and mix thoroughly. If an internal standard is to be used, add at this stage.

Solvation:

Solvate the column with isopropanol (5 mL) at a flow rate of 5 mL/min.

Equilibration:

Equilibrate the column with deionized water containing 2% (v/v) isopropanol (5 mL) at a flow rate of 5 mL/min.

Sample Application:

Load the sample (1 L) at a flow rate of 20 mL/min. Rinse the sample bottle with acetone (10 mL) dilute this to 100 mL with deionized water, apply to the extraction column. Repeat this step if necessary to prevent losses of the analytes to the walls of the sample bottle.

Interference Elution:

Rinse the column with 90:10 (v/v) deionized water:isopropanol (5 mL). Dry the column thoroughly by vacuum aspiration for 15–20 mins (use -20”Hg vacuum). If losses of the more volatile PAHs are found at this stage, consider drying the column by centrifugation.

Analyte Elution:

Elute analytes with 2 x 3 mL THF/hexane (50:50, v/v). Use a soak step of up to 2 min between aliquots to maximize recoveries of the analytes.

Evaporate the extract at low temperature under nitrogen to concentrate further if required to improve detection limit. Care should be taken at this stage to prevent losses of the more volatile analytes, the extract should NOT be taken to dryness. A keeper solvent may be of use in this situation.

General Comments:

1. Dichloromethane is a suitable alternative elution solvent. The THF/hexane elution solvent was developed and optimized in order to replace dichloromethane as necessary.
2. These compounds are very non-polar and the smaller PAHs can be extremely volatile. These properties of the analytes can lead to problems during extraction and analysis. Losses of the analytes can occur during loading of the sample, as analytes can remain on the walls of the sample container, and during concentration of the sample during evaporation.
3. This method can be scaled down for use with a lower sample volume, using ISOLUTE PAH 750 mg/ 3 mL SPE columns, p/n 927-0075-B.

Ordering Information

Part Number	Description	Quantity
927-0075-B	ISOLUTE® PAH 750 mg/3 mL	50
927-0150-C	ISOLUTE® PAH 1.5 g/6 mL	30

EUROPE

Main Office: +46 18 565900
 Toll Free: +800 18 565710
 Fax: +46 18 591922
 Order Tel: +46 18 565710
 Order Fax: +46 18 565705
 order@biotage.com
 Support Tel: +46 18 56 59 11
 Support Fax: + 46 18 56 57 11
 eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
 Toll Free: +1 800 446 4752
 Fax: +1 704 654 4917
 Order Tel: +1 704 654 4900
 Order Fax: +1 434 296 8217
 ordermailbox@biotage.com
 Support Tel: +1 800 446 4752
 Outside US: +1 704 654 4900
 us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123
 Fax: +81 3 5627 3121
 jp_order@biotage.com
 jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 2898 6655
 Fax: +86 21 2898 6153
 cn_order@biotage.com
 cn-1-pointsupport@biotage.com

To locate a distributor,
 please visit our website
 www.biotage.com

Part Number: IST1025A.V.1

© 2016 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.