

Simultaneous determination of known and unknown Extractables in pharmaceuticals packaging material by tandem mass spectrometry using Scan-MRM mode

Hemant Kesarkar¹; Sanket Anand Chiplunkar¹; Prashant Hase¹; Durvesh Sawant¹; Aseem Wagle¹; Rahul Dwivedi¹; Satyendra Singh²; Mohit Sharma²; Jitendra Kelkar¹; Pratap Rasam¹, Dheeraj Handique¹

¹Shimadzu Analytical (India) Pvt. Ltd., 1 A/B Rushabh Chambers, Makwana Road, Marol, Andheri (E), Mumbai-400059, Maharashtra, India. ²Shimadzu Analytical (India) Pvt. Ltd., Plot No-2, Sethi Chamber, M O R Land, New Rajender Nagar, New Delhi, Delhi 110060.

1. Introduction

The purpose of Extractable & Leachable (E&L) studies in the regulatory references, is to identify the traces of potential chemical substances. Because these substances may be harmful to patient due to their toxicity, or the activity of drug product may impact because of their presence. Hence, to ensure the safety and efficacy of the drug till its shelf-life period, E&L study plays an important role.

Extractables are organic or inorganic entities which are released from container closure system (CCS) into an extraction medium under laboratory condition. Whereas, leachables are the foreign organic or inorganic chemical substances which are present in the packaged drug product during normal as well as accelerated stability conditions.

Samples of CCS will be used for identification of extractables. However, packaged drug product sample will be used for the evaluation of leachable. Ophthalmic drug product packaging material was used for this study (Figure 1) and the analysis was carried out on Shimadzu triple quadrupole system GCMS-TQ8040 NX equipment (Figure 2).



Figure 1: Ophthalmic packaging system



Figure 2: GCMS-TQ8040 NX

2. Method development

Total 33 commonly identified potential extractable compounds were selected as targeted extractable which includes 9 phthalates, 2 siloxanes, 16 polyaromatic hydrocarbons (PAH) and 6 other additives. Two methods were developed to quantitate the above targeted extractables. Method-1 is used for quantification of targeted components (except PAH's) and unknowns. Method-2 is used for quantification of PAHs only.

2.1 Analytical conditions for method-1:

Instrument used	: GCMS-TQ8040 NX with AOC-20i
Column	: Rxi-5Sil MS 0.25mm I.D. X 30m d.f.=0.25µm (P/N: 221-75954-30)
Injection Mode	: Split (1:10)
Injection Temp.	: 250 °C
Flow control mode	: Linear velocity @ 36.5 cm/sec
Carrier gas	: Helium
Temp. Program	: 40 °C (1 min), 10 °C/min to 160 °C (7 min), 25 °C/min to 280 °C (5.20 min)
Ionization Mode	: Electron Impact
Ion source Temp.	: 250 °C
Acquisition mode	: Scan (from 40 m/z to 500 m/z) & MRM

2.2 Analytical conditions for method-2:

Instrument used	: GCMS-TQ8040 NX with AOC-20i
Column	: SH-I-PAH 0.25 mm I.D. X 30m d.f.=0.1 µm (P/N: 227-36074-01)
Injection Mode	: Split (1:10)
Injection Temp.	: 320 °C
Flow control mode	: Linear velocity @ 52 cm/sec
Carrier gas	: Helium
Temp. Program	: 50 °C (2 min), 12 °C/min to 270 °C (2 min), 30 °C/min to 300 °C (6.67 min)
Ionization Mode	: Electron Impact
Ion source Temp.	: 250 °C
Interface Temp.	: 300 °C

2.3 Standard preparation:

Standard solutions of linear concentration as mentioned in Table-1 and Table-2, were prepared using Ethyl Acetate as diluent.

2.4 Sample preparation:

The sample extraction procedure was carried out as per following flow diagram (Figure 3).




Figure 3: Flow diagram for sample preparation

3. Results and Discussion:

Results of Linearity and precision at LOQ for targeted compounds analyzed in both the methods are reported in below Table-1 and Table-2.

Table-1: Results for targeted compounds analysed by method-1

Sr. No.	Compound Name	Linearity range (ppb)	LOQ (ppb)	Correlation coefficient (r ²)	LOQ precision (%RSD, n=6)	S/N ratio at LOQ
1	Isobutyl Benzene	5-150	5	0.996	5.1	25
2	Decamethyl Cyclopentasiloxane	5-150	5	0.998	0.8	79
3	1,1'-Biphenyl, 2-fluoro-	25-200	25	0.996	3.3	16
4	Dimethyl phthalate	50-250	50	0.995	10.6	33
5	Tetradecamethyl Cycloheptasiloxane	5-150	5	0.999	11.5	23
6	3-tert-Butyl-4-hydroxyanisole (BHA)	25-200	25	0.999	2.1	30
7	Butylated Hydroxytoluene (BHT)	5-200	5	0.995	8.6	21
8	2,6-Di-tert-butyl-4-ethyl phenol (BHEB)	5-200	5	0.994	6.0	52
9	Diethyl Phthalate	25-250	25	0.991	6.6	114
10	3,5-di-tert-Butyl-4-hydroxybenzaldehyde	10-200	10	0.997	1.8	51
11	Di-isobutyl phthalate	10-250	10	0.997	6.6	178
12	Diamyl phthalate	5-200	5	0.994	1.4	43
13	Benzyl butyl phthalate	25-250	25	1.000	5.1	54
14	Dicyclohexyl phthalate	5-250	5	0.986	0.8	60
15	Diphenyl phthalate	10-200	10	0.992	6.3	28
16	Dinonyl Phthalate	5-250	5	0.991	12.7	11
17	Di-n-octyl phthalate	10-250	10	0.987	3.9	37

Table-2: Results for targeted compounds analysed by method-2

Sr. No.	Compound Name	Linearity range (ppb)	Correlation Coefficient (r ²)	LOQ precision (%RSD, n=6)	S/N ratio at LOQ
1	Naphthalene	(LOQ = 1 ppb)	0.998	4.4	66
2	Acenaphthylene		0.999	7.1	23
3	Acenaphthene		0.999	13.0	17
4	Fluorene		0.998	7.4	10
5	Phenanthrene		0.997	12.4	13
6	Anthracene		0.999	7.0	12
7	Fluoranthene		0.998	3.5	50
8	Pyrene		0.998	2.8	42
9	Benzo[a]anthracene		0.999	3.7	12
10	Chrysene		0.999	7.9	17
11	Benzo[b]fluoranthene		0.999	4.9	11
12	Benzo[k]fluoranthene		0.999	7.3	10
13	Benzo[a]pyrene		0.999	4.7	11
14	Indeno[1,2,3-cd]pyrene		0.997	4.1	13
15	Dibenz[a,h]anthracene		0.997	7.5	16
16	Benzo[g,h,i]Perylene		0.998	6.6	13

Representative calibration curve and chromatogram at LOQ level is depicted in Figure 4. Isobutyl Benzene, Decamethyl Cyclopentasiloxane, Dicyclohexyl Phthalate and 2-6-Di-tert-butyl-4-ethyl phenol were analysed in method-1. Naphthalene, Acenaphthalene, Fluorene and Fluoranthene were analysed in method-2.




Figure 4: a) Calibration curve, b) Chromatograms of LOQ level

Extraction profile from aqueous experiments

Aqueous sample extraction experiments covers wide pH range. This provides suitable medium for potential extractables to get extracted from packaging material when exposed to acidic or basic medium. In this CCS, extractables were not detected in scan mode. However, two targeted extractable i.e. Decamethyl Cyclopentasiloxane and Dinonyl Phthalate were quantified in MRM mode (method-1).




Figure 5: Overlaid total ion chromatograms for aqueous extracted experiment samples

Extraction profile from organic-reflux experiments:

Sample was refluxed with organic solvents such as Ethanol, Hexane, Isopropanol and Dichloromethane, with wide range of polarities. Extractables observed in this experiment are depicted in Figure 6.




Figure 6: TIC overlay for organic-reflux experimental samples, highlighting specific unknown extractables only detected in Experiment-B (Ethanol reflux)

Targeted compounds which were detected in aqueous and organic experiments were quantified using standard linearity plot and reported in Table-3. About 88 unknown impurities were identified using scan data and semi-quantified based on average slope method. However, these are not reported here.

Table-3: Content of targeted extractable compounds

Name of the Targeted comp. (Method-1)	Content (ppb)	Expt.	Name of the Targeted comp. (Method-2)	Content (ppb)	Expt.
Isobutyl Benzene	34	E	Naphthalene	62	B
Decamethyl Cyclopentasiloxane	3	A (pH 7)	Acenaphthylene	1	B
1,1'-Biphenyl, 2-fluoro-	12	E	Acenaphthene	10	E
Dimethyl phthalate	30	E	Fluorene	13	C
Tetradecamethyl Cycloheptasiloxane	3	C	Phenanthrene	37	E
3-tert-Butyl-4-hydroxyanisole (BHA)	24	B	Anthracene	9	B
Butylated Hydroxytoluene (BHT)	28	B	Fluoranthene	12	E
2-6-Di-tert-butyl-4-ethyl phenol (BHEB)	ND	NA	Pyrene	8	E
Diethyl Phthalate	236	E	Benzo[a]anthracene	0	B
3,5-di-tert-Butyl-4-hydroxybenzaldehyde	209	C	Chrysene	1	D
Di-isobutyl phthalate	1939	E	Benzo[b]fluoranthene	ND	NA
Diamyl phthalate	9	D	Benzo[k]fluoranthene	ND	NA
Benzyl butyl phthalate	12	B	Benzo[a]pyrene	ND	NA
Dicyclohexyl phthalate	350	B	Indeno[1,2,3-cd]pyrene	ND	NA
Diphenyl phthalate	6	D	Dibenz[a,h]anthracene	ND	NA
Dinonyl Phthalate	14	A (pH 7)	Benzo[g,h,i]Perylene	ND	NA
Di-n-octyl phthalate	9	B			

ND= Not Detected, NA= Not Applicable

4. Conclusion

- The migration of extractables from the CCS includes a wide range of compounds. Highly sensitive and robust techniques is required for detection at trace level. Also, Data need to be acquired in both SCAN and SIM/MRM modes to detect known and unknown extractables simultaneously without affecting sensitivity.
- Shimadzu's GCMS TQ-8040 NX triple quadrupole ultra fast mass spectrometer offers a SCAN speed of 20,000 amu/sec, 800 MRM transitions per second along with Advanced Scan Speed Protocol (ASSP) function, facilitates detection and identification extractables even at trace levels.