




Maximizing Resolution and Selectivity: Superficially Porous Column Chromatography Options

Anne Mack
Application Scientist

Stephen Luke
*LC Columns Product
Manager*

Agenda

- Superficially porous particle (SPP) specifications and benefits
- Method development challenges and objectives
- Method development with selectivity
 - Bonded phase
 - Mobile phase pH
- How to ensure the best performance from your SPP column
- When you need more help choosing a column...



SUPERFICIALLY POROUS PARTICLES SPECIFICATIONS AND BENEFITS

Current Status of Superficially Porous Particles

Status in 2000 **Status in 2010**

Status in 2015

	# of Vendors	# of Vendors	# of Vendors	Particle size range	Pore size range	# of Phase chemistries
Small molecules	0	3	16	1.3 μm to 5 μm	80Å to 120Å	>12 chemistries
Large molecules	1 (Agilent)	1 (Agilent)	9	2.6 μm to 5 μm	160Å to 450Å	>8 chemistries





Bell, LC-GC, 2015 June
Majors, LC-GC, 2014 Nov

Superficially Porous Column Technology

Poroshell 120 2.7 μm

- $d_p = 2.7 \mu\text{m}$
- Particles
 - 1.7 μm solid core
 - 0.5 μm diffusion path
 - 2.7 μm total diameter
- Efficiency (N) \approx 90% of sub-2 μm
- $N \approx 2X$ 3.5 μm (totally porous)
- Pressure \approx 40-50% of sub-2 μm
- 2 μm frit to reduce clogging
- $P_{\text{limit}} = 600$ bar for HPLC or UHPLC






Comparing Efficiency and Pressure with Different Types of Columns

Particle Size/Type	Pressure	Efficiency	LC Compatibility
3.5 μm Totally Porous	123 bar	7,800	All 400 bar instruments
2.7 μm Poroshell 120	180 bar	12,000	All LCs/UHPLCs (up to 600 bar)
1.8 μm Totally Porous	285 bar	12,500	All LCs/UHPLCs (up to 1200 bar)


Columns: 4.6 x 50mm, Mobile Phase: 60% ACN:40% Water Flow Rate: 2 mL/min

USP Method for Naproxen Tablets

Method Requirement $N > 4000$, R_s better than 11.5




Larger Diameter SPP



Poroshell 120 2.7 µm

SA = 120 m²/g

Pore size = 120-140Å



Poroshell 120 4 µm


SA = 120 m²/g

Pore size = 120-140Å

- Offers nearly 2X the performance of traditional 5 µm columns with the easy drop-in replacement for current methods


2.7 and 4 μm Poroshell 120 Have Similar Selectivity for Easy Method Transfer and Scaling

Selectivity Comparisons with 4.6 x 50 mm Columns
5-95% CH₃CN in 2 min with 0.1% formic acid, 2 mL/min



- ◆ RT on 2.7 μm Poroshell 120 EC-C18 (min)
- Linear (RT on 2.7 μm Poroshell 120 EC-C18 (min))



Comparison between Poroshell 120 2.7 μm and 4 μm for 4-quinolones in Milk



Norfloxacin
Ofloxacin
Ciprofloxacin
Pefloxacin
Lomefloxacin
Danofloxain
Enrofloxacin
Sarafloxacin
Difloxacin
Oxilinic acid
Flumequine

- $PW_{1/2}$ with 4 μm column increased by 30% compared to 2.7 μm
- Pressure on 4 μm decreased by 45% compare to 2.7 μm column. It is more suitable to use on a < 400 bar LC, while 2.7 μm column is suitable for 600 bar LC.

Long Lifetime with Poroshell 120 2.7 μm Column >1800 Injections at 550 bar - No Performance Change



A: 0.2% HCOOH in H₂O, B: 0.2% HCOOH in CH₃CN

0.833 mL/min

Time	0.00	1.25	2.50
%B	10	15	27

%B

40 °C

Agilent Poroshell 120 SB-C18, 2.1 x 100 mm, 2.7 μm

Sig=210,4nm, Ref=Off

2- μL , 3-mm micro flow cell (PN G1315-60024)

Sample: 2 μL of freshly brewed green tea

(brewed from a commercial tea bag in 6 oz of initially boiling water for six minutes)






Other Considerations when Selecting a Column

- **Robustness and batch-to-batch reproducibility**

2010



2012



Beverage Additives

METHOD DEVELOPMENT CHALLENGES AND OBJECTIVES

Challenges In Method Development

- Worldwide method transfers
 - Instruments and configurations differ from lab to lab
 - More contract labs
- Many chromatographic mode choices
 - RP, SEC, IEX, HILIC, Chiral...
 - RP most common
- Too many columns to choose from
 - Endcapped C18 is a good starting point



Defining the Objective

- How complex is the sample?
- Is high efficiency important?
- Is speed important?
- What are the instrument limitations?
- What is the skillset of the operator?



Examples of Common Separation Goals and Method Performance Criteria

Good System Suitability Parameters

- Resolution: ≥ 2
- Peak shape: USP T_f close to 1 (< 2)
- Injection Repeatability: areas, T_f , etc. (RSD 0.1 - 0.25%)
- Absolute retention factors: $1 < k < 10$
- Relative Retention: α or k_2/k_1
- Signal-to-Noise Ratio: > 10

Method Performance Criteria

- Accuracy
- Precision
 - Ruggedness
 - Robustness
- Selectivity/Specificity
- Linearity
- Range
- Quantitation Limit (LOQ, 10x S/N)
- Detection Limit (LOD, 3x S/N)


AVOID THESE for System Suitability Criteria:

*Column efficiency (theoretical plates)
& Absolute retention time*

These inhibit the ability to speed up your method in the future!

Where to Begin?

$$R_s = N^{1/2}/4 \cdot (\alpha-1)/\alpha \cdot k'/(k'+1)$$



Selectivity Impacts Resolution Most

- Change bonded phase
 - Change mobile phase
- } Typical Method Development Parameters

Plates are easiest to increase

METHOD DEVELOPMENT WITH SELECTIVITY – BONDED PHASE

Why is Changing the Bonded Phase Effective?

- Different interactions for polar and non-polar compounds.
- Exploit other interactions with bonded phase (e.g., pi-pi)
- Changing the bonded phase can improve selectivity/resolution, reduce analysis time
- Having numerous different bonded phases available on the same particle makes development easier
 - Fast SPP methods make development faster

Poroshell 120 Column Chemistries

Multiple bonded phases for flexibility in method development



Poroshell 120 EC-C18 and C8

- Robust endcapped C18 for best peak shape at pH 2-9

Poroshell 120 StableBond C18 and C8

- Robust chemistries for pH<2

Poroshell HPH-C18 and HPH-C8

- Long lifetime at high pH

Poroshell 120 Phenyl-Hexyl

- Excellent choice for pi-pi interactions
- Selectivity similar to phenyl, diphenyl, or other phenyl-hexyl columns

Poroshell 120 SB-Aq

- Proprietary bonding phase is an excellent choice for polar analytes

Poroshell 120 Bonus-RP

- Embedded polar group provides unique selectivity for polar compounds

Poroshell 120 EC-CN

- Flexible endcapped CN chemistry with Normal and Reversed Phase character

Poroshell 120 HILIC

- Bare silica HILIC for use in hydrophilic interaction chromatography of polar molecules

Poroshell 120 PFP

- Perfluorophenyl chemistry

HSM a Way to Look at Column Orthogonality

Poroshell 120	H	S*	A	B	C	k'	F
EC-C18	1.020	0.008	-0.130	-0.004	0.161	6.920	0
EC-C8	0.877	0.011	-0.232	0.023	0.127	4.840	6
EC-CN	0.421	-0.057	-0.476	0.002	0.045	0.950	17
Phenyl-Hexyl	0.752	-0.083	-0.394	0.018	0.136	3.590	13
Bonus RP	0.686	-0.030	-0.573	0.180	-0.670	3.980	75
SB-C18	0.956	-0.041	0.168	0.025	0.210	5.440	12
SB-C8	0.726	-0.087	0.068	0.044	0.087	3.560	15
SB-Aq	0.581	-0.120	-0.133	0.051	-0.014	2.150	22
PFP	0.630	-0.520	-0.520	0.430	-0.110	2.300	85

Data provide by Dwight Stoll


F_s factor describes the similarity of two columns. A small F_s indicates that two columns are very similar, while a large factor indicates that two columns are very different. Calculated according to the following equation:

$$F_s = \left\{ [12.5(H_2 - H_1)]^2 + [100(S_2^* - S_1^*)]^2 + [30(A_2 - A_1)]^2 + [143(B_2 - B_1)]^2 + [83(C_2 - C_1)]^2 \right\}^{\frac{1}{2}}$$

Further details at:

<http://www.hplccolumns.org>

HSM Data for Poroshell 120




Poroshell 120	F
EC-C18	0
EC-C8	6
EC-CN	17
Phenyl-Hexyl	13
Bonus RP	75
SB-C18	12
SB-C8	15
SB-Aq	22
PFP	85

HSM	
Hydrophobicity	H
Steric interaction	S*
Hydrogen-bond acidity	A
Hydrogen-bond basicity	B
Ion-exchange capacity	C
K ethyl benzene	EB

$$F_s = \left\{ [12.5(H_2 - H_1)]^2 + [100(S_2^* - S_1^*)]^2 + [30(A_2 - A_1)]^2 + [143(B_2 - B_1)]^2 + [83(C_2 - C_1)]^2 \right\}^{\frac{1}{2}}$$

Separation of 8 Steroids with Methanol Gradient

Best Resolution of all analytes with Poroshell 120 Phenyl-Hexyl





1. Hydrocortisone, 2. B Estradiole, 3. Androstadiene 3,17 dione, 4. Testosterone, 5. Ethyestradiene, 6. Estrone, 7. Norethindone acetate, 8. Progesterone

40-80 % Methanol/14 min, DAD 260, 80 nm 0.4 ml/min, 2.1 x 100 mm 40 C 0.1% Formic Acid in Water and Methanol, Agilent 1260 Method Development Solution

Poroshell 120 Phenyl-Hexyl vs Poroshell 120 EC-C18

Phenyl-Hexyl alternative selectivity to C18


- recommended for aromatics, especially with methanol
- compatible with highly aqueous mobile for polar compounds.



Buffer-Organic Gradients

Beta Blockers with Methanol Gradient

Best Resolution of all analytes with Poroshell Bonus-RP





1. Atenolol, 2. Pindolol, 3. Naldolol, 4. Metoprolol, 5. Acebutolol, 6. Propranolol, 7. Alprenolol

10-70 % Methanol/12 min, DAD 260 nm 0.35 ml/min, 2.1 x 100 mm 40 C 10 mM pH 3.8 Ammonium Formate Buffer and Methanol

Poroshell 120 Bonus-RP vs Poroshell 120 EC-C18


- Embedded polar group gives unique selectivity for polar compounds compared to C18.



Buffer-Organic Gradients

NSAID Separation with a Methanol Gradient

Best Resolution of all analytes with Poroshell 120 PFP




1. APAP, 2. Phenacetin, 3. Piroxicam, 4. Tolmetin, 5. Ketoprofen,
6. Naproxen, 7. Sulindac, 8. Diclofenac, 9. Diflunisal

Time	% Organic
0	8
6	100
7	100
8	8

2mL/min 254 nm

METHOD DEVELOPMENT WITH SELECTIVITY – MOBILE PHASE PH

When Does pH Affect Selectivity and Resolution? Compound Type Comparison




- **Ionizable compounds – acids and bases can change retention and selectivity most with changes in pH**

Change in Retention with pH for Ionizable Compounds is Key to Method Development

- Non-charged analytes have better retention (i.e. acids at low pH and bases at high pH)
- Silanols on silica ionize at mid-pH, increasing retention of basic analytes (i.e possible ion-exchange interactions)
- Choose mobile phase pH to optimize retention and selectivity during method development
- Ensure that your column is compatible with and stable in the mobile phase pH you select

Poroshell HPH-C18 vs Poroshell 120 EC-C18




- Acid
- ▲ Base
- × Neutral

MP-A 10mM Ammonium Formate/water
adj. to pH3 or 7 using Formic Acid
MP-B ACN or MeOH
Flow – 0.42 mL/min
Column Temp. Ambient
1 uL injection
Detection 254 nm

Gradient:	time, min	%B
	0	5
	4	95
	5	95
	6	5
	7	stop run


Use of Varied pH can Help Build Separations that are Very Different (poorly correlated)



Column: Poroshell HPH-C18 2.7 μm

Change in Retention with pH for Ionizable Compounds is Compound Dependent

**More retention for non-charged analytes
(i.e. acids at low pH and bases at high pH)**



**Column: Poroshell
HPH-C18 2.7 μ m**


Mobile Phase: 45% Methanol, 55% 20 mM Phosphate Buffer

Selectivity Can be Controlled by Changing pH

1. Procainamide
2. Caffeine
3. Acetyl Salicylic Acid
4. Hexanophenone Deg.
5. Dipyrimadole
6. Diltiazem
7. Diflunisal
8. Hexanophenone

Time	% Buffer	% MeCN
0	10	90
5	90	10
7	10	90
2 ml/min		254 nm

Poroshell HPH-C18 4.6 x 50 mm, 2.7 μ m




HOW TO ENSURE THE BEST PERFORMANCE FROM YOUR SPP COLUMN

Benefits of Installing a Guard Column

Accelerated Lifetime Test

Similac sample (milk substitute diluted 300:1) containing 2 sulfa drugs
Peak width change indicating column failure



No Guard

- Column failure @ inj. 70;
new column required

With Guard




- Guard failure @ inj. 80;
guard replaced
- Same column used throughout analysis

By installing a guard column when using dirtier samples, one can extend the life of their column, and utilize more inexpensive guard columns rather than analytical column replacements

Agilent A-Line Fittings


Importance of the Spring Loaded Feature

Most commonly used fittings in UHPLC are non-adjustable 2-piece or 3-piece metallic fittings. Since different manufacturers of column hardware have different design in column end fittings, as shown in Figure 1, a new set of tubing and fittings needs to be installed for every brand of column to guarantee that the stem length, namely the length between the bottom of the ferrule and the end of tubing, fits the column end fitting.




The spring-loaded design constantly pushes the tubing against the receiving port, delivering a reproducible connection with no dead volume for consistent chromatographic performance

Stem length is adjustable through the spring, which makes the fitting compatible with all types of LC columns.



Data Collection Rate



1 ul QC Mix, Uracil, Phenol (k=0.5), 4-Chloronitrobenzene(k=2), Napthalene(k=3.8)
55% MeCN 45 % Water 0.55 ml/min micro flow cell

WHEN YOU NEED MORE HELP CHOOSING A COLUMN...

On-Line Tool “The Navigator”

A Column and Sample Prep Selection Tool

BACK ▶

GET YOUR COLUMN SELECTION GUIDE ▶

GET HELP ▶

LC COLUMN
and SAMPLE PREP
NAVIGATOR

Let us help you find the best Agilent column and sample prep products for your application.

BEGIN ▶

The Measure of Confidence


Agilent Technologies

Contact tech support | Privacy Statement | Terms of Use | United States Home | © Agilent 2000-2012


<http://navigator.chem.agilent.com>

Literature on Poroshell 120 Columns


- There are continuous updates and additions to Poroshell 120 Columns Literature.
- Brochures, app notes, flyers and other documents are updated and added often!



5990-5951EN



5990-5951EN



5991-4893EN

Phase Selection Wall Chart

5991-6240EN

AGILENT SMALL MOLECULE LC COLUMNS OVERVIEW:

A FAMILY OF PHASE CHOICES TO PERFECT EVERY SEPARATION

Best all around—exceptional peak shape, efficiency, resolution, and lifetime



Make every LC and LC/MS in your lab work even harder



You can trust the ZORBAX name for superior reproducibility and long-term stability

Start here for method development
Poroshell 120 EC-C18
Poroshell 120 EC-C8
Poroshell 120 Phenyl-Hexyl

2.7 µm and 4 µm, 600 bar

ZORBAX Eclipse Plus
 RRHD: 1.8 µm, stable to 1200 bar
 RRHT: 1.8 µm, 600 bar

Great for Method Development
C18 (USP L1)
C8 (USP L7)
Phenyl-Hexyl (USP L11)
PAA (USP L1)



High performance and excellent peak shape with acids, bases and neutrals.

Best for low pH mobile phases—great for method development

Poroshell 120 SB-C18
Poroshell 120 SB-C8

2.7 µm, 600 bar

ZORBAX StableBond
 RRHD: 1.8 µm, stable to 1200 bar
 RRHT: 1.8 µm, 600 bar



High performance with acids, bases, and neutrals with superior lifetime at low pH.

High performance over a wide pH range

Poroshell 120 EC-CN

2.7 µm, 600 bar


ZORBAX Eclipse XDB
 RRHD: 1.8 µm, stable to 1200 bar
 RRHT: 1.8 µm, 600 bar



Good peak shape for basic, acidic, and neutral compounds with high performance over a wide pH range (pH 3-9). oXtra Dense Bonding and double endcapping help give this column a long lifetime.

High pH
 Best lifetime and peak shape at high pH

Poroshell 120 HPH-C18
Poroshell 120 HPH-C8



2.7 µm, 600 bar

ZORBAX Extend-C18
 RRHD: 1.8 µm, stable to 1200 bar
 RRHT: 1.8 µm, 600 bar



High efficiency and long life at high pH—up to pH 11.5. Improve retention, resolution and peak shape of basic compounds. High sensitivity for LC/MS separations of peptides. Unique intramolecular bonding and double endcapping provides high pH stability.

Polar Compounds
 Alternative selectivity to alkyl, phenyl, cyano phases

Poroshell 120 Bonus-RP
Poroshell 120 PFP

2.7 µm and 4 µm (PFP), 600 bar

ZORBAX Bonus-RP
 RRHD: 1.8 µm, stable to 1200 bar
 RRHT: 1.8 µm, 600 bar



Polar-embedded to improve peak shapes for basic compounds at low and mid pH.

High sensitivity for LC/MS applications and recommended for EPA 1604

Poroshell 120 SB-Aq
Poroshell 120 HILIC

2.7 µm and 4 µm (HILIC), 600 bar

ZORBAX SB-AQ
 RRHD: 1.8 µm, stable to 1200 bar
 RRHT: 1.8 µm, 600 bar



Proprietary phase best for polar compounds and high aqueous conditions. Exceptional lifetime at low pH—no endcapping.

ZORBAX HILIC
 RRHD: 1.8 µm, stable to 1200 bar
 HILIC Plus is a HILIC column based on Eclipse Plus silica for excellent peak shapes.

Learn more at
www.agilent.com/chem/discoverporoshell


Choose the right LC column for your analysis with the LC Column and Sample Prep Navigator at
www.navigator.chem.agilent.com

The Measure of Confidence



Agilent is Here to Help

Visit www.agilent.com/chem/cstechsupport



The screenshot shows the Agilent Technologies website's technical support resources page. The header includes the Agilent logo, navigation links for Register, Login, Global Sites, and CONTACT US, and a search bar. A secondary navigation bar contains links for PRODUCTS & SERVICES, SOLUTIONS, TECHNICAL SUPPORT, LIBRARY, and TRAINING & EVENTS, along with icons for myAgilent, CART (0), and QUICK ORDER. The main content area features a breadcrumb trail (Home > Promotions > Technical Support Resources) and a large blue banner with the text "AGILENT TECHNICAL SUPPORT RESOURCES" and a chromatogram graphic. Below the banner, a paragraph states Agilent's commitment to helping users succeed with their applications. A "Helpful tip" section provides advice on providing detailed information when contacting support. The "TECHNICAL SUPPORT CONTACTS" section lists contact information for various regions: US, Canada; Mexico, Brazil, Argentina and other countries in S. America; Europe, Middle East and Africa; India; Singapore, Malaysia, Australia, Korea, the Philippines; Japan; and China. An image of a technical support representative is shown on the right side of the page.

Agilent Technologies Register Login Global Sites CONTACT US Search in

PRODUCTS & SERVICES SOLUTIONS TECHNICAL SUPPORT LIBRARY TRAINING & EVENTS myAgilent CART (0) QUICK ORDER

Home > Promotions > Technical Support Resources

Technical Support Resources

AGILENT TECHNICAL SUPPORT RESOURCES

Agilent is committed to helping you succeed with your application. Send us an email with your questions and a technical support representative will contact you as soon as possible, usually within the day, provided you are emailing during business hours in your region. However, please allow up to 48 hours for a reply.

Helpful tip: Provide as much information as possible regarding the instrument you are using, your sample type, the column you are using now and what you are trying to achieve.

TECHNICAL SUPPORT CONTACTS

US, Canada – please call 800-227-9770 or email LC-column-support@agilent.com for LC columns support and sop-support@agilent.com for sample prep support

Mexico, Brazil, Argentina and other countries in S. America – please contact your preferred distributor (see www.agilent.com/chem/contactus for more information)

Europe, Middle East and Africa, please contact your distributor (see www.agilent.com/chem/contactus) or email esmc-es@agilent.com

India – please send an email to columns_helpdesk@agilent.com

Singapore, Malaysia, Australia, Korea, the Philippines: contact your local distributor (see <http://www.chem.agilent.com/en-US/Contact-Us/Pages/Singapore.aspx>) or send an email to ccc-smt@agilent.com


Japan – please call 0120-477-111 or email email_japan@agilent.com

China – please call 800-820-3278-4

Technical Support Helpdesk Representative

Conclusions

- Resolution is a common goal during method development
- Selectivity is a main driver of resolution
- Superficially porous particle columns (e.g., Poroshell 120) offer...
 - 12 chemistries, including high pH stable options
 - Faster method development
 - Higher sample throughput
 - Maintain method ruggedness
 - Compatibility with any LC system



Thank You

Questions?

