

# EPA Method 8270 Semi Volatile Organic Compounds Analysis on a New Benchtop Time-of-Flight Mass Spectrometer

Jonathan D. Byer, Todd Richards, Lorne M. Fell, Joseph E. Binkley | LECO Corporation, Saint Joseph, MI USA

## Instrument Parameters and System Performance Evaluation


### GC-TOFMS (Pegasus® BT) Conditions

Gas Chromatograph	LECO L-Pal3 Autosampler with Agilent 7890B GC
Injection:	GC inlet, split 20:1 @ 270 °C
Carrier Gas:	He @ 1.4 ml/min, Constant Flow
Column:	ZB-SemiVolatiles, 30 m x 0.25 mm ID x 0.25 µm df (Phenomenex, Torrance, CA)
Oven Program:	70 °C (1 min), to 285 °C @ 28 °C/min, to 305 °C @ 3 °C/min, to 320 °C @ 30 °C/min (1 min); Total time: 16.5 min
Transfer Line :	300 °C
Mass Spectrometer	LECO Pegasus BT
Ion Source Temperature:	250 °C
Mass Range:	30-650 m/z
Acquisition Rate:	10 spectra/s

### DFTPP tuning results for 10 replicate injections shows the Pegasus BT routinely passes

Mass	Ion Abundance Criteria	Mean	STDEV	Min	Max	Criteria
51	10-80% of Base Peak	64.5	0.25	64.1	64.9	PASS
68	< 2% of mass 69	1.93	0.04	1.87	1.98	PASS
70	< 2% of mass 69	0.38	0.01	0.37	0.40	PASS
127	10-80% of Base Peak	52.1	0.18	51.9	52.4	PASS
197	< 2% of mass 69	0.22	0.03	0.19	0.28	PASS
198	Base Peak, or > 50% of Mass 442		Base Peak			PASS
199	5-9% of mass 198	6.73	0.04	6.68	6.77	PASS
275	10-60% of Base Peak	23.0	0.08	22.9	23.1	PASS
365	> 1% of mass 198	1.90	0.02	1.86	1.92	PASS
441	Present but < 24% of mass 442	13.2	0.03	13.15	13.25	PASS
442	Base Peak, or > 50% of mass 198	76.5	0.3	76.1	76.9	PASS
443	15-24% of mass 442	19.5	0.08	19.4	19.6	PASS


Caliper - sample "GCMS Tune Mix 50 ng\_uL Split 20", 7.55214 min - 7.53208 min, Height (Abundance)






Library Hit - Similarity: 907 - Library: WR10 - Phosphine, bis(pentafluorophenyl)phenyl-, Abundance



## Quantitation in ChromaTOF® 5.0: Initial Calibration




Total ion chromatogram (TIC) of a 5 ppm standard containing 90 components.



Extracted ion chromatogram (XIC) of the quant mass for hexachlorobenzene (HCB) from 0.05 to 50 ppm and its corresponding calibration curve, with the Analytical Ion Chromatogram (AIC) of a 1 ppm standard split 20:1 indicating where HCB elutes in the chromatogram, and the profile mass spectrum of the molecular isotope cluster.

## Calibration and Quantitation on LECO Pegasus BT

Quantitation using ChromaTOF 5.0 to evaluate calibration data for 8270 compounds.



### EPA Method 8270 Calibration Curve Data for 84 Compounds from 0.05 to 50 ppm.

Name	Type	R.T. (min)	Avg RF	% RSD RF	r	Name	Type	R.T. (min)	Avg RF	% RSD RF	r
N-Nitrosodimethylamine	Analyte	1.837	0.577	8.425	-----	2,6-Dinitrotoluene	Analyte	5.819	0.171	7.38	-----
Pyridine*	Analyte	1.975	1.174	7.698	-----	1,2-Dinitrobenzene	Analyte	5.864	0.072	8.875	-----
2-Fluorophenol	Surrogate	2.475	0.367	5.14	-----	Acenaphthylene	Analyte	5.889	2.2	7.71	-----
Phenol-d6	Surrogate	3.090	0.581	3.92	-----	3-Nitroaniline	Analyte	5.965	0.173	6.112	-----
Phenol*	Analyte	3.100	0.314	7.758	-----	Acenaphthene	Analyte	6.040	1.341	2.539	-----
Aniline	Analyte	3.163	1.304	3.669	-----	Bis(2-chloroethyl) ether	Analyte	6.057	0.04	10.238	-----
Bis(2-chloroethyl) ether	Analyte	3.198	1.025	5.732	-----	4-Nitrophenol	Analyte	6.095	0.112	5.813	-----
2-Chlorophenol-D4	Surrogate	3.242	0.373	3.895	-----	2,4-Dinitrophenol†	Analyte	6.170	0.166	7.711	-----
2-Chlorophenol	Analyte	3.253	0.398	7.524	-----	4-Nitrotoluene	Analyte	6.192	2.148	1.868	-----
1,3-Dichlorobenzene	Analyte	3.380	1.032	6.147	-----	Dibenzofuran	Analyte	6.255	0.08	10.589	-----
1,4-Dichlorobenzene	Analyte	3.437	1.075	7.583	-----	2,3,5,6-Tetrachlorophenol	Analyte	6.292	0.089	11.549	-----
Benzyl alcohol	Analyte	3.513	0.234	5.162	-----	Diethyl Phthalate	Analyte	6.385	0.483	8.33	-----
1,2-Dichlorobenzene-D4	Surrogate	3.547	0.98	0.914	-----	4-Chlorophenyl phenyl ether	Analyte	6.490	0.671	3.345	-----
1,2-Dichlorobenzene	Analyte	3.558	1.01	7.825	-----	Fluorene	Analyte	6.492	1.249	6.793	-----
2-Methylphenol	Analyte	3.587	0.423	2.733	-----	4-Nitroaniline	Analyte	6.497	0.186	5.994	-----
Bis(2-chloro-1-methylethyl)ether	Analyte	3.625	1.943	5.703	-----	4,6-Dinitro-2-methylphenol	Analyte	6.525	0.052	8.015	-----
3/4-Methylphenol	Analyte	3.712	0.378	3.9	-----	Diphenylamine	Analyte	6.592	1.097	4.717	-----
N-Nitrosodi-n-propylamine	Analyte	3.729	0.508	7.917	-----	Azobenzene	Analyte	6.630	0.289	7.106	-----
Hexachloroethane	Analyte	3.837	0.429	5.399	-----	2,4,6-Tribromophenol	Surrogate	6.700	0.079	6.93	-----
Nitrobenzene-D5	Surrogate	3.867	0.522	3.019	-----	4-Bromophenyl phenyl ether	Analyte	6.920	0.438	8.076	-----
Nitrobenzene	Analyte	3.882	0.216	6.85	-----	Hexachlorobenzene	Analyte	6.970	0.336	2.276	-----
Isophorone	Analyte	4.079	1.03	6.196	-----	Pentachlorophenol*	Analyte	7.140	0.038	11.186	-----
2-Nitrophenol	Analyte	4.150	0.098	5.204	-----	Phenanthrene	Analyte	7.337	1.335	4.781	-----
2,4-Dimethylphenol	Analyte	4.167	0.289	5.908	-----	Anthracene	Analyte	7.382	1.828	3.421	-----
Bis(2-chloroethoxy)methane	Analyte	4.257	0.8	3.017	-----	Carbazole	Analyte	7.519	1.436	3.42	-----
2,4-Dichlorophenol	Analyte	4.345	0.19	5.907	-----	Dibutyl phthalate	Analyte	7.820	0.357	10.856	-----
1,2,4-Trichlorobenzene	Analyte	4.427	0.335	2.007	-----	Fluoranthene	Analyte	8.379	1.479	5.63	-----
Naphthalene	Analyte	4.499	1.146	4.684	-----	Pyrene	Analyte	8.579	1.542	3.832	-----
4-Chloroaniline	Analyte	4.537	0.772	7.899	-----	p-Terphenyl-d14	Surrogate	8.712	0.52	8.99	-----
Hexachlorobutadiene	Analyte	4.600	0.156	1.946	-----	Benzyl butyl phthalate	Analyte	9.164	0.141	9.953	-----
4-Chloro-3-methylphenol	Analyte	4.944	0.251	5.033	-----	Bis(2-ethylhexyl) adipate	Analyte	9.234	0.131	9.617	-----
2-Methylnaphthalene	Analyte	5.097	0.588	5.282	-----	Benz[a]anthracene	Analyte	9.767		0.9987	-----
1-Methylnaphthalene	Analyte	5.184	0.578	5.149	-----	Bis(2-ethylhexyl)phthalate*	Analyte	9.791	0.239	10.425	-----
Hexachlorocyclopentadiene	Analyte	5.232	0.326	2.02	-----	Chrysene	Analyte	9.811	1.379	4.165	-----
2,4,6-Trichlorophenol	Analyte	5.337	0.089	6.534	-----	Di-n-octyl phthalate	Analyte	10.676	0.279	7.858	-----
2,4,5-Trichlorophenol	Analyte	5.364	0.101	9.49	-----	Benz[b]fluoranthene	Analyte	11.264	1.041	6.505	-----
2-Fluorobiphenyl	Surrogate	5.415	1.771	3.367	-----	Benz[k]fluoranthene	Analyte	11.309	1.114	7.535	-----
2-Chloronaphthalene	Analyte	5.524	1.554	3.336	-----	Benz[a]pyrene	Analyte	11.821	0.896	7.607	-----
2-Nitroaniline	Analyte	5.607	0.127	7.753	-----	Indeno[1,2,3-cd]pyrene	Analyte	14.094	0.633	8.146	-----
1,4-Dinitrobenzene	Analyte	5.729	0.056	7.153	-----	Dimethyl phthalate	Analyte	14.164	0.704	7.091	-----
Dimethyl phthalate	Analyte	5.767	0.552	6.982	-----	1,3-Dinitrobenzene	Analyte	14.691	0.926	7.95	-----
1,3-Dinitrobenzene	Analyte	5.790	0.104	7.962	-----						

\*The 0.05 ppm standard was not included in the calibration curve because the %D was greater than 20%. † The detection limit was 0.2 ppm for 2,4-Dinitrophenol.

## Summary

- The Pegasus BT is a powerful instrument for target and non-target analysis in the same data acquisition.
- Pegasus BT meets or exceeds the EPA Method 8270 requirements for instrument performance.
- With TOFMS, all of the masses are collected all the time without compromising sensitivity and without mass spectral skewing, resulting in unsurpassed spectral fidelity for library database searching and identifying analytes with confidence.

