

Comprehensive Characterization of Natural Rubber Samples using Thermal Field-Flow-Fractionation coupled with MALS and Triple Detection

Dr. Gerhard Heinzmann, Postnova Analytics GmbH, Germany

**International Symposium on GPC/SEC and Related Techniques,
Amsterdam, The Netherlands, September 26th-29th, 2016**

- **FFF Separation Platform**
- **Thermal FFF (TF3)**
- **TF3-MALS and TF3-Triple/Tetra/Penta Detection**
- **Application Examples**

Separation and Detection Platform



Light Scattering
Detector
PN3621 MALS

Particle
Size Rg /
Molar Mass



AF2000 Flow FFF

Size Separation



CF2000 Centri FFF



SC2000 GPC



TF2000 Thermal FFF

Particle
Size Rh



Malvern
Zetasizer

Intrinsic
Viscosity /
Branching



Viscometer detector
PN3310

Concentration



RI detector PN3150

UV detector PN3211

DAD detector PN3241


Fluori detector PN3411

ELS detector PN3510



Agilent 7900 ICP-MS

Separation Mechanism



- Separation in a narrow ribbon-like channel
- Laminar flow inside the channel
- External field perpendicular to the solvent flow


Thermal FFF – Principle

TF2000: Separation

$$V_E \sim D\Delta T/D$$

Thermal Diffusion Coefficient
-> Depends on chemical composition

Temperature Field



Applications

- Rubbers Polymers
- Gels / Latexes
- Cross-linked Polymers

- Thermal gradient up to $\Delta 120^\circ\text{C}$
- Separation kDa up to several MDa
- Analysis time, 10 – 120min (no upper limit)
- Separation depends on **Size** and **Chemical Composition** (“2 Dimensional”?)



Thermal FFF vs. GPC/SEC




Advantages of Thermal FFF compared to GPC/SEC:

No shear degradation in porous columns / frits

No filtration of gel material, no clogging of columns

No interaction with stationary phase

Excellent separation for high MW → extended 'Exclusion Limit'



FFF – Triple Detection

- MALS
- Viscometer
- Concentration Detector(s)

Setup of TF2000 Thermal FFF with Triple Detection (Penta Detection) MALS / Viscometer / 3x Concentration Detector



Detector Setup: UV => MALS => Visc => RI => ELSD

Conditions: Solvent = THF, Flow rate = 0.3 mL/min

Alternative Setup: UV => MALS => Split: Visc => Waste and RI => ELSD => Waste

$$\mathbf{LS}_{\text{ signal}} = \mathbf{K}_{\text{LS}} * (\mathbf{dn/dc})^2 * \mathbf{Conc} * \mathbf{Mw}$$

$$\mathbf{RI}_{\text{ signal}} = \mathbf{K}_{\text{RI}} * (\mathbf{dn/dc}) * \mathbf{Conc}$$


$$\mathbf{UV}_{\text{ signal}} = \mathbf{K}_{\text{UV}} * \boldsymbol{\varepsilon} * \mathbf{Conc}$$

$$\mathbf{Visco}_{\text{ signal}} = \mathbf{K}_{\text{Visco}} * [\boldsymbol{\eta}] * \mathbf{Conc}$$

- **MW** from MALS + RI
- **Rg** from MALS
- **IV** from Viscometer + RI
- **Rh** from Viscometer + RI + MALS
- **Mark-Houwink** plot (log IV vs. log Mw) shows structure and degree of branching

Results

Separation of PS standards




System

- TF2000 FFF System
- PN3150 RI Detector

Conditions

- Injection Volume: 20 μ L
- Concentration: 2 mg/mL
- Temp. grad. $\Delta T = 90^\circ\text{K}$ to 0°K

DOW PS 1683 / $M_w = 250$ kDa, PD = 2.5




System

- TF2000 FFF System
- PN3621 MALS
- PN3310 Viscometer
- PN3211 UV
- PN3150 RI
- PN3510 ELSD

Conditions

- Injection Volume: 50 μ L
- Concentration: 5 mg/mL
- Temp. grad. $\Delta T = 90^\circ\text{K}$ to 0°K

DOW PS 1683 / $M_w = 250$ kDa, PD = 2.5




System



- TF2000 FFF System
- PN3621 MALS
- PN3310 Viscometer
- PN3211 UV
- PN3150 RI
- PN3510 ELSD

Conditions

- Injection Volume: 50 μ L
- Concentration: 5 mg/mL
- Temp. grad. $\Delta T = 90^\circ\text{K}$ to 0°K




Separation of PS and PMMA-Standards with same R_h SEC vs. Thermal FFF (TF3)



TF3 enables separation of molecules with the same hydrodynamic volume according to their chemical composition!


TFFF Applications – SBR Rubber with Pentadetection




- MALS
- Viscometer
- RI
- ELSD
- UV

Application Example: Natural Rubber

Molar Mass




Radius of Gyration



- **Nanoscale impurities with same size separated according to different chemical composition**

Natural Rubber Sample (Polyisoprene)



Stiff chain structure in low MW area


System

- TF2000 FFF System
- PN3621 MALS
- PN3310 Viscometer
- PN3211 UV
- PN3150 RI
- PN3510 ELSD

Conditions

- Injection Volume: 50 μL
- Concentration: 5 mg/mL
- Temp. grad. $\Delta T = 90^\circ\text{K}$ to 0°K

Natural Rubber Sample



System


- TF2000 FFF System
- PN3621 MALS
- PN3310 Viscometer
- PN3211 UV
- PN3150 RI
- PN3510 ELSD

Conditions

- Injection Volume: 50 μL
- Concentration: 5 mg/mL
- Temp. grad. $\Delta T = 90^\circ\text{K}$ to 0°K

More compact structure in high MW area (branching?)

Raw Data Fractogram of TF3 - MALS and ELSD



System

- PN5300 Auto Injector
- TF2000 Thermal FFF
- ELSD
- PN3621 MALS (92° SLS)


Conditions

- Injection Volume: 100 μ L
- Concentration: 2.0 mg/mL
- LS 92° (red trace)
- ELSD (blue trace)

- The fractogram shows a system start peak at 11 min.
- The 1st peak – the main peak - was detected between 15 - 50 min by light scattering and ELSD detection.
- A 2nd peak was detected at 54 min predominantly by light scattering.

Overlay: Molar Mass and LS Signal

- The molar mass was calculated from MALS and ELSD data
- Literature value for $dn/dc = 0.124 \text{ mL/g}$
- In the 1st peak the sample contains rubber material with a molar mass of appr. $4.7 \times 10^5 \text{ g/mol}$ (w-average) and in the 2nd peak of $3.6 \times 10^8 \text{ g/mol}$.




Conditions

- Molar Mass (red dots)
- LS Signal (blue trace)
- Fitting by Random Coil Model

		M_w [g/mol]
1. Peak 13.0 – 48.5 min	<i>n</i> -Average	2.7×10^5
	w-Average	4.7×10^5
	<i>z</i> -Average	8.1×10^5
2. Peak 48.5 – 56.0 min	<i>n</i> -Average	8.1×10^6
	w-Average	3.6×10^8
	<i>z</i> -Average	1.5×10^9

Overlay: Radius of Gyration and LS Signal

- The Radius of Gyration was calculated from MALS angular data
- The main/1st peak shows a Radius of Gyration of 42 nm (z-average) and the 2nd peak a Radius of Gyration of 218 nm.




Conditions


- Radius of Gyration (red dots)
- LS Signal (blue trace)
- Fitting by Random Coil Model

		R_g [nm]
1. Peak 13.0 – 48.5 min	<i>n</i> -Average	23
	<i>w</i> -Average	31
	z-Average	42
2. Peak 48.5 – 56.0 min	<i>n</i> -Average	65
	<i>w</i> -Average	107
	z-Average	218

Molar Mass Distribution




Zoom in Molar Mass Distribution



Differential Molar Mass Distribution (blue trace), Cumulative Molar Mass Distribution (red trace)

The sample shows a multimodal distribution. For the 1st fraction (13 – 48.5 min) the distribution is in the range of $9.2 \times 10^4 - 2.8 \times 10^6$ g/mol and for the 2nd fraction (48.5 – 56 min) in the range of $2.7 \times 10^6 - 2.3 \times 10^9$ g/mol). The sample contains 1.3 % of high molar mass material (gel). Calculation based on concentration detector signal.

Fractal Dimension




System

- PN5300 Auto Injector
- AF2000 FFF System
- PN3621 MALS Detector
- PN3150 RI Detector
- PN3510 ELSD


Conditions

- Injection Volume: 20 μ L
- Concentration: 2 mg/mL

Fractal Dimension 17.5 - 44.6 min



Fractal Dimension 51.6 – 54.8 min



The Fractal Dimension was evaluated for 2 different regions.

For the 1st region of the fractogram (17.5 – 44.6 min) the Fractal Dimension is 1.78 and for the 2nd region (51.6 – 54.8 min) 2.92.

The Fractal Dimension change indicates that in the 2nd part of the fractogram there is material with a higher density, indicating cross linking or branches.

Thermal FFF is a Powerful Method for Polymer and Biopolymer Characterization

- Increased resolution for high molar mass species
- Huge flexibility in choice of fractionation power (gradient)
- No interaction with stationary phase / no artifacts
- No degradation due to shear stress during the separation
- No filtration by stationary phase or frits
- Separation according to molecular weight and size as well as according to chemical composition offers new horizons for further applications

Thank you for your Attention