

# Thinking Outside the C18 Box



# Modes of HPLC


- Reversed-phase
    - C18, C8, other non-polar phases
    - Phenyl, PFP, other more polar phases
  - Ion Pairing
  - HILIC
    - Silica
    - Amide
    - Zwitterionic
    - Hydroxyl-based
  - Chiral
- 
- Normal phase
    - Silica (Rx-SIL)
    - Aminopropyl
  - Ion Exclusion (Hi-Plex)
    - Sulfonated polystyrene/divinylbenzene
  - Ion Exchange (Bio IEX)
    - Anion exchange
    - Cation exchange
  - Size Exclusion/Gel Permeation (AdvanceBio SEC) (PLgel)

# Retention Factor

$$k = \frac{(t_R - t_0)}{t_0}$$

$t_R$  = retention time for sample peak

$t_0$  = retention time for unretained peak



The retention factor measures the period of time that the sample component resides in the stationary phase relative to the time it resides in the mobile phase. It is calculated from the retention time divided by the time for an unretained peak.


# Pore Dewetting or Phase Collapse

- Alkyl phases such as C8 or C18 can exhibit poor retention or reproducibility of retention in low organic mobile phases
- Phenomenon known as pore dewetting or phase collapse
- Onset can be unpredictable
- A method robustness issue often mistaken as a column or lot issue
- See Przybyciel and Majors, *LCGC* **20**(6), 516-523 (2002).

# Agilent InfinityLab Poroshell Phases


Best all around	Best for low pH mobile phases	Best for high and mid pH mobile phases	Best for alternative selectivity	Best for more polar compounds	HILIC for polar compounds
Poroshell 120 <b>EC-C18</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>SB-C18</b> 2.7 µm	Poroshell <b>HPH-C18</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>Phenyl-Hexyl</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>SB-Aq</b> 2.7 µm	Poroshell 120 <b>HILIC</b> 1.9 µm, 2.7 µm, 4 µm
Poroshell 120 <b>EC-C8</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>SB-C8</b> 2.7 µm	Poroshell <b>HPH-C8</b> 2.7 µm, 4 µm	Poroshell 120 <b>Bonus-RP</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>EC-CN</b> 2.7 µm	Poroshell 120 <b>HILIC-Z</b> 2.7 µm
			Poroshell 120 <b>PFP</b> 2.7 µm		Poroshell 120 <b>HILIC-OH5</b> 2.7 µm

# Poroshell 120 Particle Sizes




InfinityLab Poroshell 120  
**1.9 µm**

Highest UHPLC  
performance



InfinityLab Poroshell 120  
**2.7 µm**


UHPLC performance at  
lower pressure




InfinityLab Poroshell 120  
**4 µm**

Improved HPLC  
performance


# Water Soluble Vitamins




Thiamine, B1




Riboflavin, B2




Nicotinic Acid, B3




Pantothenic Acid, B5




Pyridoxine, B6




Folic Acid, B9




PABA, B10



Biotin, B7




Ascorbic Acid, C



Cyanocobalamin, B12

# Water Soluble Vitamins C18 at low pH

A: 20 mM  $\text{NaH}_2\text{PO}_4$  pH 2.5  
B:  $\text{CH}_3\text{CN}$ , 10% B isocratic  
0.5 mL/min, 30 C, 210 nm  
2.1 x 50 mm, 2.7  $\mu\text{m}$   
Poroshell 120 EC-C18



- Poor retention for 4 analytes
- Coelution
- For 6 compounds,  $k' < 2$




## Now What?

- **Adjust mobile phase pH**
- Ion-pair chromatography
- Alternate column choice
- HILIC

# Why try higher pH?


## Poroshell HPH-C8 or C18



Agilent Poroshell HPH C-18, 4.6 × 50 mm, 2.7 μm


Time	% Buffer	% MeCN
0	10	90
5	90	10
7	10	90
2 mL/min		254 nm

- Peak ID
1. Procainamide
  2. Caffeine
  3. Acetyl salicylic acid
  4. Hexanophenone deg.
  5. Dipyrimadole
  6. Diltiazem
  7. Diflunisal
  8. Hexanophenone



# Water Soluble Vitamins C18 at pH 7.5

A: 20 mM  $\text{Na}_2\text{HPO}_4$  pH 7.5  
B:  $\text{CH}_3\text{CN}$ , 10% B isocratic  
0.5 mL/min, 30 C, 210 nm  
2.1 x 50 mm, 2.7  $\mu\text{m}$   
Poroshell 120 HPH-C18




- Acid groups now charged
- Loss of retention

# Now What?

- Adjust mobile phase pH
- **Ion-pair chromatography**
- Alternate column choice
- HILIC

# Ion-Pair Chromatography


Similar to reversed-phase, but an ion-pairing reagent is added to the mobile phase



- *Non-polar alkyl chain will adsorb into the non-polar stationary phase*
- *Polar part of the ion-pairing reagent will “stick-out” into the mobile phase*

# Ion-Pair Chromatography


Similar to reversed-phase, but an ion-pairing reagent is added to the mobile phase




- *Non-polar alkyl chain will adsorb into the non-polar stationary phase*
- *Polar part of the ion-pairing reagent will “stick-out” into the mobile phase*

# Some Common Ion-Pairing Reagents

## Pairs with Cations




Alkyl sulfonates



Heptafluorobutyric acid  
(HFBA)

## Pairs with Anions



Quaternary amines

# Ion-Pair Chromatography

## Suggested Experimental Conditions

Column: C8 or C18

Mobile Phase:

- Organic – often methanol
- Aqueous - Buffered with appropriate IP reagent
- Temperature controlled between 35° and 60°C

### Cations – bases

Buffer: 25 – 50 mM phosphate, pH 2- 3

IP reagent: 10-100 mM heptane sulfonate


### Anions – acids

Buffer: 25 – 50 mM phosphate, pH 6 – 7


IP reagent: 10-40 mM tetrabutyl ammonium phosphate




# Water Soluble Vitamins




Thiamine, B1




Riboflavin, B2




Nicotinic Acid, B3




Pantothenic Acid, B5




Pyridoxine, B6




Folic Acid, B9




PABA, B10



Biotin, B7




Ascorbic Acid, C



Cyanocobalamin, B12

# Water Soluble Vitamins Ion Pair Conditions



EC-C18

A: 1.5 g sodium 1-heptanesulfonate + 0.2 mL triethylamine +  
7.5 mL acetic acid + 992.5 mL water

B: CH<sub>3</sub>CN

0.5 mL/min, 10% B for 1 minute, then 10-40% B in 1 minute  
injection volume: varies according to signal strength

TCC: 30 C

260, 8 nm Ref Off, 8 nm slit, 80 Hz

The ion pairing reagent increased retention for most compounds

- 6 compounds have  $k' > 2$
- B5 and B7 could not be detected due to low signal and high background noise at 210 nm (not detectable at 260 nm)

# Ion-Pair Chromatography Limitations

- Higher level of complexity than RP, so generally chosen only if needed
- Requires careful control of IP reagent, pH, temperature
- Gradient methods are more difficult than RP
- Equilibration is much slower than RP
- Column dedicated to IP
- IP reagent in the injection solvent
- IP reagents not desirable for MS detection

## Now What?

- Adjust method conditions
- Ion-pair chromatography
- **Alternate column choice**
- HILIC


# Phase Choices

Best all around	Best for low pH mobile phases	Best for high and mid pH mobile phases	Best for alternative selectivity	Best for more polar compounds	HILIC for polar compounds
Poroshell 120 <b>EC-C18</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>SB-C18</b> 2.7 µm	Poroshell <b>HPH-C18</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>Phenyl-Hexyl</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>SB-Aq</b> 2.7 µm	Poroshell 120 <b>HILIC</b> 1.9 µm, 2.7 µm, 4 µm
Poroshell 120 <b>EC-C8</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>SB-C8</b> 2.7 µm	Poroshell <b>HPH-C8</b> 2.7 µm, 4 µm	Poroshell 120 <b>Bonus-RP</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>EC-CN</b> 2.7 µm	Poroshell 120 <b>HILIC-Z</b> 2.7 µm
			Poroshell 120 <b>PFP</b> 2.7 µm		Poroshell 120 <b>HILIC-OH5</b> 2.7 µm

These phases can be used with high aqueous mobile phases to improve retention of highly polar analytes in RPLC mode

# Water Soluble Vitamins Alternative Phases



Retention Times (min) for Water Soluble Vitamins  
 A: 20 mM sodium phosphate pH 2.5, B: acetonitrile, 0.5 mL/min, 30 C,  
 210 nm



# Water Soluble Vitamins

## Phenyl-Hexyl


A: 20 mM NaH<sub>2</sub>PO<sub>4</sub> pH 2.5  
B: CH<sub>3</sub>CN, 0% B for 2 min,  
then 0-30% B in 2.5 min  
0.5 mL/min, 30 C, 210 nm  
2.1 x 50 mm, 2.7 μm  
Poroshell 120 Phenyl-Hexyl



Phenyl-Hexyl,  
100% pH 2.5  
Buffer



- Phenyl-Hexyl has the best retention
- 7 compounds have  $k' > 2$ ;
- C18 analysis had only 4 compounds with  $k' > 2$

# Aliphatic Acids SB-Aq



Peak ID

1. Tartaric acid
2. Quinic acid
3. Malic acid
4. Citric acid
5. Fumaric acid



Column: Agilent Poroshell 120 SB-Aq, 3 × 100 mm, 2.7 μm  
(p/n 685975-314)

Eluent: 100 mM Potassium phosphate buffer, pH 2.5

Injection volume: 5 μL

Flow rate: 0.5 mL/min

Temperature: 50 °C

Detector: DAD, at 226 nm

5991-1992EN






# What Now?

- Adjust method conditions
- Ion-pair chromatography
- Alternate column choice
- **HILIC**

# What is HILIC?


## Hydrophilic Interaction Liquid Chromatography

- Retains hydrophilic compounds
- Polar stationary phase: silica, amino, diol/hydroxyl-based, zwitterionic, etc.
- Uses an organic and water mobile phase with a buffer



# HILIC Mechanism

- Retains moderate to highly polar analytes
  - A water layer is adsorbed onto the polar silica surface, creating a liquid/liquid extraction system
  - Polar analytes can **partition into and out of the water layer**, with more polar analytes having a stronger interaction (1)
  - Charged polar analytes can also undergo **ion exchange with the silica surface** (2)




1. Partitioning in and out of adsorbed water layer
2. Ion exchange with silanols


## Key points to note:

- Water is the strong solvent
- Elution is least polar to most polar, opposite of RPLC.
- Gradients run from **high organic to high aqueous** (10% aqueous to 50% aqueous is a common scouting gradient)
- Reversed-phase solvents (ACN/Water)
  - MeOH, EtOH, IPA can also be used
- Typically uses a buffer like ammonium formate or ammonium acetate
- Higher buffer concentration increases solvent strength, improves peak shape, and can change selectivity slightly


# InfinityLab Poroshell 120 HILIC Columns - Specifications




InfinityLab Poroshell Column	Bonded phase	Particle	Particle Size	Pressure Limit	Pore Size	Temp Limit	pH Range	Dimensions (ID in mm)	Dimensions (Length in mm)
<b>HILIC-Z</b>	Proprietary zwitterionic chemistry	Hybrid Poroshell superficially porous particle	2.7µm	600 bar	100 Å	80°C	2-12	2.1	5 (guard)
								3.0	50
								4.6	100
									150
									PEEK-lined version available
<b>HILIC-OH5</b>	Poly-hydroxy fructan chemistry	Poroshell superficially porous particle	2.7µm	400 bar	120 Å	45°C	1-7	2.1	50
								3.0	100
								4.6	150
<b>HILIC</b>	Bare-silica (unbonded)	Poroshell superficially porous particle	1.9µm 2.7µm 4µm	1300 bar (1.9) 600 bar (2.7) 600 bar (4)	120 Å	60°C	0-8	2.1	5 (guard)
								3.0	50
								4.6	100
									150



**HILIC-Z**  
Proprietary Zwitterionic chemistry



**HILIC-OH5**  
Poly-hydroxy fructan chemistry



**HILIC**  
Bare Silica chemistry

# HILIC Method Development: Common LC Parameters

## Type of stationary phase

- Vary retention mechanism and selectivity
- 3 phases on Agilent InfinityLab Poroshell 120 2.7  $\mu\text{m}$  particles
  - HILIC-Z, HILIC-OH5, HILIC

## Mobile phase pH

- Controls ionization of silica and analytes
- Compounds are more retained in their charged state
  - Acids should be run at high pH, bases at low pH

## Temperature

- Increasing temperature will decrease retention
- Increasing temperature will increase column efficiency
- Decreasing temperature can improve selectivity

# HILIC Method Development: Mobile Phase Considerations


## Organic solvent concentration

- Solvent strength in HILIC mode:
  - $THF < Acetone < CH_3CN < IPA < EtOH < MeOH < H_2O$
- $H_2O$  must be present — *need*  $> 3\% H_2O$
- **Ionic strength of buffer**
- Concentration of (salt) buffer increases strength
- Different anions and cations may can also affect analyte retention

## Type of buffer









- Acetates, formates good, soluble in  $CH_3CN$ —also MS friendly
- *Phosphate salts* have low  $CH_3CN$  solubility

# Less CH<sub>3</sub>CN Makes a HILIC Mobile Phase Stronger, Causing Less Retention




Column used was 2.1 x 150 mm, 2.7  $\mu$ m Agilent InfinityLab Poroshell 120 HILIC-Z (PEEK lined); A: 100 mM pH 3 Ammonium Formate in Water, B: Acetonitrile, x % B, isocratic elution, 0.25 mL/min, 30  $^{\circ}$ C, 1  $\mu$ L injection of toluene, cytosine, uracil QC mixture, 254 nm

# Water Soluble Vitamins






# Water Soluble Vitamins by HILIC



Columns used were 2.1 x 100 mm, 2.7  $\mu$ m; A: 100 mM Ammonium Acetate + 0.5% Acetic Acid (pH ~4.6) in H<sub>2</sub>O, B: CH<sub>3</sub>CN, 0.5 mL/min, 87% B for 1 min, 87-50% B in 4 min, 3 min re-equilibration, 1  $\mu$ L injection of individual vitamin standards (0.1-0.4 mg/mL each), 40  $^{\circ}$ C, 260 nm, 80 Hz

# Temperature Used to Optimize a HILIC Separation


## Water Soluble Vitamins on Agilent InfinityLab Poroshell 120 HILIC-OH5



Agilent InfinityLab Poroshell 120 HILIC-OH5 2.1 x 100 mm, 2.7  $\mu$ m; A: 100 mM Ammonium Acetate (no pH adjustment) in H<sub>2</sub>O, B: CH<sub>3</sub>CN, 0.5 mL/min, 95-60%B in 10 min, 3 min re-equilibration, 1  $\mu$ L injection of individual vitamin standards (0.1-0.4 mg/mL each), 20/40/60 °C, 260 nm, 80 Hz

# Excellent retention, peak shape and sensitivity with HILIC-Z

## Underivatized Amino Acids by LC/MS



### InfinityLab Poroshell HILIC-Z 2.1 x 100mm, 2.7µm

A: 20mM ammonium formate in H<sub>2</sub>O, pH3  
B: 9/1 ACN/H<sub>2</sub>O with 20mM ammonium formate, pH3

Gradient: 100%B to 70% B over 10 min, return to 100% B

Flow rate: 0.8 ml/min

Temp: 30deg

MS Detection: Agilent MS-QQQ, MS2 SIM mode

# HILIC Analyses Perform Best with Weak Injection Solvents

## B Vitamins on HILIC with Isocratic Elution


1  $\mu\text{L}$  injection in  $\text{H}_2\text{O}$

1  $\mu\text{L}$  injection in  $\text{H}_2\text{O}/\text{CH}_3\text{CN}$  (3:1)

1  $\mu\text{L}$  injection in  $\text{H}_2\text{O}/\text{CH}_3\text{CN}$  (1:1)

1  $\mu\text{L}$  injection in  $\text{H}_2\text{O}/\text{CH}_3\text{CN}$  (1:3)

1  $\mu\text{L}$  injection in  $\text{CH}_3\text{CN}$



Agilent ZORBAX RRHD HILIC Plus 2.1 x 50 mm, 1.8  $\mu\text{m}$ ; Mobile Phase: acetonitrile / 100 mM Ammonium Formate pH 3.2 in water (9:1), isocratic elution, 0.4 mL/min, 1  $\mu\text{L}$  injection of 5.7  $\mu\text{g}/\text{mL}$  each of 4-aminobenzoic acid, nicotinamide, riboflavin, nicotinic acid; 25  $^\circ\text{C}$ , MS Source: ESI+, 200  $^\circ\text{C}$ , 10 L/min., 30 psi, 4000 V; SIM: 138, 123, 377, 124

# Advantages of HILIC

- Retains polar analytes where reversed-phase methods may not
- Offers alternative selectivity to RPLC mode
- Can retain cations, anions, and polar neutrals in a single run
- Can improve peak shape for basic compounds
- Uses a standard LC system and common reversed-phase solvents
- Uses low viscosity mobile phases with high organic content
  - Fast methods with high flow rates
  - Longer columns for higher efficiency at lower pressures
- Enhanced detection sensitivity with MS compatible methods, in both positive and negative modes
- Can directly inject ACN extracts from C18 SPE cartridges




# Agilent InfinityLab Poroshell Phases


Best all around	Best for low pH mobile phases	Best for high and mid pH mobile phases	Best for alternative selectivity	Best for more polar compounds	HILIC for polar compounds	Chiral phases
Poroshell 120 <b>EC-C18</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>SB-C18</b> 2.7 µm	Poroshell <b>HPH-C18</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>Phenyl-Hexyl</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>SB-Aq</b> 2.7 µm	Poroshell 120 <b>HILIC</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>Chiral-CF</b> 2.7 µm
Poroshell 120 <b>EC-C8</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>SB-C8</b> 2.7 µm	Poroshell <b>HPH-C8</b> 2.7 µm, 4 µm	Poroshell 120 <b>Bonus-RP</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>EC-CN</b> 2.7 µm	Poroshell 120 <b>HILIC-Z</b> 2.7 µm	Poroshell 120 <b>Chiral-CD</b> 2.7 µm
			Poroshell 120 <b>PFP</b> 2.7 µm		Poroshell 120 <b>HILIC-OH5</b> 2.7 µm	Poroshell 120 <b>Chiral-V</b> 2.7 µm
						Poroshell 120 <b>Chiral-T</b> 2.7 µm

# Poroshell 120 Chiral Chemistries


InfinityLab Poroshell 120 Chiral-CF  
(Cyclofructan CF-6)




InfinityLab Poroshell 120 Chiral-V  
(Vancomycin)



InfinityLab Poroshell 120 Chiral-CD  
(Hydroxypropylated beta-cyclodextrin)



InfinityLab Poroshell 120 Chiral-T  
(Teicoplanin)



# Poroshell 120 Chiral Chemistries

Column Chemistry	Chiral Selector (bonded chemistry)	Typical LC Mode	Typical Applications
InfinityLab Poroshell 120 Chiral-CF	Derivatized cyclofructan (CF6)	Polar Organic (PO)	Primary amines
		Normal Phase (NP)	Primary amines
InfinityLab Poroshell 120 Chiral-CD	Hydroxypropylated- $\beta$ -cyclodextrin	Reversed Phase (RP)	Stimulants, fungicides, t-boc amino acids
		Polar Organic (PO)	Complex molecules
InfinityLab Poroshell 120 Chiral-V	Vancomycin (macrolide antibiotic)	Polar Ionic (PI)	Basic pharmaceuticals (various)
		Reversed Phase (RP)	Amines, profens
		Polar Organic (PO)	Complex neutral molecules
InfinityLab Poroshell 120 Chiral-T	Teicoplanin (macrolide antibiotic)	Polar Ionic (PI)	Beta blockers, hydroxyl acids
		Reversed Phase (RP)	Amino acids, hydroxyl acids, profens
		Polar Organic (PO)	Hydantoins, benzodiazepines



# Modes of Separation used with Infinity Lab Chiral Columns

## Polar Ionic Mode

Methanol with acid or base or volatile salt < 0.2 % wt. (MeOH + HOAc + TEA )

Non-aqueous mobile phase; fast, MS detection; for ionizable molecules – any acid or base

Dominant interactions: Ionic interaction, hydrogen bonding

Example: MeOH with 0.2 wt% ammonium formate

## Reversed Phase Mode

Methanol/Water/Buffer,

MS compatible, ideal for manufacturing QC, bioanalysis for all types of molecules

Example: 30/70 MeOH/20 mM ammonium formate (pH 4)

## Polar Organic Mode

Acetonitrile/Methanol/Ethanol/Isopropanol+ HOAc + TEA

Dominant interactions: Hydrogen bonding, dipole-dipole

Example: 60/40/0.3/0.2 ACN/MeOH/acetic acid/TEA

## Normal Phase


Heptane (or Hexane)/methanol or ethanol

Example: 60/40/0.3/0.2 ACN/MeOH/acetic acid/TEA

# Fast, High Efficiency Chiral Separations

## Traditional Chiral Separation— totally porous particle


Chirobiotic V2 (250 x 4.6 mm, 5 µm)



1. D-(+)-Amphetamine, 2. L(-)-Amphetamine, 3. D-(+)-Methamphetamine  
4. L(-)-Methamphetamine 100/0.1/0.02 MeOH/HOAc/NH<sub>4</sub>OH with a  
1.0 mL/min flow rate at room temperature and UV at 220 nm

## Agilent InfinityLab Poroshell 120 Chiral Separation— superficially porous particle

InfinityLab Poroshell 120 Chiral-V (100 x 4.6 mm, 2.7 µm)




1. D-(+)-Amphetamine, 2. L(-)-Amphetamine, 3. D-(+)-Methamphetamine  
4. L(-)-Methamphetamine 100/0.1/0.02 MeOH/HOAc/NH<sub>4</sub>OH with a  
1.0 mL/min flow rate at room temperature and UV at 220 nm

# Chiral Applications Compendium

## Publication number: 5991-8450EN

### Amines

#### 1,1-diphenyl-1-fluoro-2-aminopropane




**Retention (min):** 1.48/1.67 (Peak1/Peak 2)  
**Resolution:** 3.20

##### Method Conditions

**Column:** InfinityLab Poroshell 120 Chiral-V (10 cm x 4.6 mm, 2.7 µm)  
**Mobile phase:** 100/0.1 wt %: Methanol/Ammonium Trifluoroacetate  
**Flow Rate:** 1.0 mL/min  
**Temperature:** Ambient (23 °C)  
**Injection Volume:** 1.0 µL  
**Detection:** UV 230 nm



#### 1,1-diphenyl-2-aminopropane




**Retention (min):** 1.87/2.13 (Peak1/Peak 2)  
**Resolution:** 3.33

##### Method Conditions

**Column:** InfinityLab Poroshell 120 Chiral-V (10 cm x 4.6 mm, 2.7 µm)  
**Mobile phase:** 100/0.1 wt %: Methanol/Ammonium Trifluoroacetate  
**Flow Rate:** 1.0 mL/min  
**Temperature:** Ambient (23 °C)  
**Injection Volume:** 1.0 µL  
**Detection:** UV 230 nm




#### 1,2,2-triphenylethylamine




**Retention (min):** 1.54/1.93 (Peak1/Peak 2)  
**Resolution:** 3.88

##### Method Conditions

**Column:** InfinityLab Poroshell 120 Chiral-V (10 cm x 4.6 mm, 2.7 µm)  
**Mobile phase:** 100/0.1 wt %: Methanol/Ammonium Trifluoroacetate  
**Flow Rate:** 1.0 mL/min  
**Temperature:** Ambient (23 °C)  
**Injection Volume:** 1.0 µL  
**Detection:** UV 220 nm



#### 2-amino-4-methyl-1,1-diphenylpentane



**Retention (min):** 0.66/0.80 (Peak1/Peak 2)  
**Resolution:** 4.49


##### Method Conditions

**Column:** InfinityLab Poroshell 120 Chiral-V (10 cm x 4.6 mm, 2.7 µm)  
**Mobile phase:** 100/0.1 wt %: Methanol/Ammonium Trifluoroacetate  
**Flow Rate:** 1.0 mL/min  
**Temperature:** Ambient (23 °C)  
**Injection Volume:** 1.0 µL  
**Detection:** UV 220 nm



### Anesthetics

#### Hexobarbital




**Retention (min):** 4.63/4.99 (Peak1/Peak 2)  
**Resolution:** 1.54

##### Method Conditions

**Column:** InfinityLab Poroshell 120 Chiral-T (10 cm x 4.6 mm, 2.7 µm)  
**Mobile phase:** 30/70: Methanol/15 mM Ammonium Formate (pH 3.6)  
**Flow Rate:** 0.5 mL/min  
**Temperature:** Ambient (23 °C)  
**Injection Volume:** 1.0 µL  
**Detection:** UV 220 nm



#### Hexobarbital




**Retention (min):** 4.70/5.29 (Peak1/Peak 2)  
**Resolution:** 2.92

##### Method Conditions

**Column:** InfinityLab Poroshell 120 Chiral-V (10 cm x 4.6 mm, 2.7 µm)  
**Mobile phase:** 30/70: Methanol/15 mM Ammonium Formate (pH 3.6)  
**Flow Rate:** 0.5 mL/min  
**Temperature:** Ambient (23 °C)  
**Injection Volume:** 1.0 µL  
**Detection:** UV 220 nm



#### Mepivacaine




**Retention (min):** 4.01/4.36 (Peak1/Peak 2)  
**Resolution:** 1.42

##### Method Conditions

**Column:** InfinityLab Poroshell 120 Chiral-T (10 cm x 4.6 mm, 2.7 µm)  
**Mobile phase:** 100/0.2/0.05: Methanol/Acetic Acid/Ammonium Hydroxide  
**Flow Rate:** 1.0 mL/min  
**Temperature:** Ambient (23 °C)  
**Injection Volume:** 1.0 µL  
**Detection:** UV 220 nm



#### Mepivacaine



**Retention (min):** 4.75/5.14 (Peak1/Peak 2)  
**Resolution:** 2.25


##### Method Conditions

**Column:** InfinityLab Poroshell 120 Chiral-V (10 cm x 4.6 mm, 2.7 µm)  
**Mobile phase:** 90/10: Methanol/15 mM Ammonium Formate (pH 3.6)  
**Flow Rate:** 0.5 mL/min  
**Temperature:** 30 °C  
**Injection Volume:** 1.0 µL  
**Detection:** UV 220 nm



### Beta blockers


#### Isoproterenol




**Retention (min):** 7.76/8.56 (Peak1/Peak 2)  
**Resolution:** 2.70

##### Method Conditions

**Column:** InfinityLab Poroshell 120 Chiral-T (10 cm x 4.6 mm, 2.7 µm)  
**Mobile phase:** 100/0.2/0.1: Methanol/Trifluoroacetate/Ammonium Hydroxide  
**Flow Rate:** 0.7 mL/min  
**Temperature:** 45 °C  
**Injection Volume:** 1.0 µL  
**Detection:** UV 230 nm




#### Isoproterenol




**Retention (min):** 4.60/4.96 (Peak1/Peak 2)  
**Resolution:** 1.0

##### Method Conditions

**Column:** InfinityLab Poroshell 120 Chiral-V (10 cm x 4.6 mm, 2.7 µm)  
**Mobile phase:** 100/0.2/0.05: Methanol/Acetic Acid/Ammonium Hydroxide  
**Flow Rate:** 0.5 mL/min  
**Temperature:** Ambient (23 °C)  
**Injection Volume:** 1.0 µL  
**Detection:** UV 220 nm




#### Metaproterenol




**Retention (min):** 1.62/1.89 (Peak1/Peak 2)  
**Resolution:** 3.01

##### Method Conditions

**Column:** InfinityLab Poroshell 120 Chiral-T (10 cm x 4.6 mm, 2.7 µm)  
**Mobile phase:** 60/40/0.3/0.2: Acetonitrile/Methanol/Acetic Acid/TEA  
**Flow Rate:** 3.5 mL/min  
**Temperature:** 45 °C  
**Injection Volume:** 1.0 µL  
**Detection:** UV 230 nm




#### Metaproterenol



**Retention (min):** 4.63/4.98 (Peak1/Peak 2)  
**Resolution:** 1.89

##### Method Conditions

**Column:** InfinityLab Poroshell 120 Chiral-V (10 cm x 4.6 mm, 2.7 µm)  
**Mobile phase:** 100/0.2/0.05: Methanol/Acetic Acid/Ammonium Hydroxide  
**Flow Rate:** 0.5 mL/min  
**Temperature:** Ambient (23 °C)  
**Injection Volume:** 1.0 µL  
**Detection:** UV 230 nm



# Summary

- What do you do when C18 does not work?
- Stick with reversed-phase but--
  - Adjust pH of mobile phase
  - Try a more polar bonded phase
- Consider HILIC
- Chiral phases
- **Fast Poroshell methods make method development faster**

# Contact Agilent Chemistries and Supplies Technical Support



1-800-227-9770 Option 3, Option 3:

Option 1 for GC/GCMS Columns and Supplies

Option 2 for LC/LCMS Columns and Supplies

Option 3 for Sample Prep Products, Filtration and QuEChERS

Option 4 for Spectroscopy Supplies

**Available in the USA 8-5 all time zones**



[gc-column-support@agilent.com](mailto:gc-column-support@agilent.com)

[lc-column-support@agilent.com](mailto:lc-column-support@agilent.com)

[spp-support@agilent.com](mailto:spp-support@agilent.com)

[spectro-supplies-support@agilent.com](mailto:spectro-supplies-support@agilent.com)

# Equilibrate from high aqueous to low


## Critical factor when changing mobile phases

ZORBAX Rx-Sil, 2.1 x 150, 5um


A: 25mM ammonium acetate with 2.5 mM ammonium formate  
B: acetonitrile  
10:90 A:B  
Flow : 0.1 mL/min  
Temp: 25 °C  
Detection: UV 254nm,4, ref=360,100  
Inj. Vol 0.5 µL  
Instrument Quat 1200  
Samples: ca. 0.4 mg/ml each in ACN/water

1. thymine
2. Uracil
3. adenine
4. guanine
5. cytosine

Equilibrating from 80 to 90% ACN works



Equilibrating from 95 to 90% ACN does *not* work well



# Water Soluble Vitamins HILIC

Poroshell 120 HILIC

A: 100 mM Ammonium Formate in H<sub>2</sub>O pH 3.0 with Formic Acid


B: CH<sub>3</sub>CN

0.5 mL/min, 97% B for 2.5 minutes, then 97-60%D in 5 minutes

injection volume: varies according to signal strength

TCC: 30 C

260/210, 8 nm Ref Off, 8 nm slit, 80 Hz



# Sugar Analysis on Agilent InfinityLab Poroshell 120 HILIC-Z using ELSD

