

Polyurethane Foam as an Effective Adsorbent for Pesticides in Air (GC Analysis)

Polyurethane foam (PUF) is a lightweight, easy-to-use material for sampling pesticides in air. Low back pressure allows sampling rates of 1-5 liters/minute, ensuring short sampling times and representative samples. A PUF filter cartridge designed to attach to the front of an ORBO-1000 PUF sampling cartridge traps aerosols and particulate matter, is easily and conveniently assembled, and can be removed from the ORBO-1000 cartridge and capped separately to maintain sample integrity. Both the foam plug and the glass holder that comprise an ORBO-1000 cartridge are reusable when cleaned as described here.

Key Words:


- pesticides • chlorinated pesticides • air monitoring

Increasing concern over the presence of pesticides in indoor and ambient air has led to the use of polyurethane foam (PUF) as a lightweight, easy-to-use sampling material for these semivolatile compounds. Cylindrical PUF plugs of appropriate density effectively trap contaminants without creating excessive back pressure. Sampling rates of 1-5 liters/minute for pesticides allow shorter sampling times and ensure more representative samples. These relatively high sampling rates also can be required to attain needed detection limits. PUF plugs sometimes are used with an additional adsorbent, such as Amberlite® XAD®-2 resin, to collect more volatile contaminants, or with a filter, to trap particles.

Our PUF filter cartridge contains a replaceable 32mm diameter quartz microfiber filter with a stainless steel support which traps aerosol and particulate forms of semivolatiles, such as pesticides and PCBs. It is designed to attach to the front of an ORBO™-1000 PUF sampling cartridge (Figure A). The low pressure drop of the PUF cartridge/filter cartridge allows high sampling volume (1-5 liters/minute). The filter cartridge does not contribute contaminants to the sample, is easily and conveniently assembled, and can be removed from the ORBO-1000 cartridge and capped separately to maintain sample integrity.

Applications citing a 22mm diameter PUF plug include American Society for Testing and Materials (ASTM) Method D4861 (*Sampling and Analysis of Pesticides and Polychlorinated Biphenyls in Air*), ASTM D4947 (*Chlordane and Heptachlor Residues in Indoor Air*), US EPA Method IP-8 (*Determination of Organochlorine Pesticides in Indoor Air*), and US EPA Method TO10 (*Determination of Organochlorine Pesticides in Ambient Air*).

Figure A. Polyurethane Foam (PUF)/ORBO-1000 Pesticide Sampling Device



994-0053

The low detection limits required for these methods necessitate extensive cleaning of the PUF material, and meticulous handling during sampling and analysis. Table 1 summarizes the cleaning, sampling, and chromatographic procedures we use to evaluate the PUF plugs used with ORBO-1000 cartridges. We store the clean PUF plugs in sealed containers, protected from light, to maintain their purity.

We compared background and pesticide recovery from cleaned and uncleaned PUF plugs, to verify that our cleaning procedures do not alter the performance of the plugs. Figure B clearly shows clean plugs are essential to achieving the low background levels and detection limits required in pesticides methodology. A mix of 0.5µg of each pesticide in hexane was spiked onto a cleaned or uncleaned PUF plug, then extracted, concentrated, and analyzed as described in Table 1. Octachloronaphthalene was used as a surrogate, to monitor extraction efficiency. The high pesticide spiking concentration masked any interference by contaminants in the uncleaned foam plugs. Recovery rates for pesticides spiked onto cleaned plugs and uncleaned plugs were not significantly different (Table 2), showing that our cleaning process had no adverse effect on adsorbent performance. SPB-608 fused silica capillary columns, used in the chromatographic analyses of these analytes, are specifically tested for analyses of chlorinated pesticides.


Our investigations confirm that cleaned polyurethane foam is an effective adsorbent for pesticides in air. ORBO-1000 PUF sampling cartridges are lightweight, portable, and easy to handle during sample collection. Each consists of a cleaned PUF plug in a glass holder, capped at both ends and ready to use. The 1/4 inch (6mm) neck on the glass holder makes ORBO-1000 cartridges compatible with personal sampling pumps having sampling rates of 1 liter/minute or greater. Both foam plug and glass holder are reusable when cleaned as described in Table 1.

Table 1. Cleaning, Sampling, and Analysis for Organochlorine Pesticides (US EPA IP-8)

Sampling Medium
ORBO-1000 cartridge (7.6cm x 22mm polyurethane foam plug; density = 0.022g/cm ³)
PUF Plug Preparation*
Using Soxhlet apparatus, wash plug with acetone for 14-24 hr, followed by ether/hexane (5:95) for 16 hr. Vacuum dry 2-4 hr at room temperature. Place in glass sampling cartridge and seal until sample collection.
Sample Collection
1-5 liters/min for 4-24 hr. Pesticide concentration: 0.01-50µg/m ³ .
Sample Preparation
Using Soxhlet apparatus, extract plug with ether/hexane (5:95) for 16 hr. Concentrate extract to 10.0 mL.
GC Analysis
Column: SPB™-608, 30m x 0.53mm ID, 0.5µm phase film
Oven: 150°C (1 min) to 270°C at 4°C/min
Carrier: helium, 7mL/min
Det.: ECD, 320°C
Inj.: 1µL ether/hexane extract, direct injection (packed inlet sleeve), 220°C.

*Omit acetone wash when cleaning PUF cartridge for reuse.

Figure B. Background in PUF Cartridges for Pesticide Analysis (GC-ECD)



GC column and conditions described in Table 1.

Table 2. Pesticide Recovery from Cleaned vs. Uncleaned PUF Plugs

Pesticide	% Recovery (0.5µg Spike)	
	Cleaned	Uncleaned
Dichlorvos	68	74
2,4,5-Trichlorophenol	68	91
Pentachlorobenzene	67	56
2,4-D methyl ester	95	91
Hexachlorobenzene	71	72
α-BHC	67	73
γ-BHC (Lindane)	77	84
β-BHC	95	80
Heptachlor	79	76
Chlorothalonil	87	72
Aldrin	73	76
Ronnel	87	94
Chlorpyrifos	96	88
Heptachlor epoxide	86	84
trans-Nonachlor	90	86
Dieldrin	94	95
4,4'-DDE/Captan	99	95
4,4'-DDT	104	85
Mirex	102	85
Methoxychlor	111	82

Ordering Information:

Description	Cat. No.
ORBO-1000 PUF Cartridges, pk. of 3	20557
Glass Holder for ORBO-1000 Cartridge	20556
Cleaned PUF Plug, pk. of 3	20600-U
PUF Filter Cartridge Assembly cartridge, quartz filter, stainless steel screen, 2 endcaps (replacement parts available)	21031
SPB-608 Capillary Column 30m x 0.53mm ID, 0.5µm film	25312

Reference

1. *Compendium of Methods for the Determination of Air Pollutants in Indoor Air* US EPA, Research Triangle Park, North Carolina, USA (April 1990).
Chapter IP-8 - Determination of Organochlorine Pesticides in Indoor Air
Request as publication PB90-200288/AS from: NTIS, 5285 Port Royal Road, Springfield, VA 22161 USA (In USA, call 800-553-NTIS).

Trademarks

Amberlite, XAD – Rohm and Haas Company
ORBO, SPB – Sigma-Aldrich Co.

Contact our Technical Service Department
(phone 800-359-3041 or 814-359-3041, FAX 800-359-3044
or 814-359-5468) for expert answers to your questions.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA · Sigma-Aldrich de Argentina, S.A. · Buenos Aires 1119 AUSTRALIA · Sigma-Aldrich Pty. Ltd. · Castle Hill NSW 2154 AUSTRIA · Sigma-Aldrich Handels GmbH · A-1110 Wien
BELGIUM · Sigma-Aldrich N.V./S.A. · B-2880 Bornem BRAZIL · Sigma-Aldrich Quimica Brasil Ltda. · 01239-010 São Paulo, SP CANADA · Sigma-Aldrich Canada, Ltd. · 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC · Sigma-Aldrich s.r.o. · 186 00 Praha 8 DENMARK · Sigma-Aldrich Denmark A/S · DK-2665 Vallensbaek Strand FINLAND · Sigma-Aldrich Finland/YA-Kemia Oy · FIN-00700 Helsinki
FRANCE · Sigma-Aldrich Chimie · 38297 Saint-Quentin-Fallavier Cedex GERMANY · Sigma-Aldrich Chemie GmbH · D-82041 Deisenhofen GREECE · Sigma-Aldrich (o.m.) Ltd. · Ilioupoli 16346, Athens
HUNGARY · Sigma-Aldrich Kft. · H-1067 Budapest INDIA · Sigma-Aldrich Co. · Bangalore 560 048 IRELAND · Sigma-Aldrich Ireland Ltd. · Dublin 24 ISRAEL · Sigma Israel Chemicals Ltd. · Rehovot 76100
ITALY · Sigma-Aldrich s.r.l. · 20151 Milano JAPAN · Sigma-Aldrich Japan K.K. · Chuo-ku, Tokyo 103 KOREA · Sigma-Aldrich Korea · Seoul MALAYSIA · Sigma-Aldrich (M) Sdn. Bhd. · Selangor
MEXICO · Sigma-Aldrich Quimica S.A. de C.V. · 50200 Toluca NETHERLANDS · Sigma-Aldrich Chemie BV · 3330 AA Zwijndrecht NORWAY · Sigma-Aldrich Norway · Torshov · N-0401 Oslo
POLAND · Sigma-Aldrich Sp. z o.o. · 61-663 Poznań PORTUGAL · Sigma-Aldrich Quimica, S.A. · Sintra 2710 RUSSIA · Sigma-Aldrich Russia · Moscow 103062 SINGAPORE · Sigma-Aldrich Pte. Ltd.
SOUTH AFRICA · Sigma-Aldrich (pty) Ltd. · Jet Park 1459 SPAIN · Sigma-Aldrich Quimica, S.A. · 28100 Alcobendas, Madrid SWEDEN · Sigma-Aldrich Sweden AB · 135 70 Stockholm
SWITZERLAND · Supelco · CH-9471 Buchs UNITED KINGDOM · Sigma-Aldrich Company Ltd. · Poole, Dorset BH12 4QH
UNITED STATES · Supelco · Supelco Park · Bellefonte, PA 16823-0048 · Phone 800-247-6628 or 814-359-3441 · Fax 800-447-3044 or 814-359-3044 · email: supelco@sial.com

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.