

# Analysis of Pesticide Residues in Rice with Bond Elut QuEChERS Extraction Kits and Agilent J&W HP-5ms Ultra Inert GC Column

## Application Note

Food Safety

### Author

Andy Zhai and Yun Zou  
Agilent Technologies Co., Ltd  
412 YingLun Road  
Waigaoqiao Free Trade Zone  
Shanghai 200131 P.R.China

### Abstract

This application note describes the use of a quick, easy, effective, rugged, and safe (QuEChERS) sample preparation approach described in the European Committee (EN) for extraction and cleanup of 57 GC-amenable multiple pesticide class residues in rice. The method involves initial extraction in an aqueous/acetonitrile system, an extraction/partitioning step after the addition of salt, and a cleanup step using dispersive solid phase extraction (dispersive SPE). The target pesticides in the rice extracts were then separated by an Agilent J&W HP-5ms Ultra Inert Capillary GC column and analyzed by gas chromatography/mass spectrometry (GC/MS) operating in selective ion monitoring (SIM) mode. With Agilent RTL pesticides library, the GC/MS files were quickly screened and analyzed for pesticide residues identification. The method was validated in terms of recovery and reproducibility. The spiked levels for the recovery experiments were 50, 100, and 200 ng/g. Most of the recoveries ranged between 80 and 110%, with an average RSD of 5.53%.



**Agilent Technologies**

## Introduction

The QuEChERS method for pesticide analysis was first introduced by USDA scientists in 2003 [1]. The method was modified to address some problematic pesticides by including a buffered extraction system [2]. The EN method 15662:2007 is an European variation to the QuEChERS method [3,4]. The method uses acetonitrile extraction, followed by the salting out of water from the sample using anhydrous magnesium sulfate ( $MgSO_4$ ), NaCl and buffering citrate salts to induce liquid-liquid partitioning. A dispersive SPE is conducted for cleanup using a combination of primary secondary amine (PSA) to remove fatty acids among other components and anhydrous  $MgSO_4$  to reduce the remaining water in the extract. After mixing and centrifugation, the upper layer is ready for analysis.

GC/MS is a good tool for reliable detection and separation of pesticides. Retention Time Locking (RTL) is a simple technique for GC and GC/MS that allows a chromatographer to reproduce retention times on any Agilent GC instrument [5]. Agilent has developed GC and GC/MS RTL databases (p/n G1672AA) that include nine hundred sixty-two pesticides, metabolites, and suspected endocrine disruptors [6]. It is much easier to identify pesticides in complex food extracts by locking the method to this database.

## GC Conditions

|                        |   |
|------------------------|---|
| Column                 | Agilent J&W HP-5ms Ultra Inert, 30 m × 0.25 mm, 0.25 $\mu$ m (p/n 19091S-433UI)                                       |
| Inlet temperature      | 250 °C  |
| Carrier gas            | Helium, constant pressure mode  |
| Retention time locking | Chlorpyrifos-methyl locked to 16.596 min  |
| Injection mode         | Splitless, purge flow 50 mL/min at 0.75 min   |
| Injection volume       | 1 $\mu$ L   |
| Oven                   | 70 °C (2 min), 25 °C/min to 150 °C (0 min), 3 °C/min to 200 °C, 8 °C/min to 280 °C (10 min), post-run: 320 °C (5 min) |

## MS Conditions

|               |  |
|---------------|--|
| Solvent delay | 4 min  |
| MS temp       | 230 °C (Source); 150 °C (Quad)                                       |
| Transfer line | 280 °C   |
| MS libraries  | Agilent RTL Pesticide Library (G1672AA) and NIST08 Mass Spec Library |
| MS            | El, SIM/Scan   |
| Scan mode     | mass range (50–550 amu)  |

For other parameters, see Table 1

Table 1. Pesticides, CAS Number, Molecular Form and Target Ion

| Compound                | CAS no.    | Mol form      | Target ion |
|-------------------------|------------|---------------|------------|
| Phenanthrene-d10*       | 1517-22-2  | C14D10        | 188        |
| Triphenyl phosphate*    | 115-86-6   | C18H15O4P     | 326        |
| Methamidophos           | 10265-9-26 | C2H8N02PS     | 94         |
| Dichlorvos              | 62-73-7    | C4H7Cl2O4P    | 109        |
| Omethoate               | 1113-02-6  | C5H12N04PS    | 156        |
| Monocrotophos           | 6923-22-4  | C7H14N05P     | 127        |
| BHC alpha isomer        | 319-84-6   | C6H6Cl6       | 181        |
| Hexachlorobenzene       | 118-74-1   | C6Cl6         | 284        |
| Dimethoate              | 60-51-5    | C5H12N03PS2   | 87         |
| BHC beta isomer         | 319-85-7   | C6H6Cl6       | 219        |
| Lindane                 | 58-89-9    | C6H6Cl6       | 181        |
| Pentachloronitrobenzene | 82-68-8    | C6Cl5NO2      | 237        |
| Diazinon                | 333-41-5   | C12H21N2O3PS  | 179        |
| BHC Delta isomer        | 319-86-8   | C6H6Cl6       | 181        |
| Propanil                | 709-98-8   | C9H9Cl2NO     | 161        |
| Methyl parathion        | 298-00-0   | C8H10N05PS    | 263        |
| Chlorpyrifos methyl     | 5598-13-0  | C7H7Cl3N03PS  | 286        |
| Vinclozolin             | 50471-44-8 | C12H9Cl2N03   | 212        |
| Heptachlor              | 76-44-8    | C10H5Cl7      | 272        |
| Metalaxyl               | 57837-19-1 | C15H21N04     | 206        |
| Fenitrothion            | 122-14-5   | C9H12N05PS    | 277        |
| Pirimiphos-methyl       | 29232-93-7 | C11H20N3O3PS  | 290        |
| Aldrin                  | 309-00-2   | C12H8Cl6      | 263        |
| Malathion               | 121-75-5   | C10H19O6PS2   | 173        |
| Fenthion                | 55-38-9    | C10H15O3PS2   | 278        |
| Chlorpyrifos            | 2921-88-2  | C9H11Cl3N03PS | 197        |
| Parathion               | 56-38-2    | C10H14N05PS   | 291        |
| Triadimefon             | 43121-43-3 | C14H16ClN3O2  | 57         |
| Heptachlor exo-epoxide  | 1024-57-3  | C10H5Cl7O     | 353        |
| Pendimethalin           | 40487-42-1 | C13H19N3O4    | 252        |
| Isofenphos              | 25311-71-1 | C15H24N04PS   | 213        |
| Quinalphos              | 13593-03-8 | C12H15N2O3PS  | 146        |
| Triadimenol             | 55219-65-3 | C14H18ClN3O2  | 112        |
| Methidathion            | 950-37-8   | C6H11N2O4PS3  | 145        |
| Butachlor               | 23184-66-9 | C17H26ClN02   | 176        |
| Dieldrin                | 60-57-1    | C12H8Cl6O     | 79         |
| Isoprothiolane          | 50512-35-1 | C12H18O4S2    | 118        |
| p,p'-DDE                | 72-55-9    | C14H8Cl4      | 246        |
| Endrin                  | 72-20-8    | C12H8Cl6O     | 263        |
| Myclobutanil            | 88671-89-0 | C15H17ClN4    | 179        |
| p,p'-DDD                | 72-54-8    | C14H10Cl4     | 235        |
| o,p'-DDT                | 789-02-6   | C14H9Cl5      | 235        |
| Ethion                  | 563-12-2   | C9H22O4P2S4   | 231        |
| p,p'-DDT                | 50-29-3    | C14H9Cl5      | 235        |

|                          |             |                |     |
|--------------------------|-------------|----------------|-----|
| Iprodione                | 36734-19-7  | C13H13Cl2N3O3  | 187 |
| Phosmet                  | 732-11-6    | C11H12NO4PS2   | 160 |
| Bifenthrin               | 82657-04-3  | C23H22ClF3O2   | 181 |
| Fenpropathrin            | 64257-84-7  | C22H23NO3      | 97  |
| Tetradifon               | 116-29-0    | C12H6Cl4O2S    | 159 |
| Phosalone                | 2310-17-0   | C12H15ClNO4PS2 | 182 |
| Cyhalothrin( $\lambda$ ) | 68085-85-8  | C23H19ClF3NO3  | 181 |
| Permethrin I             | 52645-53-1  | C21H20Cl2O3    | 183 |
| Permethrin II            | 999046-03-6 | C21H20Cl2O3    | 183 |
| Cypermethrin I           | 52315-07-8  | C22H19Cl2NO3   | 181 |
| Cypermethrin II          | 65731-84-2  | C22H19Cl2NO3   | 181 |
| Fenvalerate              | 51630-58-1  | C25H22ClNO3    | 167 |
| Deltamethrin             | 52918-63-5  | C22H19Br2NO3   | 181 |

\*Internal standard

## Experimental

The experiments were performed on an Agilent 7890 Gas Chromatograph equipped with 5975C inert MSD, and Agilent 7683 Automatic Liquid Sampler (ALS). The split/splitless inlets were fitted with long-lifetime septa (p/n 5183-4761) and splitless deactivated liner (p/n 5181-3316). Separation of the compounds was achieved on an Agilent J&W HP-5ms ultra inert GC column (30 m  $\times$  0.25 mm, 0.25  $\mu$ m). Extraction and cleanup were achieved with an Agilent Bond Elut QuEChERS EN Extraction kit (p/n 5982-5650) and a Bond Elut QuEChERS EN Dispersive SPE kit (p/n 5982-5156). Injections were made using 5  $\mu$ L syringe (p/n 5181-1273). The instrument conditions are listed below.

## Chemicals and standards

All pesticides standards were purchased from Sigma-Aldrich (St. Louis, MO, USA). All reagents and solvents were HPLC or analytical grade.

## Sample preparation

Organically grown, pesticide-free rice was purchased from a local supermarket. The rice was placed into a clean plastic bag and frozen at  $-20$   $^{\circ}$ C overnight. The following day, the required amount of frozen rice was removed and thoroughly blended. Dry rice was added while comminuting, when possible; samples were comminuted thoroughly to get the best sample homogeneity. It was verified that no pieces of rice were visible in the final sample.

A 5 g ( $\pm$  0.1g) amount of precisely homogenized sample was placed into a 50 mL centrifuge tube. The QC sample was forti-

fied with 100  $\mu$ L of appropriate QC spiking solution. 100  $\mu$ L of IS spiking solution was added to all the samples except the control blank. Tubes were capped and vortexed for 1 min. Five mL of water were added to each tube using the dispenser and let sit for 30 min. Tubes were capped and vortexed for 1 min. A 10 mL aliquot of ACN was added to each tube using the dispenser. Tubes were capped and shaken by hand for 1 min. An Agilent Bond Elut QuEChERS EN extraction salt packet, containing 4 g anhydrous  $MgSO_4$ , 1 g NaCl, 1 g  $Na_3$ Citrate, and 0.5 g  $Na_2H$  Citrate sesquihydrate, was added directly to each tube. The salt bag was massaged carefully to loosen any clumped salts before pouring. No powders were left in the threads or rims of the tubes. Tubes were sealed tightly and shaken vigorously for 1 min by hand to ensure that the solvent interacted well with the entire sample and crystalline agglomerates were broken up sufficiently. Sample tubes were centrifuged at 4,000 rpm for 5 min.

A 6 mL aliquot of upper ACN layer was transferred into Bond Elut QuEChERS EN dispersive SPE 15 mL tube (p/n 5982-5156). The tubes were capped tightly and vortexed for 1 min. The tubes were centrifuged with a standard centrifuge at 4,000 rpm for 5 min. A 500  $\mu$ L aliquot from the extract was transferred into an autosampler vial, and analyzed by GC/MS.

Table 2 shows the QuEChERS EN sample preparation procedure.

Table 2. The Agilent Bond Elut QuEChERS EN Extraction Procedure

- 1) Weigh 5 g comminuted rice in a 50 mL centrifuge tube and add 5 mL water.
- 2) Add IS solution, and GC spike solution if necessary, then vortex 1 min.
- 3) Add 10 mL of ACN, vortex 1 min..
- 4) Add Bond Elut EN QuEChERS Extraction salt packet.
- 5) Cap and shake vigorously for 1 min.
- 6) Centrifuge at 4,000 rpm for 5 min.
- 7) Transfer 6 mL of upper ACN layer to Bond Elut EN Dispersive-SPE 15 mL tube.
- 8) Vortex 1 min, then centrifuge at 4,000 rpm for 5 min.
- 9) Transfer 0.5 mL of extract to a sample vial.
- 10) Analyze sample with GC-MS.

## Results and Discussion

The total ion chromatogram (TIC) of the 57 pesticides at 200 ng/mL is shown in Figure 1. These pesticides include organochlorine, organophosphate and pyrethroid pesticides. Figure 2 shows that all target pesticides can be well separated by Agilent J&W HP-5MS Ultra Inert Capillary GC column.

The GC/MS system was retention time locked (RTL) to Chropyrifos-methyl. Retention time of each compound is listed in Table 3. Agilent 30 m × 0.25 mm × 0.25 μm HP-5ms was specified in G1672AA GC and GC/MS RTL Pesticide and Disruptor database, but HP-5ms Ultra Inert Capillary GC column exhibited the same selection selectivity of HP5 ms with same dimensions (see Table 3).




Figure 1. Total ion chromatogram (TIC) of pesticides at 200 ng/mL.

## Linearity and Recovery Tests for Target pesticides in Rice

Linearity was determined using calibration curves spiked into the rice matrix. Calibration curves were constructed from data obtained by 1  $\mu$ L injections of standards at, 50, 100, 200, 300, 400 ng/mL. Each standard solution contains 200 ng/mL of internal standards (ISTDs). All the pesticides of interest have excellent linearity with calibration coefficients ( $R^2$ ) greater than 0.991, with an average  $R^2$  of 0.9990 across 57 pesticides.

Different levels of target pesticides were spiked into rice before sample preparation to evaluate purifying effect of Bond Elut QuEChERS Extraction Kits. The GC/MS TIC for the rice extract and matrix spiked extract is illustrated in Figure 2. The spiked samples were treated according to the procedure described in the sample preparation. The recovery data for spiked samples are listed in Table 3. All data were based on five replicates of matrix spikes at each level.

The levels were 0.05, 0.1, 0.2  $\mu$ g/mL. Good recoveries were achieved for most of the compounds, ranged from 80-110%, and an average % relative standard deviation (RSD) of 5.35%.

## Conclusion

Agilent Bond Elut QuEChERS EN extraction and dispersive SPE kits provide a simple, fast and effective method for the purification and enrichment of representative volatile to semi-volatile pesticide residues in rice. J&W HP-5ms ultra inert column with GC/MS provide a good method for separation and detection. With Agilent DRS software and RTL pesticides library, fifty-five pesticides were identified and analyzed quickly and accurately. The impurities and matrix effects from rice did not interfere with the quantitation of target compounds. Since the selected pesticides represented a broad variety of different classes and properties, the Agilent Bond Elut QuEChERS kits is an excellent choice for other pesticides in similar food matrices.




Figure 2. TIC of the rice extract (A) and matrix spiked extract 100 ng/mL (B) using Agilent GC/MS System and Agilent J&W HP-5MS Ultra inert 30 m  $\times$  0.25 mm  $\times$  0.25  $\mu$ m column.

Table 3. Recoveries of 57 Pesticides in Real Sample

| Compound                | Retention time (min) | 0.05 µg/mL (n = 5) |      | 0.10 µg/mL (n = 5) |      | 0.20 µg/mL (n = 5) |      |
|-------------------------|----------------------|--------------------|------|--------------------|------|--------------------|------|
|                         |                      | Recovery%          | RSD% | Recovery%          | RSD% | Recovery%          | RSD% |
| Methamidophos           | 5.70                 | 70.8               | 7.6  | 76.3               | 5.3  | 78.9               | 7.9  |
| Dichlorvos              | 5.83                 | 112.3              | 5.1  | 103.3              | 7.9  | 109.8              | 5.2  |
| Omethoate               | 10.04                | 78.4               | 3.5  | 97.6               | 4.5  | 97.1               | 6.6  |
| Monocrotophos           | 11.75                | 81.6               | 5.9  | 96.8               | 4.5  | 105.5              | 3.9  |
| BHC alpha isomer        | 12.08                | 103.4              | 6.5  | 94.7               | 9.3  | 112.0              | 3.3  |
| Hexachlorobenzene       | 12.37                | 100.5              | 2.1  | 99.3               | 3.4  | 92.3               | 2.0  |
| Dimethoate              | 12.68                | 83.4               | 6.1  | 92.7               | 5.0  | 101.9              | 3.2  |
| BHC beta isomer         | 13.21                | 95.0               | 6.0  | 100.9              | 7.2  | 87.5               | 6.2  |
| Lindane                 | 13.45                | 102.4              | 7.3  | 110.8              | 7.1  | 93.5               | 7.0  |
| Pentachloronitrobenzene | 13.68                | 86.5               | 3.1  | 93.1               | 2.9  | 101.9              | 2.4  |
| Diazinon                | 14.48                | 92.4               | 3.3  | 96.0               | 3.4  | 104.9              | 2.3  |
| BHC Delta isomer        | 14.55                | 89.6               | 4.0  | 105.5              | 4.3  | 99.4               | 3.8  |
| Propanil                | 16.12                | 102.9              | 4.5  | 103.8              | 4.9  | 103.6              | 5.3  |
| Methyl parathion        | 16.59                | 87.5               | 7.7  | 102.8              | 7.4  | 103.2              | 3.2  |
| Chlorpyrifos methyl     | 16.60                | 105.3              | 3.0  | 103.9              | 3.3  | 104.8              | 2.4  |
| Vinclozolin             | 16.63                | 116.4              | 4.8  | 115.5              | 4.2  | 115.1              | 4.4  |
| Heptachlor              | 16.79                | 102.8              | 1.6  | 96.6               | 3.9  | 96.9               | 3.2  |
| Metalaxyl               | 17.35                | 105.8              | 3.6  | 104.5              | 5.4  | 108.7              | 5.0  |
| Fenitrothion            | 18.08                | 81.9               | 2.7  | 92.7               | 6.7  | 106.9              | 4.2  |
| Pirimiphos-methyl       | 18.31                | 108.9              | 2.7  | 107.3              | 3.2  | 107.8              | 3.4  |
| Aldrin                  | 18.51                | 106.0              | 2.3  | 105.0              | 3.5  | 99.9               | 1.4  |
| Malathion               | 18.81                | 97.6               | 5.0  | 95.1               | 7.4  | 111.8              | 3.3  |
| Fenthion                | 19.12                | 105.2              | 2.5  | 102.8              | 4.1  | 104.4              | 2.2  |
| Chlorpyrifos            | 19.24                | 105.3              | 2.5  | 104.7              | 3.1  | 102.4              | 4.0  |
| Parathion               | 19.27                | 72.8               | 4.6  | 77.6               | 7.3  | 112.9              | 3.3  |
| Triadimefon             | 19.40                | 99.9               | 6.2  | 93.3               | 7.0  | 102.5              | 4.0  |
| Heptachlor exo-epoxide  | 20.71                | 106.6              | 4.7  | 97.7               | 4.0  | 100.3              | 2.1  |
| Pendimethalin           | 20.99                | 117.8              | 4.6  | 92.0               | 4.6  | 102.5              | 2.6  |
| Isofenphos              | 21.60                | 109.2              | 5.8  | 102.2              | 8.2  | 106.9              | 3.9  |
| Quinalphos              | 21.64                | 111.0              | 6.0  | 110.3              | 7.8  | 105.7              | 4.1  |
| Triadimenol             | 21.71                | 96.6               | 6.2  | 107.3              | 6.4  | 92.0               | 9.4  |
| Methidathion            | 22.29                | 94.9               | 8.9  | 95.3               | 9.5  | 111.9              | 4.2  |
| Butachlor               | 23.23                | 117.6              | 5.8  | 110.2              | 8.5  | 102.9              | 3.6  |
| Dieldrin                | 23.85                | 112.3              | 5.0  | 112.4              | 6.7  | 105.5              | 4.7  |
| Isoprothiolane          | 23.88                | 123.7              | 3.9  | 100.5              | 6.6  | 108.5              | 2.4  |
| p,p'-DDE                | 24.02                | 100.7              | 4.4  | 100.8              | 5.3  | 104.5              | 2.8  |
| Myclobutanil            | 24.44                | 100.5              | 5.1  | 89.8               | 7.9  | 106.2              | 3.7  |
| Endrin                  | 24.74                | 108.6              | 5.2  | 97.2               | 5.9  | 107.4              | 5.6  |
| p,p'-DDD                | 25.67                | 98.2               | 5.9  | 99.6               | 6.6  | 102.1              | 5.9  |

|                      |       |       |     |       |     |       |     |
|----------------------|-------|-------|-----|-------|-----|-------|-----|
| o,p'-DDT             | 25.76 | 90.9  | 5.2 | 90.6  | 6.3 | 108.9 | 3.0 |
| Ethion               | 25.99 | 89.9  | 5.4 | 86.9  | 8.3 | 111.6 | 5.5 |
| p,p'-DDT             | 26.98 | 87.0  | 5.8 | 88.6  | 7.4 | 106.0 | 3.4 |
| Iprodione            | 28.39 | 99.6  | 5.5 | 108.6 | 5.3 | 98.9  | 6.2 |
| Phosmet              | 28.49 | 91.7  | 6.0 | 96.0  | 8.8 | 99.6  | 5.5 |
| Bifenthrin           | 28.84 | 104.8 | 5.0 | 104.2 | 7.9 | 108.5 | 7.1 |
| Fenpropathrin        | 28.98 | 106.2 | 6.0 | 109.0 | 8.1 | 104.7 | 4.6 |
| Tetradifon           | 29.36 | 114.3 | 5.1 | 107.5 | 7.1 | 109.2 | 6.7 |
| Phosalone            | 29.66 | 97.2  | 6.7 | 91.4  | 8.7 | 105.2 | 4.0 |
| Cyhalothrin (lambda) | 30.37 | 84.0  | 6.5 | 91.0  | 8.1 | 105.9 | 5.3 |
| Permethrin I         | 31.37 | 98.2  | 8.1 | 97.6  | 7.6 | 103.3 | 4.3 |
| Permethrin II        | 31.55 | 100.1 | 5.7 | 96.7  | 8.1 | 101.0 | 5.2 |
| Cypermethrin I       | 32.70 | 89.3  | 7.7 | 102.9 | 9.3 | 104.4 | 5.2 |
| Fenvalerate          | 34.30 | 75.0  | 9.1 | 93.3  | 9.0 | 107.5 | 5.3 |
| Cypermethrin II      | 34.71 | 86.4  | 8.5 | 97.2  | 8.6 | 100.1 | 7.7 |
| Deltamethrin         | 35.88 | 67.3  | 8.4 | 87.6  | 8.6 | 95.4  | 4.9 |

## References

1. M. Anastassiades, S. J. Lehotay, Fast and Easy Multiresidue Method Employment Acetonitrile Extraction/Partitioning and "Dispersive Solid-Phase Extraction" for the Determination of Pesticide Residues in Produce, *J. AOAC Int.*, 2003, 86, 412-431.
2. S. J. Lehotay, Use of Buffering and Other Means to Improve Results of Problematic Pesticides in a Fast and Easy Method for Residue Analysis of Fruits and Vegetables, *J. AOAC Int.*, 2005, 88, 615-629.
3. European Committee for Standardization/Technical committee CEN/TC 275 (2007), Foods of plant origin: Determination of pesticide residues using GC-MS and/or LC-MS/MS following acetonitrile extraction/partitioning and cleanup by dispersive SPE-QuEChERS method. European Committee for Standardization, Brussels.
4. P. Payá, M. Anastassiades, Analysis of pesticide residues using the Quick Easy Cheap Effective Rugged and Safe (QuEChERS) pesticide multiresidue method in combination with gas and liquid chromatography and tandem mass spectrometric detection, *Anal Bioanal Chem.*, 2007, 389, 1697-1714.
5. V. Giarocco, B. Quimby, and M. Klee, "Retention Time Locking: Concepts and Applications," Agilent Technologies, publication 5966-2469EN, [www.agilent.com/chem](http://www.agilent.com/chem).
6. Philip L. Wylie, "Screening for 926 Pesticides and Endocrine Disruptor by GC/MS with Deconvolution Reporting Software and a New Pesticide Library", Agilent Technologies, publication 5989-5076EN, [www.agilent.com/chem](http://www.agilent.com/chem).

## For More Information

These data represent typical results. For more information on our products and services, visit our Web site at [www.agilent.com](http://www.agilent.com).

[www.agilent.com](http://www.agilent.com)

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2011  
Printed in the USA  
May 16, 2011  
5990-8108EN



**Agilent Technologies**