

TheReporter

Reprinted from Volume 14, No. 3, 1995

T295023

© 1999 Sigma-Aldrich Co.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA • Sigma-Aldrich de Argentina, S.A. • Av. Pueyrredon 2446/50 • Piso 5-B • Buenos Aires 1119
AUSTRALIA • Sigma-Aldrich Pty. Ltd. • Unit #2, 14 Anella Avenue • Castle Hill NSW 2154
AUSTRIA • Sigma-Aldrich Handels GmbH • Hebbelplatz 7 • A-1110 Wien
BELGIUM • Sigma-Aldrich N.V./S.A. • K. Cardijnplein 8 • B-2880 Bornem
BRAZIL • Sigma-Aldrich Quimica Brasil Ltda. • Rua Sabara, 566-Conj. 53 • 01239-010 São Paulo, SP
CANADA • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC • Sigma-Aldrich s.r.o. • Pobrezni 46 • 186 21 Praha 8
DENMARK • Sigma-Aldrich Denmark A/S • Vejlegaardsvej 65B • DK-2665 Vallengbaek Strand
FINLAND • Sigma-Aldrich Finland/YA-Kemia Oy • Teerisuonkuja 4 • FIN-00700 Helsinki
FRANCE • Sigma-Aldrich Chimie • Chromatographie Supelco • L'Isle d'Abeau Chesnes - B.P. 701 • 38297 Saint-Quentin Fallavier Cedex
GERMANY • Sigma-Aldrich Chemie GmbH • Geschäftsbereich Supelco • Grünwalder Weg 30 • D-82041 Deisenhofen
GREECE • Sigma-Aldrich (o.m.) Ltd. • 72 Argonafton Str. • 16346 Ilioupoli, Athens
HUNGARY • Sigma-Aldrich Kft. • Nagy Diófa u. 7., IV fl. • H-1067 Budapest
INDIA • Sigma-Aldrich Co. • Survey No. 31/1, Sitharamapalaya • Mahadevapura P.O. • Bangalore 560 048
IRELAND • Sigma-Aldrich Ireland Ltd. • Airton Road • Tallaght • Dublin 24
ISRAEL • Sigma Israel Chemicals Ltd. • Park Rabin • Rohovot 76100
ITALY • Sigma-Aldrich s.r.l. • Via Gallarate, 154 • 20151 Milano
JAPAN • Sigma-Aldrich Japan K.K. • Division Supelco • JL Nihonbashi Building • 1-10-15 Nihonbashi Horidome-cho, Chuo-ku • Tokyo 103
KOREA • Sigma-Aldrich Korea • Samhan Camus Annex, 10th Floor • 17-26 Yoido-dong, Yungdeungpo-ku • Seoul
MALAYSIA • Sigma-Aldrich (M) Sdn. Bhd. • 9-2, Jalan 2/128, Taman Gembira • Off Jalan Kuchai Lama • 58200 Kuala Lumpur • Selangor
MEXICO • Sigma-Aldrich Quimica S.A. de C.V. • Calle 6 North No. 107 • Parque Industrial Toluca 2000 • 50200 Toluca
NETHERLANDS • Sigma-Aldrich Chemie BV • Postbus 27 • 3330 AA Zwijndrecht
NORWAY • Sigma-Aldrich Norway • Sandakerveien 102 • N-0483 Oslo
POLAND • Sigma-Aldrich Sp. z o.o. • Szelagowska 30 • 61-626 Poznań
PORTUGAL • Sigma-Aldrich Quimica, S.A. • P.O. Box 131 • Sintra 2710
RUSSIA • Sigma-Aldrich Russia • TOO Techmedbiochem • Makarenko Str. 2/21 • Building 1, Flat 22 • Moscow 103062
SINGAPORE • Sigma-Aldrich Pte. Ltd. • 102E Pasir Panjang Road • #08-01 Citilink Warehouse • Singapore 118529
SOUTH AFRICA • Sigma-Aldrich (pty) Ltd. • CNR Kelly & Ackerman Streets • Southern Life Industrial Park Unit • Unit 16/17 • Jet Park 1459
SPAIN • Sigma-Aldrich Quimica, S.A. • Apt. Correos 161 • 28100 Alcobendas, Madrid
SWEDEN • Sigma-Aldrich Sweden AB • Solkraftsvägen 14C • 135 70 Stockholm
SWITZERLAND • Supelco Switzerland • Industriestrasse 25 • P.O. Box 260 • CH-9471 Buchs
UNITED KINGDOM • Sigma-Aldrich Company Ltd. • Supelco UK • Fancy Road, Poole • Dorset BH12 4QH
UNITED STATES • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email:supelco@sial.com

H

This article is archived from a past issue of The Supelco Reporter. Information in the article was appropriate at the time of publication, but product specifications, catalog numbers, and availability may have changed over time.

If you have questions about applying methodology described in this article to a current application, please contact our technical service chemists.



Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

Monitoring Antioxidant Additives in Foods, Using HPLC or Capillary GC and a New Reference Standards Kit

E. Doughty, K. Kiefer

Supelco now offers an antioxidant standards kit containing all 9 antioxidants listed in Association of Official Analytical Chemists (AOAC) Official Method 983.15: Phenolic Antioxidants in Oils, Fats, and Butter Oil, plus ethoxyquin, an additional commonly used antioxidant. These compounds are separated well, in less than 13 minutes, on a SUPELCOSIL LC-18 HPLC column, under the conditions described in the AOAC method. Several of the antioxidants also can be resolved on a SAC-5 capillary GC column.


Rancidity caused by the oxidation of unsaturated fats in foods is a major problem for the food industry. With the trace metals present in foods acting as catalysts, the unsaturated fats react with oxygen, producing free radicals and peroxides. These compounds can destroy many vitamins, reducing the nutritional value of the food, destroy pigments, bleaching the food, and cause off-odors and off-flavors. While the portion of fat undergoing such changes can be small, the resulting off-odors are strong, and very much out of proportion to the amount of oil or fat involved.

Antioxidants are added to foods and other products to prevent rancidity. Although the exact mechanism(s) is (are) not clear, antioxidants react with free radicals and peroxides, slowing the development of rancidity. Certain other additives greatly enhance the effectiveness of antioxidants. Metal scavenger and chelating agents, such as citric acid and citrates, tie up the trace metals and greatly reduce their catalytic activity. Synergism between antioxidants also has been noted, and many commercial antioxidant mixes are formulated to contain mixtures of antioxidants. The most common of these mixtures contain both butyl-hydroxyanisole (BHA) and butyl-hydroxytoluene (BHT).

In order to assure consistent product quality, these additives must be monitored. In the past, food analysts had to go to several vendors to obtain good quality standards of all of the antioxidants they might be required to monitor. Now, Supelco offers an antioxidant standards kit, containing all 9 of the antioxidants listed in Association of Official Analytical Chemists (AOAC) Official Method 983.15: *Phenolic Antioxidants in Oils, Fats, and Butter Oil*, (Figure A) plus ethoxyquin, an additional antioxidant commonly used in spices and other food products, and in cosmetics.

The AOAC has granted final approval of Method 983.15, and has included this method in the 1995 (16th) edition of the *Official Methods of Analysis* of AOAC International (1). The method describes a linear gradient HPLC analysis for resolving the nine antioxidants shown in Figure A. Figure B shows a chromatogram generated on a 15cm x 4.6mm SUPELCOSIL™ LC-18 HPLC column, using the conditions described in the AOAC method. All 9 analytes were well

Figure A. Phenolic Antioxidants Commonly Used to Preserve Food Quality



795-0404-0413


resolved in less than 13 minutes. Ethoxyquin is not detected at the 280nm UV wavelength indicated by the AOAC method, and thus is not included here. Figure B also shows the analysis on a 25cm SUPELCOSIL LC-18 column. The analytes again are well resolved, and analysis time is only slightly prolonged.

In most food products, only 1 or 2 antioxidants are used. In some cases, gas chromatography can be used to monitor the antioxidants. (Many of the antioxidants in Figure A give a poor response to FID detection, however.) Figure C shows the chromatogram generated from a mix of 7 antioxidants, including ethoxyquin, separated on a 30m x 0.25mm x 0.25µm phase film SAC™-5 capillary GC column under isothermal (200°C) conditions. Resolution is good. Ethoxyquin was eluted at approximately 7 minutes. Note that at

Figure B. Antioxidants in AOAC Method 983.15 by HPLC

Column: **SUPELCOSIL LC-18, 15cm or 25cm x 4.6mm ID, 5µm particles**
 Cat. No.: **58230-U** (15cm column), **5-8298** (25cm column)
 Mobile Phase: A = 5% acetic acid in DI water
 B = acetonitrile:methanol, 50:50
 30% B to 100% B in 10 min, hold 10 min
 Flow Rate: 2mL/min
 Temp.: ambient
 Det.: UV, 280nm
 Inj.: 10µL, 20µg/mL each analyte

1. Propyl gallate
2. 2,4,5-Trihydroxybutyrophenone
3. tert-Butylhydroquinone
4. Nordihydroguaiaretic acid
5. 2- & 3-tert-Butyl-4-hydroxyanisole
6. 2,6-Di-tert-butyl-4-hydroxymethylphenol (lonox 100)
7. Octyl gallate
8. 3,5-Di-tert-butyl-4-hydroxytoluene
9. Lauryl gallate




200°C, an analysis time of more than 15 minutes was needed to elute propyl gallate. Lauryl gallate and octyl gallate cannot be monitored under these conditions – a very high oven temperature is required to elute these high boiling/low volatility compounds from the capillary column.

If your analyses include monitoring antioxidants, this new standards kit should meet all your routine needs, and will prepare you for identifying unanticipated additional antioxidants in your samples. Note that the ethoxyquin standard is qualitative. Because of its highly reactive nature, we cannot assure a purity level of more than 70% for this compound. Each antioxidant is packaged neat, under nitrogen.

Figure C. Antioxidants with Lower Boiling Points Can Be Monitored by GC

Column: **SAC-5, 30m x 0.25mm, 0.25µm film**
 Cat. No.: **24156**
 Oven: 200°C
 Carrier: helium, 30cm/sec (set at 200°C)
 Det.: FID (250°C)
 Inj.: 2µL, 200µg/mL each analyte, split 100:1 (250°C)



Ordering Information:

Description	Cat. No.
Phenolic Antioxidants Kit	47192
Individually packaged under nitrogen, neat, in quantities listed below.	

Component	CAS No.	Qty. (mg)
Propyl gallate (PG)	121-79-9	500
2,4,5-Trihydroxybutyrophenone (THBP)	1421-63-2	500
tert-Butylhydroquinone (TBHQ)	1948-33-0	500
Nordihydroguaiaretic acid (NDGA)	500-38-9	100
2- & 3-tert-Butyl-4-hydroxyanisole (BHA)	25013-16-5	500
2,6-Di-tert-butyl-4-hydroxymethylphenol (lonox 100)	88-26-6	100
3,5-Di-tert-butyl-4-hydroxytoluene (BHT)	128-37-0	500
Lauryl gallate (Dodecyl gallate)	1166-52-5	500
Octyl gallate	1034-01-1	500
Ethoxyquin*	91-53-2	500

*Qualitative standard. Purity >70% cannot be assured.

Description	Cat. No.
SUPELCOSIL LC-18 HPLC Columns	
15cm x 4.6mm, 5µm particles	58230-U
25cm x 4.6mm, 5µm particles	58298
SAC-5 Capillary GC Column	
30m x 0.25mm, 0.25µm film	24156

Reference

1. *Official Methods of Analysis* (16th ed.), Method 983.15 Association of Official Analytical Chemists, Arlington, VA USA (1995).

Reference not available from Supelco.

SAC and SUPELCOSIL are trademarks of Sigma-Aldrich Co.