

# TheReporter

Reprinted from Volume 15, No. 2, 1996

T296032

© 1999 Sigma-Aldrich Co.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website ([www.sigma-aldrich.com](http://www.sigma-aldrich.com)), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

**ARGENTINA** • Sigma-Aldrich de Argentina, S.A. • Av. Pueyrredon 2446/50 • Piso 5-B • Buenos Aires 1119  
**AUSTRALIA** • Sigma-Aldrich Pty. Ltd. • Unit #2, 14 Anella Avenue • Castle Hill NSW 2154  
**AUSTRIA** • Sigma-Aldrich Handels GmbH • Hebbelplatz 7 • A-1110 Wien  
**BELGIUM** • Sigma-Aldrich N.V./S.A. • K. Cardijnplein 8 • B-2880 Bornem  
**BRAZIL** • Sigma-Aldrich Quimica Brasil Ltda. • Rua Sabara, 566-Conj. 53 • 01239-010 São Paulo, SP  
**CANADA** • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8  
**CZECH REPUBLIC** • Sigma-Aldrich s.r.o. • Pobrezni 46 • 186 21 Praha 8  
**DENMARK** • Sigma-Aldrich Denmark A/S • Vejlegaardsvej 65B • DK-2665 Vallensbaek Strand  
**FINLAND** • Sigma-Aldrich Finland/YA-Kemia Oy • Teerisuonkuja 4 • FIN-00700 Helsinki  
**FRANCE** • Sigma-Aldrich Chimie • Chromatographie Supelco • L'Isle d'Abeau Chesnes - B.P. 701 • 38297 Saint-Quentin Fallavier Cedex  
**GERMANY** • Sigma-Aldrich Chemie GmbH • Geschäftsbereich Supelco • Grünwalder Weg 30 • D-82041 Deisenhofen  
**GREECE** • Sigma-Aldrich (o.m.) Ltd. • 72 Argonafton Str. • 16346 Ilioupoli, Athens  
**HUNGARY** • Sigma-Aldrich Kft. • Nagy Diófa u. 7., IV fl. • H-1067 Budapest  
**INDIA** • Sigma-Aldrich Co. • Survey No. 31/1, Sitharamapalaya • Mahadevapura P.O. • Bangalore 560 048  
**IRELAND** • Sigma-Aldrich Ireland Ltd. • Airton Road • Tallaght • Dublin 24  
**ISRAEL** • Sigma Israel Chemicals Ltd. • Park Rabin • Rohovot 76100  
**ITALY** • Sigma-Aldrich s.r.l. • Via Gallarate, 154 • 20151 Milano  
**JAPAN** • Sigma-Aldrich Japan K.K. • Division Supelco • JL Nihonbashi Building • 1-10-15 Nihonbashi Horidome-cho, Chuo-ku • Tokyo 103  
**KOREA** • Sigma-Aldrich Korea • Samhan Camus Annex, 10<sup>th</sup> Floor • 17-26 Yoido-dong, Yungdeungpo-ku • Seoul  
**MALAYSIA** • Sigma-Aldrich (M) Sdn. Bhd. • 9-2, Jalan 2/128, Taman Gembira • Off Jalan Kuchai Lama • 58200 Kuala Lumpur • Selangor  
**MEXICO** • Sigma-Aldrich Quimica S.A. de C.V. • Calle 6 North No. 107 • Parque Industrial Toluca 2000 • 50200 Toluca  
**NETHERLANDS** • Sigma-Aldrich Chemie BV • Postbus 27 • 3330 AA Zwijndrecht  
**NORWAY** • Sigma-Aldrich Norway • Sandakerveien 102 • N-0483 Oslo  
**POLAND** • Sigma-Aldrich Sp. z o.o. • Szelagowska 30 • 61-626 Poznań  
**PORTUGAL** • Sigma-Aldrich Quimica, S.A. • P.O. Box 131 • Sintra 2710  
**RUSSIA** • Sigma-Aldrich Russia • TOO Techmedbiochem • Makarenko Str. 2/21 • Building 1, Flat 22 • Moscow 103062  
**SINGAPORE** • Sigma-Aldrich Pte. Ltd. • 102E Pasir Panjang Road • #08-01 Citilink Warehouse • Singapore 118529  
**SOUTH AFRICA** • Sigma-Aldrich (pty) Ltd. • CNR Kelly & Ackerman Streets • Southern Life Industrial Park Unit • Unit 16/17 • Jet Park 1459  
**SPAIN** • Sigma-Aldrich Quimica, S.A. • Apt. Correos 161 • 28100 Alcobendas, Madrid  
**SWEDEN** • Sigma-Aldrich Sweden AB • Solkraftsvägen 14C • 135 70 Stockholm  
**SWITZERLAND** • Supelco Switzerland • Industriestrasse 25 • P.O. Box 260 • CH-9471 Buchs  
**UNITED KINGDOM** • Sigma-Aldrich Company Ltd. • Supelco UK • Fancy Road, Poole • Dorset BH12 4QH  
**UNITED STATES** • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email:supelco@sial.com

H

This article is archived from a past issue of The Supelco Reporter. Information in the article was appropriate at the time of publication, but product specifications, catalog numbers, and availability may have changed over time.

If you have questions about applying methodology described in this article to a current application, please contact our technical service chemists.



Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

# New Supelcarb™ Split Vent Trap Exhibits Longer-Lasting Trapping of Toxic Compounds

W. Betz, O. Cheesman, R. Bartram

*Indoor pollution in the industrial workplace is a continuing concern. Because of this concern, Supelco has designed a cartridge containing carbonaceous media to trap pollutants from the split vent outlet of a gas chromatograph. In our laboratory studies, this Supelcarb trap outperformed conventional activated charcoal traps.*

Over the past several decades, federal regulations have facilitated the reduction of airborne contamination. Recently, adsorbent traps have been developed for eliminating air pollution caused by the capillary split vent injection systems of modern gas chromatographs.

Typical split vent flow rates (10 to 100mL/min) carry most of an injected sample into the laboratory atmosphere. Many GC manufacturers recommend venting the system into a fume hood. However, this is often impossible due to space constrictions. Another approach is to place an adsorbent trap onto the split vent outlet.


Supelco has designed the Supelcarb split vent trap, a replaceable cartridge filled with a high-capacity carbonaceous adsorbent, to trap a broad range of organic compounds. The adsorbent's narrow particle size distribution and spherical shape allow for tighter packing and less gas channeling than the irregular-shaped activated charcoal particles.

To determine the longevity of the Supelcarb trap and how it compared to three competitors' split vent traps, we tested for and compared the breakthrough volume for each trap. A breakthrough volume is defined as the amount of gas required to pass through the adsorbent bed to elute one molecule of a compound introduced at the inlet of the bed. Once the compound is injected, it will move through the bed of the trap as long as flow is maintained. When an adsorbent bed is saturated, the concentration of the toxic compound at the outlet is temporarily higher than that at the inlet. Therefore a spent trap is more harmful than no trap at all. Timely replacement of traps cannot be overemphasized.

The study simulated a GC operation in a typical lab: a 30-minute GC run and split vent flow of 65mL/min, with an automated analytical system comprised of two GCs. A trap was connected to the split vent outlet of the first GC, and an autosampler injected one microliter of dichloromethane into the GC every 30 minutes. The second GC, with an automated six-port sampling valve, was programmed to sample effluent from the trap and deliver it to the column every 30 minutes. Dichloromethane was chosen as the analyte because it is one of the fastest eluting compounds, it is a commonly used solvent, and it is a suspected carcinogen.

The collected data show breakthrough on the Supelcarb trap occurs after 1458L of helium has traveled through the trap — approximately 2 weeks at 65mL/min. The other traps did not show a comparable capacity (Figure A). In fact the breakthrough times for the other traps are less than the manufacturer's recommended replacement times. We recommend using the Supelcarb trap for longer-lasting trapping and replacing it every two weeks.

**Figure A. Breakthrough on Split Vent Traps**



	Brand A	Brand B	Brand C	Supelcarb
Time before breakthrough:	22.5 hr.	49.5 hr.	249 hr.	374 hr.
Number of injections:	45	99	498	748
Manufacturer recommended replacement time:	168 hr.	300 injections	720 hr.	336 hr.
Breakthrough volume:	88 liters	193 liters	971 liters	1458 liters

## Ordering Information:

Description	Cat. No.
<b>Supelcarb Split Vent Trap Starter Kit</b> one trap and attachment kit	<b>22536</b>
<b>Supelcarb Replacement Traps</b>	
One trap	22535,01
Two traps	22535,02
Five traps	22535,05
<b>Polypropylene Reducing Unions</b> pk. of 10	<b>21999</b>

For more information about the Supelcarb split vent trap, request publication T396096. Supelcarb is a trademark of Sigma-Aldrich Co.