

TheReporter

Reprinted from Volume 16, No. 5, 1997

T297035

© 1999 Sigma-Aldrich Co.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA • Sigma-Aldrich de Argentina, S.A. • Av. Pueyrredon 2446/50 • Piso 5-B • Buenos Aires 1119
AUSTRALIA • Sigma-Aldrich Pty. Ltd. • Unit #2, 14 Anella Avenue • Castle Hill NSW 2154
AUSTRIA • Sigma-Aldrich Handels GmbH • Hebbelplatz 7 • A-1110 Wien
BELGIUM • Sigma-Aldrich N.V./S.A. • K. Cardijnplein 8 • B-2880 Bornem
BRAZIL • Sigma-Aldrich Quimica Brasil Ltda. • Rua Sabara, 566-Conj. 53 • 01239-010 São Paulo, SP
CANADA • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC • Sigma-Aldrich s.r.o. • Pobrezni 46 • 186 21 Praha 8
DENMARK • Sigma-Aldrich Denmark A/S • Vejlegaardsvej 65B • DK-2665 Vallensbaek Strand
FINLAND • Sigma-Aldrich Finland/YA-Kemia Oy • Teerisuonkuja 4 • FIN-00700 Helsinki
FRANCE • Sigma-Aldrich Chimie • Chromatographie Supelco • L'Isle d'Abeau Chesnes - B.P. 701 • 38297 Saint-Quentin Fallavier Cedex
GERMANY • Sigma-Aldrich Chemie GmbH • Geschäftsbereich Supelco • Grünwalder Weg 30 • D-82041 Deisenhofen
GREECE • Sigma-Aldrich (o.m.) Ltd. • 72 Argonafton Str. • 16346 Ilioupoli, Athens
HUNGARY • Sigma-Aldrich Kft. • Nagy Diófa u. 7., IV fl. • H-1067 Budapest
INDIA • Sigma-Aldrich Co. • Survey No. 31/1, Sitharamapalaya • Mahadevapura P.O. • Bangalore 560 048
IRELAND • Sigma-Aldrich Ireland Ltd. • Airton Road • Tallaght • Dublin 24
ISRAEL • Sigma Israel Chemicals Ltd. • Park Rabin • Rohovot 76100
ITALY • Sigma-Aldrich s.r.l. • Via Gallarate, 154 • 20151 Milano
JAPAN • Sigma-Aldrich Japan K.K. • Division Supelco • JL Nihonbashi Building • 1-10-15 Nihonbashi Horidome-cho, Chuo-ku • Tokyo 103
KOREA • Sigma-Aldrich Korea • Samhan Camus Annex, 10th Floor • 17-26 Yoido-dong, Yungdeungpo-ku • Seoul
MALAYSIA • Sigma-Aldrich (M) Sdn. Bhd. • 9-2, Jalan 2/128, Taman Gembira • Off Jalan Kuchai Lama • 58200 Kuala Lumpur • Selangor
MEXICO • Sigma-Aldrich Quimica S.A. de C.V. • Calle 6 North No. 107 • Parque Industrial Toluca 2000 • 50200 Toluca
NETHERLANDS • Sigma-Aldrich Chemie BV • Postbus 27 • 3330 AA Zwijndrecht
NORWAY • Sigma-Aldrich Norway • Sandakerveien 102 • N-0483 Oslo
POLAND • Sigma-Aldrich Sp. z o.o. • Szlagowska 30 • 61-626 Poznań
PORTUGAL • Sigma-Aldrich Quimica, S.A. • P.O. Box 131 • Sintra 2710
RUSSIA • Sigma-Aldrich Russia • TOO Techmedbiochem • Makarenko Str. 2/21 • Building 1, Flat 22 • Moscow 103062
SINGAPORE • Sigma-Aldrich Pte. Ltd. • 102E Pasir Panjang Road • #08-01 Citilink Warehouse • Singapore 118529
SOUTH AFRICA • Sigma-Aldrich (pty) Ltd. • CNR Kelly & Ackerman Streets • Southern Life Industrial Park Unit • Unit 16/17 • Jet Park 1459
SPAIN • Sigma-Aldrich Quimica, S.A. • Apt. Correos 161 • 28100 Alcobendas, Madrid
SWEDEN • Sigma-Aldrich Sweden AB • Solkraftsvägen 14C • 135 70 Stockholm
SWITZERLAND • Supelco Switzerland • Industriestrasse 25 • P.O. Box 260 • CH-9471 Buchs
UNITED KINGDOM • Sigma-Aldrich Company Ltd. • Supelco UK • Fancy Road, Poole • Dorset BH12 4QH
UNITED STATES • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email:supelco@sial.com

H

This article is archived from a past issue of The Supelco Reporter. Information in the article was appropriate at the time of publication, but product specifications, catalog numbers, and availability may have changed over time.

If you have questions about applying methodology described in this article to a current application, please contact our technical service chemists.



Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

Solid Phase Extraction of 200 Pesticides from Fruits and Vegetables, for Analysis by GC or HPLC

L. Nolan, Sample Handling, Supelco, Bellefonte, PA, USA

Recovery rates for polar analytes are higher, and less variable, with ENVI-Carb SPE tubes than with silica-based packings. The unique carbon-based packing in ENVI-Carb tubes is another product of our carbons program.

Investigators at Agriculture and Agri-Food Canada (Ottawa, Ontario) have developed a mult-residue cleanup and analysis for monitoring more than 200 organochlorine, organophosphorus, nitrogen-containing, and carbamate pesticides in fruits and vegetables. Their system allows quick screening of "rush" samples and thorough cleanup of complex samples. Initially, the extraction procedure required preparing and using minicolumns containing a mixture of charcoal and Celite® (1), but ENVI™-Carb carbon-based SPE tubes enabled these investigators to improve recovery rates and reproducibility for several pesticides. Subsequently, they replaced the charcoal/Celite minicolumns with ENVI-Carb tubes, and eliminated the labor-intensive minicolumn preparation process. Typical results are summarized in

Table 1 – data for many additional pesticides are presented in free technical Bulletin 900, available on request.

For additional cleanup of polar and nonpolar interferences in complex samples, the investigators routinely include a C18 tube and an amino (NH₂) - phase tube in the procedure. In the Canadian group's sample-screening applications, the analysis for organochlorine, organophosphorus, and nitrogen-containing pesticides (approximately 210 compounds) is by capillary gas chromatography with mass-specific detection; residues are identified by retention time and ion ratios. Analysis for 10 carbamates is by HPLC with UV detection. Analytes monitored in each GC run (approximately 120 compounds and 90 compounds, respectively) are listed, along with recovery rates and limits of detection, in (1). Postcolumn derivatization and HPLC columns and conditions are described in (2).

Relative to liquid-liquid extraction or SPE on silica-based packings with octyl (C8) or

octadecyl (C18) phases, ENVI-Carb packing provides more uniform recovery of a wide variety of polar nonvolatile analytes, while maintaining comparable results for less polar compounds. Because ENVI-Carb carbon is nonporous, samples can be processed rapidly—adsorption does not require dispersion of the analytes into porous regions. And, although the surface area of the carbon is smaller than that of porous silica (100m²/g, versus 400-600m²/g), the bed weight required for pesticide extractions typically is only half that needed with silica-based packings. The primary surface interaction mechanism is similar for ENVI-Carb packing and silica-based packings, but ENVI-Carb carbon also acts as a weak ion exchanger (3). This expands the range of analytes these tubes can be used to extract.

For details about the pesticide extraction procedure, request Bulletin 900.

Ordering Information:

Description	Cat. No.
ENVI-Carb Solid Phase Extraction Tubes	
3mL/250mg packing, pk. of 54	57088
6mL/250mg packing, pk. of 30	57092
6mL/500mg packing, pk. of 30	57094

References

1. Fillion, J, R. Hindle, M. Lacroix, and J. Selwyn, *Journal AOAC International*, **78**: 1252-1266 (1995).
2. Chaput, D., *J. Assoc. Off. Anal. Chem.*, **71**: 542-546 (1988).
3. Andreolini, F., *et al.*, *Anal. Chem.* **59**: 1720 (1987).
References not available from Supelco.

Trademarks

Celite – Celite Corp.
ENVI - Sigma-Aldrich Co.



Table 1. High Recovery of Pesticides from Fruits and Vegetables, Using ENVI-Carb SPE Tubes

Analyte	Mean % Recov.	Std. Dev.	% C.V.	Commodity / % Recovery			
				Kiwi Fruit	Beet	Squash	Green Peas
<i>Organochlorine Pesticides</i>							
p,p'-DDT	96.6	9.1	9.4	102	99	102	79
Endosulfan sulfate	97.4	9.1	9.4	106	93	101	84
Endrin	95.0	5.5	5.8	88	98	100	91
Heptachlor	92.6	5.0	5.4	90	92	96	86
Lindane (γ-BHC)	96.1	5.0	5.2	98	90	99	92
Mirex	96.4	8.9	9.2	102	98	102	78
Methoxychlor	99.1	5.9	5.9	104	98	103	91
<i>Nitrogen-Containing Pesticides / Triazines</i>							
Atrazine	94.7	14.7	15.5	67	103	104	100
Cyanazine	96.7	6.5	6.7*	—	93	98	89
Cyprazine	95.7	6.8	7.1*	—	87	94	101
Prometon	89.6	10.4	11.6	75	87	96	99
Simazine	95.1	15.8	16.6	67	102	107	103
<i>Organophosphorus Pesticides</i>							
Chlorpyrifos	97.9	5.9	6.0	99	94	104	93
Diazinon	96.9	5.6	5.8	90	99	99	96
Dimethoate	99.3	4.6	4.6	105	96	99	100
Ethoprophos	92.0	9.8	10.7	85	83	91	91
Malathion	97.6	6.3	6.5	99	92	98	93
Methamidophos	64.9	8.3	12.7	57	59	57	65
Terbufos	89.0	8.7	9.8	76	84	97	89

*n = 6 extractions, all others n = 7.

Data provided by J. Fillion and J.C. Selwyn, Agriculture and Agri-Food Canada, Ottawa, Ontario, Canada.