



Amines

Separation of several amines and alkanes

Application Note

Materials Testing & Research

Authors

Agilent Technologies, Inc.

Introduction

Gas chromatography using an Agilent CP-Wax 51 for Amines column separates six amines and alkanes in 20 minutes.




Agilent Technologies

Conditions

Technique : GC-capillary
Column : Agilent Cp-Wax 51, 0.22 mm x 25 m fused silica
WCOT CP-Wax 51 for amines (0.2 μ m)
(Part no. CP7405)
Temperature : 70 °C \rightarrow 200 °C, 5 °C/min
Carrier Gas : N₂, 60 kPa (0.6 bar, 8.7 psi), 19 cm/s
Injector : Splitter, 60 mL/min
T = 225 °C
Detector : FID, 16 x 10⁻¹² Afs
T = 225 °C
Sample Size : 0.1 μ L
Solvent Sample : hexane

Peak identification

1. hexylamine
2. n-tridecane
3. n-tetradecane
4. morpholine
5. n-pentadecane
6. aniline



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A00093



Agilent Technologies