



Solvents

Analysis of volatile amines in DMA

Application Note

Materials Testing & Research

Authors

Agilent Technologies, Inc.

Introduction

Gas chromatography with an Agilent PorabOND U column separates 18 solvents in less than 11 minutes.



Agilent Technologies


Conditions

Technique : GC-capillary
Column : Agilent PoraBOND U, 0.32 mm x 25 m fused silica
PLOT (df = 7 µm) (Part no. CP7381)
Temperature : 100 °C (1 min) Æ 200 °C, 10 °C/min
Carrier Gas : He, 50 kPa (0.5 bar, 7 psi)
Injector : Split,
T = 250 °C
Detector : POD in HID mode;
T = 250 °C
Concentration Range : % level

Courtesy : C. Duvekot, Agilent application laboratory,
Middelburg, The Netherlands

Peak identification

1. methane
2. ethane + ethylene
3. propane
4. methylchloride
5. methanol
6. vinylchloride
7. isobutane
8. butane
9. acetaldehyde
10. ethanol
11. isopentane
12. pentane
13. 2-propanol
14. acetone
15. 1-propanol
16. tert-butanol
17. vinylacetate
18. ethylacetate



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01563



Agilent Technologies