



C₁ – C₄ hydrocarbons

Analysis of acetylenes mixture

Application Note

Energy & Fuels

Authors

Agilent Technologies, Inc.

Introduction

Often ppm levels of hydrocarbon impurities must be measured and the response for such low levels must be accurate and reproducible over time. The Agilent Select Al₂O₃ MAPD is extensively deactivated, which results in highest response for traces of polar hydrocarbons including acetylenes and dienes. Selectivity of this Al₂O₃ PLOT column is very high and it separates all C₁ - C₅ hydrocarbons. Temperature stability is 200 °C. Also, for long term monitoring of impurities in hydrocarbon streams this column is the best choice.



Agilent Technologies


Conditions

Technique : GC
Column : Agilent Select Al₂O₃ MAPD, 0.53 mm x 50 m fused silica (Part no. CP7432)
Temperature : 40 °C, 5min → 160 °C, 10 °C/min → 200 °C, 20 °C/min, hold 1 min
Carrier Gas : He, 4 psig, 4 min → 11 psig, 0.5 psig/min, 2 min
Injector : Split 60 mL/min
Detector : FID
Concentration Range : approx 100 ppm in nitrogen, synthetic standard

Courtesy : J. Luong, Dow Chemical Canada

Peak identification

1. methane
2. ethane
3. ethylene
4. n-butane
5. propadiene
6. 1-butene
7. iso-butene
8. 1,2-butadiene
9. 1,3-butadiene
10. ethyl acetylene



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A02041



Agilent Technologies