



Mineral oil in soil and water according to DIN EN ISO 9377-2

Fast analysis of diesel contamination according to DIN EN ISO 9377-2

Application Note

Environmental

Authors

Agilent Technologies, Inc.

Introduction

The analysis of mineral oil can be done highly efficiently using GC and the Agilent Select Mineral Oil column . This column was optimized for mineral oil analysis to generate the shortest analysis time. The method used is DIN-EN ISO 9377-2 which replaces the DIN H53. The Select Mineral Oil stationary phase was tuned for separation and stabilized for high temperature operation. Upper temperature limit of this column is 400 °C.





Agilent Technologies

Conditions

Technique : GC
Column : Agilent Select Mineral Oil, 0.32 mm x 15 m fused silica (optimized film thickness) (Part no. CP7491)
6 m x 0.53 mm, methyl deactivated
Temperature : 55 °C, 1.9 min → 320 °C, 80 °C/min
Carrier Gas : Nitrogen, 80 kPa
Injector : On-column
Detector : FID
Sample Size : 2 µL
Sample Size : Diesel, 1.1 mg/L and 0.11 mg/L in petroleum ether

Courtesy : Thomas Karle, Chemisches Labor; Dr. Vogt,
Karlsruhe, Germany



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A02047



Agilent Technologies