



## Quantitative Analysis of Trace Amount of Butadiene Rubber in High Impact Polystyrene (HIPS)

High impact polystyrene (HIPS) is a polystyrene (PS) copolymerized with a few percents of butadiene rubber in order to increase impact resistance of PS. Pyrolysis GC is used as a tool for quality control of HIPS. Quantitative analysis of a trace amount of butadiene present in HIPS is described here using Double-Shot Pyrolyzer®.

Fig. 1 shows a pyrogram of HIPS obtained at 550°C. Butadiene (B) and styrene (S), monomer components of HIPS, have been detected. Table 1 shows peak area ratios of B and S (B/S) in the pyrograms obtained repeatedly. The excellent reproducibility of 0.62% was obtained. In the actual quantitative analysis, various mixing ratios of samples are used to produce a calibration curve.

Table 1. Peak Area Ratios of Butadiene (B) and Styrene (S) and Reproducibility

n	B/S (%)
1	1.249
2	1.244
3	1.246
4	1.254
5	1.233
Average	1.245
RSD (%)	0.62




Fig. 1. Pyrogram of HIPS

Pyrolysis temp: 550°C, Split ratio: 1/50, GC oven temp: 40–320°C at 20°C/min  
Separation column: Ultra ALLOY+5 (5% diphenylpolysiloxane), Length 30m, Id 0.25mm, Film thickness 0.25µm  
Sample amount : 100µg, Detector : Hydrogen flame detector (FID)

Keyword : Impact Resistance, Polystyrene, Butadiene, Quantitative Analysis, Reproducibility

Applications : Plastics and Rubber Industry, General Polymer Analysis

Please forward your inquiries via our web page at: (<http://www.frontier-lab.com/>), or send us a fax message.

R&D and manufactured by:  
**Frontier Laboratories Ltd.**

1-8-14, Saikon, Koriyama,  
Fukushima, 963-8862 Japan

Phone: 81-24-935-5100 Fax: 81-24-935-5102

© : Registered trademark of Frontier Laboratories Ltd.

Your dealer: