


Sample Prep for Chromatographic Analysis of Difficult Matrixes




Real World & Real Samples


Urine Sample without sample prep



Urine Sample with sample prep




Sources of Chromatographic Errors



(R.E. Majors, LC/GC Magazine, 1991, 1997, 2002)

Time Spend on Analytical Process



(R.E. Majors, LCGC Magazine, 1991, 1997, 2002)

Sample Prep Innovations

- **Solid Phase Microextraction (SPME)**
- High specificity SPE
- Dispersive SPE
- Silver Ion SPE for FAMES
- Carbonaceous adsorbents
- Flash chromatography



Solid Phase Microextraction (SPME)

“Sample Prep Made Easy”

Enrichment technique mainly for trace analysis

Developed in collaboration with Janusz Pawliszyn, Univ. of Waterloo

Unique and proprietary to Supelco

Users are...

- GC and GC-MS analysts (**HPLC & LC-MS**)
- Analyzing compounds in gases, liquids or solids.

Interested in...

- Sample enrichment
- Solventless extraction
- Using existing GC & HPLC systems
- Economical sample prep
- **Reducing lab animal sacrifice**


Users can expect...

- Highly consistent, quantifiable results from low concentrations of analytes

Odor-Causing Compounds in Water at 2 ppt (GC/MS)

Sensitive


1. 2-Isopropyl-3-methoxypyrazine (IPMP)
2. 2-Isobutyl-3-methoxypyrazine (IBMP)
3. 2- Methylisoborneol (MIB)
4. 2,4,6-Trichloroanisole (I.S. 8ppt)
5. (\pm) Geosmin



7

Linearity of Odor-Causing Compounds from Water at ppt Levels (SPME-GC/MS)

Quantitative



SPME Overview

- Solvent-free extraction technique for nearly any sample or matrix
- Alternative to head-space GC and solid phase extraction (SPE) techniques
- Directly interfaced with GC analysis
- Non-destructive to sample
- Reusable (100+ times)
- Inexpensive
- Fast



Assembled SPME fiber and holder with fiber immersed in a liquid sample.



Manual SPME holder and inlet guide.

The SPME Concept



Sample Adsorption

Please click on the numbered steps below for an animated sequence of the instruction.

- 1 Drill down septum piercing needle to avoid breakage
- 2 Insert needle into container
- 3 Adjust needle depth for aqueous sampling or headspace sampling
- 4 Extend plunger to expose fiber
- 5 Retract fiber before removing to avoid damaging the fiber.
- 6 Drill down septum piercing needle to avoid breakage.
- 7 Remove SPME Device

©2003, Sigma-Aldrich, All Rights Reserved


Click [here](#) for animation

SPME Fiber Coating: The Business End

An equilibrium is set up between analytes dissolved in the sample (solution or gas phase) and in the liquid coating on the fiber.

The fiber coating consists of:

- GC-type phases
- Particles



Distribution Constant

Concentration of analyte in stationary phase compared to concentration of analyte in solution:

$$K = n_s / V_1 C_2^\circ$$


K = Distribution constant

n_s = Moles of analyte in stationary phase

V_1 = Volume of stationary phase

C_2° = Final analyte concentration in water

Adsorption Mechanism for SPME



Absorbent vs. Adsorbent Fibers

Absorbent-type fibers (Film-type fibers)

Analytes are extracted by partitioning

- Liquid phase
- Retains by thickness of coating

Analytes do not compete for sites

Fibers can have high capacity



Adsorbent-type fibers (Particle-type fibers)

Physically traps or interacts with analytes

- Porous particles
- High surface area

Analytes may compete for sites

Fibers have limited capacity



Types of SPME Fiber Coatings

Films – Ab sorption:


Coating	Type	Polarity
7 μm Polydimethylsiloxane (PDMS)	Absorbent	Nonpolar
30 μm PDMS	Absorbent	Nonpolar
100 μm PDMS	Absorbent	Nonpolar
85 μm Polyacrylate (PA)	Absorbent	Polar
60 μm PEG (Carbowax)	Absorbent	Polar

Particles – Ad sorption:

Coating	Type	Polarity
85 μm Carboxen-PDMS	Adsorbent	Bipolar
65 μm PDMS-DVB	Adsorbent	Bipolar
55 μm /30 μm DVB/Carboxen-PDMS	Adsorbent	Bipolar
15 μm Carbowax Z-PDMS	Adsorbent	Bipolar


PDMS-DVB Fiber SEM

Cross section of the PDMS-DVB fiber. The center is a fused silica core, surrounded by a Stableflex core. The 3-5 μ m DVB particles are suspended in PDMS and layered over the cores. 275x magnification.



PDMS-Carboxen Fiber SEM

3000X magnification of the Carboxen PDMS coating. The 3-5 μ m Carboxen-PDMS particles are suspended in PDMS.




Carboxen™ Particle – Volume Contribution

Contribution of pore types to total Carboxen pore volume:

micropores (2-20Å) = 0.29 mL/g

mesopores (20-500Å) = 0.26 mL/g

macropores (>500Å) = 0.23 mL/g



Physical Properties of Divinylbenzene and Carboxen 1006

Material	Surface Area (m ² /g)	Porosity (mL/g)*		
		macro	meso	micro
Divinylbenzene	750	0.58	0.85	0.11
Carboxen™ 1006	720	0.23	0.26	0.29

*Macropore = >500Å

Mesopore = 20-500Å

Micropore = 2-20Å


Comparison of SPME Fibers for the Extraction of Small Hydrocarbons

(Analytes at 1 ppm in air, extracted for 10 min.)


Analyte	Absorbent	Adsorbent	
	100 μ m PDMS	PDMS/DVB	Carboxen/PDMS
Ethane	0	0	750
Propane	0	0	20000
Butane	0	340	72100
Pentane	230	2150	108000
Hexane	460	9280	105000

(Absolute area responses)

Molecular Weight Range for SPME Fibers



Area Response vs. Fiber Type





Effects of Fiber Polarity & Coating Thickness

Fiber Polarity

- Analyte selectivity
- Better recovery of polar analytes
- PEG
- Polyacrylate

Coating Thickness

- Analyte selectivity
- Extraction time
- Sample capacity
- Desorption time and carryover

Effects of Phase Coating Thickness of PDMS on Analyte Recovery Relative to Chrysene*

Analyte	%Relative Recovery		
	100 μ m	30 μ m	7 μ m
Benzene	2	1	<1
Toluene	5	1	<1
Naphthalene	13	4	1
Phenanthrene	37	27	16
Anthracene	49	38	32
Pyrene	69	54	47
Benzo(a)anthracene	105	91	96
Chrysene	100	100	100
Benzo(a)pyrene	119	127	131
Indeno(1,2,3-cd)pyrene	61	140	148
Benzo(g,h,i)perylene	61	117	122

*Absolute response of chrysene set to 100%

Factors Affecting Extraction Recovery

Salts and pH

Headspace vs. direct extraction

Inlet liner volume

Stirring (sample) & agitation (fiber)

- Increases precision
- Reduces time to reach equilibrium
- Must be consistent for all analyses
- Required for analytes with high distribution constants
- Sonication may increase temperature



Effects of Salt and pH

Salt usually increases analyte uptake

Use 25-30% NaCl to salt-out samples

Salt is not necessary for large non-polar analytes, such as PAHs and large hydrocarbons, and may reduce recovery

Lower pH to extract acidic compounds

Raise pH to extract basic compounds

Beware of stability of analytes at different pH levels


The Effect of Salt and pH on Extraction of Phenols by SPME

	No Salt Neutral	No Salt pH = 2	Salt Neutral	Salt pH = 2
Phenol	810	1003	6425	6150
Methylphenol	761	882	5485	7434
2-Nitrophenol	422	474	311	2315
2,4-Dimethylphenol	1344	1476	15000	20710
2,4-Dichlorophenol	5396	8138	19803	61664
2,4,5-Trichlorophenol	3115	11097	24270	96333
2,4-Dinitrophenol	0	11	765	1182
4-Nitrophenol	626	730	6536	11438
2,3,4,6-Tetrachlorophenol	3108	27683	33938	70440
2-Methyl-4,6-dinitrophenol	55	47	920	1685
Pentachlorophenol	2305	40582	22056	143905

Headspace vs. Direct Immersion


Analytical considerations:

- Volatility of sample
- Extraction time concerns
- Sample matrix
- Selectivity of analytes




Inlet Liner Volume: Comparison for Analysis of Gaseous VOCs by SPME

**Standard Splitless Liner,
2mm ID
50ppb**



**0.75mm ID Inlet Liner
50ppb**



794-0050, 0051

SPME Automation

Compatible with common GC autosamplers (Gerstel MPS, CTC Combi PAL, *etc.*)


Improves reproducibility by automating important variables:

- Heating
- Agitation
- Equilibration time



SPME automation video (~2 mins.)


Peppermint Oil in Chocolate Cookie Bar



Milk Sample Off-Flavors by SPME-GC/MS



Prior to Exposure to Sunlight




SPME Fiber: 75 μ m PDMS/Carboxen

Sample: 3g of 2% milk + 10 μ L internal standard solution, (20 μ g/mL 4-methyl-2-pentanone) (9mL GC vial)

Column: Supel-Q™ PLOT, 30m x 0.32mm ID

Det.: GC/MS ion trap, m/z = 33-300

After 1-Hour Exposure to Sunlight




- 1. Acetone
- 2. 2-Butanone
- 3. 3-Methylpentane
- 4. Pentanal
- 5. Dimethyldisulfide
- 6. Hexanal
- IS. 4-Methyl-2-pentanone

Chromatogram provided by Ray Marsili, Dean Foods Technical Center, Rockford, IL, USA.

G00507, 508
98-0385

Residual Solvents in Commercial Ibuprofen



1. Acetaldehyde
2. Ethanol
3. Acetonitrile
4. Acetone
5. 2-Propanol
6. 2-Methylpentane
7. 3-Methyl pentane
8. Hexane
9. Ethyl acetate
10. 2,2-Dimethylpentane
11. 2,4-Dimethylpentane
12. Methylcyclopentane

10ppb Nitrosamines in Water: SPME-GC/MS

Sample: analytes in (water + 25% KCl, pH 10)

SPME Fiber: 65µm PDMS-DVB

Extraction: immersion, 15 min (rapid stirring)

Desorption: 270°C, 1 min


**Column: PTA-5 (amine deactivated,
30m x 0.32mm ID, 0.5µm film)**

Oven: 50°C (1 min) to 250°C at 10°C/min,
hold 2 min

Carrier: helium, 30cm/sec

Det.: GC/MS (quadrupole, SIM)

Inj.: splitless, 250°C (0.75mm ID liner)



34
Chromatogram courtesy of J. Clark, Liggett Group, Inc.


96-0142




New Development: Biocompatible Fiber Pipette Tips for Solvent Extraction




Single Use Biocompatible Fiber Probes for *in vivo* Analysis



Comparison of SPME *in-vivo* PK Study of Carbamazepine from Mice Whole Blood to Extracts of Plasma Removed from Mice



SPME fiber Holder with Automated DESI-1D Source



Courtesy of
Joseph Kennedy
of Prosolia

Solid Phase Microextraction (SPME) Products

Fibers

Holders

- Manual
- For autosamplers

Accessories

Instructions

Applications on CD

sigma-aldrich.com/spme



SUPELCO

Supelco Park
Belafonte, PA 16823-0048 USA
Telephone 800-247-6629 • 814-353-3441
Fax 800-447-3044 • 814-353-3044
Email: supelco@sig.com
<http://www.sigma-aldrich.com/supelco>

Bulletin 929

A Practical Guide to Quantitation with Solid Phase Microextraction

Solid Phase Microextraction (SPME) is an innovative, solvent free technology that is fast, economical, and versatile. SPME has gained wide spread acceptance as the technique of preference for many applications. This guide presents a practical introduction to quantitation using the technique based on your type of sample. We present the factors that will influence your accuracy and precision and the different quantitation approaches that you can use. To help you further, we provide specific examples for each of the different approaches discussed and suggested references for additional reading.*

Introduction 2
Quantitation Guide Table
Approaches to Quantitation
Tips to Improve Accuracy
Conclusion
Helpful Hints
Flavors & Fragrances • Food & Beverage
Environmental • Pharmaceutical • Forensics

7th Edition
Solid Phase Microextraction CD
SPME Application Reference Guide
Plus Additional SPME Literature

SUPELCO

SPME

Includes Video Demos

SUPELCO
Analytical

sigma-aldrich.com/spme

SIGMA-ALDRICH

Sample Prep Innovations

- Solid Phase Microextraction (SPME)
- **High specificity SPE (SupelMIPs)**
- Dispersive SPE
- Silver Ion SPE for FAMES
- Carbonaceous adsorbent
- Flash chromatography

Users are...

- Analytical chemists (LC, LC-MS, GC...)

Interested in...

- Very selective extraction
- Analysis at extremely low concentrations (ppb, ppt)
- Increasing specificity of sample prep from complex matrixes

Users can expect...

- More rigorous washing to remove matrix
- Detect at lower levels

High Specificity Sample Prep


The specific innovation we will describe:
SupelMIP Molecularly Imprinted Polymers

- SPE tubes
- 96-well plates




The Molecular Imprinting Process

Molecularly imprinted polymers (MIPs) are polymers that have been prepared by polymerizing either pre-formed or self-assembled monomer-template complexes together with a cross-linking monomer. After removal of the template molecule, a polymer with binding sites for the template is obtained.




The MIP Binding Site

Graphical representation of the MIP binding site, which contains a cavity of the right size and attractive molecular features that can bind to the target molecule(s).




SupelMIP Chloramphenicol: Analysis in Honey

Chloramphenicol is an antibiotic that is monitored in honey.




Background from honey sample cleaned by SupelMIP SPE and LLE for Chloramphenicol analysis



Comparison of ion suppression effect between different clean-up methods for honey. Samples were post-spiked with CAP prior to analysis.

Overview of a Typical SupelMIP SPE Procedure

Very simple methods.
Full protocols are included with each MIP product.
Protocols may require optimization depending on the sample matrix.



SupelMIP Products

- **PAHs** in edible oils
- **Nitroimidazoles** in milk, eggs and other foods
- Nonsteroidal anti-inflammatory drugs (**NSAIDS**) in wastewater and other matrices
- **Fluoroquinolones** in bovine kidney, honey and milk
- **Amphetamines** and related compounds in urine
- **Chloramphenicol** in plasma, urine, milk, honey and shrimp
- **NNAL** - nitroso compound in urine
- **TSNAs** - tobacco specific nitrosamines in urine and tobacco
- **β -agonists** and **β -blockers** in tissue, urine and wastewater
- **Clenbuterol** in urine
- **Triazines** in water
- **Riboflavin** in milk

sigma-aldrich.com/supelmip

SupelMIP™ Solid Phase Extraction
Molecularly Imprinted Polymers for the Highly Selective Extraction of Trace Analytes from Complex Matrices

SupelMIP™ SPE - Beta-receptors (beta-agonists and beta-blockers)

Product Description:
Molecularly imprinted polymers (MIPs) are a class of highly cross-linked polymer-based molecular recognition elements engineered to bind one target compound or a class of structurally related target compounds with high selectivity. Selectivity is introduced during MIP synthesis in which a template molecule, designed to mimic the analyte, guides the formation of specific cavities or imprints that are sterically and chemically complementary to the target analyte(s). It is therefore critical for analysts to use the methodology described below when using this phase. Conventional generic methodologies employed with conventional SPE chemistries (e.g., reverse-phase C18) will yield sub-optimal results when employed with this phase.

The following methods have been developed for the selective extraction of both beta-agonists and beta-blockers from human urine and waste water. The method is highly reproducible and offers extraction recoveries of 80-90% from 1 mL urine with minimal ion suppression when analyzing via LC-MS/MS conditions.

Extraction Procedure: A flow rate of ~0.5 mL/min. is recommended. For analyte elution a flow rate of ~0.2 mL/min. is recommended.

Application Name:	Extraction of beta-agonists and beta-blockers from urine or waste water
Analyte:	Altenox®, carazolol, clenbuterol, metoprolol, propranolol, ribotrine, salbutamol, terbutaline, timolol, terbutaline
Sample Matrix:	Urine or waste water
General Comments:	The method is optimized for the simultaneous extraction of both beta-agonists and beta-blockers from urine and waste water. The protocol can be optimized for specific analytes. For more information regarding method optimization for specific beta-agonists and/or beta-blockers, please contact Supelco and/or MIP Technologies All technical service. Note that extraction recoveries, especially for more polar compounds such as salbutamol, terbutaline and albuterol may be reduced when loading sample volumes larger than recommended.
SupelMIP SPE - Full beta-receptor:	25 mg/10 mL (LRC) (Cat. No. 53223-U) or 25 mg/5 mL (Cat. No. 53224-U)
Sample Pre-treatment:	For particulate laden urine samples, centrifuge at 2000 g for 10 min. Dilute 1:1 with DI Water. Add deuterium labeled internal standards as necessary. Adjust to pH 7. For particulate laden water samples, filter water sample through 1 µm filter paper.
1. Condition/eqilibrate cartridge with:	<ul style="list-style-type: none"> • 1 mL acetonitrile • 1 mL DI water
2. Load sample: Note: recommended flow rate ~0.5 mL/min.	Apply up to 2 mL diluted urine sample or 10 mL waste water sample to the cartridge
3. Wash (Interference elution): Note: Apply gentle vacuum between each wash step.	<ul style="list-style-type: none"> • 3 x 1 mL DI water (removal of salt and matrix interferences) • Apply 2 min. of full vacuum to dry the tube • 1 mL acetonitrile (selective removal of hydrophobic interferences) • 1 mL 0.1% acetic acid/DI Water (selective removal of hydrophilic interferences) • Apply 2 min. of full vacuum to dry the tube
4. Analyte elution: Note: recommended flow rate ~0.2 mL/min.	Elute beta-agonists and beta-blockers with 2 x 1 mL 1% formic acid in acetonitrile. Evaporate and reconstitute with LC mobile phase prior to analysis.
Recommended Analytical Technique:	LC-MS/MS

Topics: Sample Prep Innovations

- SPME (solid phase microextraction)
- High specificity SPE (SupelMIP)
- **Dispersive SPE**
- Silver Ion SPE for FAMES
- Carbonaceous adsorbents
- Flash chromatography

Users are...

- **Food safety analysts**

Interested in...

- **Multi-residue pesticide analysis in food and agricultural products**

Users can expect...

- **Quick, easy, inexpensive extraction method**

Dispersive SPE (dSPE or QuEChERS) Multiresidue Pesticide Method

Multi-residue (100's) pesticide analysis
Retains/removes key interferences in food samples
Analytes are un-retained

Quick (~30 min./6 samples)

Easy (no laborious steps)

Cheap

Effective (wide scope, low consumption)

Rugged (minimal sources of errors)

Safe (solvents and techniques)



Dispersive SPE Procedure



Procedure:

1. Food initially extracted with aq. miscible solvent (e.g. ACN)
2. High amounts of salts (NaCl, Mg-sulfate) and buffering agents added to induce phase separation and stabilize acid/base labile pesticides
3. Shake/centrifuge. Isolate aliquot of sup for SPE clean-up.
4. Transfer supernatant to centrifuge tube. Add bulk SPE phase(s) and salts. Shake/vortex. Centrifuge and analyze supernatant.

Standard dSPE product line configured for:

- CEN Standard Method EN – 15662
- AOAC Method 2007.01


Full details of the simple protocol is included with the product.

GC-MS of Pesticides from Oranges Following Extraction with dSPE



column: SLB-5ms, 30 m x 0.25 mm I.D, 0.25 μ m (28471-U)
 oven: 100 °C (1 min.), 10 °C/min. to 300 °C (5 min.)
 inj.: 250 °C
 MSD interface: 300 °C
 scan range: selected ion monitoring (SIM), 7 monitoring groups used
 carrier gas: helium, 1 mL/min constant
 injection: 1 μ L, pulsed (20 psi until 0.20 min.), splitless (1.0 min.)
 liner: 4 mm I.D., single taper

- | | | |
|-------------------------|----------------------------|------------------------|
| 1. Methamidophos | 11. Carbaryl | 21. Folpet |
| 2. Dichlorvos | 12. Dichlofluanid | 22. cis-Chlordane |
| 3. Acephate | 13. Chlorpyrifos | 23. Imazalil |
| 4. Propoxur | 14. p-Dichlorobenzophenone | 24. 4,4'-DDE |
| 5. Ethoprophos (I.S.) | 15. Cyprodinil | 25. Dieldrin |
| 6. Hexachlorobenzene | 16. Pencanazole | 26. Endosulfan sulfate |
| 7. γ -BHC | 17. Tolyfluanid | 27. Dicofol |
| 8. Diazinon | 18. Heptachlor epoxide | 28. cis-Permethrin |
| 9. Chlorothalonil | 19. Captan | 29. trans-Permethrin |
| 10. Methyl chlorpyrifos | 20. Thiabendazole | 30. Coumaphos |



Dispersive SPE Products

Centrifuge tubes containing pre-determined amounts of salts and SPE sorbents to support the most common method configurations used today

Product #	Description
55227-U	Dispersive SPE (dSPE) Citrate Extraction Tube, pk of 50
55237-U	Dispersive SPE (dSPE) Citrate/Sodium Bicarbonate Extraction Tube, pk of 50
55234-U	Dispersive SPE (dSPE) MgSO ₄ Extraction Tube, pk of 50
55228-U	Dispersive SPE (dSPE) PSA SPE Clean Up Tube 1, pk of 50
55229-U	Dispersive SPE (dSPE) PSA/C18 SPE Clean Up Tube 1, pk of 50
55230-U	Dispersive SPE (dSPE) PSA/ENVI-Carb SPE Clean Up Tube 1, pk of 50
55233-U	Dispersive SPE (dSPE) PSA/ENVI-Carb SPE Clean Up Tube 2, pk of 50

Also available:

- Sample packs
- Custom tubes and packing materials

sigma-aldrich.com/spe

The Extraction and Analysis of Multi-Residue Pesticides in Orange Using Dispersive SPE and GC-MS

An Trinh, Katherine Stenerson, Robbie Wolford, Olga Shimelis, and Craig Aurand

QuEChERS Multiresidue-Method

Advantages:

- > **Quick** (~30 min/batch of 6 samples)
- > **Easy** (no laborious steps)
- > **Cheap**
- > **Effective** (wide scope, low consumption)
- > **Rugged** (minimal sources of errors)
- > **Safe**

Extract in acetonitrile
→ amenable to GC-, LC- and (D)-SPE

www.quechers.com

Courtesy of M. Stenerson, CVUA Stuttgart, Germany

Topics: Sample Prep Innovations

- SPME (solid phase microextraction)
- High specificity SPE (SupelMIP)
- Dispersive SPE
- **Silver Ion SPE for FAMEs (Discovery Ag-Ion)**
- Carbonaceous adsorbents
- Flash chromatography

Users are...

- **Food analysts**

Interested in...

- **Measuring cis/trans fats or degree of unsaturation**

Users can expect...

- **To fractionate FAME samples prior to GC analysis, simplifying analytical chromatography and improving method accuracy**

Discovery Ag-Ion SPE for FAME Fractionation


Silver ion anchored onto SCX SPE support

Ag⁺ forms a charge transfer complex with unsaturated FAME double bond

- Ag⁺ = electron acceptor; double bond = electron donor

Cis configuration offers greater steric accessibility = stronger retention

Strength of interaction increases with no. of double bonds




Overview Discovery Ag-Ion SPE Procedure

- 1) Fatty acids (FA) extracted from food sample
- 2) FA converted to FAMES using BF_3
- 3) FAMES are extracted into hexane
- 4) Hexane sample applied to Discovery Ag-Ion SPE cartridge
- 5) FAMES separated using different mixtures of hexane:acetone to extract from cartridge
 - Increasing % acetone disrupts retention of strongly retained FAMES (cis and higher number of double bonds)
- 6) Fractions analyzed by GC







Cis/Trans Fractionation of FAMES from Potato Chips with and without Ag-Ion SPE



column: SP-2560, 75 m x 0.18 mm I.D., 0.14 μ m (23348-U)
oven: 180 °C, isothermal
inj.: 220 °C
det.: FID, 220 °C
carrier gas: hydrogen, 40 cm/sec. at 180 °C
injection: 0.5 μ L, 100:1 split
liner: 4 mm I.D., split, cup design



Fr. 1: hexane:acetone (96:4)
• Saturated fatty acids
• Trans monoenes




Fr. 2: hexane:acetone (90:10)
• Cis monoenes

Discovery Ag-Ion SPE Products

Description	Qty.	Cat. No.
Discovery Ag-Ion SPE		
750 mg/6 mL SPE Tube	30	54225-U
750 mg/1 mL Rezorian™ Cartridge	10	54226-U

<http://tinyurl.com/agionspe>

105 North Hartman Road
 Bellefonte, PA 16823-0942 USA
 Telephone: 800-247-6226 • 814-353-3441
 Fax: 800-447-3044 • 814-353-3044
 e-mail: supelco@supelco.com
www.sigmaaldrich.com/supelco



Technical Report




Discovery® Ag-Ion SPE for FAME Fractionation and Cis/Trans Separation

When silver ions are loaded and immobilized on to an SCX phase as counter-ions, they have the ability to form polar complexes with unsaturated fatty acid double bonds under normal-phase conditions. Discovery Ag-Ion SPE was developed for the fractionation of cis/trans isomers, and can also resolve FAMES by degree of unsaturation in which retention strength increases with increasing number of double bonds. As a result, Discovery Ag-Ion SPE allows users to fractionate FAME samples prior to GC analysis, thereby simplifying analytical chromatography and improving method accuracy. In this report we demonstrate the utility of this technology against a variety of sample matrices including potato chips, butter cookies, popcorn oil, peanut butter, and poppy seed muffins.


Although small amounts of trans fat are produced in the GI tract of cattle and are found in dairy and beef fat, trans fats are predominantly produced commercially in large quantities through a process called partial hydrogenation used to protect foods from spoilage (1). Concerns have been raised for several decades that the consumption of trans fatty acids might have contributed to the 20th century epidemic of coronary heart disease (the raising of LDL or "bad" cholesterol). As a result, on July 9, 2003, the United States Food and Drug Administration (FDA) issued a regulation requiring manufacturers to list trans fat on the Nutrition Facts panel of foods and some dietary supplements. With this rule, consumers will have more information to make healthier food choices that could lower their consumption of trans fat as a part of a heart-healthy diet. As of January 1, 2006, food manufacturers are required to list trans fat on the nutrition label (2).

Trans fats (trans unsaturated fatty acids) are fatty acids that contain double bonds that cause carbon atoms to bond in a straight configuration. As a result, they remain in a solid state at room temperature. Most naturally occurring unsaturated fatty acids are in the cis-orientation (bent), which allow for a liquid state at room temperature (Table 1).

Table 1. Types of Fatty Acids

Structure	Common Sources	Health Effects
Saturated Fatty Acids (no double bonds)		
	Palm kernel, Palm oil, Coconut (tropical oils) Butter, Hydrogenated Oils and Shortenings	Raise LDL cholesterol and increase risk of cardiovascular disease
Mono and Polyunsaturated Fatty Acids (≥ 1 cis double bond)		
	Fluid, liquid oils such as Soybean, Canola, Olive, Sunflower, and Corn Oils	Lower LDL cholesterol, associated with reduced risk of cardiovascular disease
Trans Fatty Acids (≥ 1 trans double bond)		
	Partially Hydrogenated Oils, Shortenings, Margarine, and Oils	Raise LDL cholesterol. Like saturated fat, may also lower HDL. Associated with increased cardiovascular disease and possible type II diabetes

sigmaaldrich.com/supelco



Topics: Sample Prep Innovations

- SPME (solid phase microextraction)
- High specificity SPE (SupelMIP)
- Dispersive SPE
- Silver Ion SPE for FAMES
- **Carbonaceous adsorbents (ENVI-Carbs)**
- Flash chromatography

Users are...

- Analytical chemists, HPLC, GC doing sample prep

Interested in...



- Extraction of highly polar compounds from water samples, and many others...

Users can expect...

- High extraction efficiency

Structural Classification of Carbons

Carbon Class	C–C Distance (nm)	Layer Distance (nm)
Amorphous (hexagonal)	0.139	—
Turbostratic	0.142	0.365
Graphitic	0.142	0.335
Diamond (cubic)	0.155	—




Carbon Sorbents for Sample Prep

Packed SPE tubes:


- **Supelclean™ ENVI-Carb PLUS** – spherical carbon molecular sieve for extraction of highly polar compounds from water samples
- **Supelclean ENVI-Carb-II/PSA SPE** – multilayer SPE tubes for multiresidue pesticide analysis in foods
- **Supelclean ENVI-Carb-II SPE** – isolation/removal of pigments (e.g., chlorophyll and carotenoids) and sterols commonly present in fruits, vegetables, and other natural products
- **Supelclean ENVI-Carb-II/SAX/PSA SPE** – additional ion exchange capability
- **Supelclean PSA SPE** – polymerically bonded, ethylenediamine-N-propyl phase that contains both primary and secondary amines

<http://tinyurl.com/carbonspe>



Dual Layer Supelclean ENVI-Carb-II/PSA SPE Tube

Supelclean ENVI-Carb PLUS




Spherical Carbon Molecular Sieve

- Extraction of highly polar compounds from water samples
- > 70% Abs Recovery from 0.5 L drinking water (1-100 ng/mL)


Procedure:

1. Condition w/ 10 mL MeOH & 10 mL DI water
2. Load up to 1 L sample
3. Reverse tube & elute w/ 4-5 mL MeOH in opposite direction


Examples of polar compounds:




Acephate
log P -0.85



Acrylamide
log P -0.67




1,4-Dioxane
log P -0.27



Oxamyl
log P -1.2

GC-MS Analysis of 1,4-dioxane in water extracted using ENVI-Carb Plus



- THF-d8 (Int. Std.)
- 1,4-dioxane-d8 (surrogate)
- 1,4-dioxane

Column: SPB-624, 30 m x 0.25 mm I.D. , 1.4 μ m
Oven: 30 °C (1 min.), 7 °C/min. to 90 °C, 20 °C/min. to 200 °C (3 min.)
Inj: 200 °C
Carrier: helium, 1 mL/min constant flow
Injection: 2 μ L, splitless
MS interface: 220 °C
Scan range: SIM

Topics: Sample Prep Innovations

- SPME (solid phase microextraction)
- High specificity SPE (SupelMIP)
- Dispersive SPE
- Silver Ion SPE for FAMES
- Carbonaceous adsorbents (ENVI-Carbs)
- **Flash Chromatography**

Users are...

- **Synthetic, organic chemists**
- **Medicinal chemists**

Interested in...

- **Purification of relatively large samples from reaction mixtures or other samples**


Users can expect...

- **Fast, simple, inexpensive purifications**
- **High N (spherical particles)**

Performance of Spherical vs. Irregular Silicas in Flash Application

Higher efficiency of spherical particles translates to narrower bands for more concentrated fractions and faster isolations.

samples: 5-hydroxy-DL-tryptophan and DL-tryptophan
 cartridges: 53 mm x 23 mm I.D.
 mobile phase: methanol:water (90:10)
 detection: UV 254 nm
 flow rate: 20 mL/min.



Particle Type	Total Volume	Fraction 1 Volume	Fraction 2 Volume
Spherical silica	120 mL	25 mL	40 mL
Irregular silica	400 mL	80 mL	260 mL

VersaFlash Support Literature & Products

Silica & C18 Cartridges

- Particle size options
- Cartridge size options
- Cartridges can be coupled
- Reversible
- Compatible with other systems

All system components



Summary

Solid Phase Microextraction (SPME) <http://www.sigma-aldrich.com/spme>

High specificity SPE (SupelMIPs) <http://www.sigma-aldrich.com/supelmip>

Dispersive SPE for pesticide extraction <http://www.sigma-aldrich.com/spe>

Silver Ion SPE for FAMES (Discovery Ag-Ion SPE)

<http://tinyurl.com/agionspe>

Carbon adsorbents for polar compounds (ENVI-Carbs)

<http://tinyurl.com/carbonspe>

Flash chromatography <http://www.sigma-aldrich.com/versaflash>

Acknowledgements/Collaborators

Prof. Janusz Pawliszyn, U. Waterloo, Canada

Ines de Lannoy, NoAb Biodiscovery (*in vivo* applications)

Joseph Kennedy, Prosolia (DESI)

Supelco and Fluka R&D Teams

For more information on the subjects presented here, please contact techservice@sial.com or your regional sales team.