

## GC/MS Mass Spectra Library

# Polymer Additives Library Ver.2

The analysis of additives contained in polymer materials is essential for managing and improving the quality of these materials, and for complying with regulations dealing with chemical substances. The Polymer Additives Library Ver.2 is a GC/MS mass spectral library containing information on a wide range of additives utilized in polymer materials.

In addition to mass spectra and retention indices, it contains information on the classification of additives to enable detailed analysis of additives.



GCMS-QP2050

### Contains Approximately 5,800 Mass Spectra

This library consists of a library\*, which contains 5,775 mass spectra for 590 types of polymer additives and GC/MS pyrolysates from the pyrolysis GC/MS analysis of additives, and a library containing the mass spectra for 65 compounds often targeted for analysis, selected based on their usage in the commercial market, and chemical substance regulatory information. In total, this product contains 5,840 mass spectra, which provides strong support for the analysis of polymer additives. In Ver.2, the compound registration has been newly enhanced, for example vulcanizing accelerators utilized in the rubber field. Approximately 1,000 spectra have been added.

\*ADD-MS22B F-Search Additives Library from Frontier Laboratories

### Ready-to-Use Methods

This library provides GC/MS and PY methods; both contain all necessary conditions and measurement parameters. Anyone can easily start analysis with optimal analysis conditions.

### Library Specifications

Registered compounds\*: 5,775+65

Registered information: Mass spectrum, retention index for each analytical condition, compound name, molecular weight, compositional formula, structural formula, classification of additives

Applicable models: GCMS-QP™ series, GCMS-TQ™ series (This library does not include MS/MS spectra.)

Pyrolyzer (Pyrolysis analysis system): PY-2020D, PY-2020iD, and EGA/PY-3030D

Workstations: GCMSsolution™ Ver. 4.6 and later, LabSolutions™ GCMS™ Ver. 5.131 and later, LabSolutions DB/CS GCMS Ver. 6.131 and later

### Includes Retention Indices and Classification Information of Additives

Retention indices are registered for all compounds. Compounds can be identified with a high degree of accuracy by reducing the library search results using retention indices. In addition, information on the classification of additives is included, so even without detailed knowledge of additives, users can see which additive is associated with a compound found in the library.

### Can Be Used with a Variety of GC/MS Systems

The library can be used for a variety of GC/MS applications, including pyrolysis GC/MS, which is widely used for the analysis of additives in polymer materials, and liquid sample injection GC/MS. It supports a wide range of additive analysis for customers.

\*: This consists of two types of libraries. The library with 5,775 compounds consists of the ADD-MS22B F-Search Additives Library developed by Frontier Laboratories, which has been converted to the Shimadzu mass spectral library format, and prepared for use with GCMSolution and LabSolutions GCMS. The library with 65 compounds registered is an additive library developed exclusively by Shimadzu. It contains additives often targeted for analysis, selected based on their usage in the commercial market, and chemical substance regulatory information. Using these two libraries in combination enables the analysis of an even wider range of additives.

# Polymer Additives Library Ver.2

## Filtering with Retention Index

Multiple compounds with similar mass spectra are listed as candidates when performing a library search using only the mass spectrum. Filtering with the retention index sorts the candidates by retention index, thereby providing highly accurate identification results.

Results of similarity search using mass spectrum

Hit#	Similarity	Register	Ret. Index	Compound Name	Mol Wt	Formula	Library
1	95	<input checked="" type="checkbox"/>	1729	Ethylhexyl benzoate [Original Additive: Bi	234	C <sub>15</sub> H <sub>22</sub> O <sub>2</sub>	FLAB_ADD1
2	81	<input type="checkbox"/>	2042	Undecyl benzoate [Original Additive: Dlu	276	C <sub>18</sub> H <sub>28</sub> O <sub>2</sub>	FLAB_ADD1
3	81	<input type="checkbox"/>	1796	Octyl benzoate [Original Additive: Tri-n-	234	C <sub>15</sub> H <sub>22</sub> O <sub>2</sub>	FLAB_ADD1
4	79	<input type="checkbox"/>	2108	Tridecyl benzoate [Original Additive: Diun	304	C <sub>20</sub> H <sub>32</sub> O <sub>2</sub>	FLAB_ADD1
5	66	<input type="checkbox"/>	793	1-Octene [Original Additive: Alkyl phosp	112	C <sub>8</sub> H <sub>16</sub>	FLAB_ADD1
6	65	<input type="checkbox"/>	791	Isomer of C <sub>8</sub> H <sub>16</sub> [Original Additive: Alkyl	112	C <sub>8</sub> H <sub>16</sub>	FLAB_ADD1




Results sorted using retention index filtering

Hit#	Similarity	Register	Ret. Index	Compound Name	Mol Wt	Formula	Library
1	95	<input checked="" type="checkbox"/>	1729	Ethylhexyl benzoate [Original Additive: Bi	234	C <sub>15</sub> H <sub>22</sub> O <sub>2</sub>	FLAB_ADD1


## Confirming the Additive Classification Information

The additive classification information (such as plasticizers and flame retardants) registered in the library plays a role in confirming the type of additive associated with the compounds included in search results.



## Deconvolution Analysis with LabSolutions Insight Explore™ GCMS

Even more accurate additive analysis is achieved by combining deconvolution analysis, which is capable of separating out impurities and detecting trace compounds, with this library.



## Example of System Configuration\*

Pyrolysis Analysis System: EGA/PY-3030D, GCMS-QP2050

Liquid Injection Analysis System: AOC-30i, GCMS-QP2050

\*: The accuracy of the retention indexes differs depending on the sample introduction method. Some of the pyrolysates may only be seen using a pyrolysis GC/MS analysis.

GCMS, GCMS-QP, GCMS-TQ, GCMSsolution, LabSolutions, and LabSolutions Insight Explore are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.