

Glass Insert Selection Guide for GC and GC-MS

GC/GC-MS Glass Insert/Liner Selection Guide



Contents

Chapter

1

Introduction

| | |
|---|---|
| 1-1. Glass Inserts/Liners | 3 |
| 1-2. Sample Injection Unit/Injection Method | 4 |
| 1-3. Selection Guide | 6 |

Chapter

2

Glass Liner Selection

| | |
|--|----|
| 2-1. Glass Liner Selection by GC Model/Sample Injection Unit | 10 |
| •GC-2060 series | 10 |
| •GC-2030/GC-2050 | 13 |
| •GC-2010 series | 15 |
| •GC-2014 series | 16 |
| •HS-10 | 17 |
| 2-2. Glass Liner Selection by Application | 18 |
| •SPME/SPME Arrow Analysis | 18 |
| •Aqueous Solvent Analysis | 19 |

Chapter

3

All Glass Liner List

| | |
|---|----|
| 3-1. Glass Liner Compatible Instruments | 20 |
| 3-2. Glass Liner Information Table | 23 |

Chapter

4

| | |
|--------------------------|-----------|
| Accessories | 26 |
|--------------------------|-----------|



1-1. Glass Inserts/Liners


A gas chromatograph (GC) introduces a sample into a column via an injection unit. Using a glass liner (glass insert) allows the sample to be vaporized efficiently and the vaporized sample to be mixed as uniformly as possible with the carrier gas before introduction into the column. It also reduces the extent to which contaminants such as high boiling point compounds from the sample contaminate the Injection Unit. Accordingly, selecting a glass liner appropriate for the objective improves the accuracy of analysis and leads to more reliable results.

The selection criteria for glass liners are as follows:

- Sample form
- Type of injection unit used
- Sample injection method during analysis
- Shape of the glass liner, presence of wool, and presence of surface treatment on the wool or glass liner

Glass liners with a variety of specifications are commercially available to suit these varied selection criteria, so selecting a glass liner can be challenging.

This selection guide describes how to select an appropriate glass liner for the objective/application as well as points to consider in accordance with elements of the glass liner configuration. Additionally, this guide introduces glass liners appropriate for your Shimadzu GC system.



Schematic Diagram of the Glass Liner and O-Ring



| Element | Parameter | Explanation | Role |
|-------------------|---|---|---|
| Wool | Position | Wool position within the glass liner | Controls the vaporization point |
| | Amount of filler | Weight of the wool | Controls the vaporization efficiency |
| Inner Diameter | Inner diameter | Inner diameter of the glass liner | Impact on the delivery rate of the vaporized components |
| | Capacity | Internal volume of the glass liner | Impact on the upper limit of the injection volume |
| Shape | Wool position offset prevention mechanism | Protrusion within the glass liner on the order of 1 mm | Prevents offset of the wool position |
| | Tapered | A narrowing of the inner diameter at the bottom of the glass liner on the order of 1 mm | Focuses the vaporized sample toward the column tip |
| Surface Treatment | | Chemical treatment of the surface of the glass liner | Impact on adsorption and degradation of compounds |

1-2. Sample Injection Units/Injection Methods

The glass liner is designed in accordance with the sample injection method. This section introduces typical sample injection units and injection methods.

Select a glass liner suited to the sample injection unit and injection method being used.

| Sample Injection Unit | Injection Method | Glass Liner | Target Sample | Sample Injection | Compatible Columns |
|-----------------------|--------------------------------|---|--|---|--|
| SPL | Split | Split glass liner, Split/Splitless glass liner | Medium to high-concentration samples | The sample is instantaneously vaporized, and some of it is injected into the column. | 0.1 mm to 0.53 mm capillary columns |
| | Splitless | Splitless glass liner, Split/Splitless glass liner | Low-concentration samples | The sample is vaporized, and a large portion of it is injected into the column. | 0.1 mm to 0.53 mm capillary columns (Wide diameter columns are more suitable) |
| MMI | Split | Split glass liner(78.5 mm length), Split/Splitless glass liner (78.5 mm length) | Medium to high-concentration samples | The sample is instantaneously vaporized, and some of it is injected into the column. | 0.1 mm to 0.53 mm capillary columns |
| | Splitless | Splitless glass liner(78.5 mm length) ,Split/Splitless glass liner (78.5 mm length) | Low-concentration samples | The sample is vaporized, and a large portion of it is injected into the column. | 0.1 mm to 0.53 mm capillary columns(Wide diameter columns are more suitable) |
| | Cold split | Split glass liner (78.5 mm length), Split/Splitless glass liner (78.5 mm length) | Medium to high concentration samples prone to thermal decomposition, and samples with a wide range of boiling points | A heating program is used to vaporize the compounds in the sample in order of boiling point. A portion is then injected into the column. | 0.1 mm to 0.53 mm capillary columns |
| | Cold splitless | Splitless glass liner (78.5 mm length) , Split/Splitless glass liner (78.5 mm length) | Low-concentration samples prone to thermal decomposition, and samples with a wide range of boiling points | A heating program is used to vaporize the compounds in the sample in order of boiling point. A large portion is then injected into the column. | 0.1 mm to 0.53 mm capillary columns (Wide diameter columns are more suitable) |
| | Cold Direct | Direct glass liner (78.5 mm length) | Samples contained wide range of boiling points | Almost all sample introduced into the column before vaporization via specialized shape liner. | 0.32 mm to 0.53 mm capillary columns (Wide diameter columns are more suitable) |
| | Large Volume Injection | Large-volume injection liner (78.5 mm length) | Extreme low-concentration samples | First, the injection unit will selectively vent the solvent from the sample. Then, a heating program will vaporize the solute and introduce it into the column. | 0.1 mm to 0.53 mm capillary columns |
| | Thermal desorption | ¼inch O.D TD tube, Splitless glass liner (78.5 mm length, monotrapp analysis only) | Trace amounts of VOCs/ SVOCs contained in an adsorbent. | A heating program is used to vaporize the compounds in adsorbent. | 0.1 mm to 0.53 mm capillary columns |
| WBI | Direct | Splitless glass liner/WBI glass liner | No concentration restriction | The sample is instantaneously vaporized, and the full amount is injected into the column. | 0.53 mm capillary columns |
| OCI | On-Column | OCI glass liner | Samples with a wide range of boiling points, and samples with ultra-high boiling points | The entire sample solution is injected into the column before vaporization. | 0.25 mm to 0.53 mm capillary columns |
| PTV | Cold split | PTV glass liner | Medium to high-concentration samples prone to thermal decomposition, and samples with a wide range of boiling points | A heating program is used to vaporize the compounds in the sample in order of boiling point. A portion is then injected into the column. | 0.1 mm to 0.53 mm capillary columns |
| | Cold splitless | PTV glass liner | Low-concentration samples prone to thermal decomposition, and samples with a wide range of boiling points | A heating program is used to vaporize the compounds in the sample in order of boiling point. A large portion is then injected into the column. | 0.1 mm to 0.53 mm capillary columns (Wide diameter columns are more suitable) |
| SINJ/DINJ | Full volume injection (Direct) | Packed glass liner | No concentration restriction | The sample is instantaneously vaporized, and the full amount is injected into the column. | Packed Columns |

Refer to the following for details on the types of sample injection unit.

[Details](#) 

1-3. Selection Guide


This section describes glass liner selection methods.

Selecting the Glass Liner Wool

Wool is used inside the glass liner when performing SPL, WBI, and PTV analysis. The wool traps the liquid sample injected into the glass liner, promoting efficient vaporization. The two main wool substances are hard glass and quartz glass. This section describes the role of the wool, and cautions when filling the glass liner with wool.

Wool Filling Position

Liquid samples are often vaporized on the wool surface or inside the wool plug. In other words, the wool position is the site at which the sample vaporizes, so controlling the position of the wool controls the timing of the sample vaporization.



For splitless or direct analysis with low carrier gas flow rates, it is common to fill the wool at the bottom of the glass liner.

If the vaporization site is too far away from the column inlet, it will take time for the vaporized sample to reach the column inlet, which can cause reduced sensitivity and peak shape problems. Additionally, there is a risk that the vaporized sample will spill over the top of the glass liner, leading to reversed flow or deviation from the septum purge flowrate. Caution is necessary when injecting 3 μL or more of the sample.

For split analysis with high total carrier gas flow rates, it is common to fill the wool at the top of the glass liner.

If the top is filled with wool, the syringe needle tip and the top edge of the wool tend to be in close proximity or contact. As a result, when the liquid sample is injected, the direction of scatter of the liquid from the syringe needle tip is uniform, leading to improved repeatability of area values.

Amount of Wool

The amount of wool filling differs depending on the injection unit and the injection method, but is generally between 2 mg and 20 mg.

The liquid sample is vaporized by the heat provided not only by the carrier gas but also by the filled wool. In particular, if the inside of the glass liner is not filled with wool, the liquid sample will not vaporize sufficiently, and repeatability will drop. The amount of heat provided by the wool changes depending on the amount of wool filler, so adjust the amount of wool filler if you are concerned about repeatability during liquid sample analysis.

In split analysis, the amount of wool must be increased to improve the sample vaporization efficiency.

The greater the amount of wool, the more efficient the sample vaporization, but this can also cause adsorption and degradation depending on the compound.

In splitless analysis, in which the sample retention time within the glass liner is long, if the amount of wool filler is large, excessive heat might be applied to the sample, which can cause sample adsorption or degradation. When increasing the amount of wool filler, be careful of adsorption and degradation of compounds in the sample. **With splitless analysis and direct analysis in pursuit of sensitivity, when analyzing target compounds that are comparatively prone to adsorption and degradation, it is recommended that you reduce the amount of wool filling for the analysis.**

Recommended Filling Position and Amount of Filler for Each Model

If purchasing a glass liner with no wool, and then filling the wool yourself, pay attention while filling to the amount and position of the wool. For details on the recommended wool amount and position, refer to the instruction manual and maintenance help for the applicable instrument.

| GC Model | Sample Injection Unit | Injection Method | Recommended Amount of Wool | Recommended Wool Position |
|---------------------------------------|-----------------------|---|----------------------------|---|
| GC-2060 | SPL | Split | 10 mg | 22 mm from the top surface |
| | | Splitless | 4 mg | 67 mm from the top surface (Filled to the tapered bottom) |
| | MMI | All Xtra Inert liners are pre-packed with wool. | | |
| GC-2030 GC-2050 | SPL | Split | 10 mg | 22 mm from the top surface |
| | | Splitless/Full volume injection | 4 mg | 67 mm from the top surface (Filled to the tapered bottom) |
| | PTV | Cold split/Cold splitless | 1 to 2 mg | 22 mm from the top surface |
| GC-2014 | SPL | Split | 10 mg | 20 mm from the top surface |
| | | Splitless/Full volume injection | 2 mg | 67 mm from the top surface (Filled to the tapered bottom) |
| GC-2025 (SPL) GC-2010 series (SPL) | SPL | Split | 10 mg | 25 mm from the top surface |
| | | Splitless/Full volume injection | 2 mg | 67 mm from the top surface (Filled to the tapered bottom) |
| GC-2010 series | PTV | Cold split/ Cold splitless | 1 to 2 mg | 25 mm from the top surface |

Note: The optimal wool position and amount of wool may differ depending on the sample, the measured components, and the injection volume. The above-mentioned are generally recommended values.

In split analysis, it is advisable to pack the upper surface of the wool as flat as possible, as this improves reproducibility. When injecting samples, having the upper surface of the wool as flat as possible is beneficial because it ensures stable sample vaporization, especially when the wool is in close proximity to the syringe tip.

Other Roles of the Wool


Trapping Nonvolatile Components

During analysis, as the number of injections increases, so too does the possibility that various nonvolatile components will be injected and accumulate within the sample Injection Unit.

- Nonvolatile components from the sample
- Contaminants and specks of silicon from the septum
- Metal powder from the syringe

When a large amount of nonvolatile components accumulate, this may have an impact on the analysis results.

Trapping as many of these nonvolatile components as possible with the wool, and reducing column contamination from nonvolatile components, leads to more reliable analysis results. Additionally, trapping nonvolatile components with the wool helps reduce the frequency of replacing columns and the frequency of maintenance for the GC system as a whole.



Mixing

Filling the glass liner with wool promotes mixing of the sample after vaporization with the carrier gas, reducing discrimination* during sample vaporization.


* Note: Discrimination refers to a compositional change that occurs when using a chromatograph to analyze mixed samples containing components with a range of boiling points from low to high.

Selecting the Shape of the Glass Liner

This section describes the role of the various shapes of glass liners, as well as points to consider when selecting the shape.

Wool Position Offset Prevention Mechanism


When there is a sudden pressure change in sample injection unit, such as during septum replacement, the wool position is sometimes offset, worsening the repeatability of the analysis results. Some glass liners are shaped with a protrusion to prevent wool position offset. Using a glass liner with a wool position offset prevention mechanism should provide analysis results with more stable repeatability.



Tapered Structure

A tapered part is one in which the inner diameter of the glass liner quickly narrows from the top down, from a position at the bottom of the glass liner (in proximity to the column inlet).

A tapered structure plays a role in focusing the injection of the vaporized sample toward the column. This structure is often used in glass liners for splitless and direct analysis. In these analyses, in which the total flowrate of the carrier gas is small, the delivery rate of the vaporized sample within the glass liner is slowed. For this reason, using a tapered structure helps the vaporized sample to collect at the column inlet more efficiently.



Selecting the Inner Diameter of the Glass Liner

This section describes how to select the inner diameter of the glass liner in accordance with the analysis conditions.

Inner Diameter


The inner diameter of the glass liner has an impact on the sample delivery rate and the capacity of the glass liner.

This section focuses on describing the impact on the sample delivery rate.

(Capacity is described in the next section.)

In almost all liquid sample analysis, a thicker inner diameter of glass liner is appropriate. When the carrier gas passes through the glass liner, if the glass liner has a narrow inner diameter, the linear velocity will increase. Liquid samples take some time to vaporize, and the volume expands considerably, so the use of glass liners with extremely narrow inner diameters should be avoided.

However, when analyzing some samples containing low boiling point components, the peak bandwidth tends to spread, so a fast delivery rate is required. For this sort of analysis, a glass liner with a comparatively narrow inner diameter is suitable.



Application



Capacity


The internal capacity of a glass liner is calculated as $(\text{Inner diameter of the glass liner}/2)^2 \times \pi \times \text{Length of the glass liner}$. The capacity of a GC glass liner is generally about 100 to 900 μL . Select a suitable glass liner in accordance with the vaporization volume of the sample for analysis. The following table shows the vaporization volume (at 250 °C and 140 kPa) of typical sample solvents utilized in analysis.

| Solvent Type | 1 μL Injection Volume | 2 μL Injection Volume |
|------------------|----------------------------------|----------------------------------|
| Isooctane | 110 | 220 |
| n-Hexane | 140 | 280 |
| Toluene | 170 | 340 |
| Ethyl Acetate | 185 | 370 |
| Acetone | 245 | 490 |
| Dichloromethane | 285 | 570 |
| Carbon Disulfide | 300 | 600 |
| Acetonitrile | 350 | 700 |
| Methanol | 450 | 900 |
| Water | 1010 | 2020 |

Surface Treatment of Glass Liners

Samples make direct contact with the inner walls of glass liners, so if a sample contains compounds prone to degradation, problems such as peak tailing and drops in peak intensity sometimes occur. In this case, the use of a glass liner and wool treated to deactivate the surface is recommended. Note that deactivated glass liners are rated for use up to approximately 350 °C. The non-deactivated glass liners are rated for use up to approximately 450 °C.

Example



Glass Liner Selection



2-1. Glass Liners by GC Model/Sample Injection Unit

This section introduces frequently used, general-purpose glass liners for each GC model, sample injection unit, and injection method. If you are having difficulty making a selection, it is recommended that you start by selecting from the following lineup. The full glass liner list is provided in Chapter 4, so use that in combination with this section.

GC-2060



GC-2060

SPL (Split)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature |
|----|---|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|------------------------|---------------------|
| 1 |  | 227-35007-01 | 5 | ✓ | 95 | 4.9 | 3.4 | 863 | Deactivation treatment | 350°C |

Standard accessory

SPL (Splitless)



| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature |
|----|---|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|------------------------|---------------------|
| 19 |  | 227-35008-01 | 5 | ✓ | 95 | 5 | 3.4 | 654 | Deactivation treatment | 350°C |

Standard accessory



SPL (Split/Splitless)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature |
|----|---|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|------------------------|---------------------|
| 8 |  | 227-35016-01 | 5 | ✓ | 95 | 5 | 3.5 | 847 | Deactivation treatment | 350°C |

MMI (Split)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature* |
|----|---|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|-----------------------------------|----------------------|
| 44 |  | 227-36700-01 | 5 | ✓ | 78.5 | 6.3 | 4.0 | 986 | Xtra Inert Deactivation treatment | 350°C |
| 50 |  | 227-36707-01 | 5 | --- | 78.5 | 6.3 | 4.0 | 672 | Xtra Inert Deactivation treatment | 350°C |

MMI (Splitless)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature* |
|----|---|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|-----------------------------------|----------------------|
| 45 |  | 227-36701-01 | 5 | ✓ | 78.5 | 6.5 | 4.0 | 924 | Xtra Inert Deactivation treatment | 350°C |
| 46 |  | 227-36702-01 | 5 | --- | 78.5 | 6.5 | 4.0 | 496 | Xtra Inert Deactivation treatment | 350°C |


MMI (Large-volume injection, Split/Splitless)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature* |
|----|---|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|-----------------------------------|----------------------|
| 48 |  | 227-36704-01 | 5 | ✓ | 78.5 | 6.3 | 4.0 | 924 | Xtra Inert Deactivation treatment | 350°C |

MMI (Split/Splitless)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature* |
|----|---|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|-----------------------------------|----------------------|
| 54 |  | 227-36711-01 | 5 | ✓ | 78.5 | 6.3 | 4.0 | 924 | Xtra Inert Deactivation treatment | 350°C |

MMI (Direct)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature* |
|----|---|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|-----------------------------------|----------------------|
| 47 |  | 227-36703-01 | 5 | --- | 78.5 | 6.3 | 4.0 | 697 | Xtra Inert Deactivation treatment | 350°C |

MMI (Split, Low Inner Volume)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature* |
|----|---|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|-----------------------------------|----------------------|
| 53 |  | 227-36710-01 | 5 | ✓ | 78.5 | 6.3 | 2.0 | 247 | Xtra Inert Deactivation treatment | 350°C |

MMI (Splitless, Middle Inner Volume)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature* |
|----|---|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|-----------------------------------|----------------------|
| 52 |  | 227-36709-01 | 5 | ✓ | 78.5 | 6.5 | 3.0 | 520 | Xtra Inert Deactivation treatment | 350°C |

MMI (Gas analysis)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature* |
|----|---|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|-----------------------------------|----------------------|
| 51 |  | 227-36708-01 | 5 | --- | 78.5 | 6.3 | 1.0 | 62 | Xtra Inert Deactivation treatment | 350°C |

Note: The liners for MMI inlet in the Xtra Inert series are rated for use up to 350 °C to maintain the performance of the deactivation treatment. For applications where deactivation is not critical, they can be used at temperatures up to 450 °C.

WBI (Full Volume Injection)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature |
|----|--|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|------------------------|---------------------|
| 19 |  | 227-35008-01 | 5 | ✓ | 95 | 5 | 3.4 | 654 | Deactivation treatment | 350°C |

Standard accessory

SINJ+WBC Attachment (WBC Full Volume Injection)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature |
|----|---|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 34 |  | 221-38107-02 | 1 | --- | 126 | 4.8 | 3.4 | 870 | None | 450°C |

OCI-2030 NX Liner (Metal)


| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature |
|----|---|-----------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 28 |  | 221-85694 | 1 | --- | 11.5 | - | - | - | None | 450°C |

Note: This can only be used with the OCI-2030 NX. This cannot be used with OCI-2030.

SINJ (Liner for UG-type columns)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature |
|----|---|-----------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 43 |  | 222-30935 | 1 | --- | 71 | 5 | 3.4 | 549 | None | |

SINJ (Liner for UG-type columns)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature |
|----|---|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 32 |  | 221-80902-84 | 1 | --- | 93 | 5 | 3.4 | 726 | None | |

Nexis GC-2030

GCMS-QP2020 NX
GCMS-TQ8040 NX
GCMS-TQ8050 NX

Brevis GC-2050



GC-2030

GC-2050

SPL (Split)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature | |
|----|--------|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|------------------------|---------------------|--------------------|
| 1 | | 227-35007-01 | 5 | ✓ | 95 | 4.9 | 3.4 | 863 | Deactivation treatment | 350°C | Standard accessory |

SPL (Splitless)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature | |
|----|--------|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|------------------------|---------------------|--------------------|
| 19 | | 227-35008-01 | 5 | ✓ | 95 | 5 | 3.4 | 654 | Deactivation treatment | 350°C | Standard accessory |

SPL (Split/Splitless)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature |
|----|--------|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|------------------------|---------------------|
| 8 | | 227-35016-01 | 5 | ✓ | 95 | 5 | 3.5 | 847 | Deactivation treatment | 350°C |

WBI (Full Volume Injection)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature | |
|----|--------|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|------------------------|---------------------|--------------------|
| 19 | | 227-35008-01 | 5 | ✓ | 95 | 5 | 3.4 | 654 | Deactivation treatment | 350°C | Standard accessory |

SINJ+WBC Attachment (WBC Full Volume Injection)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature |
|----|--------|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 34 | | 221-38107-02 | 1 | --- | 126 | 4.8 | 3.4 | 870 | None | 450°C |

OCI-2030 Liner (Metal)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature |
|----|--------|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 25 | | 221-49298-91 | 1 | --- | 103 | 2 | 1 | 81 | None | 450°C |

Note: Cannot be used with the OCI-2030 NX.

Simple OCI Liner

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature |
|----|---|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|---|---------------------|
| 27 |  | 221-49381-02 | 1 | --- | 95 | 3.5 | 0.8 | 48 | Deactivation treatment (Silanizing treatment) | 350°C |

Note: Cannot be used with the OCI-2030 NX.

OCI-2030 NX Liner (Metal)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature |
|----|---|-----------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 28 |  | 221-85694 | 1 | --- | 11.5 | - | - | - | None | 450°C |

Note: This can only be used with the OCI-2030 NX. This cannot be used with OCI-2030.

PTV (Cold Split/Cold Splitless)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature |
|----|---|-----------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 30 |  | 221-49300 | 1 | --- | 95 | 3.5 | 1.5 | 168 | None | 450°C |

Standard accessory

SINJ (Liner for UG-type columns)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature |
|----|---|-----------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 43 |  | 222-30935 | 1 | --- | 71 | 5 | 3.4 | 549 | None | |

SINJ (Liner for UG-type columns)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature |
|----|---|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 32 |  | 221-80902-84 | 1 | --- | 93 | 5 | 3.4 | 726 | None | |

GC-2010 series

GCMS-QP2010 SE
GCMS-QP2020
GCMS-TQ8050/8040/8030



GC-2010

SPL (Split)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature |
|----|--------|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|------------------------|---------------------|
| 1 | | 227-35007-01 | 5 | ✓ | 95 | 4.9 | 3.4 | 863 | Deactivation treatment | 400°C |
| 2 | | 221-41444-01 | 1 | --- | 95 | 4.9 | 3.4 | 863 | None | 450°C |
| 21 | | 225-20803-01 | 5 | ✓ | 95 | 5 | 3.5 | 914 | Deactivation treatment | 350°C |

Standard accessory

- 227-35007-01 can be used with the GC-2010 series. The wool filling position differs from the recommended position for the 2010 series. If it is used in the default position, the form is such that as with the GC-2030, the tip of the syringe needle will enter the wool.
- 221-41444-01 is a standard accessory for the GC-2010 unit.
- 225-20803-01 has a history of being recommended for GCMS with the GC-2010 series.

SPL (Splitless)/WBI (Full Volume Injection)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature |
|----|--------|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|------------------------|---------------------|
| 19 | | 227-35008-01 | 5 | ✓ | 95 | 5 | 3.4 | 654 | Deactivation treatment | 400°C |
| 20 | | 221-48335-01 | 1 | --- | 95 | 5 | 3.4 | 654 | None | 450°C |
| 6 | | 221-48876-03 | 5 | ✓ | 95 | 5 | 3.4 | 654 | Deactivation treatment | 350°C |

Standard accessory

- 227-35008-01 can be used with the GC-2010 series.
- 221-48335-01 is a standard accessory for the GC-2010 unit.
- 221-48876-03 has a history of being recommended for GCMS with the GC-2010 series.

SPL (Split/Splitless)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature |
|----|--------|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|------------------------|---------------------|
| 8 | | 227-35016-01 | 5 | ✓ | 95 | 5 | 3.5 | 847 | Deactivation treatment | 350°C |

OCI Liner (Metal)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature |
|----|--------|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 25 | | 221-49298-91 | 1 | --- | 103 | 2 | 1 | 81 | None | 450°C |

Simple OCI Liner

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature |
|----|--------|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|---|---------------------|
| 27 | | 221-49381-02 | 1 | --- | 95 | 3.5 | 0.8 | 48 | Deactivation treatment (Silanizing treatment) | 350°C |

PTV (Cold Split/Cold Splitless)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature |
|----|--------|-----------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 30 | | 221-49300 | 1 | --- | 95 | 3.5 | 1.5 | 168 | None | 450°C |

Standard accessory

GC-2014 series



GC-2014/GC-2014s/GC-2014c

SPL (Split)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature |
|----|--------|---------------------------|----------|--------------|-------------|---------------------|---------------------|---------------|------------------------|---------------------|
| 1 | | 227-35007-01 | 5 | 8.5 | 95 | 4.9 | 3.4 | 863 | Deactivation treatment | 350°C |
| 3 | | 221-41444 221-41444-84 | 1 5 | --- | 95 | 5 | 3.4 | 863 | None | 450°C |

Standard accessory

- 227-35007-01 can be used with the GC-2014 series. The wool filling position differs from the recommended position for the 2014 series.
- 221-41444 is a standard accessory for the GC-2014 unit. It is not wool filled.

SPL (Splitless)/WBI (Full Volume Injection)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature |
|----|--------|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|------------------------|---------------------|
| 19 | | 227-35008-01 | 5 | ✓ | 95 | 5 | 3.4 | 654 | Deactivation treatment | 400°C |
| 20 | | 221-48335-01 | 1 | --- | 95 | 5 | 3.4 | 654 | None | 450°C |

Standard accessory

- 227-35008-01 can be used with the GC-2014 series.
- 221-48335-01 is a standard accessory for the GC-2014 unit. It is not wool filled.

SINJ/DINJ+WBC Attachment (WBC Full Volume Injection)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature |
|----|--------|-----------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 33 | | 221-38107 | 1 | --- | 139 | 4.8 | 3.4 | 988 | None | 450°C |

SINJ/DINJ (Dia. 3.2 mm Packed Column Full Volume Injection Method)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature |
|----|--------|---------------------------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 37 | | 221-14093 221-14093-84 | 1 5 | --- | 139 | 4.4 (Tip 2.9) | 3.5 | 1104 | None | 450°C |

Standard accessory

- 221-14093 and 14093-84 are used for glass-packed or SUS-packed columns with a diameter of 3.0 to 3.4mm.

SINJ/DINJ (Dia. 2.6 mm Packed Column Full Volume Injection Method)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature |
|----|--------|---------------------------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 38 | | 221-14094 221-14094-84 | 1 5 | --- | 139 | 4.4 (Tip 2.3) | 3.4 | 825 | None | 450°C |

Standard accessory

- 221-14094 and 14094-84 are used for glass-packed columns with a diameter of 2.4 to 2.8mm.

HS-10

With the HS-10 headspace sampler, the transfer line is connected to the GC SPL injection port. The sample is injected from the HS to the GC.
The HS-10 glass liner is used for the SPL.



HS-10

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature |
|----|---|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|------------------------|---------------------|
| 40 |  | 227-35014-01 | 1 | --- | 95 | 5 | 1.2 | 107 | Deactivation treatment | 350°C |

Shimadzu HS-20 series Headspace sampler has a short transfer line design that reduces carryover and enables highly sensitive analysis.

A glass liner is not required for HS-20 series because a column can be connected to them without going through a injection units such as SPL.

2-2. Glass Liner Selection by Application

SPME/SPME Arrow Analysis

The solid phase micro extraction (SPME) method concentrates volatile components in SPME fibers, for sample loading into GC and GCMS systems. Analysis can be performed by the AOC-6000 Plus multi-functional autosampler system.

Refer also to C146-E424 "Smart SPME Fibers and Arrow Selection Guide".

[Brochure](#)



AOC-6000 Plus

Favorably shaped peaks can be obtained by selecting the optimal GC glass liner for the Smart SPME fibers or Smart SPME Arrow. If an ordinary glass liner for liquid injection is used when using Smart SPME fibers, the bandwidth for low boiling point compounds widens, leading to wider peak shapes. For this reason, with the SPME method, use glass liners with as narrow an inner diameter as possible to obtain sharp peaks.

SPME (SPL)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature |
|----|--------|-----------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 18 | | 221-75196 | 5 | --- | 95 | 5 | 0.8 | 48 | None | 350°C |

When using Smart SPME fibers, select a glass liner with a narrow 0.8 mm inner diameter for SPME.

SPME Arrow (SPL)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature |
|----|--------|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|-------------------|---------------------|
| 41 | | 227-35327-03 | 3 | --- | 95 | 5 | 1.3 | 126 | None | 350°C |
| 42 | | 227-35328-03 | 3 | --- | 95 | 5 | 1.7 | 216 | None | 350°C |

· When using Smart SPME Arrow with an outer diameter of 1.1 mm, use glass liners for SPME Arrow with an inner diameter of 1.3 or 1.7 mm.

· When using Smart SPME Arrow with an outer diameter of 1.5 mm, use glass liners for SPME Arrow with an inner diameter of 1.7 mm.

MMI (SPME)

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (µL) | Surface Treatment | Maximum Temperature |
|----|--------|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|-----------------------------------|---------------------|
| 49 | | 227-36705-01 | 5 | --- | 78.5 | 6.3 | 0.75 | 35 | Xtra Inert Deactivation treatment | 350°C |

Aqueous Solvent Analysis (SPL)




GC-2030/GC-2050/GC-2010 series/GC-2014 series/GC-2025

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature |
|----|--------|--------------|----------|--------------|-------------|---------------------|---------------------|---------------|------------------------|---------------------|
| 39 | | 227-35015-01 | 5 | ✓ | 95 | 5 | 3.4 | 863 | Deactivation treatment | 400°C |

227-35015-01 is a glass liner for aqueous solvent analysis. The shape, wool amount, and wool filling position are optimized, enabling stable sample vaporization and good repeatability.

Application News (01-00328) introduces examples of the use of glass liners for aqueous solvent analysis.

Application



Measurement Results for a Standard Solution for Alcohol Concentration Measurements (TCD)

All Glass Liner List



This is a list of the glass liners provided by Shimadzu. Information on compatible GC instruments is also noted. The same ID numbers used in the section list in Chapter 2 are used in this table.

3-1. List of Glass Liners and Compatible Instruments

★...Standard accessory ●...Described as a recommended item in Chapter 2 ✓...Can be used

*1 Special order for glass column + pTCD

*2 Special order for glass column + packed detector other than a TCD

*3 It can be used by remodeling the GC unit SPL to WBI with WBI Modification Kit 221-74660-41.

*4 Connections to systems other than the GC-2010, GC-2014, or GC-2030 series are handled by special order.

| ID | Figure | P/N | Quantity | Wool Filling | Category | GC-2026 | GC-2030 | GC-2050 | GC-2010 series | GC-2014 series |
|----|--------|---------------------------|----------|--------------|-------------------------|---------|---------|---------|----------------|----------------|
| 1 | | 227-35007-01 | 5 | ✓ | SPL (Split) | ★● | ★● | ★● | ● | ● |
| 2 | | 221-41444-01 | 1 | | SPL (Split) | ✓ | ✓ | ✓ | ★● | ✓ |
| 3 | | 221-41444 221-41444-84 | 1 5 | | SPL (Split) | ✓ | ✓ | ✓ | ✓ | ★ |
| 4 | | 221-75193 | 5 | ✓ | SPL (Split) | ✓ | ✓ | ✓ | ✓ | ✓ |
| 5 | | 221-48876-02 | 5 | | SPL (Splitless) | ✓ | ✓ | ✓ | ✓ | ✓ |
| 6 | | 221-48876-03 | 5 | ✓ | SPL (Splitless) | ✓ | ✓ | ✓ | ★● | ✓ |
| 7 | | 221-48876-05 | 5 | | SPL (Splitless) | ✓ | ✓ | ✓ | ✓ | ✓ |
| 8 | | 227-35016-01 | 5 | ✓ | SPL (Split / Splitless) | ★● | ★● | ★● | ● | ● |
| 9 | | 221-41544 221-41544-84 | 1 5 | | SPL (Split / Splitless) | ✓ | ✓ | ✓ | ✓ | ✓ |
| 10 | | 221-41544-05 | 5 | | SPL (Split / Splitless) | ✓ | ✓ | ✓ | ✓ | ✓ |
| 11 | | 221-75187 | 5 | ✓ | SPL (Split / Splitless) | ✓ | ✓ | ✓ | ✓ | ✓ |
| 12 | | 221-75188 | 5 | ✓ | SPL (Split / Splitless) | ✓ | ✓ | ✓ | ✓ | ✓ |
| 13 | | 221-75189 | 5 | ✓ | SPL (Split / Splitless) | ✓ | ✓ | ✓ | ✓ | ✓ |
| 14 | | 221-75190 | 5 | | SPL (Split / Splitless) | ✓ | ✓ | ✓ | ✓ | ✓ |
| 15 | | 221-75192 | 5 | | SPL (Split / Splitless) | ✓ | ✓ | ✓ | ✓ | ✓ |
| 16 | | 221-75194 | 5 | | SPL (Split / Splitless) | ✓ | ✓ | ✓ | ✓ | ✓ |



























★...Standard accessory ●...Described as a recommended item in Chapter 2 ✓...Can be used

*1 Special order for glass column + pTCD

*2 Special order for glass column + packed detector other than a TCD

*3 It can be used by remodeling the GC unit SPL to WBI with 221-74660-41 WBI Modification Kit.

*4 Connections to systems other than the GC-2010, GC-2014, or GC-2030 series are handled by special order.

| ID | Figure | P/N | Quantity | Wool Filling | Category | GC-2026 | GC-2030 | GC-2050 | GC-2010 series | GC-2014 series |
|----|---|---------------------------|----------|--------------|--|---------|---------|---------|----------------|----------------|
| 17 |  | 221-75195 | 5 | ✓ | SPL (Split / Splitless) | ✓ | ✓ | ✓ | ✓ | ✓ |
| 18 |  | 221-75196 | 5 | | SPME | ✓ | ✓ | ✓ | ✓ | ✓ |
| 19 |  | 227-35008-01 | 5 | ✓ | SPL/WBI (Splitless, Full Volume Injection) | ★● | ★● | ★● | ✓ | ✓ |
| 20 |  | 221-48335-01 | 1 | | SPL/WBI (Splitless, Full Volume Injection) | ✓ | ✓ | ✓ | ★● | ★● |
| 21 |  | 225-20803-01 | 5 | ✓ | SPL/WBI (Splitless, Full Volume Injection) | ✓ | ✓ | ✓ | ★● | ✓ |
| 22 |  | 221-75197 | 5 | ✓ | SPL/WBI (Splitless, Full Volume Injection) | ✓ | ✓ | ✓ | ✓ | ✓ |
| 23 |  | 221-41599 221-41599-84 | 1 5 | | WBI (Full Volume Injection) | ✓ | ✓ | ✓ | ✓ | ✓ |
| 24 |  | 221-41599-05 | 5 | | WBI (Full Volume Injection) | ✓ | ✓ | ✓ | ✓ | ✓ |
| 25 |  | 221-49298-91 | 1 | | OCI | ★● | ★● | | ★● | |
| 26 |  | 221-49381-01 | 1 | | Simple OCI | ✓ | ✓ | | ✓ | |
| 27 |  | 221-49381-02 | 1 | | Simple OCI | ● | ● | | ● | |
| 28 |  | 221-85694 | 1 | | OCI | ★● | ★● | | | |
| 29 |  | 225-08184 | 1 | | PTV/OCI | ✓ | ✓ | ✓ | ★ | |
| 30 |  | 221-49300 | 1 | | PTV | ★● | ★● | ★● | ★● | |
| 31 |  | 221-74830-09 | 5 | ✓ | PTV | ✓ | ✓ | ✓ | ✓ | |
| 32 |  | 221-80902 | 1 | | SINJ (Packed Column Full Volume Injection) | ● | ● | ● | | |
| 33 |  | 221-38107 | 1 | | SINJ/DINJ (WBC Full Volume Injection) | | | | | ★● |
| 34 |  | 221-38107-02 | 1 | | SINJ/DINJ (WBC Full Volume Injection) | ★● | ★● | ★● | | |
| 35 |  | 221-48993 | 1 | | SPL/WBI (Packed Column Full Volume Injection Method) | | | | ✓*1 | |
| 36 |  | 221-48886 | 1 | | SPL/WBI (Packed Column Full Volume Injection Method) | | | | ✓*2 | |
| 37 |  | 221-14093 221-14093-84 | 1 5 | | SINJ/DINJ (Packed Column Full Volume Injection Method) | | | | | ● |
| 38 |  | 221-14094 221-14094-84 | 1 5 | | SINJ/DINJ (Packed Column Full Volume Injection Method) | | | | | ● |
| 39 |  | 227-35015-01 | 5 | ✓ | Aqueous Solvent Analysis | ✓ | ✓ | ✓ | ✓ | ✓ |
| 40 |  | 227-35014-01 | 1 | | HS-10 | ● | ● | ● | ● | ● |
| 41 |  | 227-35327-03 | 3 | | AOC-6000 Plus SPME arrow | ● | ● | ● | ● | |
| 42 |  | 227-35328-03 | 3 | | AOC-6000 Plus SPME arrow | ● | ● | ● | ● | |

★...Standard accessory ●...Described as a recommended item in Chapter 2 ✓...Can be used

*1 Special order for glass column + pTCD

*2 Special order for glass column + packed detector other than a TCD

*3 It can be used by remodeling the GC unit SPL to WBI with 221-74660-41 WBI Modification Kit.

*4 Connections to systems other than the GC-2010, GC-2014, or GC-2030 series are handled by special order.

| ID | Figure | P/N | Quantity | Wool Filling | Category | GC-2026 | GC-2030 | GC-2050 | GC-2010 series | GC-2014 series |
|----|--------|--------------|----------|--------------|---|---------|---------|---------|----------------|----------------|
| 43 | | 222-30935 | 1 | | SINJ (liner for UG-type columns) | ★● | ★● | ★● | | |
| 44 | | 227-36700-01 | 5 | ✓ | MMI (Split) | ● | | | | |
| 45 | | 227-36701-01 | 5 | ✓ | MMI (Splitless) | ● | | | | |
| 46 | | 227-36702-01 | 5 | | MMI (Splitless with Cyclo) | ● | | | | |
| 47 | | 227-36703-01 | 5 | | MMI (Direct) | ● | | | | |
| 48 | | 227-36704-01 | 5 | ✓ | MMI (Large-volume injection, Split/Splitless) | ● | | | | |
| 49 | | 227-36705-01 | 5 | | MMI (SPME) | ● | | | | |
| 50 | | 227-36707-01 | 5 | | MMI (Split with Cyclo) | ● | | | | |
| 51 | | 227-36708-01 | 5 | | MMI (Gas analysis) | ● | | | | |
| 52 | | 227-36709-01 | 5 | ✓ | MMI (Splitless, Middle Inner Volume) | ● | | | | |
| 53 | | 227-36710-01 | 5 | ✓ | MMI (Split, Low Inner Volume) | ● | | | | |
| 54 | | 227-36711-01 | 5 | ✓ | MMI (Split/Splitless) | ● | | | | |

Standard Accessories for the GC-17A ver. 1-3, the GC-1700, and the GC-18A

| ID | P/N | Category |
|-----|----------------------------------|--|
| 3 | 221-41444 221-41444-84 (5 pc) | SPL (Split) |
| 9 | 221-41544 221-41544-84 (5 pc) | SPL (Splitless) |
| 23 | 221-41599 221-41599-84 (5 pc) | WBI (Full Volume Injection) |
| 100 | 221-41484 | Packed Column Full Volume Injection |

Standard Accessories for the GC-14A/B

| ID | P/N | Category |
|-----|----------------------------------|--|
| 101 | 221-32574-01 | SPL (Split) without taper |
| 102 | 221-32544-01 | SPL (Split) with taper |
| 103 | 221-32544 | SPL (Splitless) |
| 104 | 221-38151-04 | Septum purge unit (Full Volume Injection) |
| 33 | 221-38107 | WBC attachment (Full Volume Injection) |
| 105 | 221-32998-01 | CLH (Injection unit) |
| 106 | 221-33000 | CLH (Detector) |
| 37 | 221-14093 221-14093-84 (5 pc) | Dia. 3.2 mm packed Full Volume Injection |
| 38 | 221-14094 221-14094-84 (5 pc) | Dia. 2.6 mm packed Full Volume Injection |

Standard Accessories for the GC-8A

| ID | P/N | Category |
|-----|--------------|---|
| 107 | 221-25822-03 | SPL (Split) |
| 108 | 221-25944-03 | SPL (Splitless) |
| 109 | 221-39148 | WBC attachment (Full Volume Injection) |
| 110 | 221-18384-04 | CLH (Injection unit) |
| 111 | 221-18756-02 | CLH (Detector) |

3-2. Glass Liner Information

| ID | Figure | P/N | Quantity | Wool Filling | Length (mm) | Outer Diameter (mm) | Inner Diameter (mm) | Capacity (μL) | Surface Treatment | Maximum Temperature* |
|----|--------|---------------------------|----------|--------------|-------------|---------------------|---------------------|---------------|------------------------|----------------------|
| 1 | | 227-35007-01 | 5 | ✓ | 95 | 4.9 | 3.4 | 863 | Deactivation treatment | 350°C |
| 2 | | 221-41444-01 | 1 | | 95 | 4.9 | 3.4 | 863 | None | 450°C |
| 3 | | 221-41444 221-41444-84 | 1 5 | | 95 | 5 | 3.4 | 863 | None | 450°C |
| 4 | | 221-75193 | 5 | ✓ | 95 | 5 | 3.4 | 863 | Deactivation treatment | 350°C |
| 5 | | 221-48876-02 | 5 | ✓ | 95 | 5 | 3.4 | 654 | Deactivation treatment | 350°C |
| 6 | | 221-48876-03 | 5 | ✓ | 95 | 5 | 3.4 | 654 | Deactivation treatment | 350°C |
| 7 | | 221-48876-05 | 5 | | 95 | 5 | 3.4 | 654 | Deactivation treatment | 350°C |
| 8 | | 227-35016-01 | 5 | ✓ | 95 | 5 | 3.5 | 847 | Deactivation treatment | 350°C |
| 9 | | 221-41544 221-41544-84 | 1 5 | | 95 | 5 | 2.6 | 504 | None | 450°C |
| 10 | | 221-41544-05 | 5 | | 95 | 5 | 2.6 | 504 | Deactivation treatment | 350°C |
| 11 | | 221-75187 | 5 | ✓ | 95 | 5 | 3.4 | 799 | Deactivation treatment | 350°C |
| 12 | | 221-75188 | 5 | ✓ | 94.5 | 5 | 3.4 | 858 | Deactivation treatment | 350°C |
| 13 | | 221-75189 | 5 | ✓ | 95 | 5 | 3.4 | 863 | Deactivation treatment | 350°C |
| 14 | | 221-75190 | 5 | | 95 | 5 | 3.3 | 813 | Deactivation treatment | 350°C |
| 15 | | 221-75192 | 5 | | 95 | 5 | 3.5 | 847 | Deactivation treatment | 350°C |
| 16 | | 221-75194 | 5 | | 95 | 5 | 3.4 | 863 | Deactivation treatment | 350°C |
| 17 | | 221-75195 | 5 | ✓ | 95 | 5 | 3.4 | 863 | Deactivation treatment | 350°C |
| 18 | | 221-75196 | 5 | | 95 | 5 | 0.8 | 48 | Deactivation treatment | 350°C |
| 19 | | 227-35008-01 | 5 | ✓ | 95 | 5 | 3.4 | 654 | Deactivation treatment | 350°C |
| 20 | | 221-48335-01 | 1 | | 95 | 5 | 3.4 | 654 | None | 450°C |
| 21 | | 225-20803-01 | 5 | ✓ | 95 | 5 | 3.5 | 914 | Deactivation treatment | 350°C |
| 22 | | 221-75197 | 5 | ✓ | 95 | 5 | 3.3 | 599 | Deactivation treatment | 350°C |
| 23 | | 221-41599 221-41599-84 | 1 5 | | 95 | 5 | 2.6 | 504 | None | 450°C |



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.