

High Performance Liquid Chromatograph

Method Scouting System




Maximizing Efficiency for Method Development

Fully Automated Method Scouting

Method Scouting System is a method development system based on Shimadzu's UHPLC technology.

The combination with Method Scouting Solution dedicated control software achieves a fast and accurate method scouting workflow, offering excellent support for method development.



Improving Method Development Productivity

Method Scouting System assists in analyses involving a large number of parameters for columns, mobile phases, etc., such as method scouting for new compounds or chiral compounds in the pharmaceutical and chemical fields. It can automatically determine optimum analysis conditions, improving method development efficiency.

› Flexible System Options

Method scouting can be carried out with the various LC models provided by Shimadzu. Choose from more compact instruments with excellent ease-of-use, systems designed for SFC, etc., depending on your analysis needs.

› Fully Automated Method Scouting

Method Scouting System is able to perform highly efficient and robust method scouting by automating entirely the switching between different combinations of mobile phases and columns.

› Optimization and Robustness Evaluation

In addition to column and mobile phase selection, the system supports gradient optimization and robustness evaluations. Analysis is run with varying parameters such as gradient, oven temperature, and flow rate. It is also compatible with continuous analysis based on a user-determined experimental plan.




Flexible System Options

Method scouting can now be performed with Shimadzu's diverse LC product line-up. Choose the system suited to your laboratory environment, whether for SFC investigations or simple investigations on compact systems.


Nexera™ Method Scouting System

- Investigate analytical conditions with automated switching of columns and mobile phases
- Compatible with systems including 130 MPa UHPLC models or UHPLC-like models with method transfer.




Supercritical Fluid Chromatograph Nexera/UCs UHPLC/SFC Switching System

- Optimized for method investigation of chiral compounds
- Also compatible with LC/SFC switching, and can investigate various separation modes



i-Series Method Scouting System


- Saves costs in transferring to routine testing
- Space-saving, low-cost design




Fully Automated Method Scouting

Method Scouting Solution dedicated software automates the switching of column and mobile phase combinations, leading to efficient and reliable method development.

Previous Method




Method Scouting System + Method Scouting Solution



Leave Your Method and Analysis Scheduling to Method Scouting System

Usually when method scouting, a method has to be created every time the column or mobile phase is switched. For example, 100 combinations of analytical conditions leads to 100 method files, not only creating a huge amount of labor but also increasing the risk of operational mistakes.


Method Scouting System is intuitive to operate and can automate these changes in analytical conditions, allowing you to use your resources more efficiently.



Labor-Saving Mobile Phase Blending

The mobile phase blending function in Method Scouting System allows you to quickly adjust the salt concentration, pH, and organic solvent ratio on-line using only a few pre-prepared solvents, and quickly search for suitable conditions.

This significantly reduces the time required to prepare and evaluate mobile phases. The integrated i-Series LC system can also perform different types of mobile phase blending, two for water-based and two for organic solvents, using a low-pressure gradient.



Acquire data Analyze data


Column 3 → Column 4 → Overnight → Column 5 → Column 6



Improved Operating Rates with Automation of Analysis Sequences

Purge conditions and equilibration during automatic switching of mobile phases and columns are key issues in method scouting.

With Method Scouting Solution, these and all laboratory operations can be performed automatically based on conditions set in advance, from automated control to system checks and system shutdown.



Seamless Evaluation of Results in a Multi-Data Report

The data can be evaluated quantitatively with a multi-data report. Shimadzu's evaluation method, based on the degree of separation and the number of peaks detected, can quantitatively analyze chromatograms to derive the optimum method. Register templates provided in Method Scouting Solution and the report will be produced as soon as the analysis is completed, so the optimum conditions can be checked right away.



Ordered by evaluated score
A browser also allows for visual evaluation of multiple chromatograms.

Optimization and Robustness Evaluation

In addition to column and mobile phase selection, the system supports gradient optimization and robustness evaluations. Analysis is run with varying parameters such as gradient, oven temperature, and flow rate. It is also compatible with continuous analysis based on a user-determined experimental plan.

Simple to Set Conditions for Optimization

After choosing the columns and mobile phases, gradient conditions need to be optimized. Using Method Scouting Solution, it is simple to create an analysis schedule which changes these conditions automatically.

1. Pattern Settings

Gradient Mode: Linear

Pattern Generation Mode: Time

B. Conc. (%)

Time (min)

2. Generated Patterns

#	Use	I1	I2	I3	I4	Flow Rate	Oven Temp.
1	<input checked="" type="checkbox"/>	0.00	5.00	3.00	2.00	1.0000	40
2	<input checked="" type="checkbox"/>	0.00	6.00	3.00	2.00	1.0000	40
3	<input checked="" type="checkbox"/>	0.00	7.00	3.00	2.00	1.0000	40
4	<input checked="" type="checkbox"/>	0.00	8.00	3.00	2.00	1.0000	40
5	<input checked="" type="checkbox"/>	0.00	9.00	3.00	2.00	1.0000	40
6	<input checked="" type="checkbox"/>	0.00	1.00	3.00	2.00	1.5000	40
7	<input checked="" type="checkbox"/>	0.00	6.00	3.00	2.00	1.5000	40
8	<input checked="" type="checkbox"/>	0.00	7.00	3.00	2.00	1.5000	40
9	<input checked="" type="checkbox"/>	0.00	8.00	3.00	2.00	1.5000	40
10	<input checked="" type="checkbox"/>	0.00	9.00	3.00	2.00	1.5000	40
11	<input checked="" type="checkbox"/>	0.00	1.00	3.00	2.00	1.2000	40
12	<input checked="" type="checkbox"/>	0.00	6.00	3.00	2.00	1.2000	40
13	<input checked="" type="checkbox"/>	0.00	6.00	3.00	2.00	1.3000	40
14	<input checked="" type="checkbox"/>	0.00	7.00	3.00	2.00	1.3000	40

Parameter

Parameter	Initial Value	# of Steps	Step Size
Initial Conc. Holding Time	0.00	1	1.000
Gradient Time 1	5.00	5	1.000
Gradient Time 2	3.00	1	1.000
Equilibration Time	2.00	1	1.000
Flow Rate	1	3	0.1
Oven Temp.	40	5	2

Generate

100
80
60
40
20
0

B. Conc. (%)

Time (min)

OK Cancel

Easily create an analysis schedule for changing parameters:

- Gradient conditions (concentration/time)
- Column oven temperature
- Flow rate
- etc.

Set number of steps and step size for changing parameters

Efficient Robustness Evaluations

Method Scouting Solution can import any experimental design (Excel® format) and execute the analysis schedule. It can finely adjust analytical parameters such as flow rate and oven temperature to streamline robustness assessments that require multiple runs.

B. Conc.	pH	Oven Temp.	Column	Flow Rate
39	2.5	40	Scepter-OOS	1
40	2.5	40	Scepter-OOS	1
41	2.5	40	Scepter-OOS	1
40	2.4	40	Scepter-OOS	1

Import data

Experimental Design (Excel)

Import Settings

Batch creation


Batch schedule including column equilibration is automatically created

Link experimental plan and HPLC conditions

Improve Efficiency Even With Standard LCs

Operational efficiency can be improved even with a standard hardware configuration in the optimization and robustness evaluation stages that do not require automatic switching of mobile phases and columns.


Method Scouting Solution is compatible with LCs that do not include mobile phase or column switching valves, so you can use your current LC for more efficient method development.



Interface for a simple configuration

Seamless Method Migration with LabSolutions™ CS

Once the optimal conditions are determined, they may be applied as general-purpose analysis conditions and transferred to other working departments. The data needs to be managed over long time periods and retrieved smoothly when needed. The combination of LabSolutions CS and Method Scouting Solution not only makes the method development process more efficient, but also streamlines a series of workflow tasks from reporting results to other departments to transferring methods.



Application

Productivity Improvements for Method Scouting Process

When looking for the optimum mobile phase, column, and gradient conditions for a simultaneous analysis of multiple components, it is more efficient to break down the process into a number of steps. In the example shown below, simultaneous analysis conditions for a 13-component mixture of cephem antibiotics are determined by first narrowing down the mobile phase and column to items that are suited to the type of analysis (STEP 1) then optimizing the gradient conditions (STEP 2). This sequence of steps allows for efficient method scouting with reduced mobile phase consumption.

13-Component Mixture of Cephem Antibiotics

1. Cefsulodin
2. Cefadroxil
3. Cephapirin
4. Cefaclor
5. Cephalixin
6. Cephadrine
7. Cefotaxime
8. Cefazolin
9. Cefuroxime
10. Cefmetazole
11. Cefoxitin
12. Cefoperazone
13. Cephalothin


Preparation of Mobile Phases and Columns

- Mobile Phase: (A) (a) Phosphoric acid (sodium) buffer solution (pH 2.6)
 (b) Citric acid (sodium) buffer solution (pH 3.1)
 (c) Acetic acid (ammonium) buffer solution (pH 4.7)
 (d) Acetic acid (ammonium) aqueous solution (pH 6.7)
- (B) (a) Acetonitrile
 (b) Methanol
 (c) Acetonitrile/Methanol = 1/1 (v/v)
- Column: (1) Reversed phase column A (50 mm L. x 3.0 mm i.d., 2.2 mm)
 (2) Reversed phase column B (50 mm L. x 3.0 mm i.d., 2.2 mm)
 (3) Reversed phase column C (50 mm L. x 3.0 mm i.d., 2.2 mm)
 (4) Reversed phase column D (50 mm L. x 3.0 mm i.d., 2.6 mm)
 (5) Reversed phase column E (50 mm L. x 3.0 mm i.d., 2.6 mm)
 (6) Reversed phase column F (50 mm L. x 3.0 mm i.d., 2.6 mm)
- 4 aqueous mobile phases
 3 organic mobile phases
 6 columns

STEP 1 Select mobile phase and column.

STEP 2 Investigate gradient conditions.

Determine conditions for simultaneous analysis.




Mobile Phase : (A) Phosphoric acid (sodium) buffer solution (pH 2.6)
 (B) Acetonitrile
 Column : Reversed phase column E (50 mm L. x 3.0 mm i.d., 2.6 μm)
 Time Program : B. Conc. 5% (0 min) → 65% (5 min) → 90% (5.01–7 min)
 Flowrate : 1.0 mL/min
 Injection Volume : 5 μL
 Column Temp : 40°C
 Detection Wavelength: 260 nm (SPD-M20A)

Method Scouting Example for 13-Component Mixture of Cephem Antibiotics


Reducing the Labor Involved in Mobile Phase Blending

Method development for HPLC requires the selection of mobile phases and columns that are suited to the type of analysis. Mobile phase optimization has a particularly large influence on separation, and is therefore treated as important. However, mobile phase optimization requires the examination of a wide variety of parameters, including salt concentration, pH, and proportions of organic solvents. Because of this, mobile phase optimization commonly requires substantial time, effort, and skill. Nexera Method Scouting System has a mobile phase blending function that adjusts salt concentration, pH, and proportions of organic solvent on-line from multiple solvents that have been prepared in advance, and is quickly able to find conditions suited to a given type of analysis. This function can not only reduce considerably the time required to prepare and investigate mobile phases, it can also reduce solvent consumption.




Water, aqueous phosphoric acid solution, and aqueous disodium hydrogen phosphate solution are connected to solvent delivery unit A, and the pH and salt concentration can be adjusted at will by changing the mixing ratio.

Pump A Mixing Ratio (a/b/c)		Buffer Concentration (mmol/L)		
		2	5	10
pH	3.01	90/6/4	75/15/10	50/30/20
	5.67	90/5/5	75/12.5/12.5	50/25/25
	6.50	90/4/6	75/10/15	50/20/30



In reversed-phase chromatography, the type of organic solvent used affects retention behavior and peak shape. When a single organic solvent is used and a component cannot be separated or there is a desire to improve peak shape, mixing more than one organic solvent can resolve the problem.



Nexera, LabSolutions and ACQ are trademarks of Shimadzu Corporation.
Excel is either a registered trademark or a trademark of Microsoft Corporation in the United States and/or other countries.



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.