



# Pesticides Analysis Using the Agilent 5977A Series GC/MSD

## Application Note

Food Testing and Agriculture

### Authors

Jia-jia Wu and Wen-wen Wang  
Agilent Technologies Co. Ltd (China)  
Beijing,  
People's Republic of China

### Abstract

A method has been developed using the Agilent 5977A Series GC/MSD that can deliver excellent linearity ( $R^2 \geq 0.991$ ), reproducibility (RSD 11% in solvent and apple matrix), and sensitivity, most minimum detection limits (MDLs)  $\leq 6.7$  ppb of quantitation for 42 pesticides commonly analyzed in China.

### Introduction

With the increasing globalization of the food industry, there is greater scrutiny on food safety, resulting in major changes in the number of pesticides that are being regulated and monitored, as well as the allowable levels of those pesticides in food. There are more than 1,000 registered pesticides worldwide, and many countries employ strict regulations on pesticide residues in food and animal feed. Currently, the lowest maximum residue level (MRL) for most pesticides in China is 0.01 mg/kg, parts per million (ppm), and can range as high as 5 mg/kg, depending on the food matrix being tested.

This application note demonstrates the ability of the 5977A Series GC/MSD to provide sensitive, accurate, and reproducible analysis for 42 commonly tested pesticides in China, with MDLs well within the required MRLs.



**Agilent Technologies**

## Experimental

### Standards and Reagents

Standards for 42 pesticides were commercially obtained (Table 1). The working calibration and matrix spiking standards were prepared in hexane using the Agilent 7696A Sample Prep WorkBench.

Table 1. 42 Pesticide Standards

| No. | Classification   | CAS        | Name                    |
|-----|------------------|------------|-------------------------|
| 1   | Conazole         | 43121-43-3 | Triadimefon             |
| 2   | Dicarboximide    | 32809-16-8 | Procymidone             |
| 3   | Dicarboximide    | 50471-44-8 | Vinclozolin             |
| 4   | Dicarboximide    | 36734-19-7 | Iprodione               |
| 5   | Organochlorine   | 319-84-6   | BHC- <i>alpha</i>       |
| 6   | Organochlorine   | 319-85-7   | BHC- <i>beta</i>        |
| 7   | Organochlorine   | 58-89-9    | BHC- <i>gamma</i>       |
| 8   | Organochlorine   | 319-86-8   | BHC- <i>delta</i>       |
| 9   | Organochlorine   | 72-55-9    | DDE-p,p'                |
| 10  | Organochlorine   | 72-54-8    | DDD-p,p'                |
| 11  | Organochlorine   | 789-02-6   | DDT-o,p'                |
| 12  | Organochlorine   | 50-29-3    | DDT-p,p'                |
| 13  | Organochlorine   | 82-68-8    | Pentachloronitrobenzene |
| 14  | Organophosphorus | 62-73-7    | Dichlorvos              |
| 15  | Organophosphorus | 298-02-2   | Phorate                 |
| 16  | Organophosphorus | 122-14-5   | Fenitrothion            |
| 17  | Organophosphorus | 55-38-9    | Fenthion                |
| 18  | Organophosphorus | 2921-88-2  | Chlorpyrifos            |
| 19  | Organophosphorus | 56-38-2    | Parathion               |
| 20  | Organophosphorus | 24353-61-5 | Isocarbophos            |
| 21  | Organophosphorus | 732-11-6   | Phosmet                 |

Table 2. Agilent 7890B GC and Agilent 5977A Series GC/MSD Conditions

#### GC run conditions

|                           |  |
|---------------------------|--|
| Analytical column         | HP-5ms Ultra Inert 30 m × 250 µm, 0.25 µm (p/n 19091S-433UI)   |
| Injection volume          | 1 µL   |
| Injection mode            | Splitless  |
| Inlet temperature         | 280 °C   |
| Liner                     | UI, splitless, single taper, glass wool (p/n 5190-2293)  |
| Plated seal kit           | Gold Seal, Ultra Inert, with washer (p/n 5190-6144)  |
| Carrier gas               | Helium, constant flow, 1 mL/min  |
| Oven program              | 60 °C for 1 minute,<br>then 40 °C/min to 170 °C,<br>then 10 °C/min to 310 °C,<br>then hold for 2 minutes |
| Transfer line temperature | 280 °C   |

## Instruments

The study was performed on an Agilent 7890B gas chromatograph equipped with a Split/Splitless Inlet and coupled to an Agilent 5977A Series GC/MSD, using SIM and Electron Ionization (EI) acquisition modes. Table 2 lists the instrument conditions.

| No. | Classification      | CAS         | Name                    |
|-----|---------------------|-------------|-------------------------|
| 22  | Organophosphorus    | 2310-17-0   | Phosalone               |
| 23  | Organophosphorus    | 60-51-5     | Dimethoate              |
| 24  | Organophosphorus    | 333-41-5    | Diazinon                |
| 25  | Organophosphorus    | 298-00-0    | Parathion-methyl        |
| 26  | Organophosphorus    | 121-75-5    | Malathion               |
| 27  | Organophosphorus    | 41198-08-7  | Profenofos              |
| 28  | Organophosphorus    | 24017-47-8  | Triazophos              |
| 29  | Organophosphorus    | 99675-03-3  | Isofenphos-methyl       |
| 30  | Pyrazole            | 120068-37-3 | Fipronil                |
| 31  | Pyrethroid          | 39515-41-8  | Fenpropathrin           |
| 32  | Pyrethroid          | 91465-08-6  | Cyhalothrin             |
| 33  | Pyrethroid          | 68359-37-5  | Cyfluthrin              |
| 34  | Pyrethroid          | 70124-77-5  | Flucythrinate           |
| 35  | Pyrethroid          | 102851-06-9 | Fluvalinate- <i>tau</i> |
| 36  | Pyrethroid          | 82657-04-3  | Bifenthrin              |
| 37  | Pyrethroid          | 51877-74-8  | Permethrin              |
| 38  | Pyrethroid          | 52315-07-8  | Cypermethrin            |
| 39  | Pyrethroid          | 51630-58-1  | Fenvalerate             |
| 40  | Pyrethroid          | 52918-63-5  | Deltamethrin            |
| 41  | Pyrimidine          | 53112-28-0  | Pyrimethanil            |
| 42  | Substituted Benzene | 1897-45-6   | Chlorothalonil          |

#### MS conditions

|                        |             |
|------------------------|-------------|
| Solvent delay          | 3.5 minutes |
| Acquisition mode       | SIM         |
| Tune                   | Etune.u     |
| Gain factor            | 5.00        |
| Source temperature     | 250 °C      |
| Quadrupole temperature | 150 °C      |
| TID                    | On          |

## Sample Preparation

A 20 g apple sample was homogenized with 40 mL of acetonitrile and mixed vigorously for 1 minute in a vortex mixer. A 5-g amount of NaCl was added, and the sample was mixed on the vortex mixer for another minute. The sample was then centrifuged at 4,200 rpm for 5 minutes. A 20-mL aliquot of the supernatant, equivalent to 10 g of sample, was removed and concentrated to approximately 1 mL. The extract was then cleaned up using a Bond Elut carbon/NH<sub>2</sub> column (500 mg/500 mg, 6 mL, p/n 12252202). The effluent was

collected and evaporated to near dryness with a gentle stream of nitrogen. The residue was dissolved in hexane to a volume of 10 mL. One milliliter of the resulting blank extract, corresponding to 1 g of sample, was used to prepare a 50 ng/mL parts per billion (ppb) fortified matrix-matching standard (pesticide standards added to blank extracts).

## Acquisition Parameters

Table 3 shows the ions used for acquisition.

Table 3. Retention Times and Acquisition Parameters

| Target compound         | RT    | Q <sub>0</sub> | Q <sub>1</sub> | Q <sub>2</sub> | Q <sub>3</sub> |
|-------------------------|-------|----------------|----------------|----------------|----------------|
| Dichlorvos              | 4.57  | 109            | 185            | 79             | 187            |
| Phorate                 | 7.33  | 75             | 121            | 260            | 97             |
| BHC- <i>alpha</i>       | 7.47  | 181            | 219            | 183            | 217            |
| Dimethoate              | 7.61  | 87             | 93             | 125            | 143            |
| BHC- <i>beta</i>        | 7.84  | 217            | 181            | 183            | 219            |
| BHC- <i>gamma</i>       | 7.96  | 217            | 183            | 219            | 111            |
| Pentachloronitrobenzene | 8.04  | 237            | 249            | 295            | 214            |
| Pyrimethanil            | 8.10  | 198            | 199            | 200            | 77             |
| Diazinon                | 8.11  | 179            | 137            | 152            | 199            |
| BHC- <i>delta</i>       | 8.31  | 181            | 219            | 183            | 217            |
| Chlorothalonil          | 8.39  | 266            | 264            | 268            | 270            |
| Vinclozolin             | 8.91  | 212            | 285            | 198            | 187            |
| Parathion-methyl        | 8.94  | 263            | 109            | 125            | 79             |
| Fenitrothion            | 9.38  | 277            | 125            | 109            | 260            |
| Malathion               | 9.53  | 173            | 127            | 125            | 93             |
| Fenthion                | 9.71  | 278            | 125            | 109            | 169            |
| Chlorpyrifos            | 9.75  | 197            | 199            | 314            | 97             |
| Parathion               | 9.76  | 291            | 109            | 97             | 139            |
| Triadimefon             | 9.79  | 57             | 208            | 85             | 210            |
| Isocarbophos            | 9.87  | 136            | 121            | 120            | 110            |
| Isofenphos-methyl       | 10.19 | 199            | 58             | 121            | 231            |

| Target compound         | RT    | Q <sub>0</sub> | Q <sub>1</sub> | Q <sub>2</sub> | Q <sub>3</sub> |
|-------------------------|-------|----------------|----------------|----------------|----------------|
| Fipronil                | 10.43 | 367            | 369            | 213            | 351            |
| Procymidone             | 10.62 | 96             | 283            | 285            | 67             |
| Profenofos              | 11.31 | 339            | 139            | 206            | 208            |
| DDE- <i>p,p'</i>        | 11.39 | 246            | 318            | 316            | 248            |
| DDD- <i>p,p'</i>        | 12.13 | 235            | 237            | 165            | 236            |
| DDT- <i>o,p'</i>        | 12.20 | 235            | 237            | 165            | 236            |
| Triazophos              | 12.41 | 161            | 162            | 172            | 77             |
| DDT- <i>p,p'</i>        | 12.79 | 235            | 237            | 165            | 236            |
| Iprodione               | 13.45 | 314            | 187            | 189            | 244            |
| Phosmet                 | 13.65 | 160            | 161            | 77             | 93             |
| Bifenthrin              | 13.69 | 181            | 165            | 166            | 182            |
| Fenpropathrin           | 13.81 | 181            | 97             | 125            | 265            |
| Phosalone               | 14.33 | 182            | 121            | 184            | 367            |
| Cyhalothrin             | 14.63 | 181            | 197            | 208            | 209            |
| Permethrin              | 15.40 | 183            | 163            | 165            | 184            |
| Cyfluthrin              | 16.12 | 163            | 206            | 165            | 227            |
| Cypermethrin            | 16.24 | 181            | 163            | 165            | 77             |
| Flucythrinate           | 16.44 | 199            | 157            | 44             | 207            |
| Fenvalerate             | 17.14 | 167            | 125            | 181            | 152            |
| Fluvalinate- <i>tau</i> | 17.32 | 250            | 252            | 209            | 181            |
| Deltamethrin            | 17.85 | 181            | 253            | 251            | 255            |

RT – retention time in minutes

Q<sub>0</sub> – quantifier ion

Q<sub>1</sub>, Q<sub>2</sub>, Q<sub>3</sub> – qualifier ions

## Results and Discussion

### Linearity

Calibration curves were constructed from 20 to 200 ng/mL (ppb) in hexane for most target compounds. All 42 compounds had calibration coefficient values  $\geq 0.991$  (Table 4).

Table 4. Calibration Coefficients ( $R^2$ ) for the 42 Pesticides

| Target compound         | $R^2$ | Target compound         | $R^2$ |
|-------------------------|-------|-------------------------|-------|
| Dichlorvos              | 0.998 | Fipronil                | 0.997 |
| Phorate                 | 0.998 | Procymidone             | 0.998 |
| BHC- <i>alpha</i>       | 0.998 | Profenofos              | 0.997 |
| Dimethoate              | 0.995 | DDE-p,p'                | 0.998 |
| BHC- <i>beta</i>        | 0.995 | DDD-P,P'                | 0.998 |
| BHC- <i>gamma</i>       | 0.998 | DDT-o,p'                | 0.998 |
| Pentachloronitrobenzene | 0.997 | Triazophos              | 0.995 |
| Pyrimethanil            | 0.998 | DDT-p, p'               | 0.996 |
| Diazinon                | 0.998 | Iprodione               | 0.995 |
| BHC- <i>delta</i>       | 0.998 | Phosmet                 | 0.995 |
| Chlorothalonil          | 0.999 | Bifenthrin              | 0.996 |
| Vinclozolin             | 0.999 | Fenpropathrin           | 0.996 |
| Parathion-methyl        | 0.993 | Phosalone               | 0.995 |
| Fenitrothion            | 0.991 | Cyhalothrin             | 0.995 |
| Malathion               | 0.998 | Permethrin              | 0.992 |
| Fenthion                | 0.998 | Cyfluthrin              | 0.995 |
| Chlorpyrifos            | 0.998 | Cypermethrin            | 0.995 |
| Parathion               | 0.989 | Flucythrinate           | 0.994 |
| Triadimefon             | 0.994 | Fenvalerate             | 0.996 |
| Isocarbophos            | 0.991 | Fluvalinate- <i>tau</i> | 0.995 |
| Isofenphos-methyl       | 0.998 | Deltamethrin            | 0.992 |

### Reproducibility and Minimum Detection Limits in Solvent

Table 5 illustrates the good reproducibility obtainable with the 5977A Series GC/MSD in solvent. The majority of the relative standard deviation values (RSDs) across all 42 pesticide compounds were  $\leq 5\%$ , with none exceeding 11%, across eight injections.

The corresponding MDLs, calculated as the product of the RSD multiplied by the Student t test value, ranged from 3 to 14.3 ppb, with the majority being  $\leq 6.7$  ppb.

Table 5. Reproducibility RSDs and Calculated MDLs for a 50 ppb Standards Mix Sample in Solvent\*

| Compound                | RSD (%) | MDL (ppb) | Target compound         | RSD (%) | MDL (ppb) |
|-------------------------|---------|-----------|-------------------------|---------|-----------|
| Dichlorvos              | 5       | 6.4       | Fipronil                | 6       | 8.3       |
| Phorate                 | 5       | 6.3       | Procymidone             | 4       | 4.7       |
| BHC- <i>alpha</i>       | 3       | 4.3       | Profenofos              | 10      | 12.0      |
| Dimethoate              | 7       | 8.3       | DDE-p,p'                | 3       | 4.5       |
| BHC- <i>beta</i>        | 3       | 3.6       | DDD-p,p'                | 5       | 6.7       |
| BHC- <i>gamma</i>       | 2       | 3.0       | DDT-o,p'                | 3       | 3.6       |
| Pentachloronitrobenzene | 3       | 3.8       | Triazophos              | 8       | 9.9       |
| Pyrimethanil            | 4       | 5.0       | DDT-p,p'                | 3       | 4.0       |
| Diazinon                | 4       | 4.7       | Iprodione               | 10      | 12.2      |
| BHC- <i>delta</i>       | 3       | 3.5       | Phosmet                 | 8       | 9.7       |
| Chlorothalonil          | 3       | 4.1       | Bifenthrin              | 6       | 8.2       |
| Vinclozolin             | 3       | 4.2       | Fenpropathrin           | 6       | 8.6       |
| Parathion-methyl        | 5       | 5.9       | Phosalone               | 7       | 9.1       |
| Fenitrothion            | 5       | 6.6       | Cyhalothrin             | 8       | 10.6      |
| Malathion               | 5       | 6.2       | Permethrin              | 8       | 9.6       |
| Fenthion                | 5       | 5.9       | Cyfluthrin              | 9       | 11.3      |
| Chlorpyrifos            | 4       | 5.7       | Cypermethrin            | 9       | 11.0      |
| Parathion               | 5       | 5.6       | Flucythrinate           | 11      | 13.1      |
| Triadimefon             | 4       | 5.5       | Fenvalerate             | 9       | 12.0      |
| Isocarbophos            | 8       | 9.1       | Fluvalinate- <i>tau</i> | 11      | 14.3      |
| Isofenphos-methyl       | 5       | 6.7       | Deltamethrin            | 9       | 11.2      |

\*Eight consecutive injections were used to calculate the RSDs.

### Reproducibility and Minimum Detection Limits in Apple Matrix

Reproducibility was quite good in apple matrix, with many of the RSDs being lower than those calculated for analysis in solvent (Table 6 versus Table 5). All but two of the MDLs were  $\leq 7.5$  ppb, none was higher than 10.1 ppb, and all were well below the Chinese MRL requirements (Table 6).

Table 6. Reproducibility RSDs and Calculated MDLs for a 50 ppb Standards Mix Sample Spiked into Apple\*

| Target compound         | RSD (%) | Calculated MDL (ppb) | Required MRL† (ppb) | Target compound         | RSD (%) | Calculated MDL (ppb) | Required MRL† (ppb) |
|-------------------------|---------|----------------------|---------------------|-------------------------|---------|----------------------|---------------------|
| Dichlorvos              | 2       | 2.6                  | 200                 | Fipronil                | 2       | 3.4                  | 20 (cereal)         |
| Phorate                 | 3       | 4.2                  | 10                  | Procymidone             | 3       | 4.5                  | 5,000 (grape)       |
| BHC- <i>alpha</i>       | 3       | 4.1                  | 50**                | Profenofos              | 4       | 5.8                  | 50                  |
| Dimethoate              | 3       | 3.7                  | 1,000               | DDE-p,p'                | 3       | 3.5                  | 50***               |
| BHC- <i>beta</i>        | 4       | 5.9                  | 50**                | DDD-p,p'                | 1       | 1.8                  | 50***               |
| BHC- <i>gamma</i>       | 3       | 5.2                  | 50**                | DDT-o,p'                | 2       | 2.4                  | 50***               |
| Pentachloronitrobenzene | 3       | 4.0                  | 20 (watermelon)     | Triazophos              | 2       | 3.7                  | 200                 |
| Pyrimethanil            | 3       | 3.9                  | 1,000 (pear)        | DDT-p,p'                | 2       | 2.9                  | 50***               |
| Diazinon                | 3       | 4.1                  | 100 (cereal)        | Iprodione               | 3       | 5.3                  | 5,000               |
| BHC- <i>delta</i>       | 3       | 4.3                  | 50**                | Phosmet                 | 3       | 4.1                  | 5,000 (orange)      |
| Chlorothalonil          | 2       | 2.9                  | 1,000               | Bifenthrin              | 3       | 5.0                  | 500                 |
| Vinclozolin             | 3       | 4.0                  | 1,000 (cucumber)    | Fenpropathrin           | 3       | 3.8                  | 5,000               |
| Parathion-methyl        | 3       | 4.1                  | 10                  | Phosalone               | 2       | 2.6                  | 1,000 (spinach)     |
| Fenitrothion            | 2       | 3.2                  | 500                 | Cyhalothrin             | 4       | 5.1                  | 200                 |
| Malathion               | 3       | 3.7                  | 2,000               | Permethrin              | 4       | 5.5                  | 2000                |
| Fenthion                | 3       | 5.0                  | 50                  | Cyfluthrin              | 7       | 10.1                 | 500                 |
| Chlorpyrifos            | 3       | 4.5                  | 1,000               | Cypermethrin            | 5       | 7.3                  | 2,000               |
| Parathion               | 3       | 4.5                  | 10                  | Flucythrinate           | 5       | 7.5                  | 500                 |
| Triadimefon             | 3       | 3.7                  | 1,000               | Fenvalerate             | 4       | 5.7                  | 1,000               |
| Isocarbophos            | 3       | 6.0                  | 10                  | Fluvalinate- <i>tau</i> | 3       | 4.1                  | 500 (spinach)       |
| Isafenphos-methyl       | 3       | 3.7                  | 10                  | Deltamethrin            | 2       | 3.1                  | 100                 |

\* Eight consecutive injections were used to calculate the RSDs.

† Chinese regulation GB 2763-2012, MRL: Maximum Residue Level

\*\* Four isomers in total

\*\*\* Four DDT, DDE, DDD in total

## Conclusion

The Agilent 5977A Series GC/MSD can deliver sensitive, accurate and reproducible results for the analysis for 42 commonly tested pesticides in China, including MDLs in some food matrices that can be well within the required limits of Chinese regulation GB 2763-2012. The use of automation such as the Agilent 7696A Sample Prep WorkBench to prepare standards also results in more linear calibration curves.

## For More Information

These data represent typical results. For more information on our products and services, visit our Web site at [www.agilent.com/chem](http://www.agilent.com/chem).

[www.agilent.com/chem](http://www.agilent.com/chem)

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2013  
Printed in the USA  
April 12, 2013  
5991-2212EN



**Agilent Technologies**