

Poster Reprint

ASMS 2019

TP305

An Optimization Tool for MS Signal Acquisition in GC Triple Quadrupole Mass Spectrometry

Anastasia Andrianova¹, Melissa Churley²

¹Agilent Technologies Inc., Wilmington, DE 19808 USA; ²Agilent Technologies Inc., Santa Clara, CA 95051 USA

Introduction

Development of GC/MS/MS MRM transitions is a challenging and time-consuming multi-step process, which may be further complicated by analyte co-elution and matrix interferences.

A new modular optimization tool allows for fully automated end-to-end MRM development. Alternatively, each of the modules, i.e., MRM development steps, can be completed individually. These steps include:

- ✓ Precursor ion identification
- ✓ Product ion scan at various collision energies
- ✓ Selection of product ions
- ✓ Collision energy optimization

The optimization tool uses spectral deconvolution to reliably identify precursor ions, even in the presence of chromatographic interferences such as column bleed or co-eluting peaks. The tool significantly reduces the need for user intervention and generates the complete GC/MS/MS data acquisition method, containing the optimized MRM transitions.

Experimental

GC coupled to a triple quadrupole MS was used for MRM development. The tool is fully integrated with GC/MS data acquisition software.

It enables:

- ✓ Identification of the analytes using library search of deconvoluted spectra
- ✓ Selection of precursor ions
- ✓ For each selected product ion, a product ion scan at various collision energies
- ✓ Selection of the best product ions
- ✓ Collision energy optimization
- ✓ Generation of time-segment MRM or dynamic MRM acquisition methods


The tool allows the user to perform only those selected steps that are necessary, depending on the information already available. For example, precursor ions can be imported from an existing GC/MS SIM acquisition method.

The software is capable of generating the sequences necessary for gathering all the data required for optimization. All generated data files, methods, and log files are stored within a Project folder, and may be re-opened at a later date for further optimization.

Experimental


Optimizer for GC Triple Quad

Installed automatically with MassHunter GC/MS Data Acquisition Version 10.0. Supported for use with Agilent 7000 series and 7010 series TQ GC/MS.



Who will benefit from using the Optimizer?

- **New TQ users**
 - ✓ Convert existing SQ SIM or scan methods to MRM
- **Existing TQ users**
 - ✓ Re-optimize collision energies for current MRM method
 - ✓ Create MRMs for new targets using scan data
 - ✓ Update retention times for a new column configuration or oven program



The type of method (Scan, SIM, or MRM) supplied during setup will determine the optimization steps and the options available during setup. The Optimization Tool can automate portions of the optimization procedure, or the user may choose to manually review the results of each optimization step and optimize only the ions of interest.

MRM Optimization for an Environmental Contaminant: Fenarimol

Retention Time Determination and Precursor Ion Selection

Scan TIC of a pesticide mixture.


Deconvoluted mass spectrum of Fenarimol.

Available and selected precursor ions.

Fenarimol				Selected Precursor Ions			
Select	Mass	Abundance	%	Compound Name	RT (min)	Mass	
<input checked="" type="checkbox"/>	139	2,620,894.00	1.00	1	Fenarimol	15.062	295
<input checked="" type="checkbox"/>	107	1,864,054.00	0.71	2	Fenarimol	15.062	251
<input checked="" type="checkbox"/>	111	984,514.00	0.38	3	Fenarimol	15.062	219
<input checked="" type="checkbox"/>	219	983,875.00	0.38	4	Fenarimol	15.062	139
<input type="checkbox"/>	141	858,237.00	0.33	5	Fenarimol	15.062	111
<input checked="" type="checkbox"/>	251	851,637.00	0.32	6	Fenarimol	15.062	107
<input checked="" type="checkbox"/>	79	719,715.00	0.27	7	Fenarimol	15.062	79
<input type="checkbox"/>	75	650,970.00	0.25				
<input type="checkbox"/>	253	570,180.00	0.22				
<input type="checkbox"/>	191	433,197.00	0.17				
<input type="checkbox"/>	112	432,854.00	0.17				
<input type="checkbox"/>	113	431,769.00	0.16				
<input type="checkbox"/>	221	323,040.00	0.12				
<input type="checkbox"/>	189	319,126.00	0.12				
<input type="checkbox"/>	152	279,740.00	0.11				
<input type="checkbox"/>	140	231,718.00	0.09				
<input checked="" type="checkbox"/>	295	231,046.00	0.09				


Product Ion Selection

In a product ion scan experiment, a maximum of four collision energies may be used. In this experiment, the collision energies used were 5, 15, 25, and 35 eV.



Collision Energy Optimization

Once the MRM transitions are selected, the next step is to determine the optimal collision energy. The optimization can be performed over the defined range or within several steps around the optimal collision energy from the previous step.



Conclusions

- A highly-automated optimization tool for MS signal acquisition allowed for efficient development of method-ready MRM transitions.
- The Optimizer may be used as either a complete workflow with automation, or the user may review the results of each optimization step before continuing.
- The optimized results can be saved as MRM or dMRM method.