



Featured Application: Extractable Semivolatile Organic Compounds on Rxi-5Sil MS

Analyze More Semivolatile Samples per Shift Using Split Injection and an Rxi-5Sil MS Column

- Faster oven cycle increases sample throughput.
- Split injection reduces maintenance frequency.
- Reliably meets or exceeds method requirements for sensitivity and linearity.

Semivolatiles are typically analyzed using splitless injection, but this approach results in slow analysis times, frequent maintenance due to matrix accumulation, and injection-to-injection variability. Combined, these factors reduce the number of samples that can be analyzed before quality control criteria are no longer met. Using split injection is a simple way to increase sample throughput for semivolatiles analysis because it provides faster oven cycle times due to the use of higher initial oven temperatures. In addition, because less matrix is injected onto the column, more samples can be analyzed before downtime is required for maintenance.


As an alternative to splitless injection, using split injection with a robust, low-bleed Rxi-5Sil MS column produces highly repeatable results for trace-level analysis of semivolatiles, such as those listed in EPA Method 8270D. The selectivity of this column ensures reliable separation of problematic compounds, like benzo[b]fluoranthene and benzo[k]fluoranthene. An Rxi-5Sil MS column was selected for this work because it provides consistent chromatographic performance for semivolatiles analysis, making it easier to meet Method 8270 requirements while benefiting from the faster sample throughput obtained with split injection.

Peaks

1. 1,4-Dioxane-d8 (IS)	25. 2-Nitrophenol	49. Acenaphthylene	73. Carbazole
2. 1,4-Dioxane	26. 2,4-Dimethylphenol	50. 3-Nitroaniline	74. di-n-Butyl phthalate
3. N-Nitrosodimethylamine	27. Benzoic acid	51. Acenaphthene-d10 (IS)	75. Fluoranthene
4. Pyridine	28. Bis(2-chloroethoxy)methane	52. Acenaphthene	76. Benzidine
5. 2-Fluorophenol (SS)	29. 2,4-Dichlorophenol	53. 2,4-Dinitrophenol	77. Pyrene-D10 (SS)
6. Phenol-d6 (SS)	30. 1,2,4-Trichlorobenzene	54. 4-Nitrophenol	78. Pyrene
7. Phenol	31. Naphthalene-D8 (IS)	55. 2,4-Dinitrotoluene	79. p-Terphenyl-d14 (SS)
8. Aniline	32. Naphthalene	56. Dibenzofuran	80. 3,3'-Dimethylbenzidine
9. Bis(2-chloroethyl)ether	33. 4-Chloroaniline	57. 2,3,5,6-Tetrachlorophenol	81. Butyl benzyl phthalate
10. 2-Chlorophenol	34. Hexachlorobutadiene	58. 2,3,4,6-Tetrachlorophenol	82. Bis(2-ethylhexyl)adipate
11. 1,3-Dichlorobenzene	35. 4-Chloro-3-methylphenol	59. Diethyl phthalate	83. 3,3'-Dichlorobenzidine
12. 1,4-Dichlorobenzene-D4 (IS)	36. 2-Methylnaphthalene	60. 4-Chlorophenyl phenyl ether	84. Benz[a]anthracene
13. 1,4-Dichlorobenzene	37. 1-Methylnaphthalene	61. Fluorene	85. Chrysene-D12 (IS)
14. Benzyl alcohol	38. Hexachlorocyclopentadiene	62. 4-Nitroaniline	86. Bis(2-ethylhexyl)phthalate
15. 1,2-Dichlorobenzene	39. 2,4,6-Trichlorophenol	63. 4,6-Dinitro-2-methylphenol	87. Chrysene
16. 2-Methylphenol	40. 2,4,5-Trichlorophenol	64. N-nitrosodiphenylamine	88. Di-n-octyl phthalate
17. Bis(2-chloroisopropyl)ether	41. 2-Fluorobiphenyl (SS)	65. 1,2-Diphenylhydrazine	89. Benzo[b]fluoranthene
18. 4-Methylphenol	42. 2-Chloronaphthalene	66. 2,4,6-Tribromophenol (SS)	90. Benzo[k]fluoranthene
19. 3-Methylphenol	43. 2-Nitroaniline	67. 4-Bromophenyl phenyl ether	91. Benzo[a]pyrene
20. N-nitroso-di-n-propylamine	44. 1,4-Dinitrobenzene	68. Hexachlorobenzene	92. Perylene-D12 (IS)
21. Hexachloroethane	45. Dimethyl phthalate	69. Pentachlorophenol	93. Indeno[1,2,3-cd]pyrene
22. Nitrobenzene-D5 (SS)	46. 1,3-Dinitrobenzene	70. Phenanthrene-D10 (IS)	94. Dibenz[a,h]anthracene
23. Nitrobenzene	47. 2,6-Dinitrotoluene	71. Phenanthrene	95. Benzo[ghi]perylene
24. Isophorone	48. 1,2-Dinitrobenzene	72. Anthracene	

Chromatogram and conditions on page 2

See peak list
on page 1




Column Rxi-5Sil MS, 30 m, 0.25 mm ID, 0.25 μ m (cat.# 13623)
Sample 8270 MegaMix (cat.# 31850)
 8270 Benzidines mix (cat.# 31852)
 Benzoic acid (cat.# 31879)
 1,4-Dioxane (cat.# 31853)
 Revised B/N surrogate mix (cat.# 31888)
 Acid surrogate mix (4/89 SOW) (cat.# 31063)
 Revised SV internal standard mix (cat.# 31886)
Diluent: Methylene chloride
Conc.: 40 μ g/mL (IS/SS 20 μ g/mL)
Injection
Inj. Vol.: 1 μ L split (split ratio 10:1)
Liner: Premium 4 mm Precision liner w/wool (cat.# 23305.5)
Inj. Temp.: 270 $^{\circ}$ C
Split Vent Flow Rate: 12 mL/min
Oven
Oven Temp.: 70 $^{\circ}$ C (hold 1 min) to 285 $^{\circ}$ C at 28 $^{\circ}$ C/min to 305 $^{\circ}$ C at 3 $^{\circ}$ C/min to 320 $^{\circ}$ C at 30 $^{\circ}$ C/min (hold 1 min)
Carrier Gas He, constant flow
Flow Rate: 1.2 mL/min
Detector MS
Mode: Scan
Transfer Line Temp.: 280 $^{\circ}$ C
Analyzer Type: Quadrupole
Source Temp.: 270 $^{\circ}$ C
Quad Temp.: 150 $^{\circ}$ C
Electron Energy: 70 eV
Solvent Delay Time: 1.3 min
Tune Type: DFTPP
Ionization Mode: EI
Scan Range: 35-550 amu
Scan Rate: 5.36 scans/sec
Instrument Agilent 7890A GC & 5975C MSD

Rxi-5Siil MS Columns (fused silica)

low-polarity phase; Crossbond 1,4-bis(dimethylsiloxy)phenylene dimethyl polysiloxane

Description	temp. limits	qty.	cat.#
30 m, 0.25 mm ID, 0.25 µm	-60 to 320/350 °C	ea.	13623



8270 MegaMix Standard (76 components)

Acenaphthene (83-32-9)	1,2-Dichlorobenzene (95-50-1)	1-Methylnaphthalene (90-12-0)
Acenaphthylene (208-96-8)	1,3-Dichlorobenzene (541-73-1)	2-Methylnaphthalene (91-57-6)
Aniline (62-53-3)	1,4-Dichlorobenzene (106-46-7)	2-Methylphenol (<i>o</i> -cresol) (95-48-7)
Anthracene (120-12-7)	2,4-Dichlorophenol (120-83-2)	3-Methylphenol (<i>m</i> -cresol) (108-39-4)
Azobenzene (103-33-3)*	Diethylphthalate (84-66-2)	4-Methylphenol (<i>p</i> -cresol) (106-44-5)
Benz(a)anthracene (56-55-3)	2,4-Dimethylphenol (105-67-9)	Naphthalene (91-20-3)
Benzo(a)pyrene (50-32-8)	Dimethylphthalate (131-11-3)	2-Nitroaniline (88-74-4)
Benzo(b)fluoranthene (205-99-2)	Di- <i>n</i> -butyl phthalate (84-74-2)	3-Nitroaniline (99-09-2)
Benzo(ghi)perylene (191-24-2)	1,2-Dinitrobenzene (528-29-0)	4-Nitroaniline (100-01-6)
Benzo(k)fluoranthene (207-08-9)	1,3-Dinitrobenzene (99-65-0)	Nitrobenzene (98-95-3)
Benzyl alcohol (100-51-6)	1,4-Dinitrobenzene (100-25-4)	2-Nitrophenol (88-75-5)
Benzyl butyl phthalate (85-68-7)	4,6-Dinitro-2-methylphenol (Dinitro- <i>o</i> -cresol) (534-52-1)	4-Nitrophenol (100-02-7)
Bis(2-chloroethoxy)methane (111-91-1)	2,4-Dinitrophenol (51-28-5)	N-Nitrosodimethylamine (62-75-9)
Bis(2-chloroethyl)ether (111-44-4)	2,4-Dinitrotoluene (121-14-2)	N-Nitroso-di- <i>n</i> -propylamine (621-64-7)
Bis(2-ethylhexyl)adipate (103-23-1)	2,6-Dinitrotoluene (606-20-2)	2,2'-Oxybis(1-chloropropane) (108-60-1)
Bis(2-ethylhexyl)phthalate (117-81-7)	Di- <i>n</i> -octyl phthalate (117-84-0)	Pentachlorophenol (87-86-5)
4-Bromophenyl phenyl ether (101-55-3)	Diphenylamine (122-39-4)**	Phenanthrene (85-01-8)
Carbazole (86-74-8)	Fluoranthene (206-44-0)	Phenol (108-95-2)
4-Chloroaniline (106-47-8)	Fluorene (86-73-7)	Pyrene (129-00-0)
4-Chloro-3-methylphenol (59-50-7)	Hexachlorobenzene (118-74-1)	Pyridine (110-86-1)
2-Chloronaphthalene (91-58-7)	Hexachlorobutadiene (87-68-3)	2,3,4,6-Tetrachlorophenol (58-90-2)
2-Chlorophenol (95-57-8)	Hexachlorocyclopentadiene (77-47-4)	2,3,5,6-Tetrachlorophenol (935-95-5)
4-Chlorophenyl phenyl ether (7005-72-3)	Hexachloroethane (67-72-1)	1,2,4-Trichlorobenzene (120-82-1)
Chrysene (218-01-9)	Indeno(1,2,3- <i>cd</i>)pyrene (193-39-5)	2,4,5-Trichlorophenol (95-95-4)
Dibenz(a,h)anthracene (53-70-3)	Isophorone (78-59-1)	2,4,6-Trichlorophenol (88-06-2)
Dibenzofuran (132-64-9)		

1,000 µg/mL each in methylene chloride (3-methylphenol and 4-methylphenol at 500 µg/mL),
1 mL/ampul cat.# 31850 (ea.)

*1,2-diphenylhydrazine (8270-listed analyte) decomposes to azobenzene (mix component) in the injector.

**N-nitrosodiphenylamine (8270-listed analyte) decomposes to diphenylamine (mix component) in the injector.

Easier calibration!

8270 MegaMix and 8270 Matrix Spike Mix include 3-methylphenol and 4-methylphenol at 1/2 x concentration of other components.

Benzoic Acid

Benzoic acid (65-85-0)

2,000 µg/mL in methylene chloride, 1 mL/ampul	cat.# 31879 (ea.)
1,000 µg/mL in methanol, 1 mL/ampul	cat.# 31415 (ea.)

8270 Benzidines Mix (3 components)

Benzidine (92-87-5)

3,3'-Dichlorobenzidine (91-94-1)

3,3'-Dimethylbenzidine (*o*-tolidine) (119-93-7)

2,000 µg/mL each in methanol, 1 mL/ampul	cat.# 31688 (ea.)
2,000 µg/mL each in methylene chloride, 1 mL/ampul	cat.# 31852 (ea.)



Acid Surrogate Mix (4/89 SOW) (3 components)

2-Fluorophenol (367-12-4)

Phenol-d6 (13127-88-3)

2,4,6-Tribromophenol (118-79-6)

2,000 µg/mL each in methanol, 1 mL/ampul	cat.# 31025 (ea.)
2,000 µg/mL each in methanol, 1 mL/ampul	cat.# 31025.15 (15-pk.)
2,000 µg/mL each in methanol, 1 mL/ampul	cat.# 31025.25 (25-pk.)



1,4-Dioxane

1,4-Dioxane (123-91-1)

2,000 µg/mL in P&T methanol, 1 mL/ampul	cat.# 30287 (ea.)
2,000 µg/mL in methylene chloride, 1 mL/ampul	cat.# 31853 (ea.)
1.9 mg/mL in dimethyl sulfoxide, 1 mL/ampul	cat.# 36294 (ea.)

Revised B/N Surrogate Mix (4 components)

2-Fluorobiphenyl (321-60-8) *p*-Terphenyl-d14 (1718-51-0)
 Nitrobenzene-d5 (4165-60-0) Pyrene-d10 (1718-52-1)

1,000 µg/mL each in methylene chloride, 1 mL/ampul	cat.# 31887 (ea.)
--	-------------------

SV Internal Standard Mix (6 components)

Acenaphthene-d10 (15067-26-2) 1,4-Dichlorobenzene-d4 (3855-82-1) Perylene-d12 (1520-96-3)
 Chrysene-d12 (1719-03-5) Naphthalene-d8 (1146-65-2) Phenanthrene-d10 (1517-22-2)

2,000 µg/mL each in methylene chloride, 1 mL/ampul	cat.# 31206 (ea.)
2,000 µg/mL each in methylene chloride, 1 mL/ampul	cat.# 31206.15 (15-pk.)
2,000 µg/mL each in methylene chloride, 1 mL/ampul	cat.# 31206.25 (25-pk.)



Topaz 4.0 mm ID Precision Inlet Liner w/ Wool

for Agilent GCs equipped with split/splitless inlets

ID x OD x Length	qty.	cat.#
Precision, Premium Deactivation, Borosilicate Glass with Quartz Wool		
4.0 mm x 6.3 mm x 78.5 mm	5-pk.	23305

* 100% SATISFACTION GUARANTEE: If your Topaz inlet liner does not perform to your expectations for any reason, simply contact Restek Technical Service or your local Restek representative and provide a sample chromatogram showing the problem. If our GC experts are not able to quickly and completely resolve the issue to your satisfaction, you will be given an account credit or replacement product (same cat.#) along with instructions for returning any unopened product. (Do not return product prior to receiving authorization.) For additional details about Restek's return policy, visit www.restek.com/warranty



CE Ex

22655



22658



22652



22657



22653

Restek Electronic Leak Detector

Don't let a small leak turn into a costly repair—protect your analytical column by using a Restek leak detector.

Description	qty.	cat.#
Leak Detector with Hard-Sided Carrying Case and Universal Charger Set (U.S., UK, European, Australian)	ea.	22655
Small Probe Adaptor for Leak Detector	ea.	22658
Dynamic Duo Combo Pack (Restek Leak Detector and ProFLOW 6000 Flowmeter)	kit	22654
Soft-Sided Storage Case for Leak Detector or ProFLOW 6000 Flowmeter	ea.	22657
Car Charger/Adaptor	ea.	22652
Universal AC Power Adaptor	ea.	22653

Avoid using liquid leak detectors on a GC! Liquids can be drawn into the system and/or into the leak detector.

*Caution: The Restek electronic leak detector is designed to detect trace amounts of hydrogen in a noncombustible environment. It is NOT designed for determining leaks in a combustible environment. A combustible gas detector should be used for determining combustible gas leaks under any condition. When using it to detect hydrogen, the Restek electronic leak detector may only be used for determining trace amounts in a GC environment.

Restek Recommends: If you think that your Restek electronic leak detector needs service or repair, please contact Restek Customer Service before sending your unit in (cat.# 22655-R). Leak detector service/repair will include inspection and testing of the unit.