

Poster Reprint

ASMS 2025
Poster number ThP 083

New Pollutant Screening of Environmental Samples Using High-Resolution GC/Q-TOF and Accurate Mass Spectral Library of Emerging Environmental Pollutants

Ye Kong; Zhe Cao; Jiajia Wu; Shuai Wang;
Agilent Technologies, Inc. Beijing, China

Introduction

This paper conducts qualitative and semi-quantitative analysis of complex GC/Q-TOF sample data by using high-resolution, accurate-mass gas chromatography/quadrupole time-of-flight mass spectrometry (GC/Q-TOF), an accurate mass spectral library of environmental pollutants specifically created for emerging pollutant screening, and analytical software with a screening process platform. Actual environmental water samples are used to verify the sensitivity, accuracy, and spiked recovery rate of the solution. All results confirm that this solution is simple, stable, and effective, providing technical support for emerging pollutant screening.

Experimental

Sample Preparation

The pretreatment of environmental water samples entails the extraction of organic substances via an optimized liquid-liquid extraction solvent. This procedure can be performed either manually or automatically using the PAL3 RTC automated processing platform. Following this, data acquisition is carried out using the Agilent 7250 GC/Q-TOF instrument platform.




Figure 1. 7250 GC/Q-TOF system and hardware setup

Experimental

Instrument Analysis

The samples were analyzed by Agilent 8890 GC and 7250 high resolution accurate mass GC/Q-TOF system, It was configured with a mid-column backflush setup (Figure 1). A 55.5 - min retention time locked GC oven temperature program is utilized which is listed in Table 1. The key mass spectrometry parameters of the GC/Q-TOF are as follows: It operates in the Electron Ionization (EI) mode with an energy of 70 electron volts, conducts detection in the time-of-flight scan mode within a mass-to-charge ratio range of 50 - 1000. The spectral acquisition rate is 5 Hz, and both the interface temperature and the source temperature are set at 280 °C. Data analysis is performed with the assistance of MassHunter data analysis software. The high-resolution accurate mass spectral library employed in the screening process is independently created by the laboratory.

Table 1. GC-Q/TOF Operational Conditions.

GC and MS Conditions	Value
Column (2 ea.)	HP-5MS UI, 15 m, 0.25 mm ID, 0.25 μm film
Inlet	MMI, 4-mm UI liner with wool
Injection	1 μL by automatically in splitless mode 5 μL by manually in solvent vent mode
Injection temperature	280 °C
Inlet flow (column 1)	~ 1.0 mL/min
PUU flow (column 2)	column 1 flow + 0.2 mL/min
Oven program	60 °C for 1 min 40 °C/min to 120 °C, 0 min 5 °C/min to 310 °C, 15 min
Backflushing conditions	5 min (Post-run), 310 °C (Oven) 45 psi (Aux EPC), 2 psi (Inlet)
Transfer line temperature	280 °C
Ion source	EI, 70 eV
Ion source temperature	280 °C
Quadrupole temperature	150 °C
Spectral Acquisition	50-1000 m/z, 5 spectra/sec

Data Analysis

- The data processing used Agilent Masshunter data analysis software, including Suremass.
- The target screening of new pollutants was based on an independently created by the laboratory GC/Q-TOF PCDL library¹.

Results and Discussion

Results of Compound Screening and Semi-Quantification

Analysis by GC/Q-TOF provided high-efficiency screening of new pollutants. An accurate mass pesticides library (Figure 2) was employed to perform targeted screening using a tool to search by principle ions of each compound with the identified hits verified via mass accuracy, S/N, library match factor, RT match and coelution score. After rapid and efficient screening, Agilent's MassHunter quantitative software can be used to conduct semi-quantitative analysis on the positive results of the samples (Figure 3), meeting the qualitative and quantitative analysis needs of different operators.




Figure 2. List of compound distribution in the new pollutant PCDL

Figure 3. Overview of the New Contaminants PCDL Database

Sample Screening

Targeted screening was set up to automatically extract six ions per new pollutants from the PCDL, and to require at least two of these to produce EICs with a coelution score ≥ 75 and an $S/N \geq 3$. If a compound passing these requirements had an RT within ± 0.1 minutes, it was considered identified. This workflow is very simple and efficient. The complex data collected by GC/Q-TOF can be directly imported into the Masshunter screening process software to achieve automatic spectral library retrieval, generation of screening reports, and generation of semi-quantitative reports. In this paper, actual water samples were collected from two different areas, namely rivers near farmland and ordinary surface water. Using this screening process, trace amounts of new pollutant compounds were detected in both cases (Figure 5).

Method Sensitivity and Stability Validation

To verify the reliability of the method, experiments for verifying sensitivity, stability, and spiked recovery of actual samples were designed. In terms of sensitivity and stability, this scheme can achieve five consecutive repeated injections at trace concentrations of $0.1 \mu\text{g/L}$ and $0.5 \mu\text{g/L}$, with the data stability reaching $RSD\% < 10\%$. Regarding the spiked recovery of actual samples, for the automated pretreatment scheme, solvent standards were used to prepare reference substances, and the calculated recovery rates of samples spiked at a concentration of $5 \mu\text{g/L}$ were in the range of $84.4\% - 151.7\%$. For the manual scheme, matrix-spiked reference substances were used to calculate the recovery rates of $5 \mu\text{g/L}$ spiked samples, which were in the range of $91.6\% - 137\%$ (Figure 4).





Figure 4. The recovery results obtained after spiking at a level of $5 \mu\text{g/L}$ for the two pretreatment schemes respectively

Results and Discussion

Figure 5. Based on the PCDL quantitative screening process, pollutants such as PCBs (Polychlorinated Biphenyls), PAHs (Polycyclic Aromatic Hydrocarbons), organophosphates, UV absorbers, and plasticizers were detected in a certain surface water.



Semi-Quantification Result

In response to the current requirement of sometimes needing to issue semi-quantitative analysis reports for screening analysis, this method also explored the results of automatically generating semi-quantitative analysis reports using software. By using the setting of inherited calibration reference standards available in the Masshunter software, all other compounds were treated as inherited targets of the three target substances (PAH) for semi-quantitative concentration calculation. Finally, combined with the customized screening semi-quantitative analysis report template of the software, the semi-quantitative results of samples can be presented in the form of a classic report (Figure 6).

Statu	Compound	CAS#	Formula	RT/min	RT/min	Area	Height	Final Units Conc.	Lib Match Score	Mass Accuracy (PPM)	# of Verified Ions
+	Acenaphthene	83-32-9	C12H10	10.137	0.005	683541	258443	4.8633 ng/ml	99.8	153.0699	-1.69
+	Pyrene	129-00-0	C16H10	22.125	0.005	1252422	347907	8.2978 ng/ml	99.1	202.0777	-1.26
+	Benzo[a]pyrene	50-32-8	C20H12	33.817	0.003	680707	97954	72.2926 ng/ml	99.8	252.0934	-2.39
+	Hexachlorobutadiene	87-68-3	C4Cl6	5.820	0.005	26901	9622	0.1914 ng/ml	99.9	224.8408	-3.01
+	2-Chloronaphthalene, 2z	91-58-7	C10H7Cl	8.198	0.001	738558	287299	5.2547 ng/ml	99.4	162.0231	-0.72
+	o-BHC	319-84-6	C6H6Cl6	14.293	0.009	133340	44531	0.9487 ng/ml	99.9	182.9344	-1.45
+	Hexachlorobenzene	118-74-1	C6Cl6	14.559	0.001	446883	152238	3.1795 ng/ml	99.8	283.8096	0.26
+	p-BHC	319-85-7	C6H6Cl6	15.338	0.008	122828	40414	0.8739 ng/ml	99.9	182.9344	-0.82
+	TCDF / Tris(2-chloroethyl)phosphate	115-96-8	C6H12Cl3O4P	15.530	0.005	45448	13074	0.3234 ng/ml	97.5	248.9845	-0.61
+	γ-BHC (Lindane)	58-89-9	C6H6Cl6	15.560	0.001	132349	43927	0.9416 ng/ml	99.4	180.9373	-1.49
+	Phenanthrene	85-01-8	C14H10	15.896	0.000	1191643	383778	8.4783 ng/ml	99.3	178.0777	-0.57
+	m-BHC	319-86-8	C6H6Cl6	16.501	0.007	87911	28479	0.5824 ng/ml	99.4	180.9373	-1.84
+	2,4,4'-Trichlorobiphenyl, 28z	7012-37-5	C12H7Cl3	17.762	0.005	485843	156513	3.2189 ng/ml	99.9	255.9608	-0.78
+	Phosphoric acid, dibutyl phenyl ester	2528-36-1	C14H23O4P	17.838	0.002	111576	33001	0.7392 ng/ml	95.3	175.0155	-2.40
+	Chlordane-trans (γ-Chlordan)	5103-74-2	C10H6Cl8	22.002	0.000	164388	51719	1.0891 ng/ml	99.9	372.8254	-0.39
+	2,2',4,4',5,5'-Hexachlorobiphenyl, 01z	37680-73-2	C12H5Cl5	22.354	0.001	301065	98409	1.9947 ng/ml	100.0	325.8799	-0.46
+	Chlordane-cis (α-Chlordan)	5103-71-9	C10H6Cl8	22.570	0.004	142741	46055	0.9457 ng/ml	99.9	372.8254	-0.73
+	Endrin	72-20-8	C12H8Cl6O	24.183	0.009	42793	12767	0.2835 ng/ml	99.9	262.8564	-2.75
+	2,3,4,4',5,5'-Hexachlorobiphenyl, 14z	74472-37-0	C12H5Cl5	25.034	0.005	351629	104173	2.3297 ng/ml	99.9	325.8799	-0.66
+	o,p'-DDT	789-02-6	C14H9Cl5	25.061	0.001	149773	46315	0.9923 ng/ml	98.8	235.0076	-2.10
+	Endosulfan sulfate	1031-07-8	C9H6ClO4S	26.095	0.013	18221	5345	0.1202 ng/ml	99.3	271.8096	-4.03
+	p,p'-DDT	50-29-3	C14H9Cl5	26.295	0.005	106953	34491	0.7086 ng/ml	97.4	235.0076	-1.63
+	Phenol, 2-(2H-1,2,3-benzotriazol-2-yl)-4,6-bis(1,1-dimethylethyl)-	3846-71-7	C20H25N3O	28.370	0.000	283203	77554	30.0768 ng/ml	99.1	308.1758	-0.92
+	Hlrex	2385-85-5	C10Cl12	29.657	0.001	234855	70603	24.9421 ng/ml	99.9	271.8096	-1.05
+	2,3,3',4',4',5,5'-Hexachlorobiphenyl, 89z	39035-31-9	C12H3Cl7	30.827	0.003	203719	50871	21.6354 ng/ml	99.9	393.8020	-0.73

Figure 6. semi-quantitative screening analysis report can show the semi-quantitative results <https://www.agilent.com/en/promotions/asms>

Conclusions

- Based on the currently collected list of new pollutants, a truly meaningful new pollutant PCDL high-resolution spectral library has been created.
- This work demonstrated a total workflow for targeted screening and semi-quantification in real water sample.
- The confidence in identification of new pollutants was enhanced by stable RT and excellent mass accuracy as a result of using an RTL backflush method and high-resolution accurate mass measurement.

References

- Wong, L., Black, G., Young, T., Nieto, S. Accurate Mass Library for PFAS Analysis in Environmental Samples and Workflow for Identification of Pollutants in Drinking Water Using GC/Q-TOF. *Agilent Technologies Application Note*, 5994-6966EN, 2023

This information is subject to change without notice.

DE-006381

© Agilent Technologies, Inc. 2025
Published in USA, May 15, 2025