

# Perfluoro Acid Anhydrides

## Product Specification

Trifluoroacetic acid anhydride (TFAA), pentafluoropropionic acid anhydride (PFPA), and heptafluorobutyric acid anhydride (HFBA) are perfluoroacylated acylation reagents. These reagents form stable, volatile derivatives with alcohols, amines, and phenols. The major use for perfluoro acid anhydrides is to prepare electron-capturing derivatives for GC/electron capture detection. This detector provides greatly enhanced responses for halogenated derivatives. If reaction times and temperatures are adjusted to allow for the differences in reactivity, the similarities among TFAA, PFPA, and HFBA are sufficient to allow these reagents to be interchangeable.

Although the anhydride may be used on its own, the acylation reactions go most smoothly and quickly in a solvent and with a catalyst. Bases such as triethylamine (TEA) and trimethylamine (TMA) often are added to promote reactivity. PFPA and HFBA should be used with an acid scavenger, to help drive the reaction to completion and to prevent column damage from acidic byproducts of the derivatization reaction.

**TFAA** is the most reactive and most volatile of the three anhydrides. No acid byproducts are formed in derivatization reactions with TFAA. Amino acids and steroids are the analytes most commonly derivatized with TFAA. Amino acids are derivatized primarily to obtain volatile derivatives, not for monitoring with an ECD.

**PFPA** is used to prepare volatile derivatives of alcohols, amines, and phenols for use with either an ECD or an FID. PFPA should be used with an acid scavenger to prevent column damage. The derivatives require a low analysis temperature.

**HFBA** provides the derivatives most sensitive to ECD. Like PFPA, it yields volatile derivatives of alcohols, amines, and phenols for ECD or FID, and should be used with an acid scavenger to prevent column damage.

### Features/Benefits

Produce stable, volatile derivatives of alcohols, amines, and phenols for electron capture or flame ionization detection.

Frequently used in confirmation testing for drugs of abuse by GC/MS. (TFAA is used to identify methamphetamine, PFAA is used to identify opiates and benzoyllecgonine, HFBA is used to identify amphetamines and phencyclidine.)

### Typical Procedure

This procedure is intended to be a guideline and may be adapted as necessary to meet the needs of a specific application. Always take proper safety precautions when using an acetylating reagent – consult MSDS for specific handling information.


Prepare a reagent blank (all components, solvents, etc., *except sample*), following the same procedure as used for the sample.

1. Dissolve 50µg sample (250µg for FID) in 0.5mL benzene.
2. Add 0.1 mL 0.05M trimethylamine (acid scavenger) in benzene, followed by 10µL PFPA, HFBA, or TFAA.
3. Cap the vial and heat at 50°C for 15 minutes.
4. Cool the mix and add 1mL 5% ammonia in water.
5. Shake for 5 minutes, allow the layers to separate, and inject an aliquot of the benzene (upper) layer onto the chromatograph.

### Properties

#### Trifluoroacetic Acid Anhydride

Structure:



CAS Number: 407-25-0

Molecular Formula: (CF<sub>3</sub>CO)<sub>2</sub>O

Formula Weight: 210.03

bp: 39.5-40°C


d: 1.487

n<sub>D</sub>: <1.300

Appearance: clear, colorless liquid

#### Pentafluoropropionic Acid Anhydride

Structure:



CAS Number: 356-42-3

Molecular Formula: (CF<sub>3</sub>CF<sub>2</sub>CO)<sub>2</sub>O

Formula Weight: 310.05

bp: 69-70°C


d: 1.571

n<sub>D</sub>: <1.300

Appearance: clear, colorless liquid

#### Heptafluorobutyric Acid Anhydride

Structure:



CAS Number: 336-59-4

Molecular Formula: (CF<sub>3</sub>CF<sub>2</sub>CF<sub>2</sub>CO)<sub>2</sub>O

Formula Weight: 410.06

mp: -43°


d: 1.665

n<sub>D</sub>: 1.2870

Appearance: clear, colorless liquid

797-0244,0245,0246

Derivatization times vary widely, depending upon the specific compound(s) being derivatized. If derivatization is not complete under the procedure described here, the addition of a catalyst, use of another solvent, higher reaction temperature, longer reaction time, and/or higher reagent concentration should be evaluated.



797-0399

## Mechanism (1,2)

Acylation involves the introduction of an acyl group into a molecule that has a replaceable hydrogen atom (OH, NH, or SH group). **Except for TFAA**, anhydride acylating reagents form acidic byproducts that must be removed prior to GC analysis, to prevent destructive effects on the phase in the column. Consequently, acylations with anhydride reagents normally are performed in pyridine, tetrahydrofuran, or another solvent capable of accepting the acid byproducts.

## Toxicity – Hazards – Storage – Stability

Perfluoroacyl anhydrides are corrosive, flammable, and moisture-sensitive. Store in a bottle or ampul at room temperature, in a dry, well ventilated area. Use only in a well ventilated area and keep away from ignition sources. Moisture can hinder the reaction.

Recommended storage conditions for the unopened product are stated on the label. If you store an opened container or transfer the contents to another container for later reuse, validate that your storage conditions adequately protected the reagent.

## References

1. K. Blau and J. Halket *Handbook of Derivatives for Chromatography* (2nd ed.) John Wiley & Sons, New York, 1993.
2. D.R. Knapp *Handbook of Analytical Derivatization Reactions* John Wiley & Sons, New York, 1979.

## Additional Reading

- F.F. Lawrence and J.J. Ryan, *J. Chromatogr.* 130, 97 (1977).  
 D.E. Coffin, *J. Assoc. Off. Anal. Chem.* 52, 1044 (1969).  
 N.P. Sen, *J. Food Sci.* 34, 22 (1969).  
 D.D. Clarke, *et al.*, *J. Gas Chromatogr.* 5, 307 (1967).

## Ordering Information

Description	Cat. No.
<b>TFAA</b>	
10 x 1mL	3165-U
25mL	33164
<b>PFPA</b>	
10 x 1mL	33167
25mL	33168
<b>HFBA</b>	
10 x 1mL	33170-U
<b>Micoreaction Vessels with Hole Caps and Septa</b>	
1mL, pk. of 12	33293
3mL, pk. of 12	33297
5mL, pk. of 12	33299
<b>Books</b>	
<i>Handbook of Derivatives for Chromatography</i> K. Blau and J. Halket	Z246220
<i>Handbook of Analytical Derivatization Reactions</i> D.R. Knapp	23561

Contact our Technical Service Department (phone 800-359-3041 or 814-359-3041, FAX 800-359-3044 or 814-359-5468) for expert answers to your questions.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website ([www.sigma-aldrich.com](http://www.sigma-aldrich.com)), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA · Sigma-Aldrich de Argentina, S.A. · Buenos Aires 1119 AUSTRALIA · Sigma-Aldrich Pty. Ltd. · Castle Hill NSW 2154 AUSTRIA · Sigma-Aldrich Handels GmbH · A-1110 Wien  
 BELGIUM · Sigma-Aldrich N.V./S.A. · B-2880 Bornem BRAZIL · Sigma-Aldrich Quimica Brasil Ltda. · 01239-010 São Paulo, SP CANADA · Sigma-Aldrich Canada, Ltd. · 2149 Winston Park Dr., Oakville, ON L6H 6J8  
 CZECH REPUBLIC · Sigma-Aldrich s.r.o. · 186 00 Praha 8 DENMARK · Sigma-Aldrich Denmark A/S · DK-2665 Vallensbaek Strand FINLAND · Sigma-Aldrich Finland/YA-Kemia Oy · FIN-00700 Helsinki  
 FRANCE · Sigma-Aldrich Chimie · 38297 Saint-Quentin-Fallavier Cedex GERMANY · Sigma-Aldrich Chemie GmbH · D-82041 Deisenhofen GREECE · Sigma-Aldrich (o.m.) Ltd. · Ilioupoli 16346, Athens  
 HUNGARY · Sigma-Aldrich Kft. · H-1067 Budapest INDIA · Sigma-Aldrich Co. · Bangalore 560 048 IRELAND · Sigma-Aldrich Ireland Ltd. · Dublin 24 ISRAEL · Sigma Israel Chemicals Ltd. · Rehovot 76100  
 ITALY · Sigma-Aldrich s.r.l. · 20151 Milano JAPAN · Sigma-Aldrich Japan K.K. · Chuo-ku, Tokyo 103 KOREA · Sigma-Aldrich Korea · Seoul MALAYSIA · Sigma-Aldrich (M) Sdn. Bhd. · 58200 Kuala Lumpur  
 MEXICO · Sigma-Aldrich Química S.A. de C.V. · 50200 Toluca NETHERLANDS · Sigma-Aldrich Chemie BV · 3330 AA Zwijndrecht NORWAY · Sigma-Aldrich Norway · Torshov · N-0401 Oslo  
 POLAND · Sigma-Aldrich Sp. z o.o. · 61-663 Poznań PORTUGAL · Sigma-Aldrich Química, S.A. · Sintra 2710 RUSSIA · Sigma-Aldrich Russia · Moscow 103062 SINGAPORE · Sigma-Aldrich Pte. Ltd.  
 SOUTH AFRICA · Sigma-Aldrich (pty) Ltd. · Jet Park 1459 SPAIN · Sigma-Aldrich Química, S.A. · 28100 Alcobendas, Madrid SWEDEN · Sigma-Aldrich Sweden AB · 135 70 Stockholm  
 SWITZERLAND · Supelco · CH-9471 Buchs UNITED KINGDOM · Sigma-Aldrich Company Ltd. · Poole, Dorset BH12 4QH  
 UNITED STATES · Supelco · Supelco Park · Bellefonte, PA 16823-0048 · Phone 800-247-6628 or 814-359-3441 · Fax 800-447-3044 or 814-359-3044 · email: [supelco@sial.com](mailto:supelco@sial.com)

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.