



15th Multidimensional Chromatography Workshop

Alkene Quantitation in Plastic Waste-Derived Alternative Fuels using GC×GC FID

Genesis Barzallo, Hung Gieng, Petr Vozka

Department of Chemistry and Biochemistry, California State University, Los Angeles

gbarzal@calstatela.edu

C³AL
COMPLEX CHEMICAL COMPOSITION ANALYSIS LAB

Background

- California landfills enough plastic each day to fill 285 Olympic-size swimming pools (more than 12,000 tons).¹
- In 2021, there were 2.3 billion units sold in the world's tire market, which amounted to 54.5 million tons of tires (2.7 billion units in 2027)².



Chemical Conversion Processes



**UNIVERSITY OF
CHEMISTRY AND TECHNOLOGY
PRAGUE**


- Pyrolysis (+ hydrotreating) of plastic foils and waste tires

**P PURDUE
UNIVERSITY®**

- Low Pressure Hydrothermal Processing of polyolefin plastic waste



THE FRACTIONAL DISTILLATION OF PLASTIC WASTE



Chemical Conversion Processes



- Pyrolysis (+ hydrotreating) of plastic foils and waste tires



- Low Pressure Hydrothermal Processing of polyolefin plastic waste

Statement of the Problem


- High olefin content can be detrimental to engine cleanliness.
- Causes gum formations.

There are currently no methods for detailed quantitation of olefins, especially in high concentrations.




Methodology

- LECO's QuadJet GC×GC-FID with LN₂ thermal modulator
- Reversed-phase column configuration
 - (30 m × 0.25 mm × 0.25 μm Rxi-17 Sil MS, 0.9 m × 0.25 mm × 0.10 μm DB-1HT)
- Flame ionization detector (FID)




Methodology

- He flow rate of 1.5 mL/min
- Front inlet and FID temperatures: 275 and 280 °C, respectively
- Second oven and modulator temperature offsets: 35 and 20 °C, respectively
- Modulation period: 1.5 s (0.45 s hot pulse)
- Oven temperature program: from 40 °C (hold 0.2 min) to 200 °C (hold 1 min) with 3 °C/min temperature ramp rate
- Split injection: 1:20



GC×GC-FID Classification



Sample Preparation

Sample + Internal Standard +
DMDS + I₂ solution



70 °C oven for at least an hour





Sodium Thiosulfate



Organic phase analyzed




Theory



Model Compounds


Before Derivatization



Standard Compounds			
1-Hexadecene	1-Octadecene	4-Methyl-1-Pentene	1-Tetradecene
2,4-Dimethyl-1-Heptene	1,7-Octadiene	Cycloheptene	1,5-Cyclooctadiene
Cyclo-dodecatriene	Alpha-terpinene	4,4-Dimethyl-1-Cyclohexene	1-Heptene
Cyclohexene	Cyclo-dodecene	Trans-2,2-Dimethyl-3-Hexene	3-Ethyl-5-Methyl-2-Hexene
2,4-Dimethyl-1,3-Pentadiene	1-Octene	2-Ethyl-1-Hexene	2,5-Dimethyl-2,4-Hexadiene
1-Decene	Cis-Cyclooctene	Cis-Trans-3-Heptene	

Model Compounds


1 hr Post Derivatization



Standard Compounds			
1-Hexadecene	1-Octadecene	4-Methyl-1-Pentene	1-Tetradecene
2,4-Dimethyl-1-Heptene	1,7-Octadiene	Cycloheptene	1,5-Cyclooctadiene
Cyclo-dodecatriene	Alpha-terpinene	4,4-Dimethyl-1-Cyclohexene	1-Heptene
Cyclohexene	Cyclo-dodecene	Trans-2,2-Dimethyl-3-Hexene	3-Ethyl-5-Methyl-2-Hexene
2,4-Dimethyl-1,3-Pentadiene	1-Octene	2-Ethyl-1-Hexene	2,5-Dimethyl-2,4-Hexadiene
1-Decene	Cis-Cyclooctene	Cis-Trans-3-Heptene	

Model Compounds


2 hr Post Derivatization



Standard Compounds			
1-Hexadecene	1-Octadecene	4-Methyl-1-Pentene	1-Tetradecene
2,4-Dimethyl-1-Heptene	1,7-Octadiene	Cycloheptene	1,5-Cyclooctadiene
Cyclo-dodecatriene	Alpha-terpinene	4,4-Dimethyl-1-Cyclohexene	1-Heptene
Cyclohexene	Cyclo-dodecene	Trans-2,2-Dimethyl-3-Hexene	3-Ethyl-5-Methyl-2-Hexene
2,4-Dimethyl-1,3-Pentadiene	1-Octene	2-Ethyl-1-Hexene	2,5-Dimethyl-2,4-Hexadiene
1-Decene	Cis-Cyclooctene	Cis-Trans-3-Heptene	

Model Compounds


3 hr Post Derivatization



Standard Compounds			
1-Hexadecene	1-Octadecene	4-Methyl-1-Pentene	1-Tetradecene
2,4-Dimethyl-1-Heptene	1,7-Octadiene	Cycloheptene	1,5-Cyclooctadiene
Cyclo-dodecatriene	Alpha-terpinene	4,4-Dimethyl-1-Cyclohexene	1-Heptene
Cyclohexene	Cyclo-dodecene	Trans-2,2-Dimethyl-3-Hexene	3-Ethyl-5-Methyl-2-Hexene
2,4-Dimethyl-1,3-Pentadiene	1-Octene	2-Ethyl-1-Hexene	2,5-Dimethyl-2,4-Hexadiene
1-Decene	Cis-Cyclooctene	Cis-Trans-3-Heptene	

Model Compounds


4 hr Post Derivatization



Standard Compounds			
1-Hexadecene	1-Octadecene	4-Methyl-1-Pentene	1-Tetradecene
2,4-Dimethyl-1-Heptene	1,7-Octadiene	Cycloheptene	1,5-Cyclooctadiene
Cyclo-dodecatriene	Alpha-terpinene	4,4-Dimethyl-1-Cyclohexene	1-Heptene
Cyclohexene	Cyclo-dodecene	Trans-2,2-Dimethyl-3-Hexene	3-Ethyl-5-Methyl-2-Hexene
2,4-Dimethyl-1,3-Pentadiene	1-Octene	2-Ethyl-1-Hexene	2,5-Dimethyl-2,4-Hexadiene
1-Decene	Cis-Cyclooctene	Cis-Trans-3-Heptene	

Model Compounds


5 hr Post Derivatization



Standard Compounds			
1-Hexadecene	1-Octadecene	4-Methyl-1-Pentene	1-Tetradecene
2,4-Dimethyl-1-Heptene	1,7-Octadiene	Cycloheptene	1,5-Cyclooctadiene
Cyclo-dodecatriene	Alpha-terpinene	4,4-Dimethyl-1-Cyclohexene	1-Heptene
Cyclohexene	Cyclo-dodecene	Trans-2,2-Dimethyl-3-Hexene	3-Ethyl-5-Methyl-2-Hexene
2,4-Dimethyl-1,3-Pentadiene	1-Octene	2-Ethyl-1-Hexene	2,5-Dimethyl-2,4-Hexadiene
1-Decene	Cis-Cyclooctene	Cis-Trans-3-Heptene	

Model Compounds


24 hr Post Derivatization



Standard Compounds			
1-Hexadecene	1-Octadecene	4-Methyl-1-Pentene	1-Tetradecene
2,4-Dimethyl-1-Heptene	1,7-Octadiene	Cycloheptene	1,5-Cyclooctadiene
Cyclo-dodecatriene	Alpha-terpinene	4,4-Dimethyl-1-Cyclohexene	1-Heptene
Cyclohexene	Cyclo-dodecene	Trans-2,2-Dimethyl-3-Hexene	3-Ethyl-5-Methyl-2-Hexene
2,4-Dimethyl-1,3-Pentadiene	1-Octene	2-Ethyl-1-Hexene	2,5-Dimethyl-2,4-Hexadiene
1-Decene	Cis-Cyclooctene	Cis-Trans-3-Heptene	

Model Compounds

36 hr Post Derivatization




Standard Compounds			
1-Hexadecene	1-Octadecene	4-Methyl-1-Pentene	1-Tetradecene
2,4-Dimethyl-1-Heptene	1,7-Octadiene	Cycloheptene	1,5-Cyclooctadiene
Cyclo-dodecatriene	Alpha-terpinene	4,4-Dimethyl-1-Cyclohexene	1-Heptene
Cyclohexene	Cyclo-dodecene	Trans-2,2-Dimethyl-3-Hexene	3-Ethyl-5-Methyl-2-Hexene
2,4-Dimethyl-1,3-Pentadiene	1-Octene	2-Ethyl-1-Hexene	2,5-Dimethyl-2,4-Hexadiene
1-Decene	Cis-Cyclooctene	Cis-Trans-3-Heptene	

Calculations

Olefins in sample

- Iso-alkenes
- Olefins with 1 double bond + Linear-alkenes
- Olefins with 2 double bonds
- Olefins with 3 double bonds



Equation

$$Wt. \% Olefin, C\# = P.A. Pre-Derivatization, C\# - P.A. Post-Derivatization and Normalization, C\#$$

Example

$$Wt. \% Iso-alkene, C11 = P.A. Pre-Derivatization, C11 - P.A. Post-Derivatization and Normalization, C11$$


Scrap Tire Pyrolysis Oil Results

Before	wt. %		After	wt. %
<i>n</i> -alkanes	1.54		<i>n</i> -alkanes	1.54
Iso-alkanes + Iso-alkenes	0.86		Iso-alkanes	0.64
			Iso-alkenes	0.22
Monocycloalkanes + Olefins with 1 Double Bond + Linear-alkenes	4.43		Monocycloalkanes	0.96
			Linear alkenes	1.04
			Olefins with 1 Double Bond	2.43
Dicycloalkanes + Olefins with 2 Double Bonds	27.01		Dicycloalkanes	1.80
			Olefins with 2 Double Bonds	25.20
Tricycloalkanes + Olefins with 3 Double Bonds	4.60		Tricycloalkanes	1.10
			Olefins with 3 Double Bonds	3.50
Aromatics	53.66		Aromatics	53.66
Light Hydrocarbons	7.90		Light Hydrocarbons	7.90


Results

	Waste Tire Pyrolysis, Gasoline		Post Derivatization				
Monocycloalkanes	p.a.	wt. %	p.a.	After Norm p.a.	Olefins (p.a.)	Olefin (wt.%)	Real Cyclo (wt.%)
C5	52883.75	6.00	38374.69	31454.66	21429.09	2.32	3.68
C6	81707.23	9.27	65393.26	53601.03	28106.20	3.05	6.23
C7	86175.05	9.78	64515.59	52881.63	33293.42	3.61	6.17
C8	50744.41	5.76	50393.38	41306.05	9438.36	1.02	4.74
C9	56524.36	6.41	16947.64	13891.51	42632.85	4.62	1.80
C10	1839.73	0.21	553.02	453.30	1386.43	0.15	0.06
C11	0	0.00	0	0	0	0	0.00
C12	0	0.00	0	0	0	0	0.00
C13	0	0.00	0	0	0	0	0.00
C14	0	0.00	0	0	0	0	0.00
C15	0	0.00	0	0	0	0	0.00
C16	0	0.00	0	0	0	0	0.00
Total	329874.53	37.44	236177.58	193588.18	136286.35	14.77	22.67


Validation (Compared to other ASTM methods)




Validation (Compared to other ASTM methods)




Validation (Compared to other ASTM methods)




Validation (Compared to other ASTM methods)



Validation (Compared to other ASTM methods)



Validation (Compared to other ASTM methods)




Conclusion

- We introduced a possible method for olefin quantitation in complex alternative fuels made from the conversion of plastic waste.
- Validate the results using GC×GC with time-of-flight mass spectrometry (TOFMS).




Pegasus BT 4D GC×GC-TOFMS



Chromatogram – diesel-like sample pre-derivatization



Chromatogram – diesel-like sample post-derivatization



3



Identified Compounds

e	4-Methyl-1-Pentene	1-Tetradecene
e	Cycloheptene	1,5-Cyclooctadiene
	4,4-Dimethyl-1-Cyclohexene	1-Heptene
	Trans-2,2-Dimethyl-3-Hexene	3-Ethyl-5-Methyl-2-Hexene
	2-Ethyl-1-Hexene	2,5-Dimethyl-2,4-Hexadiene
	Cis-Trans-3-Heptene	

Conclusion

- We introduced a possible method for olefin quantitation in complex alternative fuels made from the conversion of plastic waste.
- Validate the results using GC×GC with time-of-flight mass spectrometry (TOFMS).



Pegasus BT 4D GC×GC-TOFMS

References

1. Lam, M. (2022, November 3) “You Can’t Recycle Your Way Out”: California’s plastic problem and what we can do about it. *KQED*.
2. Tire Market: Global Industry Trends, Share, Size, Growth, Opportunity and Forecast 2022-2027, 2022.
3. Carlson, D. A.; Roan, C. S.; Yost, R. A.; Hector, J. Dimethyl disulfide derivatives of long chain alkenes, alkadienes, and alkatrienes for gas chromatography/mass spectrometry. *Anal. Chem.* **1989**, *61*, 1564-1571.
4. Curiale, J. A.; Frolov, E. B. Occurrence and origin of olefins in crude oils. A critical review. *Org. Geochem.* **1998**, *29*, 397-408.
5. Auersvald, M.; Staš, M.; Šimáček, P. Electrometric bromine number as a suitable method for the quantitative determination of phenols and olefins in hydrotreated pyrolysis bio-oils. *Talanta (Oxford)* **2021**, *225*, 122001
6. Lecocorp. (2022, February 17) ChromaTOF® tile analytical software. *LECO Website Main Page*. LECO Corporation.
7. Jin, K.; Vozka, P.; Gentilcore, C.; Kilaz, G.; Wang, N. L. Low-pressure hydrothermal processing of mixed polyolefin wastes into clean fuels. *Fuel (Guildford)* **2021**, *294*, 120505

Dr. Vozka's Research Lab



calstatela.edu/research/c3al

pvozka@calstatela.edu

www.linkedin.com/in/genesis-b-042363b3

Follow us on
Instagram!



@CAL_C3AL



CAL STATE LA
CALIFORNIA STATE UNIVERSITY, LOS ANGELES

