

Implementing Tile-Based Fisher Ratio Analysis of GC×GC-TOFMS Data to Obtain a Master Peak Table of All Detected Analyte Compounds in Many Petroleum-Based Samples

MDCW 2026

Rachel Halvorsen¹, Wenjing Ma¹, Caitlin Cain¹, Hep Ingham², Rachel Mohler², Rob Synovec¹

¹ Department of Chemistry, University of Washington
Seattle, WA 98195

² Chevron Technical Center
Richmond, CA 94801







Study Background: Fuel Sample Analysis




- Chevron's Technical Center contracted group to develop more advanced fuel analysis techniques
- Direct improvement by simply moving from GC to GC×GC for separation
- Increased peak capacity and resolution necessitated improved chemometrics for peak tables
- Novel method adapted from tool used to find peak differences between samples forwarded




Quick Overview of GC×GC



Modulate




Visualize





Chemometric Method Overview




Tile-Based
F-Ratio

Single Comprehensive Hit List

Hit #	¹ Dt _R (min)	² Dt _R (sec)	F-Ratio	ID	MV
1	3.97	1.66	High	Pentane	877
...
N	31	5.16	Low	Dodecane	923

p-Test



SSPT: Sample Specific Peak Table



Tile-Based F-Ratio (Through ChromaTOF Tile)

Utilize a tile-based fisher ratio method that discovers analytes that are statistically different in concentration between sample classes.

$$\text{Fisher Ratio} = \frac{\text{Between Class Variance}}{\sum(\text{Within Class Variance})}$$




Figure 1. GCxGC chromatographic plane divided into rectangular tiles.






Figure 2. Feature intensity vs. sample number for two significant (high F-ratio or low p=-value) features. Samples 1, 2, and 3 (purple) belong to one sample class and features 4, 5, and 6 (orange) belong to the other sample class.



Sample Class Linear-Scale Chromatograms




FRH: Fisher Ratio Hitlist



- Samples run on Rxi-1 column (30 m × 0.25 mm × 0.25 μm) first dimension and Rxi-17MS column (2.0 m × 0.18 mm × 0.20 μm) on second dimension.
- Sample classes are all fuels, mixture of different manufacturing methods and samples taken from different batches
- Major components all elute earlier than 15-27 minutes depending on sample





Sample Classes: Log-Scale Chromatograms



- Several samples have many more late-eluting minor components!
- Variations in these minor components also more pronounced between sample batches



Histogram of 719 hits' F-Ratios




- As expected from chromatograms showing a variety of large and small peaks across differing samples, we get a wide range of F-Ratios for our hits. While traditionally we'd discard hits below a certain F-Ratio, we take a different approach here
- Too many variables for F-Ratio alone to tell us much about each hit




Use of Blank Tiles for p-testing



- Using ChromaTOF Tile tells us which samples have relatively higher peaks for given compounds but does not tell us directly which compounds are detectable.
- To resolve this, three blank chromatograms are used to generate per-tile p-tests to be applied to all hits across nine sample classes.
- To p-test, we need a Confidence Interval – how to choose?



3 Blank Replicates



**Inferred Gaussian
Distribution of Summed
Tile Blank Signal**

99% CI
p-test threshold

Hit #	H11	H18	...
1	Pass	Pass	...
2	Pass	Pass	...
3	Pass	Fail	...
4	Fail	Pass	...


**Master Peak Table
(MPT)**




Example of Detectable Peak: 1-Methylnaphthalene @ m/z 115




Hydrobate 11




Reformat 11




Naphtha 007




Naphtha 009




FCC Gasoline




Hydrobate 18




Reformat 18




Naphtha 008



Naphtha 010




Blank






Summed m/z 115 Tile Signal with 99% CI Limit





Summed m/z 115 Tile Signal with 99% CI Limit: Zoomed





More Example Hits at 99% CI




Hit #	F-ratio	F-ratio <i>m/z</i>	Nominal Identification	MV	RT1(min)	RT2 (sec)	In H11	In H18	In R11	In R18	In N007	In N008	In N009	In N010	In FCC
1	2.85E+04	110	1,2,3-trimethylcyclopentene	877	3.97	1.66	3	3	3	3	3	3	3	3	3
17	4.51E+03	98	2-octene	882	4.57	1.91	3	3	3	3	3	3	3	3	3
73	1.29E+03	202	β-vatirenene	818	29.26	4.06	3	3	0	0	3	0	0	3	0
80	1.22E+03	172	2,6,10,14-tetramethylheptadecane	872	14.87	1.94	0	0	0	0	0	0	0	3	3
83	1.17E+03	151	1-methylene-1H-indene	948	15.87	5.44	0	0	3	3	2	3	3	3	3
95	1.07E+03	119	1,2,4-timethylbenzene	867	9.07	2.89	3	3	3	3	3	3	3	3	3
169	5.47E+02	169	2,4-diethyl-1-methylbenzene	901	15.57	3.38	0	0	0	0	0	3	3	1	3
247	2.56E+02	64	1,2-dimethyl-4-(phenylmethyl)-benzene	872	31.46	5.18	0	0	3	3	0	3	3	0	0
292	2.06E+02	79	2,6-dimethylundecane	874	19.33	2.02	0	0	3	0	0	3	3	0	3
349	1.68E+02	115	1-methylnaphthalene	923	19.7	5.26	0	0	3	3	3	3	3	3	3

- Example hits align with our qualitative analysis (heavier analytes are found in heavier samples, more aromatic analytes are found in samples with higher aromatic content, etc.)
- Out of all classifications across all 719 hits, only 6.3% are heterogeneous across replicates (2.53% 1/3 passed and 3.77% 2/3 passed), this indicates good reproducibility.

Sample-Specific Peak Tables, Total #s of Peaks

(Unshaded for Entire Chromatogram, Shaded for 99% Summed TIC)





Conclusions



- Developed method gains benefits of tiling, particularly minimization of retention time shifting, and allows recovery of trace analytes with reduced need for user input/intervention
- Method should be widely applicable to other samples
- Similar approaches could likely be applied with different chemometric metrics other than F-Ratio to better fit experimental needs, e.g. variance ranking for cases where replicates are impractical



Principal Investigator
Rob Synovec

Graduate Students

Rachel Halvorsen

Augie Dobrecevich

Wenjing Ma

Haylee Meissner

Cassy Padilla

Valencia Parker

Jungho Ahn

Undergraduate Students

Stephanie Nguyen




Synovec Lab

Gas Chromatography, Liquid Chromatography,
and Mass Spectrometry, with Multi-Dimensional
Data Analysis



TIC Signal-Based Rough Simulated Distillation Curves





Simulated Distillation Sample Class Comparisons

