


Analysis of organochlorine pesticides (USEPA 608.1)

BPX5

Column part number	054142		
Phase	BPX5	Detector	MS
Column	30 m x 0.22 mm x 0.25 μ m	Injection mode	Splitless, 0.5 min
Initial temperature	50°C, 1 min	Carrier gas	He, 30 cm/sec
Rate 1	30°C/min to 150°C	Injection volume	1 μ L
Rate 2	6°C/min to 260°C	Interface/source/quad	320°C, 250°C, 150°C
Final temperature	260°C		



For more information visit www.trajanscimed.com or contact techsupport@trajanscimed.com