


## Analysis of raspberry flavours

BPX5

Column part number	054114		
Phase	BPX5	Final temperature	250°C, 10 min
Column	50 m x 0.22 mm x 0.25 µm	Detector	FID, 320°C
Initial temperature	60°C, 10 min	Carrier gas	He, 25 psi
Rate	2°C/min	Injection mode	Split, 240°C



For more information visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com)